# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**TERMOMETRO AURICOLARE**  
**IR EAR THERMOMETER**  
**THERMOMETRE AURICOLAIRE**  
**OHRENTHERMOMETER**  
**TERMÓMETRO AURICULAR**  
**TERMÔMETRO AURICULAR**  
**ΘΕΡΜΟΜΕΤΡΟ ΑΥΤΙΟΥ**

**Model: TH889**

**REF** 25580

**Gima S.p.A.**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)  
Made in China

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| Guidance and manufacturer's declaration – electromagnetic emissions   |            |  |
|---|------------|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. |            |  |
| Emissions test  | Compliance | Electromagnetic environment – guidance   |
| RF emissions<br>CISPR 11  | Group 1    | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  |
| RF emissions<br>CISPR 11  | Class B    | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes |

| Guidance and manufacturer's declaration – electromagnetic immunity   |   |   |   |
|--|---|---|---|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.  |   |   |   |
| Immunity test  | IEC 60601 test level                          | Compliance                                    | Electromagnetic environment – guidance  |
| Radiated RF<br>IEC 61000-4-3   | 10 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1KHz | 10 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1KHz | <p><b>Recommended separation distance</b></p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.7 \text{ GHz}$ <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p><b>NOTA 1</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>NOTA 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> |   |   |   |

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

| Immunity test  | IEC 60601 test level  | Compliance level  | Electromagnetic environment - guidance  |
|--|---|---|---|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2             | contact $\pm 8$ Kv<br>air $\pm 2, \pm 4, \pm 8,$<br>$\pm 15$ kv | contact $\pm 8$ Kv<br>air $\pm 2, \pm 4, \pm 8,$<br>$\pm 15$ Kv | Floors should be wood, concrete or ceramic tile.<br>If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Power frequency<br>(50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m<br>50 Hz or 60Hz   | 30 A/m<br>50 Hz or 60Hz   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.         |

**Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter m |                                       |  |
|---|---|---------------------------------------|--|
|   | 150KHz to 80MHz<br>$d = 1.2 \sqrt{P}$                       | 80MHz to 800MHz<br>$d = 1.2 \sqrt{P}$ | 800MHz to 2.7GHz<br>$d = 2.3 \sqrt{P}$ |
| 0.01  | N/A   | 0.12                                  | 0.23                                   |
| 0.1   | N/A   | 0.38                                  | 0.73                                   |
| 1   | N/A   | 1.2                                   | 2.3                                    |
| 10  | N/A   | 3.8                                   | 7.3                                    |
| 100   | N/A   | 12                                    | 23                                     |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. **NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Manufacturer's declaration-electromagnetic immunity

#### Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

The device is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment

| Test frequency (MHz) | Band <sup>a)</sup> (MHz) | Service <sup>a)</sup>   | Modulation <sup>b)</sup>                    | Maximum power (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) | Compliance LEVEL (V/m) (for home healthcare) |
|----------------------|--------------------------|---|---|-------------------|--------------|---------------------------|--|
| 385                  | 380 – 390                | TETRA 400   | Pulse modulation b) 18 Hz                   | 1,8               | 0,3          | 27                        | 27   |
| 450                  | 430 – 470                | GMRS 460, FRS 460   | FM c) $\square$ ±5 kHz deviation 1 kHz sine | 2                 | 0,3          | 28                        | 28   |
| 710                  | 704 – 787                | LTE Band 13, 17   | Pulse modulation b) 217 Hz                  | 0,2               | 0,3          | 9                         | 9  |
| 745                  |                          |   |   |                   |              |                           |  |
| 780                  |                          |   |   |                   |              |                           |  |
| 810                  | 800 – 960                | GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5          | Pulse modulation b) 18 Hz                   | 2                 | 0,3          | 28                        | 28   |
| 870                  |                          |   |   |                   |              |                           |  |
| 930                  |                          |   |   |                   |              |                           |  |
| 1 720                | 1700 – 1990              | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation b) 217 Hz                  | 2                 | 0,3          | 28                        | 28   |
| 1 845                |                          |   |   |                   |              |                           |  |
| 1 970                |                          |   |   |                   |              |                           |  |
| 2 450                | 2400 – 2570              | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7            | Pulse modulation b) 217 Hz                  | 2                 | 0,3          | 28                        | 28   |
| 5 240                | 5100 – 5800              | WLAN 802.11 a/n   | Pulse modulation b) 217 Hz                  | 0,2               | 0,3          | 9                         | 9  |
| 5 500                |                          |   |   |                   |              |                           |  |
| 5 785                |                          |   |   |                   |              |                           |  |

**NOTE** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.