**Gebrauchsanweisung
Infrarot-Multifunktions-
Thermometer**
Instructions
**Infrared Multifunction
Thermometer**
Mode d'emploi
**Thermomètre à infrarouge
multifonctions**
Instrucciones para el uso
**Termómetro multifuncional
infrarrojo**
Инструкция по эксплуатации
**Инфракрасный
многофункциональный термометр**
Istruzioni per l'uso
**Termometro multifunzionale
a infrarossi**

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Table of Contents **Page**

1. Important information – please read before start-up	22
2. Applications	23
3. The Advantages of this Ear Thermometer	23
4. Important Safety Instructions	25
5. Product Description	26
6. How this Ear Thermometer Measures Ear Temperature	27
7. Control Displays and Symbols	28
8. How to Reload a New Probe Cover	28
9. Directions for Use	30
10. Changing from Fahrenheit to Celsius and viceversa	32
11. How to recall 12 readings in Memory Mode	33
12. Cleaning and Storage	33
13. Error Messages	34
14. Replacing the Battery	35
15. Technical Specifications	35
16. Symbol	36
17. Calibration	36
18. EMC requirements	37

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1. Important information – please read before start-up

You have acquired a high-quality Riester multi-function infrared thermometer manufactured in compliance with Directive 93/42/ EEC for medical products and subject to continuous stringent quality control, ensuring reliable results for years to come.

Clinical thermometers

Part 5: Performance of infra-red ear thermometers
(with maximum device)

German version EN 12470-5

Please read these Instructions carefully prior to start-up and keep them in a safe place. Should you have any queries, we will be pleased to assist you at any time. You will find our address on the last page of these Instructions. We will be pleased to supply you with the address of your authorized RIESTER distributor on request.

Please note that correct and safe operation of this thermometer will be guaranteed only when RIESTER accessories are used throughout.

Please note the following prior to using the thermometer:

1. Ensure that the measuring probe (1) is shielded by a disposable probe cover (2) (even when the thermometer is not being used).
2. Keep the multi-function thermometer out of reach of children.
3. Keep the multi-function thermometer away from sources of intense heat, otherwise it could be damaged. (Do not hold the thermometer in a flame).
4. The multi-function thermometer must not be exposed to the sun for any length of time, as this will result in inaccurate temperature readings.