· HEINE OMEGA 600 HEINE OMEGA 600 wired

HEINE OMEGA 600
HEINE OMEGA 600 wired

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

SVENSKA

NEDERLANDS

DANSK

NORSK

SUOMI

PORTUGUÊS

HEINE OMEGA 600

HEINE OMEGA 600 wired

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Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	<p>Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.</p> <p>The supply voltage quality should be that of a typical hospital environment.</p> <p>Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.</p>
Performance features of the ME system that have been determined to be essential to the performance	None
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Note	Flickering of the LED of the OMEGA 600 wired is possible because of radiated electromagnetic fields.
	Brief illumination of the LED during charging of the OMEGA 600 is possible because of electrostatic discharge.
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	EMC compatibility is only ensured if original HEINE spare parts, accessories and spower sources are used as described in the chapter „Accessories“. The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.
Test	Compliance
RF emissions CISPR 11	Group 1 Class B
Conducted emissions (EN 55011/CISPR 11)*	Passed
Radiated emissions (EN 55011/CISPR 11)	
Harmonic current emissions (IEC 61000-3-2)*	
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*	
Immunity	See attached immunity test levels

Immunity test levels

Test	Test level	
	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ± 8 kV Air Discharge: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m 80–2700 MHz 80 % AM at 1kHz	
Electrical fast transients / bursts (IEC 61000-4-4)*	± 2 kV 100 kHz repetition frequency	
Surges (IEC 61000-4-5)*	± 0.5 kV, ± 1 kV	
Conducted disturbances induced by RF fields (IEC 61000-4-6)*	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM: ± 5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 50Hz or 60 Hz	30 A/m; 60 Hz
Voltage dips (IEC 61000-4-11)*	0 % U_T ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles Single phase: at 0°	
Short interruptions (IEC 61000-4-11)*	0 % U_T ; 250/300 cycles	

*n/a: "Not applicable" in the internally powered mode

Technical Specification OMEGA 600 and OMEGA 600 wired

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Classification according to EN ISO 15004-2:2007 and ANSI Z80.36-2016	Group II The classification was performed together with a Ø54mm/16 Diopter HEINE ophthalmoscopy lens.

Technical Specification OMEGA 600

CB1	Li-Po cell
Input	USB 2.0 Type C: 5 V, 1.2 A
Power consumption	6 W
Protection class	Charging: class II Operating: internally powered
Charging time	typ. 1.5 h
Operating time (at maximum charge capacity)	typ. 4 h
Operating time visionBOOST)	typ. 1.5 h
Weight	475 g incl. rechargeable battery

Technical Specification OMEGA 600 wired

Input	100–240 V~ / 50–60 Hz / 160–80 mA
Power consumption	6 W
Protection class	Class II
Weight	655 g

Accessories

CW1 – Wall Charger	X-095.17.320
CC1 – Charging Case	X-000.99.091
E4-USBC (USB-C cord approx. 2 m)	X-000.99.300

Options

TM2 – Teaching Mirror	C-000.33.212
OMEGA 600 Breath Shield	C-000.33.019
Fundus charts Pad with 50 pcs.	C-000.33.208

Spare parts


CB1 – OMEGA 600 Battery	X-007.99.687
OMEGA 600 eyepiece +2D	C-000.17.116
OMEGA 600 eyepiece 0D	C-000.17.115

Erläuterung der verwendeten Symbole
Explanation of utilized symbols
Explicación des symboles utilisés
Explicación de los símbolos utilizados
Spiegazione dei simboli utilizzati
Förklaring av symboler som används
Verklaring van de gebruikte symbolen
Forklaring af de anvendte symboler
Symbolforklaring
Käyttöttyjen symbolien selitys
Explicação dos símbolos utilizados

	<p>CE-Kennzeichnung kennzeichnet die Übereinstimmung mit der Europäischen Medizinprodukteverordnung (EU) 2017/745.</p> <p>The CE mark indicates that the product complies with the medical device regulation (EU) 2017/745.</p> <p>Le marquage CE indique que le produit est conforme au règlement (UE) 2017/745 relatif aux dispositifs médicaux.</p> <p>La marca CE indica que el producto cumple el Reglamento europeo sobre productos sanitarios (UE) 2017/745.</p> <p>Il marchio CE indica la conformità con il regolamento sui dispositivi medici (UE) 2017/745.</p> <p>CE-märkning markerar en överensstämmelse med förordningen om medicinska produkter (EU) 2017/745.</p> <p>CE-märkning duidt de overeenstemming aan met de verordening betreffende medische hulpmiddelen (EU) 2017/745.</p> <p>CE-mærkningen angiver overensstemmelse med forordningen om medicinsk udstyr (EU) 2017/745.</p> <p>CE-merket angir at produktet er i samsvar med forskriften om medisinsk utstyr (EU) 2017/745.</p> <p>CE-merkintä tarkoittaa, että laite lääkinnällisiä laitteita koskevan asetuksen (EU) 2017/745 kanssa.</p> <p>O simbolo CE identifica a concordância com o regulamento de Dispositivos Médicos (EU) 2017/74.</p>
	<p>Katalog- oder Bestellnummer Catalogue- or order number Numéro de catalogue ou de commande Número de catálogo o de pedido Codice catalogo e di dell'ordine numero Katalog- eller Beställningsnummer Catalogus- of bestelnummer Katalog- eller Ordrenummer Katalog- eller bestillingsnummer Luettelo- tai viitenumero Número de catálogo ou pedido</p>
	<p>Hersteller Manufacturer Fabricant Fabricante Produttore Tillverkare Fabrikant Producent Produsent Valmistaja Fabricante</p>
	<p>Hersteldatum Date of manufacture Date de fabrication Fecha de fabricación Data di produzione Tillverkningsdatum Productiedatum Produktionsdato Produksjonsdato Valmistuspäivä Data de fabricação</p>

	<p>Getrennte Sammlung von Elektro- und Elektronikgeräten. (Europäische WEEE Richtlinie) Product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal according to local provisions. (European Waste Electrical and Electronic Equipment Directive, WEEE) Tri sélectif des appareils électriques et électroniques. (Directive européenne DEEE) Desechado separado de aparatos eléctricos y electrónicos. (Directiva Europea RAEE) Raccolta differenziata di apparecchi elettrici ed elettronici (direttiva europea RAEE). Separat insamling av elektriska och elektroniska apparater (det europeiska WEEE-direktivet). Gescheiden inzameling van elektrische en elektronische apparaten (Europese AEEA-richtlijn). Separat indsamling af elektrisk og elektronisk udstyr (det europæiske WEEE-direktiv). Produkter med dette symbolet skal ikke kasseres sammen med vanlig husholdningsavfall, men krever separat kassering i henhold til lokale bestemmelser. (European Waste Electrical and Electronic Equipment Directive, WEEE) Sähkö- ja elektroniikkalaitteille tarkoitettu erillinen keräyspiste (eurooppalainen WEEE-standardi). Coleção separada de aparelhos elétricos e eletrônicos (Diretrizes Europeias WEEE).</p>
	<p>Batterien müssen einer zentralen Sammelstelle zugeführt werden. Dispose the rechargeable batteries at your local collection point. Mise au rebut au point de collecte local. Desechar en un punto limpio de la comunidad. Smaltimento in un centro di raccolta comunale. Avfallshantera på ett kommunalt insamlingsställe. Afvoer naar een gemeentelijke inzamellocatie. Bortskaffelse på et kommunalt indsamlingssted. Avhønd hos ditt lokale avfallshåndteringsanlegg. Laitte on hävitettävä paikalliseen keräyspisteeseen. Proceda à eliminação do aparelho em um ponto de coleta municipal.</p>
	<p>Zulässiger Temperaturbereich in °C für Lagerung und Transport Temperature limits in °C for storage and transport Plage de température admise en °C pour le stockage et le transport Rango de temperatura permitida en °C para almacenar y transportar el producto Temperatura ammessa in °C per conservazione e trasporto Tillåtet temperaturintervall i °C för lagring och transport Toegestane temperaturen in °C voor opslag en transport Tillått temperaturområde i °C ved oppbevaring og transport Temperaturbegrensning i °C for oppbevaring og transport Näyttää pakkauksen sallittu säilytys- ja kuljetuslämpötilan (°C). Limite de Temperatura permitida em °C para armazenamento e transporte</p>
	<p>Zulässiger Temperaturbereich in °F für Lagerung und Transport Temperature limits in °F for storage and transport Plage de température admise en °F pour le stockage et le transport Rango de temperatura permitida en °F para almacenar y transportar el producto Temperatura ammessa in °F per conservazione e trasporto Tillåtet temperaturintervall i °F för lagring och transport Toegestane temperaturen in °F voor opslag en transport Tillått temperaturområde i °F ved oppbevaring og transport Temperaturbegrensning i °F for oppbevaring og transport Näyttää pakkauksen sallittu säilytys- ja kuljetuslämpötilan (°F) Limite de Temperatura permitida em °F para armazenamento e transporte</p>
	<p>Zulässige Luftfeuchtigkeit für Lagerung und Transport Humidity limitation for storage and transport Humidité admise pour le stockage et le transport Humedad del aire permitida para almacenar y transportar el producto Umidità atmosferica ammessa durante il trasporto e la conservazione Tillåten luftfuktighet för transport och lagring Toegestane luchtvochtigheid voor opslag en transport Tillått luftfugtighed ved oppbevaring og transport Fuktighetsbegrensning for oppbevaring og transport Sallittu ilmastokosteus kuljetuksen ja varastoinnin aikana Umidade do ar admissível para o armazenamento e transporte</p>

	<p>Zulässiger Luftdruck für Lagerung und Transport Pressure limitation for storage and transport Pression atmosphérique admise pendant le transport et le stockage Presión de aire permitida para almacenar y transportar el producto Pressione atmosferica ammessa durante il trasporto e la conservazione Tillåten lufttryck för lagring och transport Toegestane luchtdruk voor opslag en transport Tilladt lufttryk ved opbevaring og transport Trykkbegrensning for oppbevaring og transport Sallittu ilmanpaine kuljetuksen ja varastoinnin aikana Pressão do ar admissível para o armazenamento e transporte</p>
	<p>Vorsicht Bruchgefahr! Fragile, handle with care! Fragile ! Manipuler avec soin Atención. Frágil. Attenzione: pericolo di rottura! Försiktigt! Risk för brott Voorzichtig, kans op breuk! Forsigtig, risiko for brud! Ømtålig, behandles forsigtigt! Varo särkymisvaaraa! Perigo de quebra!</p>
	<p>Trocken lagern! Keep dry! Conserver au sec ! Conservar en un lugar seco! Evitare ambienti umidi! Förvaras torr! Droog bewaren! Opbevares tørt! Hold tørt! Säilytetään kuivassa paikassa! Armazenar em ambiente seco!</p>
	<p>Gebrauchsanweisung verbindlich befolgen. (Hintergrundfarbe: blau, Vordergrundfarbe: weiß) Follow instructions for use! (Background color: blue, foreground color: white.) Suivre le mode d'emploi. (Couleur de fond : bleu ; couleur du premier plan : blanc) Seguir obligatoriamente las instrucciones de uso. (Color de fondo: azul, color de primer plano: blanco) Attenersi obbligatoriamente alle istruzioni per l'uso. (Colore dello sfondo: blu, colore in primo piano: bianco) Bruksanvisningen ska alltid följas. (Bakgrundsfärg: blå, förgrundsfärg: vit) De gebruiksaanwijzing is bindend en dient gevolgd te worden. (achtergrondkleur: blauw, voorgrondkleur: wit) Følg altid brugsanvisningen. (Baggrundsfarve: Blå; forgrundsfarve: Hvid) Følg brugsanvisningen! (Bakgrunnsfarge: blå, forgrunnsfarge: hvit) Käyttöohjeita on noudatettava tarkasti. (Taufaväri: sininen, etualan väri: valkoinen) Siga as instruções de uso! (Cor de fundo: azul, cor de primeiro plano: branco)</p>
	<p>Unique Device Identification</p>
	<p>Medical Device</p>
	<p>Gebrauchsanweisung Instructions for use Mode d'emploi Manual de instrucciones Istruzioni per l'uso Bruksanvisning Gebruiksaanwijzing Brugsanvisning Bruksanvisning Käyttöohjeet Instruções de utilização</p>

	<p>Nur in geschlossenen Räumen benutzen. For indoor use only. Utiliser uniquement dans des locaux fermés. Sólo utilizar en espacios cerrados. Utilizzare solo in ambienti chiusi. Får endast användas i slutna rum. Uitsluitend in afgesloten ruimten toepassen. Må udelukkende benyttes i lukkede rum. Skal kun brukes i lukkede rom. Käyttöön sallittu ainoastaan suljetuissa tiloissa. Utilizar apenas em espaços fechados.</p>
	<p>The Regulatory Compliance Mark (RCM) The Regulatory Compliance Mark (RCM) Marque réglementaire de conformité (RCM) The Regulatory Compliance Mark (RCM) The Regulatory Compliance Mark (RCM) The Regulatory Compliance Mark (RCM) The Regulatory Compliance Mark (RCM) The Regulatory Compliance Mark (RCM) Vaatimustenmukaisuusmerkki (RCM) (RCM)- Marca de Conformidade Regulamentar</p>
	<p>Geräte der Schutzklasse II Class II equipment Appareil de classe de protection II Aparato de clase de protección II Apparecchio di classe di protezione II Apparat i skyddsklass II Apparaat van beschermingsklasse II Apparater i sikkerhedsklasse II Klasse II-utstyr Suojausluokan II laite Aparelho da classe de proteção II</p>
	<p>Wechselstrom (AC) Alternating current (AC) Courant alternatif (CA) Corriente alterna (CA) Corrente alternata (AC) Växelström (AC) Wisselstroom (AC) Vekselstrøm (AC) Vekselstrøm (AC) Vaihtovirta (AC) Corrente alternada (CA)</p>
	<p>Gleichstrom Direct current DC Tension continue Tensión continua (CC) Tensione continua Likspänning Gelijkspanning Jævnspænding DC Likestrøm dc Tasajännite Tensão contínua (CC)</p>
	<p>Ladestandanzeige Charge status indicator Indicateur de l'état de charge Indicador del estado de carga Indicatore dello stato di carica Indikator för laddningsstatus Indicatielampje oplaadstatus Ladestatusindikator Ladestatusindikator Lataustilan merkivalo Indicador do estado de carga</p>

5V ~ 1.2A