PORTABLE AUDIOMETER

PICCOLO

MULTILANGUAGE USER MANUAL

Document title:	AU1P-Piccolo User Manual DA
Code:	AU1-MA303_D
Revision:	Rev. 02
Date:	2025.03.07

PICCOLO

BÆRBART AUDIOMETER

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SIGNALS SPECIFICATIONS	
Attenuators step	1, 3, 5 dB
Presentation mode	Continuous Pulsed, with rate of 0.5, 1 or 2 Hz
Frequency accuracy	0,1 %
Intensity accuracy	±3 dB between 125 Hz and 4 kHz ±5 dB above 4 kHz
Total Harmonic Distortion (THD)	AC: less than 2,5 % BC: less than 5,5 %
Warble tone	Frequency of the modulating signal: 5 Hz Modulation waveform: sine wave Modulation range: ±12%
NBN	Band: ½ octave, i.e.: - lower cut-off frequency $f_l = f / 1.1892$ - upper cut-off frequency $f_u = f \cdot 1.1892$ where f is the centre frequency
WN	Lower cut-off frequency: 100 Hz Upper cut-off frequency: 24 kHz
SN	As specified in IEC 60645-2 §13
External signals	EXT1 and EXT2 input: max 3 Vrms

AVAILABLE OUTPUTS			
Output	Basic	Plus	Speech
Air conduction (TDH-39, DD45 or DD65 headphones)	•	•	•
Air conduction (ER-3, IP30 or ER-5 insert earphones)	•	•	•
Bone conduction (B-71 bone vibrator)	-	•	•

PURE AND WARBLE TONE AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES					
Freq. (Hz)	AC TDH39 DD45 (dB HL)	AC DD65 (dB HL)	AC ER-3 IP30 (dB HL)	AC ER-5 (dB HL)	BC B71 (dB HL)
125	80	80	90	90	-
250	100	95	105	100	45
500	110	110	110	110	65
750	115	110	115	120	70
1.000	120	110	120	120	75
1.500	120	110	120	120	80
2.000	120	110	120	115	80
3.000	120	110	120	115	75
4.000	120	110	110	110	75
6.000	105	95	95	100	55
8.000	95	90	90	90	50

(*) Levels refer to device powered by mains power adapter. In case of device powered by USB, maximum levels decrease by 10 dB.

SPEECH AUDIOMETRY MAXIMUM INTENSITIES				
AC TDH39 DD45 (dB HL)	AC DD65 (dB HL)	AC ER-3 IP30 (dB HL)	AC ER-5 (dB HL)	BC B71 (dB HL)
100	90	100	100	55

(*) Levels refer to device powered by mains power adapter. In case of device powered by USB, maximum levels decrease by 10 dB.

EXTERNAL SIGNALS LEVEL INDICATOR (only for Speech model)	
Type of indicator	VU-meter
Dynamic range	+3..-20dB
Input Voltage at 0 dB	1.5 Vrms
Speed of the volume following	Increase: 60 dB/s Decrease: 60 dB/s

MASKING AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES				
Freq. (Hz)	AC TDH39 DD45 (dB EM)	AC DD65 (dB EM)	AC ER-3 IP30 (dB EM)	AC ER-5 (dB EM)
125	60	55	70	65
250	80	75	85	85
500	95	90	95	95
750	100	90	100	100
1.000	105	95	105	100
1.500	105	95	105	100
2.000	105	95	105	100
3.000	105	95	105	100
4.000	105	95	100	100
6.000	100	85	90	95
8.000	90	85	80	80
WN	90	80	80	80
SN	90	75	80	80

(*) Levels refer to device powered by mains power adapter. In case of device powered by USB, maximum levels decrease by 10 dB.

ACOUSTICAL SAFETY OF THE DEVICE	
Alert condition	Tone intensity higher than 100 dB HL (IEC 60645-1, §5.2)
Safety measures in alert condition	1) The operator has to press “Higher dB” button to increase the intensity over 100 dB HL 2) Warning on the display 3) “Normally on” function disabled

COMPATIBLE TRANSDUCERS		
<i>Type</i>	<i>Manufacturer</i>	<i>Model</i>
Supra-aural headphones	Telephonics Corp.	TDH39
Supra-aural headphones	Radioear Corp.	DD45
Circum-aural headphones	Radioear Corp.	DD65
Insert earphones	Etymotic Research Inc.	ER-3
Insert earphones	Etymotic Research Inc.	ER-5
Insert earphones	Radioear Corp.	IP30
Bone vibrator	Radioear Corp.	B71

PATIENT – OPERATOR COMMUNICATION			
	<i>Basic</i>	<i>Plus</i>	<i>Speech</i>
Talk-over through external microphone (not included)	•	•	•
Patient response trigger	•	•	•

AUDIOMETER CONTROL	
Required Software	Inventis Maestro
Communication protocol	USB 1.1
PC communication requirements	USB 2.0 port minimum

MECHANICS	
Size (LxDxH)	160 x 160 x 35 mm / 6.3 x 6.3 x 1.2 in
Weight (device only)	300 g / 10.6 oz

SOCKETS ON THE REAR PANEL		
<i>Description</i>	<i>Type</i>	<i>Connector</i>
Power supply	In	DC plug 2.5 mm
L and R headphones	Out	2 audio jack, 1/4” mono
Bone vibrator (not in Piccolo Basic)	Out	Audio jack, 1/4” mono
Patient response trigger	In	Audio jack, 1/4” mono
External microphone for talk-over	In	Audio jack, 3.5 mm mono
USB	In - Out	USB type B
<i>Only on Piccolo Speech</i>		
LINE IN	In	Audio jack, 3.5 mm stereo

SPECS OF INPUT FACILITIES	
<i>Input</i>	<i>Electrical property</i>
Power supply	Internal pin +6V, external pin 0V
Patient response	Switches 3V to logical input (switch current: 10mA)
LINE IN.	Sensitivity: 3mV at max volume and 0Vu Impedance: 10K Ω Freq. response: 75-12000Hz +/- 3dB
Microphone	Electret or 200 Ω dynamic microphone Impedance: 47K Ω Freq. response: 100-12KHz +/- 3dB Electret Bias: 2.2V trough 2.2K Ω

SPECS OF OUTPUT FACILITIES		
<i>Output</i>	<i>Available Voltage</i>	<i>Nominal Impedance</i>
L and R phones	8V _{pp}	10 Ω
Bone vibrator	8V _{pp}	10 Ω

SOUND ATTENUATION VALUES			
Frequency	TDH 39 (*) DD45 (*)	DD65	ER-3 ER-5 IP30
[Hz]	[dB]	[dB]	[dB]
125	3.0	8.3	33.5
250	5.0	15.5	34.5
500	7.0	26.1	34.5
750	-	-	-
1000	15.0	32.4	35
1500	-	-	-
2000	26.0	43.6	33
3000	-	-	-
4000	32.0	43.8	39.5
6000	-	-	-
8000	24.0	45.6	43.5

(*) With MX41\AR or PN 51 cushions

REFERENCE EQUIVALENT THRESHOLD LEVELS						
	TDH 39	DD45	DD65	ER-3 IP30	ER-5	B71*
<i>Ref. std.</i>	ISO 389-1 (ANSI S3.6)	ISO 389-1 (ANSI S3.6)	Vendor Tech. Specificat.	ISO 389-2 (ANSI S3.6)	ISO 389-2	ISO 389-3 (ANSI S3.6)
<i>Freq.</i> <i>[Hz]</i>	<i>dB</i> <i>[re 20μPa]</i>	<i>dB</i> <i>[re 20μPa]</i>	<i>dB</i> <i>[re 20μPa]</i>	<i>dB</i> <i>[re 20μPa]</i>	<i>dB</i> <i>[re 20μPa]</i>	<i>dB</i> <i>[re 1μN]</i>
125	45	47.5	30.5	26	26	-
250	25.5	27	17.0	14	14	67
500	11.5	13	8.0	5.5	5.5	58
750	7.5 (8)	6.5	5.5	2	2	48.5
1000	7	6	4.5	0	0	42.5
1500	6.5	8	2.5	2	2	36.5
2000	9	8	2.5	3	3	31
3000	10	8	2.0	3.5	3.5	30
4000	9.5	9	9.5	5.5	5.5	35.5
6000	15.5	20.5	21.0	2	2	40
8000	13	12	21.0	0	0	40

(*) Calibration of bone vibrator (B71) refers to the placement on the mastoid.

Annex B

Electromagnetic compatibility

Piccolo has been thoroughly tested and respects the limits for electro-medical devices specified by IEC 60601-1-2 standards. These limits ensure reasonable protection against hazardous interference in typical medical installations.

The instrument generates, uses and radiates radio frequency energy. If not installed and used according to the instructions in this manual, it may interfere with other nearby devices. No guarantee is given that interference will not occur under certain conditions.

This instrument is suitable for use in professional healthcare facility environment, i.e. in hospital environments, except for near active HF surgical equipment and RF shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

Piccolo should not be used adjacent to or stacked with other equipment. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

The existence of electromagnetic interference can be verified easily by switching the instrument off and back on again. If it is proven that the device is indeed interfering with other devices, try to solve the problem by adopting one of the following solutions:

- change the orientation and/or position of the affected device;
- move the two devices further away from each other;
- contact the manufacturer or authorised service organisation for further assistance.

List of cables, transducers and accessories

Cables, transducers and accessories with which Inventis claims the compliance with the IEC 60601-1-2 standard are those ones supplied with the device, in particular the followings:

- 1) 6Vdc Medical Grade mains power adapter
- 2) Power supply cable (maximum length: 1.8 m)
- 3) TDH39, DD45 and DD65 transducers with 2 m double signal shielded cable
- 4) Insert earphone: ER-3 (Etymotic) or IP30 (RadioEar)
- 5) Bone vibrator B-71 transducer with 2m shielded signal cable

- 6) Patient response switch (manufactured by Inventis) with 2 m shielded cable
- 7) Talk over microphone with 2m shielded cable
- 8) Stereo cable 3.5mm plug to 3.5mm plug, 1.8m, shielded
- 9) USB cable, shielded, maximum length: 2 m

The use of accessories, transducers and cables other than those specified, except for transducers and cables sold by the manufacturer as spare parts for internal components, may result in increased emissions or decreased immunity of the device.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of Piccolo, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

The instrument does not have any ESSENTIAL PERFORMANCE as per the IEC 60601 1 standard.

Note: All necessary instruction for maintaining compliance with regard to electromagnetic compatibility can be found in the maintenance section in this manual. No further steps are required.

Guidance and manufacturer's declaration – electromagnetic emissions		
Piccolo is intended for use in the electromagnetic environment specified below. The customer or the user of Piccolo should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR11	Group 1	Piccolo uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	Piccolo is suitable for use in professional healthcare facility environment and directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flickers emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
Piccolo is intended for use in the electromagnetic environment specified below. The customer or the user of Piccolo should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ⁽¹⁾ ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air ⁽¹⁾	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a professional healthcare facility environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a professional healthcare facility environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. $40\% U_T$ (60% dip in U_T) for 5 cycles. $70\% U_T$ (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. $40\% U_T$ (60% dip in U_T) for 5 cycles. $70\% U_T$ (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	Mains power quality should be that of a professional healthcare facility environment. If the user of Piccolo requires continued operation during power mains interruptions, it is recommended that Piccolo be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields at power frequencies must correspond to the levels typical of professional healthcare facilities.
Note: ⁽¹⁾ A lock or reboot of the device with no permanent damage is acceptable ⁽²⁾ U_T is the a.c. main voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
Piccolo is intended for use in the electromagnetic environment specified below. The customer or the user of Piccolo should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	Portable and mobile RF communications equipment should be used no closer than 30cm (12 inches) to any part of Piccolo, including cables specified by the manufacturer. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey⁽¹⁾, should be less than the compliance level in each frequency range.⁽²⁾ Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 Mhz - 2,7 Ghz	3 V/m 80 Mhz to 2,7 Ghz	
Radiated RF (from RF Wireless communication equipment) IEC 61000-4-3	9 V/m 704-787 MHz 5100 – 5800 MHz	9 V/m 704-787 MHz 5100 – 5800 MHz	
	27 V/m 380 - 390 Mhz	27 V/m 380 - 390 Mhz	
	28 V/m 430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	28 V/m 430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	
Proximity Fields IEC 61000-4-39	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz 7.5 A/m 13.56MHz Pulse Mod 50KHz	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz 7.5 A/m 13.56MHz Pulse Mod 50KHz	
<i>Note:</i> At 80 MHz and 800 MHz, the higher frequency range applies. <i>Note:</i> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. <i>Note:</i> ⁽¹⁾ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If measured field strength in the location Triangle is used in exceeds the applicable RF compliance level above, Triangle should be observed to verify its normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Triangle. ⁽²⁾ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Guidance and manufacturer's declaration – electromagnetic immunity	
Function to verify for freedom from unacceptable risk	Acceptance pass/fail criteria
Sound generation operating correctly	No sound from transducers exceeding 80dBHL; a lock or reboot of the device is acceptable