X36-1 LETTINO DA VISITA ELETTRICO
X36-1 ELECTRIC EXAMINATION TABLE
X36-1 DIVAN D'EXAMEN ÉLECTRIQUE
X36-1 CAMILLA DE EXPLORACIÓN ELÉCTRICA
X36-1 ELEKTRISCHE UNTERSUCHUNGSLIEGE
X36-1 MESA DE EXAME ELÉTRICA
X36-1 ΗΛΕΚΤΡΙΚΟ ΕΞΕΤΑΣΤΙΚΟ ΚΡΕΒΑΤΙ
X36-1 ЕЛЕКТРИЧЕСКА МАСА ЗА ПРЕГЛЕД
X36-1 ELEKTRICKÝ VYŠETŘOVACÍ STŮL
X36-1 ELEKTRISK UNDERSØGELSESELEJE
X36-1 ELEKTRILINE LÄBIVAATUSLAUD

X36-1 SÄHKÖINEN TUTKIMUSPÖYTÄ
X36-1 ELEKTROMOS VIZSGÁLÓASZTAL
X36-1 ELEKTRINĚ APŽIŮROS KUŠETĚ
X36-1 ELEKTRISKAIS IZMEKLĒŠANAS GALDS
X36-1 ELEKTRISCHE ONDERZOEKSTAFEL
X36-1 ELEKTRYCZNY STÓŁ DO BADAŃ
X36-1 CANAPEA ELECTRICĂ DE EXAMINARE
X36-1 ELEKTRICKÝ VYŠETROVACÍ STŮL
X36-1 ELEKTRIČNA PREGLEDNA MIZA
X36-1 ELEKTRISKT UNDERSÖKNINGSBORD

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Jiangsu Saikang Medical Equipment Co., Ltd
 No. 35 Lehong Road, Modern Agriculture Demonstration Park,
 Zhangjiagang City, Jiangsu Province, China
 Made in China

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

The whole examination table can be controlled by hand operator

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Guidance and Manufacturer's Declaration – Electromagnetic Immunity		
The V8v8c is intended for use in the following electromagnetic environment. The customer or user of the V8v8c should make sure it is used in such an environment		
Emission test	Compliance	Electromagnetic Environment Guidance
RF Emission to CISPR 11	Group 1	The bed uses RF energy for its internal function
RF Emission to CISPR 11	Class B	The bed is intended for use in all types of establishment including residential and similar user that are directly connected to a public supply network that also serves buildings used for residential purposes
Harmonics according to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker acc. to IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity		
The V8v8c is intended for use within healthcare facility environment. Users or customers must use bed in this environment		
Immunity test	IEC 60601-1-2 EN 60601-1-2 Test level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)..... IEC 61000-4-2 /EN 61000-4-2	Contact discharge: +8kV Air discharge: +2kV, +4kV, +8kV, +15kV	Contact discharge: +8kV Air discharge: +2kV, +4kV, +8kV, +15kV Floor should be made of wood and concrete or be tiled with ceramic tiles. If the floor is covered with synthetic flooring material, the relative air humidity must be at least 30% Can be used when higher ESD levels are present
Electrical fast transient/burst IEC 61000-4-4 /EN 61000-4-4	+2kV Power supply line +1kV Input/Output lines	+2kV Power supply line The quality of the supply voltage should be equivalent to that of a typical business or hospital environment.
Surge IEC 61000-4-5 /EN 61000-4-5	Line - line: +0.5kV, +1.0kV Line - ground: +0.5kV, +1.0kV, +2.0kV	Line - line: +0.5kV, +1.0kV Line - ground: +0.5kV, +1.0kV, +2.0kV The quality of the supply voltage should be equivalent to that of a typical business or hospital environment.
Power Frequency Magnetic Fields IEC 61000-4-8 /EN 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz Magnetic fields with a network frequency should be equivalent to those to be found in a typical business or hospital environment.
Voltage Dips IEC 61000-4- 11 /EN 61000-4-11	0% UT: 0.5 cycles Phase: 00, 450, 900, 1350, 1800, 2250, 2700, 3150	0% UT: 0.5 cycles Phase: 00, 450, 900, 1350, 1800, 2250,
		0% UT: 1 cycle 70% UT: 25/30
	0% UT: 1 cycle 70% UT: 25/30	cycles Single phase: 00
	cycles Single phase: 00 0% UT ; 300 cycle	0% UT ; 300 cycle The quality of the supply voltage should be equivalent to that of a typical business or hospital environment. If the person using the bed requires that the bed functions must continue despite any interruptions in the energy supply, it is recommended that the bed be connected to an uninterruptible electricity supply or a battery

Short interruptions IEC 61000-4-11 /EN 61000-4-11	0% UT: 250/300 cycles	0% UT: 250/300 cycles
		The quality of the supply voltage should be equivalent to that of a typical business or hospital environment. If the person using the bed requires that the bed functions must continue despite any interruptions in the energy supply, it is recommended that the bed be connected to an uninterruptible electricity supply or a battery
Note: UT is the AC Power voltage before applying in the test level		

Guidance and Manufacturer's Declaration - Electromagnetic Immunity							
The V8v8c is intended for use in an electromagnetic environment that used wireless RF transceiver. Electromagnetic interference can be prevented by maintaining the minimum distance between RF transceiver and V8v8c depend on the maximum output and frequency of the communication device as recommended below.							
Immunity test	IEC 60601-1-2 /EN 60601-1-2 Test level	Compliance level (V/m)	Band d (MHz)	Service	Modulation	Maximum power (W)	Distance (m)
Enclosure port to RF wireless communication equipment IEC 61000-4-3 /EN 61000-4-3	27V/m 385MHz	27	380-390	TETRA 400			
	28V/m 450MHz	28	430-470	GMRS 460, FRS 460			
	9V/m 710MHz	9	704-787	LTE Band 13, 17			
	9V/m 745MHz						
	9V/m 780MHz						
	28V/m 810MHz	28	800-960	GSM 800/900, TETRA 800, iDen 820, CDMA 850, LTE Band 5			
	28V/m 870MHz						
	28V/m 930MHz						
	28V/m 1720MHz	28	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4,25; UMT S			
	28V/m 1845MHz						
	28V/m 1970MHz						
	28V/m 2450MHz	28	2400-2570	Bluetooth, WLAN, 802.11 b/g/n RFID 2450, LTE band 7			
	9V/m 5240MHz 9V/m 5500MHz 9V/m 5785MHz	9	5100-5800	WLAN 802.11 a/n			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The V8v8c is intended for use in the following electromagnetic environment. The customer or user of the V8v8c should make sure it is used in such an environment.			
Immunity test	IEC 60601-1-2 /EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment - Guidance
RF Conductivity IEC 61000-4-6 EN 61000-4-6	3V 0.15MHz-80MHz	3V 0.15MHz-80MHz	Field intensity from the fixed RF transmitter determined by a field investigation of the magnetic field must be lower than the compliance level of the respective frequency band. Near devices with the following symbol, interference could occur. 
	6V in ISM bands between 0.15MHz- 80MHz 80%AM at 1kHz	6V in ISM bands between 0.15MHz- 80MHz 80%AM at 1kHz	
Radioactive RF electromagnetic field IEC 61000-4-3 /EN 61000-4-3	3V/M 80MHZ-2.7GHz 80%AM at 1kHz	3V/M 80MHZ-2.7GHz 80%AM at 1kHz	

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