
HANDHELD TYMPANOMETER SCREENING AUDIOMETER

TIMPANI

MULTILANGUAGE USER MANUAL

TYMPANOMETER

TIMPANI

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TYMPANOMETRY 226HZ	
Probe tone	With AGC
Frequency and intensity	226 Hz $\pm$ 1%; 85 $\pm$ 1.5 dB SPL
Measurement range and accuracy	0.2 to 8.0 ml $\pm$ 0.1 ml or $\pm$ 5% (whichever is greater)

Representation	Meatus-compensated
Influence of ambient temperature	-0.003 ml/°C
Influence of atmospheric pressure	-0.0002 ml/daPa
Scan range	Standard +200 to -400 daPa Low +100 to -300 daPa ± 10 daPa or ± 10 % (whichever is greater)
Scan rate	400 daPa/s
Pressure control	automatic
Pressure safety limits	Upper limit 550 daPa Lower limit -750 daPa

TYMPANOMETRY 1000HZ – Tympanometry 1000Hz license required	
Frequency and intensity	1000 Hz ± 1%; 75 ± 1.5 dB SPL
Measurement range and accuracy	0.9 to 16 mmho ± 0.5 mmho or ± 5% (whichever is greater)
Representation	Meatus-compensated

ACOUSTIC REFLEXES – Reflex and Reflex 1000Hz license required	
Type of stimulation	Ipsilateral, pulsed (50ms ON, 70ms OFF)
Stimulus frequencies and accuracy	1kHz ± 1% (with 1000Hz-only reflex license) 500Hz, 1kHz, 2kHz, 4kHz ± 1% (with full reflex license)
Harmonic distortion	THD 2.5% maximum
Intensity and accuracy	70 to 100dB HL ± 3 dB HL
Duration of stimulus	1s
Type of test	- Fixed intensity, adjustable intensity, 5dB steps - Threshold search, attenuator steps 5dB or 10dB, initial value and final value adjustable through 5dB steps
Recognition threshold	Adjustable 0.04 or 0.06 ml ± 0.01ml There is a negligible risk of artifacts in measurements at high stimulus levels, and they do not influence the reflex identification system
Test pressure	Automatic - Tympanogram peak pressure - Peak pressure less gradient (width of pressure curve at half peak height) - Atmospheric pressure

AUDIOMETRY – audiometry license required	
Stimulus	Pure tone
Attenuator step	5 dB
Method of presentation	Continuous, Pulsed (cadence: 1Hz)

Frequency accuracy	0.1%	
Intensity accuracy	± 3 dB between 125Hz and 4kHz; ± 5 dB over 4kHz	
Total harmonic distortion (THD)	VA: less than 2.5%	
Compatible transducers		
Type	Manufacturer	Model
Supra-aural headphones	Radioear Corp.	DD45
Insert earphones	Etymotic Research Inc.	ER-3C

PURE TONE AND WARBLE TONE AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES		
Freq [Hz]	AC DD45 [dB HL]	AC ER-3C [dB HL]
125	60	70
250	80	90
500	90	95
750	100	100
1.000	100	100
1.500	100	100
2.000	100	100
3.000	100	100
4.000	100	100
6.000	90	80
8.000	85	70

EQUIVALENT THRESHOLD LEVELS OF REFERENCE FOR PURE TONE			
	DD45	ER-3C	IPSI
Ref. standard	ISO 389-1 (ANSI S3.6)	ISO 389-2 (ANSI S3.6)	ISO 389-2 (ANSI S3.6)
Coupler	IEC 60318-3	IEC 60318-5	IEC 60318-5
Freq. [Hz]	dB [re 20 μ Pa]	dB [re 20 μ Pa]	dB [re 20 μ Pa]
125	47.5	26.0	-
250	27.0	14.0	-
500	13.0	5.5	11
750	6.5	2.0	-
1000	6.0	0.0	5.5
1500	8.0	2.0	-

2000	8.0	3.0	7
3000	8.0	3.5	-
4000	9.0	5.5	2
6000	20.5	2.0	-
8000	12.0	0.0	-

SOUND ATTENUATION VALUES		
Freq. [Hz]	DD45(*) [dB]	ER-3C [dB]
125	3.0	33.5
250	5.0	34.5
500	7.0	34.5
750	-	-
1000	15.0	35.0
1500	-	-
2000	26.0	33.0
3000	-	-
4000	32.0	39.5
6000	-	-
8000	24.0	43.5

(*) With MX41\AR or PN 51 cushions

PATIENT MANAGEMENT	
Maximum number of patients	50
Details memorized	Patient details (first name, last name, date of birth, gender), test date and time, Tympanogram (RH + LH), reflex tracing (RH + LH), audiometric thresholds (RH + LH)

PHYSICAL SPECIFICATIONS	
Display	TFT LCD 2.8" RGB, 240x320 pixels Viewing area 43.2 mm x 57.6 mm
Touchscreen	Capacitive
Dimensions of Tympanometer instrument	(W x D x H) 65 x 44 x 240 mm / 2.6 x 1.8 x 9.5 in
Weight of instrument	340 g / 12 oz
Dimensions of docking station	(W x D x H) 109 x 83 x 131 mm / 4.3 x 3.3 x 5.2 in
Weight of base only	280 g / 9.9 oz

TYMPANOMETER CONNECTORS	
AC headphone	Output, 3.5mm 4-pole audio jack, 8Vpp max with 10Ω load

Patient response button	Input, 2.5 mm mono audio jack, 4Vpp max
USB	I/O, mini B type, 5.5Vdc max
Docking station contacts	I/O, target type for spring-loaded contact, +/-10Vpp

DOCKING STATION CONNECTORS	
USB	I/O, 2x mini B type, 5.5Vdc max
Docking station contacts	I/O, spring-loaded contact, +/-10Vpp

INTERFACE WITH COMPUTER	
Connection:	USB (no driver needed)
Minimum pc requirements	USB 2.0 port
Compatible software products	Inventis Maestro suite

BLUETOOTH PRINTER INTERFACE	
Module type	Bluetooth v4.2 – dual mode
Frequency	2402 - 2480 GHz
Maximum power in transmission	Class 1 +8 dBm from antenna
Sensitivity	94 dBm
Coverage	100m maximum
Conformity	EC: Essential requirements as in article 3 of EU directive 2014/53/EU; Radio Equipment Directive (RED); FCC ID: SQGBT850; Industry Canada IC: 3147A-BT850

Upon request, Inventis will make available circuit diagrams, parts lists, descriptions, calibration instructions or other information that may assist service personnel in repairing parts of the device that Inventis has deemed repairable by staff.

Electromagnetic compatibility

The instrument has been tested and respects the limits for electro-medical devices specified by IEC 60601-1-2 standards. These limits ensure reasonable protection against hazardous interference in typical medical installations.

The device generates, uses and radiates radio frequency energy. If not installed and used according to the instructions in this manual, it may interfere with other nearby devices. No guarantee is given that interference will not occur under certain conditions.

This device is suitable for use in professional healthcare facility environments, i.e. in hospitals, except for near active HF surgical equipment and RF-shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

The instrument should not be used adjacent to or stacked on other equipment. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

The existence of electromagnetic interference can be verified easily by turning the device off and back on again. If it is found that the device is indeed interfering with other equipment, try to solve the problem by adopting one of the following solutions:

- Change the orientation and/or position of the affected device.
- Move the two devices further away from each other.
- Contact the manufacturer or authorised service center for further assistance.

List of cables, transducers and accessories

Cables, transducers and accessories for which Inventis declares compliance with IEC 60601-1-2 standard are those supplied with the device, in particular:

- Medical grade USB power adapter
- USB cable: 2m length, screened.
- Docking station
- Transducers: 2m length, screened.
- Patient button: 2m length, screened.

The use of accessories, transducers and cables other than those specified, except for transducers and cables sold by the manufacturer as spare

parts for internal components, may result in increased emissions or decreased electromagnetic immunity of the device and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Anyone connecting other items of equipment must ensure that the overall system complies with IEC 60601-1-2.

The instrument does not have any ESSENTIAL PERFORMANCE as per the IEC 60601 1 standard.

Note: All instructions necessary for maintaining compliance with regard to electromagnetic compatibility can be found in the maintenance section in this manual. No further steps are required.

Manufacturer's indications and declaration - electromagnetic emissions		
Timpani is intended for use in an electromagnetic environment as specified below. The customer or the user of Timpani should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment Guidance
RF emissions CISPR11	Group 1	Timpani uses RF energy for its internal functions, and contains a Bluetooth radio module that responds to pertinent regulations. Consequently, the RF emissions generated are minimal and unlikely to interfere with other equipment operating nearby.
RF emissions CISPR11	Class B	Timpani is suitable for use in professional healthcare facilities and for connection direct to the low voltage public electricity grid.
Harmonic emissions IEC 61000-3-2	Class A	
Fluctuations in voltage/emissions (flicker) IEC 61000-3-3	Complies	

Manufacturer's indications and declaration - electromagnetic immunity			
Timpani is intended for use in the electromagnetic environment as specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity tests	Test level to IEC 60601	Level of Compliance	Electromagnetic environment – guidelines
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a professional healthcare environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a professional healthcare environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% $U_T^{(1)}$ (dip >95% in U_T) for half cycle. 40% U_T (dip >60% in U_T) for 5 cycles. 70% U_T (dip >30% in U_T) for 25 cycles. <5% U_T (dip >95% in U_T) for 5 s.	< 5% $U_T^{(1)}$ (dip >95% in U_T) for half cycle. 40% U_T (dip >60% in U_T) for 5 cycles. 70% U_T (dip >30% in U_T) for 25 cycles. <5% U_T (dip >95% in U_T) for 5 s.	Mains power quality should be that of a professional healthcare environment. If the user of Timpani requires continued operation during power mains interruptions, it is recommended that the instrument be powered from an uninterruptable power supply or its battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a professional healthcare environment.
<i>Note:</i> ⁽¹⁾ U_T is the A.C. mains voltage prior to application of the test level.			

Manufacturer's indications and declaration - electromagnetic immunity			
Timpani is intended for use in an electromagnetic environment as specified below. The customer or the user of Timpani should assure that it is used in such an environment.			
Immunity tests	Test level to IEC 60601	Level of Compliance	Electromagnetic environment – guidelines
Conducted RF IEC 61000-4-6	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	Portable and mobile RF communications equipment should be used no closer than 30 cm (12 inches) to any part of Timpani, including cables specified by the manufacturer. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ⁽¹⁾ , should be less than the compliance level in each frequency range. ⁽²⁾ Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz	3 V/m 80 MHz - 2.7 GHz	
Radiated RF (from RF Wireless communication equipment) IEC 61000-4-3	9 V/m 704-787 MHz 5100 – 5800 MHz	9 V/m 704-787 MHz 5100 – 5800 MHz	
	27 V/m 380 - 390 Mhz	27 V/m 380 - 390 Mhz	
	28 V/m 430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	28 V/m 430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	
Proximity Fields IEC 61000-4-39	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz 7.5 A/m 13.56MHz Pulse Mod 50KHz	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz 7.5 A/m 13.56MHz Pulse Mod 50KHz	
<i>Note:</i> At 80 MHz and at 800 MHz, the higher frequency range is applied. <i>Note:</i> These guidelines may not be valid in all situations. Electromagnetic propagation is affected by absorption and reflection on encountering structures, objects and persons.			
<i>Note:</i> ⁽¹⁾ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If measured field strength in the location Timpani is used in exceeds the applicable RF compliance level above, Timpani should be observed to verify its normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Timpani. ⁽²⁾ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Manufacturer's indications and declaration - electromagnetic immunity	
Function requiring verification to exclude unacceptable risks	Pass/fail acceptance criteria
Sound generator operating correctly	No unwanted sound from transducers exceeding 80 dB; lockup or restart of the device is acceptable
Tympanometry on test cavity, conducted correctly in normal operating conditions	Flat tympanometry curve, ECV indication equivalent to nominal value of cavity +/- 0.1ml ESD: presence of artifacts in tympanometry recognizable by skilled professional, HW error, lockup or restart of device — all these are acceptable

Contains transmitter module in compliance with ETSI EN 301 489-1 and ETSI EN 300 328 standards

*The device emits radio frequencies in the 2.4 GHz band, class 1
Contains transmitter module in compliance with ETSI EN 301 489-1 and ETSI EN 300 328 standards*