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REF 20775 / BP2020V3 -KIDS

MYCROCLEAN ITALIA SRL
Via Carlo Porta, 14
20064 Gorgonzola - Italy
Made in Italy

LOT Vedere etichetta / See label / Voir l'étiquette
Siehe Etikett / Ver la etiqueta / Ver rótulo
Se etikett / Δείτε ετικέτα / Patrz etykieta
Consultați eticheta / انظر الملصق

Vedere etichetta / See label / Voir l'étiquette
Siehe Etikett / Ver la etiqueta / Ver rótulo
Se etikett / Δείτε ετικέτα / Patrz etykieta
Consultați eticheta / انظر الملصق

UNI EN 14683:2019

000001

The result obtained on the BP2020 series masks is 7.5, therefore good or very breathable

4 - Instructions for use and warnings

Handle with clean and sanitised hands, making sure the mask adheres well to the face, in the right position, i.e. with the smallest part towards the nose, and the larger one over the mouth and chin.

DO NOT MODIFY THE MASK; DO NOT CUT THE ELASTIC BANDS; WEAR IT RIGHT WAY OUT (as shown in the figure) Otherwise the mask loses its protective function.

5 - Washing

By hand or in the washing machine, with a detergent, up to 60°C. **DO NOT BLEACH**

By complying with the instructions for use and maintenance, the product maintains the same physical and technical features up to approx. 500 washes.

6 - Disposal

Mask to be disposed of as SUW (unless it has been contaminated, in which case it must be disposed of as "special waste"). The outer bag is disposed of as plastic.

REF	Product code	LOT	Lot number
	Manufacturer		Date of manufacture
	Non-sterile		Consult instructions for use
	Medical Device complies with Directive 93/42/EEC		

ENGLISH

1 - Product description and intended use

Washable mask, in 99% polyester and 1% carbon fibre fabric, designed to protect the nose and mouth

2 - Classification

The BP2020V3 series masks are classified as Type II, according to the UNI EN 14683:2019 standard.

3 - Features

The device consists of one layer and is machine washable up to 60°C (the fabric withstands up to 150°C) It is not sterile
PFE values: Particle Filtration efficiency (%) 0.3 µm 0.5 µm at 95%

The fabric is **NON-CYTOTOXIC** according to tests carried out on fibroblasts in compliance with the UNI EN ISO 10993-5 standard.

BFE Filtration

BFE = 99.8% (from laboratory tests carried out in compliance with the UNI EN 14683:2019 standard)

Biological load

UFC <30 UFC (from laboratory tests carried out in compliance with the UNI EN 14683:2019 standard) **Air permeability (cm²/cm³/sec) ASTM E96** (Also known as ISO - 11092 or Ret or Hohenstein test)