CATETERE URETRALE STERILE MONOUSO
STERILE URETHRAL CATHETER FOR SINGLE USE
CATHÉTERA URÉTRAL STÉRILE À USAGE UNIQUE
CATÉTER URETRAL ESTÉRIL DE UN SOLO USO
CATETER URETRAL ESTÉRIL DE UTILIZAÇÃO ÚNICA
STERILER HARNRÖHRENKATHETER FÜR DEN EINMALGEBRAUCH
ΑΠΟΣΤΕΙΡΩΜΕΝΟΣ ΟΥΡΗΘΡΙΚΟΣ ΚΑΘΗΤΗΡΑΣ ΜΙΑΣ ΧΡΗΣΗΣ
STERYLNY CEWNIK UROLOGICZNY DO JEDNORAZOWEGO UŻYTKU
STERILNÍ URETRÁLNÍ KATÉTR NA JEDNO POUŽITÍ
STERIL URINRÖRSKATETER FÖR ENGÅNGSBRUK
STERIILI JA KERTAKÄYTTÖINEN VIRTSAPUTKIKATETRI
STERILNI URETRALNI KATETER ZA ENKRATNO UPORABO
STERILNÝ MOČOVÝ KATÉTER NA JEDNO POUŽITIE
CATETER URETRAL STERIL DE UNICĂ FOLOSINȚĂ
STERIELE URETHRALE KATHETER VOOR EENMALIG GEBRUIK
STERILNI KATETER ZA URETRU ZA JEDNOKRATNU UPORABU
STERIL URINRÖRSKATETER FÖR ENGÅNGSBRUK
STERILT URINRØRSKATETER TIL ENGANGSBRUG
СТЕРИЛЕН УРЕТРАЛЕН КАТЕТЪР ЗА ЕДНОКРАТНА УПОТРЕБА
VIENKARTINIO NAUDOJIMO STERILUS ŠLAPLĖS KATETERIS
STERILS URĪNIZVADKANĀLA KATETRS VIENREIZĒJAI LIETOŠANAI
STERIILNE UREETRA KATEETER ÜHEKORDSEKS KASUTAMISEKS

REF F5110.10 (Gima 22003) - F5110.12 (Gima 22004) - F5110.14 (Gima 22005) -
F5110.16 (Gima 22006) - F5110.18 (Gima 22007) - F5100.08 (Gima 22022)
- F5100.10 (Gima 22023) - F5100.12 (Gima 22024) - F5100.14 (Gima 22025)
- F5100.16 (Gima 22026)

Jiangsu Suyun Medical Materials Co., Ltd. No.18
Jin Qiao Road Dapu Industrial Park 222002 Lianyungang
JIANGSU PROVINCE, PEOPLE'S REPUBLIC OF CHINA
Made in China

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße, 80, 20537, Hamburg, Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

 0123

ENGLISH**Product name:**

Sterile Urethral Catheter for Single Use (Nelaton catheter)

Type:

Paediatric male, Paediatric female, Male, Female

Sizes:

Straight Tip/ Standard tip

Intended purpose:

For temporary emiction by inserting it to bladder through urethra when patients can not micturate by themselves.

Intended users:

It can be used by trained doctors, nursing staff, patients or accompanying personnel.

Structural Composition:

The product is composed of drainage funnel and tube

Indications:

Urinary obstruction and urine re-tention, etc.

Contraindications:

Those who are allergic to PVC materials should not use this product.

Patient group:

Patients with urinary retention caused by urinary tract obstruction. Rescue critically ill patients.

Patients in need of diagnosis and treatment of bladder disease. People undergoing surgery.

Patients with spinal cord injury.

Performance characteristics:

Single-use and sterile

Performance Characteristics:

- a. The tube, the head end and the hole are free of impurities.
- b. The surface is smooth, the top is round and without pitting point, the inner cavity is clean and smooth.

Method:

- a. Wash your hands and put on gloves to clean around the urethral opening. Those who have foreskin should pull the foreskin down to expose the urethral opening and clean it.
- b. After preparing for catheterization, open the individual package from the opening.
- c. Lubricate the catheter with medical liquid paraffin or medical silicone oil.
- d. Insert the urinary catheter into the bladder slowly through the urethra for normal catheterization.

Any residual risk and any undesirable side-effects:

none

Installation information:

none

Clinical benefits:

The product help patients to temporary emission.

Caution:

a. When using the product, it should strictly implement the requirements of aseptic operation specifications and relevant regulations, and can be used by trained doctors, nursing staff, patients or accompanying personnel.

However, when patients or accompanying staff conduct urinary catheterization for the first time, please use it under the guidance of professional medical staff. In the course of use, the operator or patient should be treated in a timely manner according to medical regulations.

b. If you find any possible quality doubt during use, please stop using it immediately and notify our company quickly.

c. The product should be sterilized with EO, and the product should be sterile.

d. The product is disposable and should be destroyed after use.

e. Please check the individual package before using it.

f. Don't use if package is opened or damaged.

g. Sterile and validity for 5 year. Don't use if overdue.

h. This product needs to be used in conjunction with a container for urine. When using it, it should avoid splashing urine and polluting the environment.

i. Discard after use. After use, the products shall be collected, transported, stored and disposed to ensure that the products are harmless to human health and the environment. Or the product can be disposed in accordance with the local medical regulations for waste disposal.

j. If the product is used in a household environment, after use, the product shall be disposed as household hazardous waste.

The product shall be collected, transported, stored and disposed by a dedicated establishment or an undertaking.

Or the product can be disposed in accordance with the local medical regulations for waste disposal.

k. All serious incidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

Storage:

The product should be stored in a well-ventilated room temperature.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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