

Instructions for use

ETHILON™
POLYAMIDE 6 OR POLYAMIDE 6,6
STERILE SYNTHETIC NON-ABSORBABLE SURGICAL SUTURE,
USP / Ph. Eur.

DESCRIPTION

ETHILON™ Suture is a sterile, monofilament, synthetic, non-absorbable, surgical suture composed of polyamide 6 [NH-(CH₂)₆-NH-CO-(CH₂)₆-CO]_n, or polyamide 6,6 [NH-(CH₂)₄-NH-CO-(CH₂)₆-CO]_n. ETHILON™ Suture is available undyed and dyed black with Hematine HCK (Color Index 75290), and dyed green with D&C Green No.5 (Color Index 61570) to enhance visibility in the surgical field. ETHILON™ Suture is available in a range of gauge sizes and lengths, non-needled or attached to needles of various types and sizes, and in presentations as described in the HOW SUPPLIED section. ETHILON™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Polyamide 6 or Polyamide 6,6 Suture and United States Pharmacopoeia (USP) for Non-Absorbable Surgical Sutures. The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. sizes as equivalent which is reflected on the labeling.

INDICATIONS

ETHILON™ Suture is indicated for use in general soft tissue approximation and/or ligation, including cardiovascular, ophthalmic and neurosurgical procedures.

APPLICATION

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

PERFORMANCE / ACTIONS

ETHILON™ Suture elicits an initial, minimal inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While polyamide is not absorbed, progressive hydrolysis of the polyamide *in vivo* may result in gradual loss of tensile strength over time.

en

CONTRAINDICATIONS

Adverse reactions associated with the use of this device include wound dehiscence, gradual loss of tensile strength over time, calculi formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, minimal inflammatory tissue reaction, and transient local irritation at the wound site.

WARNINGS

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing ETHILON™ Suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users. Like all foreign bodies, this product may potentiate infection.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilament sutures.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse reactions associated with the use of this device include wound dehiscence, gradual loss of tensile strength over time, calculi formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, minimal inflammatory tissue reaction, and transient local irritation at the wound site.

STERILITY

ETHILON™ Suture is sterilized by irradiation. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures.

STORAGE

No special storage conditions required. Do not use after expiry date.

HOW SUPPLIED

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability. ETHILON™ Suture is available as sterile monofilament strands in USP sizes 11-0 through 2 (metric sizes 0.1 – 5.0) in a variety of lengths, with and without permanently attached needles. The sutures are also available in presentations containing the following:
1. Lead Seal and Surgical Bolster in which the lead seal is used to maintain the position of the bolster relevant to the suture knot to maintain proper tension.
2. Retention tubing (elastomer tubing), which is used to spread the load of the suture at the surface of the skin.

ETHILON™ Suture is available in one, two or three dozen units per box.

ar

- cs Označení zákonného výrobce naleznete na etiketě výrobku.
- de Für den behördlich zugelassenen Hersteller siehe Produktetikett.
- el Για πληροφορίες σχετικά με τον αναγνωρισμένο νόμιμο κατασκευαστή, ανατρέξτε στην ετικέτα του προϊόντος.
- en For recognized legal manufacturer, refer to product label.
- es Para conocer el fabricante legal, remítase a la etiqueta del producto.
- fi Tunnistettu laillinen valmistaja käy ilmi tuotetarrasta.
- fr Fabricant légal reconnu, voir l'étiquette du produit.
- hu A jogilag elismert gyártót lásd a termék címkéjén.
- it Per il fabbricante legalmente riconosciuto, consultare l'etichetta del prodotto.
- kk Белгілі заңды өндірушісін анықтау үшін өнім жапсырмасын қараңыз.
- ko 인증된 법적 제조사에 대한 정보는 제품 라벨을 참조하십시오.
- nl Raadpleeg het productetiket voor de erkende wettelijke fabrikant.
- no For anerkjent juridisk produsent, se produktets etikett.
- pl Informacje dotyczące uprawnionego legalnego producenta podano na etykietce produktu.
- pt Para conhecer o fabricante legal reconhecido, consulte a etiqueta do produto.
- ro Pentru producătorul legal recunoscut, a se consulta eticheta produsului.
- ru См. данные об официальном производителе на ярлыке изделия.
- sk Zákonne uznaný výrobca je uvedený na etikete výrobku.
- sv För erkänd, laglig tillverkare, se produktetiketten.
- tr Yetkili yasal imalatçı için ürün etiketine bakınız.
- zh-cn 关于认定的合法制造商的信息, 请参阅产品标签。
- zh-tw 有關認定的合法製造廠的資訊, 請參閱產品標籤。

ETHICON, LLC

475 C Street
Los Frailes Industrial Park
Suite 401
Guaynabo, Puerto Rico 00969
USA
1-877-ETHICON
+1-513-337-6928

ETHICON, LLC

Highway 183 Km 8.3
San Lorenzo, Puerto Rico 00754
USA
1-877-ETHICON
+1-513-337-6928

Johnson & Johnson International

c/o European Logistics Centre
Leonardo Da Vinciiaan, 15
BE-1831 Diegem
Belgium
+1-513-337-6928

Johnson & Johnson Medical GmbH

Robert-Koch-Strasse 1
Norderstedt
22851
Germany

© Ethicon, Inc. 2015

ETHILON™

- ar خط جراحی
- cs ŠÍČÍ MATERIÁL
- da SUTUR
- de NAHTMATERIAL
- el ΠΑΜΜΑ
- en SUTURE
- es SUTURA
- fi OMMELAINE
- fr FIL DE SUTURE
- hu VARRÓANYAG
- it SUTURA
- kk ТИҢУ МАТЕРИАЛЫ
- ko 봉합사
- nl HECHTMATERIAAL
- no SUTUR
- pl NICI CHIRURGICZNE
- pt FIO DE SUTURA
- ro FIR DE SUTURĂ
- ru ШОВЫЙ МАТЕРИАЛ
- sk CHIRURGICKÁ NIŤ
- sv SUTUR
- tr SÜTÜR
- zh-cn 缝线
- zh-tw 縫合線



390323R02
LAB100601593v2
08/2019



	<p>ar - لا تعيد تعبئة وحدة هذا فقط da-Må ikke genbruges / de-Nicht zur Wiederverwendung / el-Μην αναποχρησιμοποιείτε / en-Do not re-use / es-No reutilizar / fi-Ei saa käyttää uudelleen / fr-Ne pas réutiliser / hu-Ne használja újra / it-Non riutilizzare / kk-Қайта пайдаланбаңыз / ko-재사용하지 마십시오 / nl-Niet opnieuw gebruiken / no-Skal ikke brukes flere ganger / pl-Nie używać powtórnie / pt-Não reutilizar / ro-Nu reutilizați / ru-Не использовать повторно / sk-Nepoužívať opakovane / sv-Får ej återanvändas / tr-Tekrar kullanılmayın / zh-cn-不得重复使用 / zh-tw-不得重複使用</p>
	<p>ar - الرجاء الرجوع إلى التعليمات للاستخدام خلال da-Anvendes inden / de-Verwendbar bis (Datum) / el-Χρησιμοποιήστε πριν / en-Use by date / es-Usar antes de fecha / fi-Viimeinen käyttöpäivä / fr-A utiliser avant / hu-Lejárató idő / it-Utilizzare entro il / kk-Жарамдылык мерзімі / ko-사용기한일 / nl-Uiterste gebruiksdatum / no-Utløpsdato / pl-Wykorzystać do dnia / pt-Data de validade / ro-A se utiliza până la data de / ru-Использовать до (даты) / sk-Datum spotreby / sv-Används före-datum / tr-Son kullanma tarihi / zh-cn-使用有效期 / zh-tw-使用有效日期</p>
LOT	<p>ar - رمز التمييز da-Batch kode / de-Chargenbezeichnung / el-Κωδικός παρτίδας / en-Batch code / es-Código de lote / fr-Étiquetage / fi-Koodi de lot / hu-Tételkód / it-Codice di lotto / kk-Бума коды / ko-회차 코드 / nl-Batchcode / no-Batchkode / pl-Kod partii / pt-Código do lote / ro-Codul lotului / ru-Код партии / sk-Kód šarže / sv-Satskod / tr-Parti kodu / zh-cn-批号 / zh-tw-批號</p>
STERILE	<p>ar - معقم بالاشعاع da-Steriliseret ved stråling / de-Sterilisiert durch Bestrahlung / el-Αποστειρωμένο με χημική ακτινοβολία / en-Sterilized using irradiation / es-Esterilizado mediante irradiación / fi-Steriloitu säteittäminällä / fr-Stérilisé par irradiation / hu-Besugárzással sterilizálva / it-Sterilizzato con radiazioni ionizzanti / kk-Сәулелену арқылы зарарсыздандырылған / ko-방사선조사를 사용하여 멸균되었음 /</p>

STERILE	<p>nl-Gesteriliseerd door bestraling / no-Sterilisert med stråling / ph-Wysterylizowano promieniowaniem / pt-Esterilizado por irradiação / ro-Sterilizat prin iradiere / ru-Стерилизовано облучением / sk-Sterilizované žiarením / sv-Steriliserad med strålning / tr-İşleme ile sterilize edilmiştir / zh-cn-使用辐照灭菌 / zh-tw-經放射線滅菌</p>
REF	<p>ar - رقم الكاتالوج cs-Katalogové číslo / da-Varenummer / de-Katalognummer / el-Αριθμός καταλόγου / en-Catalog Number / es-Número de catálogo / fi-Tuotenumero / fr-Numéro de référence au catalogue / hu-Katalógusszám / it-Numero di catalogo / kk-Каталог бойынша нөмөрі / ko-카탈로그 번호 / nl-Catalogusnummer / no-Katalognummer / pl-Numer katalogowy / pt-Número de catálogo / ro-Număr de catalog / ru-Номер по каталогу / sk-Katalógové číslo / sv-Katalognummer / tr-Katalog numarasi / zh-cn-目录编号 / zh-tw-目錄編號</p>
	<p>ar - تحذير cs-Upozornění / da-Bemærk / de-Achtung / el-Προσοχή / en-Caution / es-Atención / fi-Huomio / fr-Attention / hu-Vigyázat / it-Attenzione / kk-Назар аударыңыз / ko-주의 / nl-Let op / no-Forsiktig / pl-Przeostroga / pt-Atenção / ro-Atentie / ru-Внимание / sk-Upozornenie / sv-ÖBS / tr-Dikkat / zh-cn-注意 / zh-tw-注意</p>
	<p>ar - صانع da-Producent / de-Hersteller / el-Κατασκευαστής / en-Manufacturer / es-Fabricante / fi-Valmistaja / fr-Fabricant / hu-Gyártó / it-Fabbricante / kk-Дайындаушы / ko-제조사 / nl-Fabrikant / no-Produsent / pl-Producent / pt-Fabricante / ro-Producător / ru-Производитель / sk-Výrobca / sv-Tillverkare / tr-İmalatçı / zh-cn-制造商 / zh-tw-製造廠</p>

	<p>ar - لا تعيد تعبئة هذا المنتج cs-Neprováděte opětovnou sterilizaci / da-Må ikke resteriliseres / de-Nicht erneut sterilisieren / el-Μην επαναποστειρωσετε / en-Do not resterilize / es-No reesterilizar / fi-Ei saa steriloida uudelleen / fr-Ne pas restériliser / hu-Ne sterilizálja újra / it-Non ristilizzare / kk-Қайта зарарсыздандырмаңыз / ko-재멸균하지 마십시오 / nl-Niet opnieuw steriliseren / no-Skal ikke resteriliseres / pl-Nie wyjalawiać powtórnie / pt-Não reesterilizar / ro-Nu reesterilizați / ru-Не стерилизовать повторно / sk-Nesterilizujte opakovane / sv-Får ej omsteriliseras / tr-Tekrar sterilize etmeyin / zh-cn-不得重复灭菌 / zh-tw-不得重複滅菌</p>
EC REP	<p>ar - ممثل مُعتمد لدى الاتحاد الأوروبي cs-Zplnomocněný zástupce v Evropském společenství / da-Autoriseret repræsentant i EU / de-Bevollmächtigter in der Europäischen Gemeinschaft / el-Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / en-Authorized Representative in the European Community / es-Representante autorizado en la Comunidad Europea / fi-Valtuutettu edustaja Euroopan yhteisössä / fr-Mandataire agréé dans la Communauté européenne / hu-Engedélyelt rendelkező képviselő az Európai Közösségen / it-Rappresentante autorizzato nella Comunità Europea / kk-ЕР Resmi өкілі / ko-유럽연합 공식 대리인 / nl-Gemachtigd vertegenwoordiger in de Europese Unie / no-Autorisert representant i EU / pl-Autoryzowany przedstawiciel we Wspólnocie Europejskiej / pt-Representante autorizado na Comunidade Europeia / ro-Reprezentant autorizat în Uniunea Europeană / ru-Официальный представитель в ЕС / sk-Autorizovaný zástupca v Európskom spoločenstve / sv-Auktoriserad representant i Europeiska gemenskaperna / tr-Avrupa Topluluğu Yetkili Temsilcisi / zh-cn-欧洲共同体授权代理 / zh-tw-歐洲共同體授權代表</p>


	<p>ar - لا تستخدم هذا المنتج إذا وجدت الصوة تالفة da-Må ikke anvendes, hvis pakningen er beskadiget / de-Bei beschädigter Verpackung nicht verwenden / el-Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά / en-Do not use if package is damaged / es-No usar si el envase está dañado / fi-Ei saa käyttää, jos pakkaus on vaurioitunut / fr-Ne pas utiliser si l'emballage est endommagé / hu-Ne használja, esli csomagolás megsérült / it-Non utilizzare se l'imballaggio non è integro / kk-Егер қаптама зақымданған болса пайдаланбаңыз / ko-포장이 파손되어 있는 경우에는 사용하지 마십시오 / nl-Niet gebruiken als de verpakking beschadigd is / no-Skal ikke brukes hvis pakningen er skadet / pl-Nie używać, jeśli opakowanie jest uszkodzone / pt-Não utilizar se a embalagem estiver danificada / ro-Nu se utilizează dacă ambalajul este deteriorat / ru-Не использовать, esli упаковка повреждена / sk-Nepoužívať, esli obal poškodený / sv-Får inte användas om förpackningen är skadad / tr-Ambalaj hasarlıysa kullanılmayın / zh-cn-如果包装有损坏, 不得使用 / zh-tw-如果包装有破損, 不得使用</p>
	<p>ar - عدد الوحدات da-Antal enheder / de-Stückzahl / el-Αριθμός μονάδων / en-Number of Units / es-Número de unidades / fi-Yksiköiden lukumäärä / fr-Nombre d'unités / hu-Az egységek száma / it-Numero di unità / kk-Бірліктер саны / ko-제품 개수 / nl-Aantal eenheden / no-Antall enheter / pl-Liczba jednostek / pt-Número de unidades / ro-Număr de unități / ru-Число единиц / sk-Počet kusov / sv-Antal enheter / tr-Ünite sayısı / zh-cn-件数 / zh-tw-數量</p>

MS/X	<p>ar - متعدد الخيوط da-Flertåret / de-Mehrfachstrang / el-Πολυκλώβητο νήμα / en-Multi-strand / es-Multitbra / fi-Moniliikainen / fr-Multibrins / hu-Többszáls / it-Multifilo / kk-Көп түйін / ko-다중 가닥 / nl-Meerdere draden / no-Flertåret / pl-Wielopasmowy / pt-Multifios / ro-Multifilament / ru-Мультиниточный / sk-Viacvláknová / sv-Multi-strand / tr-Çok iplikli / zh-cn-多股 / zh-tw-多股</p>
	<p>ar - زون الفتحة / da-TRÆK AF/LØFT/ÅBEN HER / de-HIER ABZIEHEN/ANHEBEN/ÖFFNEN / el-ΑΠΟΚΛΙΝΕΤΕ/ΑΝΑΨΗΚΟΤΕ/ΑΝΟΙΞΤΕ ΕΔΩ / en-PEEL LIFT/OPEN HERE / es-DESPEINDE/LEVANTAR/ABRIR AQUÍ / fr-LÉVÉ/NOUVEAU TÂTE/FR-DETAACHER/SOULEVER/OUVRIR ICI / hu-ITT HÚZZA LE/EMELJE FEL/NYISSA FEL / it-STRAPPARE/SOLLEVARE/APRIRE QUI / kk-ОСЫ ЖЕРДЕЛІЗ СЫЛДИРЫҢЫЗ/КӨТЕРІҢІЗ/АШЫҢЫЗ / ko-여기를 떼어 내십시오 / nl-Hier openen / no-Her åpne / pl-Odklej/odwinić/otwócić / pt-Puxar/levantar/abrir aqui / ro-Zona de deschidere / sk-TU ODLUPNÚT/ODVIHNÚT/OTVORIŤ / sv-DRA AV/LYFT/OPPNÄ HÄR / tr-BURADAN SOYUN / zh-cn-剥开/提起/此处开 / zh-tw-剝開/提起/此處開啟</p>
	<p>ar - انظر هنا da-RIV HER / de-HIER AUFREISSEN / el-ΞΕΚΛΙΝΕ ΕΔΩ / en-TEAR HERE / es-RASGAR AQUÍ / fr-REVI TÂTE / fr-DÉCHIRER ICI / hu-ITT SZAKTÁSA LE / it-ROMPERE QUI / kk-ОСЫ ЖЕРДЕЛІЖЫПТЫҢЫЗ / ko-여기를 찢으십시오 / nl-Hier afscheuren / no-Riv her / pl-ROZERWAC W TYM MIEJSCU / pt-RASGAR AQUÍ / ro-RUPETI AICI / ru-ОТΠΑТЬ ЗДЕСЬ / sk-TU ODRHNÚT / sv-RIV HÄR / tr-BURADAN YIRTIN / zh-cn-由此撕开 / zh-tw-在此撕開</p>

THIS PANEL IS LEFT BLANK INTENTIONALLY.


IFU PRINTING SPECIFICATION SHEET

PAGE LAYOUT



Flat Size

FOLD PATTERN



Folded Size

TITLE ETHILON™		DESCRIPTION Map IFU		LAB NUMBER LAB100601593v2	SPECIAL INSTRUCTIONS/COMMENTS n/a			BINDING n/a	COLORS Black			
FLAT SIZE 27.56" x 17.28" 700 mm x 439 mm		FOLDED SIZE 2.3" x 4.32" 58.42 mm x 110 mm		RMC NUMBER 390323R02	PAGE COUNT 2		LANGUAGES ar, cs, da, de, el, en, es, fi, fr, hu, it, kk, ko, nl, no, pl, pt, ro, ru, sk, sv, tr, zh-cn, zh-tw		SELF COVER <input checked="" type="checkbox"/>	PLUS COVER <input type="checkbox"/>	SEALING METHOD	WAFER SEAL <input checked="" type="checkbox"/>
BLEED SIZE .5" (12.7 mm) <input type="checkbox"/> .125" (3.175 mm) <input type="checkbox"/>		NONE <input checked="" type="checkbox"/> BLEED ALL SIDES <input type="checkbox"/>		BLEED TOP <input type="checkbox"/>	BLEED RIGHT <input type="checkbox"/>	BLEED LEFT <input type="checkbox"/>	BLEED BOTTOM <input type="checkbox"/>					
STOCK Refer to MS159-007					<p>DRAWING IS NOT TO SCALE: DRAWINGS REFLECT INFORMATION FOR PRODUCTION OF PRINTED PIECES AND DO NOT CONTAIN ACTUAL ARTWORK. This document or data herein or herewith is not to be reproduced, used or disclosed in whole or part without the permission of Ethicon, Inc.</p> <p style="text-align: center;">ETHICON</p>							