## Instructions for use

IFU

### MONOCRYL™ Poliglicaprone 25 (MONOFILAMENT) STERILE SYNTHETIC ABSORBABLE SURGICAL SUTURE

**DESCRIPTION**  
MONOCRYL™ Suture is a sterile, monofilament, synthetic, absorbable surgical suture composed of a copolymer of glycolide and ε-pirolidone ε-caprolactone. The empirical molecular formula of the polymer is (C<sub>10</sub>H<sub>16</sub>O<sub>2</sub>)<sub>x</sub> (C<sub>8</sub>H<sub>13</sub>O<sub>2</sub>)<sub>y</sub>.

Poliglicaprone 25 copolymer has been found to be nonaromatic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

MONOCRYL™ Suture is available in a range of gauge sizes and lengths, non-needled or attached to needles of varying types and sizes as described in the HOW SUPPLIED section.

MONOCRYL™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Synthetic Absorbable Monofilament Sutures and United States Pharmacopoeia (USP) for Absorbable Surgical Sutures except for a slight oversize in diameter for all gauges as shown in the table below.

USP Suture Size	Metric Suture Size	Maximum Oversize
5-0	0.7	0.03 mm
4-0	1.5	0.05 mm
3-0	2	0.07 mm
2-0	3	0.05 mm
0	3.5	0.08 mm
0	4	0.06 mm

The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. Sizes as equivalent, which is reflected on the labeling.

**INDICATIONS**  
MONOCRYL™ Sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

**APPLICATIONS**  
Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound site.

**PERFORMANCE / ACTIONS**  
MONOCRYL™ Suture elicits a minimal initial inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption occur by means of hydrolysis, where the polymer degrades to lactic acid, which is subsequently absorbed and metabolized in the body. Absorption begins as loss of tensile strength followed by a loss of mass. Implantation studies in rats show the following profile:

DYED		UNDYED	
Post Implantation	Approx. Original Strength Remaining	Post Implantation	Approx. Original Strength Remaining
7 days	60% to 70%	7 days	50% to 60%
14 days	30% to 40%	14 days	20% to 30%

All of the original tensile strength is essentially lost by 28 days post-implantation.

Absorption is essentially complete between 91 and 119 days.

**CONTRAINDICATIONS**  
This suture (dyed and undyed), being absorbable, should not be used where extended approximation of tissues under tension is required.

**WARNINGS**  
The safety and effectiveness of MONOCRYL™ Sutures has not been established in the following areas – neural tissue, cardiovascular tissue, microsurgery and ophthalmic surgery. Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing MONOCRYL™ Sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Medical professionals should consider the in vivo performance UNDER PERFORMANCE/ACTIONS section when selecting a suture.

**PRECAUTIONS**  
In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crushing damage due to application of surgical instruments such as forceps or needle holders. Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be removed after 7 days if possible. Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

This suture may be inappropiate to elderly, malnourished or debilitated patients, or in patients suffering from conditions that may delay wound healing.

As with any suture material, adequate knot security requires the standard surgical technique of flat square throws with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilament sutures.

Care should be taken to avoid damage when handling surgical needles. Grasping the needle in an one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause failure of the needle. Grasping at the attachment end could cause bending or breakage. Resharpening needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Broken needles may result in extended or additional sutures or residual suture bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens. Discard used needles in "Sharp" containers.

**ADVERSE REACTIONS**  
Adverse reactions associated with the use of this device include wound dehiscence, callus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occur, transient inflammatory tissue reaction, and transient local induration at the wound site.

**STERILITY**  
MONOCRYL™ Suture is sterilized by ethylene oxide. Do not resterilize. Do not use if package is opened or damaged. Discard opened suture.

**STORAGE**  
No special storage conditions required. Do not use after expiry date.

**HOW SUPPLIED**  
Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

MONOCRYL™ Suture is available as sterile, monofilament strands in sizes 6-0 through 1 (metric sizes 0.7 - 4), in a variety of lengths, with and without permanently attached needles.

MONOCRYL™ Suture is also available in sizes 3-0 through 1 (metric sizes 2 - 4) as sterile sticks attached to CONTROL RELEASE™ removable needles, which enable the needles to be pulled off instead of being cut.

MONOCRYL™ Suture is available in one, two or three drain units per box.

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© Ethicon, Inc. 2017