
# Introcan Safety<sup>®</sup> 3

---


---

**B | BRAUN**

LLDorder 6029 - grunewald 221544

**en** Fig. 1: Before needle removal  
**de** Abb. 1: Vor dem Entfernen der Kanüle  
**bg** Фиг. 1: Преди изваждане на иглата  
**cs** Obr. 1: Před vytažením jehly  
**da** Fig. 1: Inden fjernelse af kanylen  
**el** Σχ. 1: Πριν από την αφαίρεση της βελόνας  
**es** Fig. 1: antes de la extracción de la aguja  
**et** Joon 1: Enne nõela eemaldamist  
**fi** Kuva 1: Ennen neulan poistamista  
**fr** Fig. 1 : Avant le retrait de l'aiguille  
**hr** Slika 1: prije vađenja igle  
**hu** 1. ábra: A tű eltávolítása előtt  
**it** Fig. 1: Prima della rimozione dell'ago  
**lt** 1 pav.: prieš nuimant adatą  
**lv** 1. att. Pirms adatas noņemšanas  
**nl** Afb. 1: Voor verwijderen van de naald  
**no** Fig. 1: Før fjerning av nål  
**pl** Ryc. 1: Przed wyjęciem igły  
**pt** Fig. 1: Antes da remoção da agulha  
**pt BR** Fig. 1: Antes da remoção da agulha  
**ro** Fig. 1: înainte de scoaterea acului  
**ru** Рис. 1. До извлечения иглы  
**sk** Obr. 1: Pred vyťahnutím ihly  
**sl** Slika 1: Pred odstranitvijo igle

**sr** Sl. 1: Pre uklanjanja igle  
**sv** Fig. 1: Innan nålen tas bort  
**tr** Şek. 1: İğneyi çıkarmadan önce  
**vi** Hình 1: Trước khi tháo kim  
**zh** 图 1: 撤出针芯之前

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





## en Instructions for use

### Description

The Introcan Safety® 3 Closed IV Catheter is an indwelling peripheral catheter to obtain intravascular and subcutaneous tissue access for short-term use. It is designed to:

- Minimize inadvertent needlestick injuries with a passive safety shield
- Improve catheter stability with an integrated stabilization platform designed to minimize movement within the vessel
- Control the flow of blood from the catheter hub after needle removal and during subsequent (dis-) connections of Luer access devices, helping to reduce blood exposure.

The Introcan Safety® 3 Closed IV Catheter can be used with Luer access devices which are in compliance with ISO 80369-7.

This device is intended for single use only. It is provided sterile and non-pyrogenic.

### Materials Used

- PUR, PP, ABS, Chrome-Nickel steel, POM, Elastomer
- FEP, PP, ABS, Chrome-Nickel steel, POM, Elastomer

Components are not made with natural rubber latex, PVC, or DEHP.

### Indication

- The Introcan Safety® 3 Closed IV Catheter is inserted into a patient's vascular system for short term use to
  - sample venous or arterial blood,
  - monitor blood pressure, or
  - administer infusion solutions, blood products or pharmaceuticals (e.g. in or without carrier solutions) according to the SPC of the drug/solution.

- The Introcan Safety® 3 Closed IV Catheter is also indicated for subcutaneous infusion therapies according to the SPC of the drug/ solution.
- The 14-24 gauge catheters may be used with power injectors set at a maximum pressure of 325 psi and a recommended maximum flow rate. Use Luer lock connection only. See table below for recommended maximum flow rates when using a power injector. Testing performed at room temperature (22°C / 72°F). Warming contrast media per manufacturers' recommendations can reduce the pressure necessary to achieve recommended flow rates.

	Contrast Media [mPa*s]	Flow Rate (ml/sec)
24 Gauge*	2.3	5.0
	27.5	2.5
22 Gauge*	2.3	8.0
	27.5	3.5
20 Gauge*	2.3	10.5
	27.5	4.0
18 Gauge*	2.3	14.0
	27.5	5.0
16 Gauge*	2.3	15.5
	27.5	5.5
14 Gauge*	2.3	16.5
	27.5	5.5

\* applicable to all lengths

The above are flow rates tested in a laboratory setting and are the maximum flow rates our IV catheters are able to achieve and do not constitute guarantees, warranties, or predictions regarding the outcome of your matter. In any case, it is the responsibility of the user to accommodate the flow rates to each patient's condition and/ or therapy needed.

### Patient Population & Intended User

The catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness for the procedure being performed, the solution being infused and duration of therapy.

The catheter is intended to be used by qualified Healthcare Professionals who are trained on the use of peripheral IV catheters according to national regulations.

### Contraindications

- The Introcan Safety® 3 Closed IV Catheter should not be used in patients with known hypersensitivity to any of the materials employed.
- The Introcan Safety® 3 Closed IV Catheter is not intended to be used to facilitate the placement of vascular access devices such as guidewires, indwelling central venous catheters (CVC), peripherally inserted central catheters (PICC) and midline catheters into the vascular system.

### Precautions & Warnings

- Observe standard precautions on all patients. Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.
- When placing or maintaining any IV Catheter, standard precautions must be adhered in accordance with Centers for Disease Control and Prevention / Occupational Safety and Health Administration (CDC / OSHA) standards for blood-borne pathogens to avoid the risk of exposure to contaminated blood.
- Use only if packaging is intact. This device is sterile unless packaging has been opened or damaged.
- Do not re-use. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

- In the case of an unsuccessful IV Catheter placement attempt, remove the needle first to activate safety mechanism, then remove catheter from patient and discard both.
- Never reinsert the needle inside the catheter once the needle has been partially or completely withdrawn as it may pierce and/or sever the catheter.
- Do not attempt to bypass the safety mechanism.
- In the unlikely event that the safety mechanism has not been triggered/activated, keep the tip of the needle away from the body and fingers at all times and immediately dispose the IV Catheter into a suitable sharp-proof container. Care must be taken to avoid needlestick injuries.
- Report needlestick injuries immediately and follow established institutional protocols
- Use of only ISO Luer slip and Luer lock connectors ensures proper functioning.
- Luer Slip connections should not be left unattended due to potential of disconnection.
- Extreme care should be taken not to damage, pierce, cut or sever the catheter. Therefore, do not bend the catheter and/or needle during insertion, advancement, or removal of the needle.
- Do not use scissors or sharp instruments at or near the insertion site.
- If blood is present, flush device according to institutional protocol.
- Always clearly mark arterial lines to avoid inadvertent injection.
- Verify adequate collateral circulation prior to arterial puncture.
- Patency of catheter must be assured immediately before power injecting.
- Measures should be taken to avoid kinking or obstructing the catheter system during power injection to avoid product failure.

- Prior to use with power injectors, ensure that secure connection exists between the catheter and power injector. Use Luer Lock connection only.
- In the event of occlusion or lack of patency, power injector pressure-limiting function may not prevent catheter failure and/or infiltration.
- Monitor tightness of all connections during the whole infusion procedure.

#### Residual Risks & Side Effects

- General risks related to IV Catheters: infiltration, extravasation, leakage, air embolism, phlebitis, thrombophlebitis, thrombosis, catheter-related bloodstream infection, local infection, inflammation, catheter occlusion, catheter fracture or clot formation.
- Risks of arterial catheterization include injury to adjacent structures during insertion, infection, vascular spasm and thrombotic or embolic complications that may result in occlusion of the artery causing ischemia.
- Non-ischemic complications of arterial catheterization, or even repeated attempts at catheterization, include bleeding, pseudo aneurysms, arteriovenous fistulae, nerve palsy, infection and injury to the tendon sheaths and adjacent structures during insertion.
- Risks related to subcutaneous infusion therapies:
  - systemic side effects: acute cardiac failure and hyponatremia
  - local side effects: edema, inflammation, bleeding, cellulitis, erythema and pain.

#### Application

Use aseptic technique. Do not rotate catheter hub prior to insertion.

1. Gather all necessary items for insertion and stabilization procedure.
2. Select and disinfect per institutional protocol.
3. Apply tourniquet.
4. Flex stabilization wings before returning them to horizontal position. Remove protective guard

in a straight outward motion. Inspect device and confirm catheter hub is properly seated on flashback chamber.

5. Anchor the vessel with gentle skin traction and access vessel. Observe blood return in flashback chamber to confirm successful vessel entry (see figure A).
6. Lower and advance the entire catheter and needle unit slightly to ensure catheter tip is in vessel (see figure B).
7. Using the push-off plate, advance catheter off the needle (approximately 1/8" or 3mm). You can observe blood return between the needle and catheter to confirm catheter is in the vessel (see figure C). After confirmation, continue advancing catheter off the needle and into vessel.
8. Release tourniquet.
9. Gently press the stabilization platform to the skin to stabilize the catheter (see figure D). Withdraw the needle straight out with a controlled and continuous motion (minimize rotation or bending of the needle). The metal safety shield will automatically attach to the needle tip as the needle tip exits the catheter hub (see figure E). Blood flow from catheter hub is restricted after needle is removed.
10. Dispose of the shielded needle immediately into a sharps-proof container.
11. Immediately connect the infusion line or accessory device and cover the puncture site with a sterile and transparent dressing (see figure F) per institutional protocol.
12. After removal, dispose catheter according to local guidelines and/ or institutional protocols.

#### Duration of Use

- Duration of use is dependent on the applied therapy according to the SPC of the drug/ solutions and should follow the national guidelines and/or hospital protocols.

- The insertion site should be checked frequently and at regular intervals. The catheter should be removed in event of local or systemic signs of infection.

**Disposal**

Disposal according to local guidelines and/ or institutional protocols.

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

<b>en</b>	Do not re-use
<b>de</b>	Nicht wiederverwenden
<b>bg</b>	Да не се използва повторно
<b>cs</b>	Nepoužívat opětvorně
<b>da</b>	Må ikke genbruges
<b>el</b>	Να μην επαναχρησιμοποιείται
<b>es</b>	No reutilizable
<b>et</b>	Mitte korduvkasutada
<b>fi</b>	Ei saa käyttää uudelleen
<b>fr</b>	Ne pas réutiliser
<b>hr</b>	Nemojte ponovno upotrebljavati
<b>hu</b>	Ne használja fel újra
<b>it</b>	Non riutilizzare
<b>lt</b>	Negalima naudoti kartotinai
<b>lv</b>	Neizmantot atkārtoti
<b>nl</b>	Niet opnieuw gebruiken
<b>no</b>	Skal ikke gjenbrukes
<b>pl</b>	Nie używać ponownie
<b>pt</b>	Não reutilizar
<b>pt</b>	Não reutilizar
<b>BR</b>	
<b>ro</b>	A nu se reutiliza
<b>ru</b>	Не использовать повторно
<b>sk</b>	Opätovne nepoužívať
<b>sl</b>	Samo za enkratno uporabo
<b>sr</b>	Ne koristiti ponovo
<b>sv</b>	Får ej återanvändas
<b>tr</b>	Yeniden kullanmayınız
<b>vi</b>	Không tái sử dụng
<b>zh</b>	不得二次使用

Caution
Achtung
Внимание
Pozor (výstraha)
Forsigtig
Προσοχή
Atención
Ettevaatust
Huomio
Attention
Oprez
Figyelem!
Attenzione
Atsargiai
UZMANĪBU!
Let op
Viktig
Uwaga
Cuidado
Cuidado
Atenție
Внимание!
Varovanie
Previdno
Oprez
Obs
Dikkat
Thận trọng
警告

Consult instruction for use
Gebrauchsanweisung beachten
Вижте инструкциите за употреба
Čtěte návod k použití
Læs brugsanvisningen
Συμβουλευτείτε τις οδηγίες χρήσης
Consulte las instrucciones de uso
Vt kasutusjuhendit
Katso käyttöohje
Consulter les instructions d'utilisation
Pogledajte upute za uporabu
Nézze meg a használati utasítást
Consultare le istruzioni per l'uso
Žr. naudojimo instrukcijas
Lasit lietošanas instrukciju
Raadpleeg gebruikersinformatie
Se i bruksanvisningen
Przed użyciem należy zapoznać się z instrukcją
Consulte as instruções de utilização
Consulte as instruções de uso
Consultați instrucțiunile de utilizare
См. руководство по эксплуатации
Pozří návod na použití
Glejte navodila za uporabo
Pogledajte uputstvo za upotrebu
Se bruksanvisningen
Kullanma kılavuzuna bakınız
Tham khảo hướng dẫn sử dụng
查阅使用说明

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**CH REP** B. Braun Medical AG Seesatz 17 6204 Sempach Switzerland

Manufacturing site: B. Braun Medical Ind. Sdn. Bhd., Bayan Lepas Industrial Zone, 11900 Penang, Malaysia

B. Braun Melsungen AG  
34209 Melsungen  
Germany  
www.bbraun.com

**CE** 0123

1222 15328211

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- CA** Distributed by: **B. Braun Medical Inc.**