GIMA 23916
GIMA 23918
GIMA 23921

Shandong Lianfa Medical Plastic Products Co., Ltd
No.1 Shuangshan Sanjian Road, Zhangqiu, Jinan City,
250200, Shandong P. R. China
Made in China

01-1328-100, 01-1330-100, 01-1332-100

Linkfar Healthcare GmbH
Niederrheinstraße 71, 40474 Düsseldorf, Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

M23916-M-Rev.0.10.24

ENGLISH

BLOOD LANCETS

Model

I, II, III

Specification:

14G-33G

Intended Use:

The Blood lancets is intended for capillary blood sampling in order to obtain a small blood sample for various tests. It is designed for use by both healthcare professionals and lay users, providing a safe and convenient method for blood specimen collection in clinical and home healthcare settings.

Contraindication:

Use with caution in individuals with coagulation mechanism disorders.

Information about the matching lancing device:

LINKFAR Lancing Devices/ Lianfa Lancing Devices/ OneTouch® UltraSoft®/ OneTouch® UltraVue/ FreeStyle Lite®/ Bayer Microlet®/ Bayer Microlet® Vaculance/ Glucolet Lancing Device/ BD Lancet Device/ Easy Touch® lancing device/ Mylife Unio (Ypsomed)/ Mylife Pura (Ypsomed)/ Beurer Lancing Device Autolet Lancing Device/ Glucolet Lancing Device/ SD Biosensor Lancing Device/ TaiDoc Lancing Device/ CareSoft Lancing Device/ Tempo Lancing Device/ Autolet/ MICROLET@ NEXT/ Medifun/ SteriLance/ AIMSCO Lancing device


Product main structure:

The Blood Lancets primarily consists of lancet body, needle and cap.


The specific structure of the product is as follows:

Model I Model II Model III


Model I

Model II

Model III


1.Cap - 2.Needle - 3.Lancet body

Instructions for use

1. Wash hands with soap and warm water. Dry thoroughly.
2. Insert a new blood lancet into the lancet holder of lancing device.
3. Remove lancet cover (steps 1 & 2 above).
4. Consult the instructions for use of your lancing device.
5. Select a fingertip site, and proceed to puncture the skin.
6. After use, discard the used lancet in appropriate container.
7. Wash hands with soap and warm water again. Wear gloves if use in healthcare facilities.

Precautions and cautions

- Please read the instructions before use.
- Check the package before use. Do not use it if the package is damaged.
- Do not reuse lancets.
- Do not store in Lancing Devices.
- For use only on a single patient.
- It is forbidden to use after the expiration date
- Discard the entire device to biohazard container after use.
- Used Lancets should be treated as sharps waste and/or medical waste and disposed of properly following

Warning:

- Not intended for more than one use. Do not use on more than one patient.
- Improper use of blood lancets can increase the risk of inadvertent transmission of blood borne pathogens, particularly in settings where multiple patients are tested.
- The expiration time is 5 years. Using the product beyond the expiration date is strictly prohibited.
- For single use only. Do not use more than one patient. The device is deactivated after a single use, cannot be used more than once.
- The product has no therapeutic or diagnostic function.
- If the lancet's protective cap is damaged or lost, please do not use it.
- Use with caution for those with coagulation mechanism disorders.
- For obese users or individuals who require larger blood volumes, it is recommended to choose products with a smaller "G" value or large puncture depth for blood collection.

Shelf life: Valid period of product is 5 years.

Storage and Transportation:

1. The device shall be stored at an environment where the temperature is between -20 ° C and 40 ° C and the relative humidity of less than 80%. The storage environment shall be free of corrosive gases, dry, away from sunlight, well-ventilated, and clean.
2. The device shall be handled with care during transportation and handling.

Disposal:

- The product must be disposed of after use.
- Products contaminated with blood must be disposed of in accordance with the valid regulations for hazardous materials.
- The product must not be used after the expiry date!
- Safe disposal must be carried out in accordance with current national regulations.

Clinical benefit:

Help medical staff and patients complete capillary blood sampling in order to obtain a small blood sample for various tests.

Adverse events:

We searched for adverse events and state of the art literature on equivalent and similar devices in the market, and found no adverse events for equivalent and similar products. The products comply with the current technical standard. The product has been sold in the market and no adverse events have been found so far. Therefore, product safety and efficacy have not been affected. If any adverse events occur during the use of the product, please contact us and we will resolve the relevant issues as soon as possible.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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