
**Etilometro digitale GIMA - GIMA Digital Alcohol
 Computer - Alcootest digital GIMA - Alcoholimetro
 digital GIMA - Alcoolimetro digital**

**MANUALE D'USO E MANUTENZIONE
 USE AND MAINTENANCE BOOK
 INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
 MANUAL DE USO Y MANTENIMIENTO
 MANUAL DE USO E MANUTENÇÃO**

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ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

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ATTENTION: The user of the item agrees not to consider the manufacturer, the distributor or the retailer liable for the consequences of any decision based on the use of the item, as far as driving/using vehicles, boats, planes or other means of transport or harming anyone and causing any possible damage to things, are concerned.

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Thanks for purchasing GIMA Breath Analyzer, a hand held Digital Alcohol Computer. It gives you both digital and graphic indication of your breath alcohol content and thus your possible concentration of intoxication.

Operation Instructions

- 1) Install 3 x 1.5V alkaline batteries (AAA size) into the battery compartment (auto-cleaning will automatically activate according to the CAUTION manual)
 - 2) To begin, wash your mouth with water before using the unit. Abstain from smoking or drinking alcohol for at least 20 minutes.
 - 3) Hold **MODE** Key for 0.5 seconds to toggle Real Time clock / Alcohol mode. If less than 0.5 seconds, turn backlight on.
 - 4) Insert the provided blowing tube (or drinking straw) into the gas inlet located on top of the LCD screen.
 - 5) In Alcohol mode, press the **ST** key once to activate an approximate 17 seconds sensor-warm-up process. "PH" flashes for about 4 seconds and then followed by 13-second count down. "READY" flashes with beep tone when count down to zero and ready to start. (or directly press the **ST** key to start the test cycle in the clock mode)
 - 6) Press & hold the **BLOW** Key, while blowing continuously through the blowing tube for approximately 4 seconds. Stop blowing when you hear a single beep tone and a reading will be displayed on the display.
 - I. A measurement of less than 0.15% BAC will show "LO"
 - II. A measurement between 0.15% & 1.50% BAC will show its digital reading with circuit bar indicators.
 - II. A measurement above 1.50% BAC, "HI" & "NO DRIVING" will flash on the display.
- Press **SET** key once to convert the alcohol unit in Alcohol mode
 Press **ADV** key once to scroll forward the 4 memory readings
 Hold **ADV** key for 2 seconds to clear the memory readings in alcohol mode.
- 7) Hold **MODE** Key for 0.5 sec to enter into Real Time clock: Hold **SET** key for 2 second to enter into time setting in real time mode. Press **SET** key once to decrease the selected digit in real time setting mode, press **BLOW** key once to shift the select digit in real time setting mode, press **ADV** key once to advance the selected digit in real time setting mode.

Alert Signals

GIMA Breath Analyzer is designed to alert user as following:

- 8) Low Level – For measurement below 0.15% BAC, "LO" will be shown on display with GREEN backlight.
 - 9) Safety Level – For measurements between Low Level (0.15% BAC) and Pre-alert Level (Warning Level), the alcohol concentration level shows on the display with GREEN backlight.
 - 10) Pre-alert Level – For measurements within 0.15% BAC below the alert point. The analyzer beeps together with "NO DRIVING" & "HI" flashes alternatively for 3 times. GREEN & ORANGE backlight exchanges alternatively as well. The "HI" will then remain on the LCD with ORANGE backlight. The alcohol concentration level displays as per the bar indicator.
- DON'T DRIVE: You are very close to the alert level or you have reached the alert level if take tolerance of the unit into account**
- 11) Alert Level – For measurements equal to or above the alert concentration level, the analyzer beeps together with "NO DRIVING" & "HI" flashes alternatively for 15 seconds. ORANGE & RED backlight exchanges alternatively. The alcohol concentration level displays and "NO DRIVING" will then remain on the LCD. If reading greater than 1.5% BAC, "NO DRIVING" & "HI" will remain on the LCD. Press **MODE** key can stop the alert sound.
- YOU ARE PROHIBITED TO DRIVE since you already exceed the alert level.**

Note

- 1) To ensure high precision readings, user should perform the auto-cleaning if unused for a long period. (especially for the first time of using). User can start the auto-cleaning process by holding **MODE** & **ST** keys together for 2 seconds or by removing the batteries and power up again. Put a few drops of alcohol solvent (~2%) inside the provided cleaning cap accessory. Place the cap onto the gas inlet when you see the animation of "CAP u" > "CAP n" for 10 seconds. Then take off the cap when "Auto" & "Clr" flashes and leave the GIMA Breath Analyzer perform the 5-minutes auto cleaning. The analyzer will then return to normal Alcohol mode automatically when the cleaning process finishes. Make sure no alcohol residues left by blowing fresh air into the gas inlet before formal testing.
- 2) Please notice that pressing the **MODE** key at any time can easily terminate the auto-cleaning process.
- 3) To achieve an accurate reading, user should ensure there is no alcohol residue inside the chamber before starting a new alcohol test. Users can pre-check the tester until the 'Lo' indication is shown on the display. This may be achieved by pressing **ST** key once. Holding the **BLOW** key and blow fresh air into gas inlet for about few seconds. Alternatively, shaking the unit forwards and backwards several times will allow any alcohol residue to leave before starting a new alcohol test.
- 4) After any 5 consecutive tests, allow the device to rest for 2 minutes.
- 5) An "Err" flashes on the display for 5 seconds for the following improper operations:
 - I. When user does not hold the **BLOW** key for at least 4 seconds during blowing until a beep sound is heard.
 - II. When user fails to start the test within 5 seconds after "READY" flashes.
 User can press **MODE** Key once to stop the alarm sound of Err.
- 6) When the product is unused, please store in a dry environment.
- 7) When low battery inductor flashes, the reading is not valid, please replace the batteries immediately.
- 8) If the difference between two successive reading >0.1%, user should test again.

The GIMA Digital Alcohol Computer should be used as a general guide to the level of alcohol in your body. Many factors (e.g. environment gas residue, smoke & various smell in human mouth, ect...) may affect the BAC level, so the Readings are For Reference Only. Neither the reseller nor manufacturer accept any liability arising from the use of this unit
DO NOT drink before you drive!

Product Specifications (specification subject to change without notice)

1. Range: 0.15 ~ 1.50 ‰ BAC	7. Dimensions: 106.5 X 45 X 24 mm
2. Power Input: 3 X 1.5V AAA batteries	8. Weight: 50g (without batteries)
3. Display: Graphic & Digital display with tri-color backlighting	9. Alert level calibrated: 0.5% BAC
4. Real time clock power save mode: 2 minutes after any key in	10. Agency Approval: CE
5. Breath sampling: About 4 seconds continuous breath sample	11. Accessories: 3 pcs blowing tubes, 1 pc Auto-cleaning cap
6. Recycle time: ~20 seconds	

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.
 For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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