

TERMOMETRO DIGITALE BL3 - AMPIO SCHERMO
BL3 DIGITAL THERMOMETER - WIDE SCREEN
DIGITAL THERMOMETRE BL3 - ÉCRAN LARGE
TERMOMETRO DIGITAL BL3 - PANTALLA ANCHA
TERMOMETRO DIGITAL BL3 - ECRĂ AMPLO
ΨΗΦΙΑΚΟ ΘΕΡΜΟΜΕΤΡΟ BL3 - ΕΥΡΕΙΑΣ ΘΟΝΗΣ

Manuale d'uso - User manual - Manuel de l'utilisateur - Guia de uso
Guia para utilização - Οδηγίες χρήσης

5. Once the degree sign °C (°F) on the display stopped flashing, the correct temperature is indicated.

The temperature reading will not change after the °C (°F) stops flashing.

6. The unit will automatically turn off in 1 minutes. However, to prolong battery life, it is best to turn the thermometer-off by pressing the ON/OFF button once the temperature has been noted.

Special Features:

Normal Temperature Alarm An alarm will sound when the peak temperature has registered and is ready-to-read. The "Normal Temperature" alarm (four quick beeps followed by a pause) will sound repeatedly for approximately 10 seconds.

Last-Temperature Memory Recall Switch the thermometer on, a short beep will sound. At the same time the thermometer runs a self-check test. During which all the digital segments appear on the LCD. After that the last reading value will appear automatically, and the 37.0°C (98.6°F), on the LCD for 2 seconds.

The reading is only over-written when a new temperature value is recorded.


For Models with Fever Alarm If the temperature registers higher than 38.0°C (100.4°F), the fever alarm (a quick beeping sound which occurs continuously for about 10 seconds) will sound when the peak temperature has registered. Also, the result will keep on the LCD until the thermometer is turned off.

For Switchable Centigrade and Fahrenheit Models Hold down ON/OFF button, until °C or °F appeared. After that the °C or °F can be switched. When the button is released, the thermometer will be in "ready-in-measure" mode.

NOTE: Always disinfect the thermometer before using.

NOTE: If the room temperature is less than 32.0°C (89.6°F), it will be displayed rather than L °C (or L°F). **NOTE:** You may put a sanitary plastic probe shield on the sensor tip to help prevent cross-infection and contamination. Use of the probe cover may result in a 0.1°C (0.2°F) difference from actual temperature. Dispose of any used probe shield and store the unit in a safe place.

BATTERY REPLACEMENT


1. When the symbol  appears, the battery needs to be replaced.









Disposal: The product must not be disposed of along with other domestic waste.
The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.














ENGLISH

SPECIFICATION

Range:	32.0-43.9°C (89.6-111.1°F) Below 32.0°C (89.6°F) displays L°C (°F) Above 43.9°C (111.1°F) displays H°C (°F)
Resolution:	0.1°C (°F)
Accuracy:	±0.1°C (34.0°C to 42.0°C); ±0.2°C (other range) ±0.2°F (93.2°F to 107.6°F); ±0.4°F (other range)
Display:	Liquid crystal display 3 ½ digits
Battery: (included)	CR2032 (3V)
Power Consumption:	0.40 milliwatts in measurement mode
Battery Life:	More than 180 hours of continuous operation
Dimensions:	Approx. 117mmx30mmx16mm
Weight:	Approx. 20grams including battery
Beeper:	Approx. 5 seconds sound signal when peak temperature reached
Memory:	Built-in memory enabling display of previously measured value
Working Conditions:	Temperature: 10-40°C (50-104°F) Relative Humidity: 15-95% non condensing
Storage Conditions:	Temperature: -25-55°C (-13-131°F) Relative Humidity: 15-95% non condensing
Guarantee of Quality:	Certification ISO 13485 Comply with ASTM E112, EN12470-3, EN60601-1
Product Classification:	Type BF equipment 

HOW TO USE

- Disinfect the probe with ethyl (rubbing) alcohol before using.
- Depress the ON/OFF button. The display will read  or  →  or  (Last Memory)
- Release the power switch and the display will show L°C (L°F) with °C (°F) flashing.
- Place the probe in the appropriate position (oral, axillary or rectal).  or 

	IT Codice prodotto GB Product code FR Code produit ES Código producto PT Código produto GR Κωδικός προϊόντος		IT Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB Caution: read instructions (warnings) carefully FR Attention: lisez attentivement les instructions (avertissements) ES Precaución: lea las instrucciones (advertencias) cuidadosamente PT Cuidado: leia as instruções (avisos) cuidadosamente GR Προσοχή: διαβάστε προσεκτικά τις οδηγίες (ειδοποιήσεις)
	IT Numero di lotto GB Lot number FR Numéro de lot ES Número de lote PT Número de lote GR Αριθμός παρτίδας		IT Seguire le istruzioni per l'uso GB Follow instructions for use FR Suivez les instructions d'utilisation ES Siga las instrucciones de uso PT Siga as instruções de uso GR Ακολουθήστε τις οδηγίες χρήσης
	IT Fabricante GB Manufacturer FR Fabricant ES Fabricante PT Fabricante GR Παραγωγός		IT Parte applicata di tipo BF GB Type BF applied part FR Appareil de type BF ES Aparato de tipo BF PT Aparelho de tipo BF GR Συσκευή τύπου BF
	IT Rappresentante autorizzato nella Comunità europea GB Authorized representative in the European community FR Représentant autorisé dans la Communauté européenne ES Representante autorizado en la Comunidad Europea PT Representante autorizado na União Europeia GR Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση		IT Smaltimento RAEE GB WEEE disposal FR Disposition DEEE ES Disposición WEEE PT Disposição REEE GR Διάθεση WEEE
	IT Conservare in luogo fresco ed asciutto GB Keep in a cool, dry place FR A conserver dans un endroit frais et sec ES Conservar en un lugar fresco y seco PT Armazenar em local fresco e seco GR Διατηρείται σε δροσερό και στεγνό περιβάλλον		IT Dispositivo medico conforme alla Direttiva 93/42/CEE GB Medical Device complies with Directive 93/42/EEC FR Dispositif médical conforme à la directive 93/42/CEE ES Dispositivo médico según la Directiva 93/42/CEE PT Dispositivo médico em conformidade com a Diretiva 93/42/CEE GR Ιατρική συσκευή σύμφωνα με την οδηγία 93/42/CEE
	IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR A conserver à l'abri de la lumière du soleil ES Conservar al amparo de la luz solar PT Guardar ao abrigo da luz solar GR Κρατήστε το μακριά από ηλιακή ακτινοβολία		
	IT Data di fabbricazione GB Date of manufacture FR Date de fabrication ES Fecha de fabricación PT Data de fabrico GR Ημερομηνία παραγωγής		