
GIMA

ARTICOLI & APPARECCHI PER MEDICINA
PROFESSIONAL MEDICAL PRODUCTS

27235

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel 199 400 401 - fax 199 400 403 - gima@gimaitaly.com
Export: phone + 39 02 953854209/221/225 - fax + 39 02 95380056
-export@gimaitaly.com - www.gimaitaly.com

1. Technical data

Maximum capacity: 150kg/330lb


Minimum capacity: 3kg/6.6lb

Weight tolerance: $\pm 1\%$


Weight graduation: 100g

Battery :2xCR2032

2. Preference SET

(1) Press "SET" button to switch on. Use button "▲" and "▼" to choose a memory unit, for example, "P-01". It will back to "0.0" within 5 seconds without operation.


(2) Press "SET" button to confirm, then Press button "▲" and "▼" to choose weight unit (Kg, lb, st), and then press "SET" button to confirm.

(3) After setting is finished, the last measurement will be cleared away.


3. Use of Bathroom Scale

Weighing Mode


(1) Knock on the plastic lightly, the screen will display "8888" then back to "0.0", and then the user can stand on the scale to measure weight. (In this mode, the scale just measure weight only and will not remember the measuer value.)

Weighing and Memory Mode

(2) Press "SET" to turn on scale. Press "▲" or "▼" to choose a memory unit. The chosen memory unit will flash for a few time, and then it display "0.0".

(3) Stand on scale steadily. It will display the weight value. If any, it will also display last weight value, weight difference and percentage of difference respectively for two rounds. Minus will appear if weight value is smaller that last weight value. However weight difference in lb does not come with minus. If weigt difference in lb is negative and over 100lb, it will display "0-Er".

(4) When the scale back to "0.0kg", you can press the swith button to change wheight unit.

4. Inserting and Changing Battery

The battery compartments of the scale are situated on the back. Before the first use, open the battery compartments of the scale and insert two new "CR2032" batteries respectively in the right polarity.

Important:

When you change your batteries, be careful!! You must discard the batteries according to the rules in force. Dump them in an appropriate place in order to assure their recycling and the protection of environment.

5. Maintenance Care

- (1).The scale is intended for household use only.
- (2).The scale should be placed on a firm and level surface (avoid carpet).
- (3).Be ware of slippery platform.
- (4).The scale should never get wet.
- (5).Do not pile sundries on it.
- (6).Only clean the scale with soft cloth and neutral cleaning products.
- (7).Batteries should be removed from the battery compartments if the scale is not to be used for a long period of time.
- (8).Do not strike the scale with hard and sharp objects and avoid strong vibration.
- (9)Do not attempt to disassemble or repair the scale by yourself.

6. Troubleshooting

Indication	Reason	Remedy
O-Ld	Maximum weighting capacity exceeded	Step off the scale immediately. Avoid overloading.
LO	Low voltage indication.	Replace batteries with new ones.

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Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel travelling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.