OTOPLUS DC

0476

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English

These instructions apply to all the versions of the OTOPLUS unit and all the related accessories. For this reason, not all the paragraphs of this manual will refer to the equipment you have bought.

EUROCLINIC will provide the authorised service center with all the information necessary to install and repair the equipment.

No part of this manual is to be reproduced, stored in a retrieval system or transmitted in any form or in any means, i.e. electronic, mechanical, photocopying, translation or otherwise, without the prior written permission of EUROCLINIC.

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1. INTRODUCTION

- 1.1 General characteristics and identification of the parts
- 1.2 Cautions
- 1.3 Symbols

2. INSTALLATION**3. USE OF THE INSTRUMENTS**

- 3.1 Front commands
- 3.2 Insufflator
- 3.3 Syringe
- 3.4 Suction
 - 3.4.1
- 3.5 Medicine nebulizer

4. USE OF THE ACCESSORIES

- 4.1 Back Socket
- 4.2 Germicide compartment
- 4.3 Clar lamp
- 4.4 Larynx mirror-warmer

5. MAINTENANCE

- 5.1 Empty condense filter
- 5.2 Filling the water tank
- 5.3 General maintenance and spare parts
- 5.4 Cleaning secretions tank
- 5.5 Cleaning of the unit
- 5.6 Discharge and scrap
- 5.7 Disinfectant products chart

6. TECHNICAL CHARACTERISTICS


- 6.1 Technical data charts – absorbed power
- 6.2 Working periods of the instruments summary – identification of the unit

7. IDENTIFICATION LABELS

- 7.1 Serial Number Label
- 7.2 Protection Mass Label
- 7.3 Disposal

8. GUIDE CHARTS AND PRODUCER'S DECLARATIONS

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

1. INTRODUCTION

1.1 General characteristic and identification of the parts

Unit for otorhinolaryngology completely autonomous and mobile on castor wheels composed of:

- A) Instrument tray holder which is able to hold up to a maximum of 4 instruments. It has articulated arms which enables one to position it in the most comfortable of ways.
- B) Suction instruments holder
- C) Duropall working bench complete with waterproof back panel: a back panel which houses different sockets is available by request.
- D) Housing for the auxiliary instruments complete with a series of drawers of various sizes and dimensions according to one's choice

1.2 CAUTIONS

- 1) Do not install the machinery in environments which are exposed to explosions or frost.
- 2) The temperature and relative humidity and environmental temperature in which the machinery is installed should be respectively included between 10°C and 40°C, 30% and 75%, 700 mBar and 1060 mBar.
- 3) During installing, check that the electrical tension on the sticker on the machine corresponds to that in use in order to avoid damaging.
- 4) The building, repairing, modifying, calibration and all the other operations that require the opening of the unit should be carried out by technical staff authorized by EUROCLINIC.
- 5) The electrical circuit in which the unit is connected to, should comply to the IEC 64.8 sec. 710 regulation (regulation regarding electrical circuits in rooms used for medical purposes).
- 6) The machinery should be used exclusively following the instructions illustrated in this manual.
- 7) EUROCLINIC refuses any responsibility for the following: damages or accidents caused by the carelessness of the user, the non compliance to the safety rules, a non correct use of the machine, a use of the machine which is different to that of its purpose. EUROCLINIC denies all responsibility which has nothing to do with faults in the material and non correct building.
- 8) Before installation and functioning check that there are not loose parts, burnt electric cables, or various anomalies and if they are evident please contact EUROCLINIC or authorised technical staff.

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

9) Mobile phones and radios can influence the good functioning of medical device OTOPLUS.

10) If during the use of unit OTOPLUS happen strange changes of performance or various anomalies we suggest to interrupt all operations, switch off the machine and contact authorised staff EUROCLINIC in order to avoid any kind of damage.

11) The medical device OTOPLUS must be used matching (electrical connection) just other medical devices **and with exclusive use destination the patient area**

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

1.3 Symbols

Classification of the "BF" type machinery.

"BF" type machinery (machinery of the type I,II,III or machinery housing an internal electric source with a protection degree against direct and indirect contact and in particular for:

- permitted current dispersion
- reliability of the ground contact)

Class I machinery.

This symbol indicates: FOR FURTHER INFORMATION CONSULT THE ATTACHED DOCUMENTATION

This symbol indicates: ATTENTION !!! CONSULT THE ATTACHED DOCUMENTATION

If the sockets situated on the back panel are used for feeding other instruments different to the unit, it is advisable to check that they comply with the IEC 62,5 law and that they need the following current: 220/240 V AC 4A

Earthing contact of the machinery

Alternated current

Symbol discharge and scrap

2. INSTALLATION

The installation of the equipment must be made respecting the instructions of the following booklet of instructions. It must be made by authorised staff Euroclinic, or thanks to the simplicity of operations also by the user themselves under their responsibility.

The user should always make sure that the electric circuit complies to IEC values (in other words: 220/240 V, alternate monophasic 50 Hz, the plant should comply with IEC 64.8 sec. 710 regulation). The operations are:

- Full the tank but not to the rim
- Insert the plug in a socket

It is advisable to insert the plug in a socket with a bipolar switch compliant to IEC and having the following technical characteristics:

230 V

50 A

contact opening of 3 mm.

To check that the installation has been properly made it is enough to check that the switch on button is properly lit up (SEE PICTURE 1).

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

Figura 1
Picture 1

00	Gianni Andraghetti	08-01-2005	Longo M.	08-01-2005
Release	Redatto	Data	Approvato	Data

3. USE OF THE INSTRUMENTS

In order to use the instruments of the machine it is necessary to proceed with the switch on of the machine by pressing the button 1 on front of the work bench. When the green light is on, this means that the machine is working. Turn it off by pressing the same switch again.

3.1 Frontal comands

The following comands are present on the front panel of the unit: it is possible to use them only after turning the general switch on (see picture 1)

- 1- general switch
- 2- water heating switch
- 3- potentiometer for Clar lamp.
By turning it clockwise the light increases by turning it the opposite way it decreases

3.2 Insufflator

For a correct use of the instrument it is advisable to work according to the following periods:
 work 5 min
 pause 10 sec

The functioning of the instrument stars when it is taken out of its housing and pushing the pedal. When the corresponding light on the instrument holder is on, this means that the instrument is working. In order to vary the intensity of air, one should adjust it with the lever behind the cannula holder.

N.B. In the "DC" version, the pressure is visualized in a pressure gauge, head-on the unit
N.B. Use only spare parts by EUROCLINIC.

3.3 Ear washings syringe

For a correct functioning of the instrument we suggest to use it according to the following periods:
 Work 5 min
 Off 10 sec

Before using the syringe it is necessary to fill up, not completely full, the tank.

The functioning of the instrument is activated by taking it out of its housing. When the corresponding indicator-light on the instruments tablet is lit up, this means that the instrument is functioning.

For obtaining warm water flux, firstly turn on the warm water switch placed in the front of the unit (see section 3.1): the led stand-by on the tablet starts blinking. When the right temperature is reached, the led is lit up and the acoustic sign is activated. Now it is possible to push one of the two buttons placed in the syringe hand-piece.

N.B. Use only spare parts by EUROCLINIC.

N.B. The First cycle of heating lasts approximately 15 minutes

N.B. During the first heating cycle or after a prolonged pause with the unit off, as soon as the right temperature is reached, it is always necessary to make overflow from the syringe a water jet for 10 second, in order to settle the unit.

00	Gianni Andraghetti	08-01-2005	Longo M.	08-01-2005
Release	Redatto	Data	Approvato	Data

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N.B. In case that the water temperature overcome the set value, a safety thermostat blocks automatically the ear washing syringe and a fix acoustic advising signal is activated, until the normal cycle is not re-activated.

WARNING!! In case of reserve please immediately switch off the equipment and fill up the tank.

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For a correct use of the instrument it is advisable to work according to the following periods:

Work 5 min
 Off 10 sec

The starting of the instrument is obtained by taking it out of its housing. When the corresponding light on the tray is on this means that the instrument is working. To vary the intensity of the aspiration it is necessary to adjust the lever placed on the end of the cannula holder. (see previous picture)

N.B. In the "DC" version, the aspiration is visualized in a vacuum gauge, head-on the unit.

Notice: If the container is full of secret operative system overflow aspiration is stopped, so to drain (see 5.2)

Notice: The useful capacity is stated on the vessel, the meeting can not work on inclined planes.

NB: Use only spare parts by EUROCLINIC.

3.4.1. Accessori x Aspiratore

The unit Otoplus is set with a pulling-out neck headed by a piece of silicone pipe (see pic. 3.4.1)); moreover together with the machine is supplied an apposite Kit for setting the funnel for secretions collecting (see pic. 3.4.2) – which contains a pipe fitting to be coupled to the suction pipe + funnel (see assembling phase).

Antimicrobial hydrophobic filter

On the Unit is planned an antimicrobial hydrophobic filter on suction (see pic. 3.4.3), in order to ensure the killing (99%) of any bacteria that come from the suction system. The replacement **is advisable** every 2 – 5 working day.

Notice: EUROCLINIC RECOMMEND THE USE OF PVC PIPES OR FROSTED SILICONE SIZE OF 8X12mm

3.5 Medicaments Atomizer (See Pic. 2)

For a correct use of the instrument it is advisable to work according to the following periods:

Work 5 min
 Off 10 sec

The functioning of the instrument is obtained by taking it out of its housing and pushing the pedal. When the corresponding light on the tray is on this means that the instrument is working. In order to obtain the nebulization of the medicine products it is necessary to adjust the handpiece of the instrument. (see previous picture).

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
**LAMPADA DI CLAR/
CLAR LAMP**


**SUPPORTO OTTICHE
*Fixed Holder for endoscope***

**VANO GERMICIDA/
*GERMICIDE COMPARTMENTS***

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4. USE OF THE ACCESSORIES

4.1 Socket back panel

When installing the machine, the plugs are automatically connected to the electric system. Check that other instruments connected to supplementary sockets can function with 220-240 V alternated monophase 50Hz.

The plug is protected by 4 A fuse that operates with charge of over 900VA.

Always take the plugs of these other instruments out after their use and during cleaning; do not connect more than one instrument to each socket.

4.2 Germicide compartments (See Pic.3)

On request, the machine can be equipped with this auxiliary drawer, that permits only the maintenance in hygienically conditions of the instruments already sterilized.

The functioning of the UV lamp is obtained by turning the general switch of the unit on (see chapter 3).

The germicide lamp can be turned on and off by a switch placed on the side of the drawers turns the lamp off automatically if the germicide compartment is accidentally opened, and it will turn the lamp on again when the drawer is closed.

WARNING! We repeat once more that the germicidal space does not constitute substitutive element of the sterilization process, but function only as co-adjuvant and supporter.

NB: do not use lamps with technical features different from those indicated by EUROCLINIC.

N.B. We suggest to substitute the lamp each 500 hours of functioning

ATTENTION! Take the plug of the Unit out before carrying out any operations on UV circuit.

4.3 Clar Lamp (Pic. 3)

The lamp is turned on automatically by taking it from its housing, to turn it off it is necessary to place it back.

The adjusting of the intensity of the light is obtained by adjusting the potentiometer placed on the front (see chapter 3.1)

4.4 Mirror-warmed larynx

The functioning of the mirror-warmed larynx is obtained by pressing the switch placed next to the fan (see figure).

Place the surgical instrument in front of the warm air flux and wait a few seconds.

WARNING !

N.B. When overused, the system has a thermal safety switch which stops the air flow until normal values of temperature are reached again; it is necessary in that case to wait for some minutes before using it again.

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

5. MAINTENANCE

Attention! The following proceeding have to be carried out when the machine isn't plugged into the mains.

5.1 Empty of the condense filter

Near the panel there is a filter for collecting the condense created by the compressor; in case it is filled it is necessary to empty it; place a glass under the filter, turn the knob and pull downside the knob, empty the condense and close again the knob!

Warning! If the knob is not closed the compressor leaks air through it!

5.2 Filling the water tank

In acoustic case of signal intermittent it means that it is intervened the reserve of the tank and that it is necessary to fill it up. In order to fill the tan please proceed in the following way:

- 1- take the plug out of the machine
- 2- open the door of the auxillary instruments compartment
- 3- take the tank out, unwind the tap and with a funnel pour the water,without reach the maximum level .
- 4- close the tap and the door
- 5- plug the machine in again

Notice: Use only distilled water

Vano alloggiamento
 Auxiliary instruments
 compartment

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

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5.3 General maintenance and spare parts



To guarantee the correct functioning of the instrument EUROCLINIC suggests to effectuate a periodical revision Of the entire machine.

EUROCLINIC does not take charge of any responsibility for damages, hitches or inefficiencies due to the substitution of non-original spare parts.

5.3.1 OIL PUMP MAINTENANCE (optional)

1. Oil Level Check

It is suggested to check the oil level of the pump's engine every 2-3 working months, according to the hours of activity sustained by the Pump:



2. Oil Refill

After approximately 500h of suction use, and in case of darkoil due to its deterioration, proceed to its replacement by following point 5.4 of instructions manual DVP (see pict. 1).

3. Filter Replacement

The pump's filter is able to reduce residual fumes for 2000 working hours, after such period of time it is necessary to replace it with a new filter; see the following pictures 2,3:

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

5.4 Cleaning the unit

The choice of materials used for building the unit has been carefully taken in order to correspond to the cleaning and disinfecting needs. Given the big differences between chemical products used in a medical study, it is possible to ruin the tops. The major damages are caused by the amount of time that these chemical products are left on the tops. It is important to immediately dry the interested parts with a damp cloth. To clean the painted or metal parts, use a specific product for these materials which doesn't damage their composition. IMPORTANT! Do not use products with alcohol, ammonia or petrol

5.4.1 Cleaning and maintenance of secretion tank

The tank which contains the aspirated liquids can be cleaned with water and detergents which are compatible with the materials (see chart) used and/or sterilized in an autoclave at 135° C. Separate the main parts of the tank before the sterilization, insert the glass upside down in the autoclave and do not place anything on top of it whilst disinfecting.

DO NOT USE SOLVENTS OR ALCOHOL FOR CLEANING AND DISINFECTING. THE USE OF THESE PRODUCTS CAN DAMAGE THE PLASTIC PIECES OF THE SYSTEM.

Always substitute the damaged or faulty parts with original spare parts supplied by EUROCLINIC.

In case of cleaning in autoclave please stick strictly to the temporary chart shown in following table 5.5

5.5 Vase secreted disposable

The unit can be setted for single-use kit of different type. **Please keep attending during the replacement of the ampoules in order to avoid to invert the suction tubes VACUUM and PATIENT, because it could damage the motor. If the intake of liquids or solids into the engine, immediately turn off the suction and call the service center for maintenance.**

5.6 Empty of the condense

Near the panel there is a filter for collecting the condense created by the compressor; in case it is filled it is necessary to empty it; place a glass under the filter, turn the knob and pull downside the knob, empty the condense and close again the knob!

Warning! If the knob is not closed the compressor leaks air through it!

5.7 Discharge or scrap

The machine does not contains dangerous elements and can therefore follow normal procedures of discharging and scraping, in harmony with rules of Country where the unit has worked.

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

5.8 Condensate drainage

Near the door there is a filter for the condensate collection, that is created by the compressor; if it will be full it is necessary to discharge it. Put a glass under the filter, turn and pull under, drain the condensate and close the knob!

Attention! If the knob will not be closed the compressor let pass the air through it.


Fig. 4
 Pic. 4

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

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5.8 Instrument sterilization chart

Instrument sterilization chart		
INSTRUMENT	TEMP	KIND OF STERILIZATION
Syringe handpiece	135°C	In autoclave 5 – 7 min
Nebulizer terminal	135°C	In autoclave 5 – 7 min
Insufflator cannula	135°C	In autoclave 5 – 7 min
Aspirator cannula	135°C	In autoclave 5 – 7 min
Cleaning and disinfecting products chart		
Disinfectants	Pharmaceutic products and cosmetics	
+ bactol 5% + clorammina - DDT + delegal 5% = dimamina 5% = iodine dye - fenic acid + lysoform 2% - TB lysoform + bactol + merfen 2% + octozone + perydrol + resorcina 1% = sagrotane 5% + alcohol to be purely burnt + mercury chloride + trosilina G + oxiginated water = zefirol	+ blood plasma + sun tan lotion + hydroplex = iodine dye + lanolyn = mantolo 90% in alcohol -nail varnish - solvent + mouth wash + periston R + vaselyn + vicks rub	
CHART SIGNS + compatible - non compatible = compatible but with caution		

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

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6. TECHNICAL FEATURURES
6.1 Technical data chart

Electrical tension	220/240 V
Frequency	50 Hz
Power absorbed by the unit (according to the configuration)	1200 W
Weight	120 Kg
Maximum number of instruments	4 pcs
Maximum air pressure	2.8 Atm
Air capacity	26 l/m
Maximum water pressure	2.2 Atm
Water capacity	0.4 l/m
Dimensions (in mm)	
Height of the working bench	850 mm
Width	850 mm
Depth	540 mm

Absorbed power chart

DESCRIPT.	ABSORB.	TENSION	SERVICE
Germicide Lamp	8 W	230 V ac	Continuos
Aspirator	70 W	230 V ac	Continuos
Compressor	70 W	230 V ac	Continuos
Pump	48 W	230 V ac	Continuos

6.2 Summary of the working periods of the instruments
Device identification summary

INSTRUMENT	WORK	OFF
Syringe	5'	10"
Insufflator	5'	10"
Aspirator	5'	10"
Nebulizer	5'	10"

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
7.3 DISPOSAL

The symbol of crossed out garbage bin (on the back of the equipment) indicates that the product at the end of its working life must be collected separately from other garbage. This refers to the following law: art. 13 of D.L: 25 /07/2005, n.151 : "Implementation of the regulations 2002/95/CE, 2002/96/CE and 2003/108/CE, about decrease of dangerous substances on electric and electronic equipments" and moreover it refers to garbage disposal.

The separate collection of rubbish of this equipment (felt into disuse) must be arranged and coordinated by the manufacturer. When the user wants to throw away the equipment, he will contact the Manufacturer Company that will indicate all procedures to follow for the disused equipment separate collection of rubbish.

Therefore, the appropriate separate collection of rubbish, aimed to actuate afterwards disused equipment recycling, treatment and to ambient garbage disposal, contributes to avoid possible bad effects on the environment and on the health. In addition, these procedures raise the reemployment and/or the recycling of the materials of which the equipment is composed.

The unauthorized garbage disposal, on the part of the holder, requires administrative law sanctions application.

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The **OTOPLUS** is intended for use in the electromagnetic environment specified below. The customer or the user of the **OTOPLUS** should assure that it is used in such an environment

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At 80 Mhz and 800 Mhz is to apply the higher frequency range.

These guidelines could not be applicable to every situation.

The electro-magnetic propagation is influenced by absorption and reflection of structures, objects and people. The field intensities for fixed transmitters as base-stations for radio-telephones (mobiles and cordless) and earth radio-mobiles, radio-amateurs sets, radio transmitters in AM or FM and TV transmitters cannot be forecasted with absolute certainty and precision. To detect an electro-magnetic environment caused by RF fixed transmitters, an electro-magnetic inquiry of the site itself should be considered.

*If the field intensity measured in the place where a **SURPLUS** is used overcome the applicable conformity level specified before, the normal functioning of **SURPLUS** unit could be placed under observation. If strange performances are noticed, it could be necessary to proceed to further measures, like a different positioning and orientation of the **SURPLUS**. The field intensity on a frequency range from 150 KHz to 80 MHz should be lower than 3V/m.*

01	Francesco Marrone	30-03-10	Longo M.	30-03-10
Release	Redatto	Data	Approvato	Data

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Recommended separation distances between portable and mobile RF communications equipment and the L'OTOPLUS

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The OTOPLUS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the OTOPLUS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the OTOPLUS as recommended below, according to the maximum output power of the communication equipment.

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*For transmitters rated at a maximum output power not listed above, the recommended distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.*

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