HAPPYNEB II

GUARANTEED

3

YEARS


PROFESSIONAL

NEBULIZZATORE
NEBULIZER
INHALATOR
NÉBULISEUR
NEBULIZADOR

HAPPYNEB II

3A HEALTH CARE S.R.L.
Via Marziale Cerutti 90/F, Loc. S. Tomaso
25017 Lonato (BS) - ITALY
Ph. +39 030 9133177
Fax +39 030 9919114
www.3-a.it

SIMBOLOGIA - SYMBOLS - SYMBOLIQUE - VERWENDETE SYMBOLE - SIMBOLOGÍA

□ Apparecchio di classe II/Class II/Appareil de classe II/Schutzklasse II/Aparato de clase II.

⚡ Apparecchio di tipo B/Type B device/Appareil de type B/Geräteyp B/Aparato de tipo B.

○ Interruttore spento/Switch "OFF"/Interrupteur éteint/Ein-Ausschalter ausgeschaltet/Interruptor apagado.

I Interruttore acceso/Switch "ON"/Interrupteur allumé Ein-Ausschalter eingeschaltet/Interruptor encendido.

~ Corrente alternata/Alternating current/Courant alternatif/Wechselstrom/Corriente alterna.

 Non utilizzare l'apparecchio mentre si fa il bagno o la doccia/Never use the unit while having a bath or shower/Ne pas utiliser l'appareil pendant le bain ou la douche/Das Gerät darf nicht während des Badens oder in der/No utilizar el aparato mientras se hace el baño o la ducha.

 Leggere attentamente le istruzioni prima di utilizzare questo dispositivo/ Carefully read the instructions before using this device/ Veuillez lire attentivement les instructions avant d'utiliser le dispositif/Vor dem Gebrauch dieser Vorrichtung die Anleitung aufmerksam lesen/ Leer detenidamente las instrucciones antes de utilizar este dispositivo.

PROFESSIONAL PISTON COMPRESSOR NEBULISER

1	Piston compressor	4	Air filter compartment	7	Mouthpiece	10	Nosepiece
2	Power cord	5	Fasterjet Nebulizer	8	Adult face mask	11	Replacing air filter
3	ON/OFF Switch	6	Air hose	9	Child face mask	12	Assembling nebuliser kit

Dear Customer, This nebuliser is a high-quality device for inhalation therapy for asthma, chronic bronchitis, and other respiratory illnesses. This medical device generates aerosols with variable characteristics to adapt to the patient's respiratory pathology and generating a medicinal mist of particles that are small enough to reach even the deepest regions of your lungs and be of maximal benefit. The operation and usage of the device is very straightforward. All common liquid medication can be used for inhalation therapy. Please read through these instructions carefully so that you understand all functions and safety information. We want you to be happy with your Happyneb II product. For any further questions you may have, please ask your chemist or usual retailer.

IMPORTANT SAFETY INSTRUCTIONS

This device may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application. This device is not suitable for anaesthesia and lung ventilation. This device should only be used with original accessories as shown in these instructions. Do not use the device if you think it is damaged or notice anything unusual. Never open the device. This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section. Protect it from: - water and moisture - extreme temperatures - impact and dropping - contamination and dust - direct sunlight - heat and cold. Only use the medication prescribed for you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the therapy.

All accessories should only be used by one person. Never bend the nebuliser over 60°.

 Ensure that children do not use the device unsupervised some parts are small enough to be swallowed.

 Use of this device is not intended as a substitute for consultation with your physician.

PREPARATION AND USAGE OF THIS DEVICE

Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning and Disinfecting».

1 Assemble the nebuliser kit 12 Ensure that all parts are complete.
2 Fill the nebuliser with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level. 3 Connect the nebuliser with the air hose 6 to the compressor 1 and plug the power cord 2 into the socket. 4. Switch the ON/OFF switch 3 to position «I» to turn on the device and place the mouthpiece 7 in the mouth or fit one of the face masks over mouth and nose. The mouthpiece gives you a better drug delivery to the lungs. Choose between adult 8 or child face mask 9 and make sure that it encloses the mouth and nose area completely. 5. Breathe in and out calmly during the therapy. Sit in a relaxed position with the upper body upright. **Do not lie down while inhaling.** Stop inhalation if you feel unwell. 6 After completing the inhalation period recommended by your doctor, switch the ON/OFF switch 3 to position «O» to turn off the device and unplug it from the socket. 7 Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning and Disinfecting».

FASTERJET NEBULIZER

This compressed air system for aerosoltherapy is equipped with the innovative Fasterjet nebuliser with a breath enhanced valve system. This enables the correct quantity of medication to be featured to the respiratory capacity of each patient.

The valve controls the flow of medication during inhalation, thus reducing medication loss during exhalation.

Fasterjet increases the quantity of inhaled medication.

Fasterjet is capable of nebulising efficiently all medications for inhalation therapy.

A synchronized functioning of the nebuliser is possible only by using the mouthpiece with expiratory valve. In order to increase nebulising speed, it is recommended to remove the inhalation valve (white support) when using a mask or a nose fork.

CLEANING AND DISINFECTING

Thoroughly clean all components to remove medication residuals and possible impurities after each treatment. The compressor 1 and the air hose 6 should be cleaned with a clean, moist cloth.

-  Always wash your hands well before cleaning and disinfecting the accessories.
-  Do not expose the compressor to water or heat.
-  Replace the air hose for each treatment with a new patient or in case of impurities.

 Unplug the device before cleaning.

Cleaning with water

Wash all nebuliser components (except air hose) under warm tap water (max 60 °C) for about 5 minutes adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer. Rinse thoroughly making sure that all deposits are washed away and leave to dry.

Disinfecting

All nebuliser components (except air hose) can be disinfected with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer. Disinfectants are usually available at pharmacies.

Sterilizing with steam

All nebuliser components (except air hose and masks) can be heat steam sterilized up to 121 °C (20 min.) or 134 °C (7 min.). EN554/ISO11134. The sterilization packaging must conform to EN868/ISO11607 and be suitable for steam sterilization. After sterilization always let all components cool down to ambient temperature before further use. Do not repeat sterilization cycle when components are still warm.

MAINTENANCE, CARE, AND SERVICE

Order all spare parts from your dealer or pharmacist, or contact Microlife-Service (see foreword).

It is advisable replace nebulizer after some 100 to 120 treatments on single patient or after about 20 sterilization cycles. Check the filter continually for cleanliness and replace it if dirty, or after a maximum of 3 months use. Spare filters are provided with the device. To replace the filter, open the air filter compartment 4 marked «Filter» and change the filter 11

MALFUNCTIONS AND ACTIONS TO TAKE

The device cannot be switched on

Ensure the power cord 2 is correctly plugged into the socket. Ensure the ON/OFF switch 3 is in the position «I». The device has automatically switched off due to overheating. Wait until it has cooled down and try again.

The nebuliser functions poorly or not at all

Ensure the air hose 6 is correctly connected at both ends. Ensure the air hose is not squashed, bent, dirty or blocked. If necessary, replace with a new one. Ensure the nebuliser 5 is correctly assembled 12. Ensure the required medication has been added.

GUARANTEE

This device is covered by a 3 year guarantee from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) conforming date of purchase or the receipt. The guarantee covers only the compressor. The replaceable components like nebuliser, masks, mouth-piece, air hose, and filters are not included. Opening or altering the device invalidates the guarantee. The guarantee does not cover damage caused by improper handling, accidents or non-compliance with the operating instructions.

TECHNICAL DATA

Piston electrocompressor with thermal cut

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WARRANTY CERTIFICATE

Validity 36 months from date of purchase

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