



GIMA

PROFESSIONAL MEDICAL PRODUCTS

ALTERNATING PRESSURE MATTRESS

REF 28578



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China



All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

INSTRUCTIONS

The Alternating Pressure Mattress consists of two parts: the mattress and the pump. The pump utilizes a small compressor to remain quiet and energy efficient. The control panel is simple and very user friendly. The mattress relieves pressure by sequentially deflating and inflating alternate air cells on approximately 12 minutes timed intervals. It is widely recognized that constant pressure to a bony prominence is the leading cause of skin breakdown. The continuous movement supplied by the unit lessens the areas of constant pressure and enhances circulation.

OPERATION














1. Place the mattress on the bed frame with the hose end at the foot section of the bed frame.
2. Using the integrated hooks, securely hang the pump on the bed end at the foot end or place on a smooth flat surface.
3. Connect the air hoses from the mattress to the pump.
4. Plug the pump into a wall outlet. Be sure the power cord is safely away from possible hazards.
5. Turn on the power switch on the control panel of the pump. The pump will begin to inflate the mattress.
6. After inflation, adjust the mattress using the dial on the pump.

CAUTION

1. Do not smoke on or near the pump.
2. Keep the pump away from heated surfaces.
3. Explosion risk if used in the presence of flammable anesthetics.
4. Replace fuse as marked: T1A 250V. Power Input: ~ 220 V - 50 Hz - 0,1 A

PRODUCT SPECIFICATIONS

WEIGHT CAPACITY Maximum patient weight:	200 kg
ELECTRICAL Operating power: Fuses:	~ 220 V 50 Hz 0,1 A max. T1A 250 V
ENVIRONMENTAL CONDITIONS OPERATING CONDITIONS Ambient Temperature: Relative Humidity:	 10°C ~ 40°C (50°F ~ 104°F) 10% ~75%
STORAGE AND SHIPPING CONDITIONS Ambient Temperature: Relative Humidity:	-18°C ~ +43°C (0°F ~ 110°F) 10% ~ 95%

Simboli - Symbols					
	IT Conservare in luogo fresco ed asciutto GB Keep in a cool, dry place		IT Conservare al riparo dalla luce solare GB Keep away from sunlight		IT Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB Caution: read instructions (warnings) carefully
	IT Seguire le istruzioni per l'uso GB Follow instructions for use		IT Parte applicata di tipo BF GB Type BF applied part		IT Apparecchio di classe II GB Class II applied
	IT Smaltimento RAEE GB WEEE disposal		IT Codice prodotto GB Product code		IT Fabbricante GB Manufacturer
	IT Dispositivo medico conforme al regolamento (UE) 2017/745 GB Medical Device compliant with Regulation (EU) 2017/745		IT Numero di serie GB Serial number		IT Data di fabbricazione GB Date of manufacture
	IT Dispositivo medico GB Medical Device				