



# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**TERMOFORO (CON COPERTURA)**  
**HEATING PAD (WITHOUT COVER)**  
**COUSSIN CHAUFFANT AVEC COUVERTURE**  
**ALMOFADA DE AQUECIMENTO COM FRONHA**  
**HEIZKISSEN MIT HÜLLE**  
**ΘΕΡΜΟΦΟΡΑ ΜΕ ΚΑΛΥΜΜΑ**  
**PODUSZKA ELEKTRYCZNA Z POSZEWKĄ**  
**VYHŘÍVACÍ DEČKA S KRYTEM**  
**SĀCULETĪ ÎNCĂLZITOR CU HUSĂ**  
**VÄRMEDYNA MED FODRAL**  
**WARMTEMAT MET AFDEKKING**  
**LÄMPÖTYYNYSUOJAN KANSSA**

**REF** SS03N(42\*63)S (GIMA 28671)

Fabbricante / Manufacturer / Fabricant / Fabricante / Hersteller / Κατασκευαστής  
/ Producent / Výrobce / Producent / Producent / Tillverkare / Valmistaja  
**Dongguan Sunbright Electric Appliance Co., Ltd.**  
No.102, Shangdong Yiheng East Road, Qishi Town,  
523455 Dongguan City, Guangdong Province, China  
Made in China

Importato da / Imported by / Importé par / Importado por / Importiert von / Εισάγεται από  
/ Importowane przez / Importováno / Importat de / Importerad av / Geïmporteerd door /  
Maahantuoja:  
**Gima S.p.A.**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)



## Instruction Manual SS03N(42\*63)S

**IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!**  
Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.  
Should you give this device to another person, it is vital that you also pass on these instructions for use.

This instruction manual belongs to this device.  
It contains important information about starting up and operation.  
Read the instruction manual thoroughly.  
Non-observance of these instructions can result in serious injury or damage to the device



### WARNING

These warning notes must be observed to prevent any injury to the user.



### CAUTION

These notes must be observed to prevent any damage to the device.



### NOTE

These notes give you useful additional information on the installation or operation.



Class II applied



Product complies with European Directive



Follow instructions for use



Product code



Caution: read instructions (warnings) carefully



Lot number



WEEE disposal



Do not puncture the heated cushion!



Do not use the heated cushion when it is folded!



Not to be used by very young children (0-3 years)!



Only use the heated cushion indoors!



Machine wash 30°



Do not chlorinate!



Do not dry the heated cushion in a tumble drier!

The cushion cover may not be put in a dryer!



Do not iron the heated cushion!

You may iron the cushion cover at a low temperature!



Do not dry clean!

### SAFETY INFORMATION

- Examine the appliance frequently for signs of wear or damage. If there are such signs, if the appliance has been misused or does not work, return it to the supplier before switching it on again.
- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Do not fold or crease the heat padding used.
- Neither pierce the pad with, nor attach safety pins or other spiky or sharp objects to the pad.
- Do not use the device on children or persons who are disabled, asleep or insensitive to heat (persons who cannot react on overheating).
- Children below the age of 3 years are not allowed to use this device, because they are not able to react to overheating.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.
- Never use the pad unsupervised.
- The appliance can be used by children older than 3 years and younger than 8 years under supervision and with the control always set to minimum temperature value.
- The cable and control unit of the appliance can give rise to risks of entanglement, strangulation, tripping or treading if not correctly arranged. The user shall make sure that excess cables and electric cords shall be arranged in a safe way.
- The device is only for private use and is not to be used in hospitals.
- Do not fall asleep while the heat pad is switched on.
- Prolonged use of the heated cushion at a high setting could lead to skin burns.
- Do not cover the pad with any other cushions. The control switch should not be

covered or placed on or below the pad while the unit is in operation.

- Do not sit on the heat pad. Instead, rest the pad on or against the part of the body to be. The intended use of the is to only back and neck.
- Never touch a heat pad that has fallen into water. Unplug the unit from the mains outlet immediately.
- Keep the mains lead away from hot surfaces.
- Never carry, pull or turn the heat pad by the mains lead and do not let the lead get tangled.
- Do not use the heat pad when wet and only use it in dry environments (not in the bathroom or similar environments).
- The controls and leads must not be exposed to any kind of moisture.
- If a fault occurs, do not attempt to repair the unit yourself. Repairs must only be carried out by an authorized specialist dealer or other suitably qualified personnel.
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
- A damaged power supply cable must be replaced with a power cable of the same type.
- Heat pad may only be operated with the control unit (SSO3) supplied.
- When storing the heat pad, allow it to cool down before folding.
- Should you have any concerns about health issues, consult your doctor before using the heat pad.
- Do not use the heat pad on parts of the body that are swollen, injured or inflamed.
- Should you experience persistent pains in the muscles and joints, please inform your doctor. Persistent pains can be a symptom of serious illness.
- If you experience any pain or discomfort while using the unit, stop using it immediately.
- Do not crease the appliance by placing items on top of it during storage.
- Do not use if wet
- When the appliance is likely to be used for an extended period, for instance by the user falling asleep, controls should be adjusted to position 1.

## OPERATION

Plug the unit into the mains outlet and move the slider control from position 0 to position 1. The power LED will light up. After a few minutes you will be able to feel the pad getting noticeably warmer. If you require a higher temperature setting, move the switch to the next positions 2, 3, 4, 5, up to the maximum heat setting, position 6. If the pad gets too hot for you, move the slider back to position 2 or 1. After a short while the heating pad cools down noticeably. To switch the unit off, move the switch back to position 0. The power LED will go out, indicating that the unit is switched off. Unplug the unit from the mains outlet. The heat pad will switch off automatically after approx. 90 minutes of continuous use. To switch it back on again, move the slider control to position 0 and then back to the desired heat setting. Unplug the unit from the mains outlet once you have finished using it.

## CLEANING AND MAINTENANCE

### Heat pad

- Before cleaning the pad, unplug it from the mains outlet and let it cool down for at least ten minutes.
- The heating pad is equipped with a detachable cable. Disconnect the plug and remove the cable from the heating pad.
- The heating pad can be machine washed at 30° in a delicate programme or by hand gently. The best way to do this is to put the heated.
- Rinse the heating pad several times to remove all the detergent.
- Store the unit in a clean and dry place, laid out flat with nothing covering or resting on top of it.
- Only use the appliance once it is completely dry

### GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies



**Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.**