# **MANICO FLYER FLYER HANDLE MANCHE PORTE-LUMIÈRE FLYER MANGO FLYER HANDLICHTQUELLE FLYER**

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**ATTENTION:** The operators must carefully read and completely understand the present manual before using the product.

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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**REF** 819-139-20 (GIMA 30794)  
819-138-01 (GIMA 30795)

Medical Devices (Pvt) Ltd  
Wazirabad Road, Ugoki,  
Sialkot - Pakistan  
Made in Pakistan

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## Product

- 30794 flyer handle(LED)
- 30795 flyer handle (halogen)

A convenient handle that provides focused and bright light. The system is portable and easy to use.

Flyer Handle contains two AAA batteries connected in series and gives 3.0 volt DC thus eliminating the need for adapters. This is the perfect solution for hectic everyday clinical and hospital life.

## Features

**Handle:** Handle contains two AAA size batteries.

**Bottom Cap:** Situated at the bottom of each handle. Rotate clockwise to turn ON and Counter clockwise to turn OFF

**Connector:** Located on the top of each handle and contains LED or halogen Lamp. Connect your instrument here.

## LED Specifications

**LED Voltage:** 2.5V

**LED Power:** 0.33 W

**Lux:** 60000

## Bulb Specifications

**Voltage:** 2.5V


**Current:** 0.7 amp

Code 31478 bulb (for handle 30795).

## Assembling

First of all connect part No. 1 with part No. 2 then insert batteries in the handle and connect part No.3 with part No.2. Now the handle is ready for use.

**Note:** Insert the batteries in the right polarity as shown.

## Replacement of batteries

Replace 2 batteries at the same time. New batteries should not be mixed with partially exhausted ones. Different batteries or different brands should not be mixed. Failure to observe these precautions may result problems and thus increase the probability of leakage.

Rechargeable batteries are not recommended. Do not short circuit the batteries terminals. Do not expose the batteries to water or fire.

## Precriptions

The product shall be used by qualified personnel only. Do not use the equipment in case it is damaged. Apply to your retailer. Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Broken LEDs can cause injuries or cut. Carefully handle and immediately replace.

## Use and Maintenance

Place the product in such a way to make sure that it is protected from dust and water to assure its hygienic conditions and protection against short circuit.

- It is not good practice to use or leave equipment exposed to extreme conditions.
- It is advantageous to remove batteries immediately from equipment which has ceased to function satisfactorily, or when not in use for a long period.
- Be sure to turn off the equipment after use.
- Store batteries in a cool, dry place and out of direct sunlight. Handle the product with great care.

Before removing the diagnostic unit, make sure that instrument is off. Any unauthorized modification/service may cause short circuit that may damage the equipment, so always contact the authorized person.

**Symbols**

	Caution: read instructions (warnings) carefully		Follow instructions for use		Medical Device compliant with Regulation (EU) 2017/745
	Keep in a cool, dry place		Keep away from sunlight		WEEE disposal
	Manufacturer		Date of manufacture		Type B applied part
	Product code		Lot number		Imported by
	Authorized representative in the European community		Medical Device		

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

**GIMA WARRANTY CONDITIONS**

The Gima 12-month standard B2B warranty applies.