GIMA

PROFESSIONAL MEDICAL PRODUCTS

CASCHETTO LED GIMA 5W - NERO
GIMA LED HEADLIGHT 5W - BLACK
LAMPE FRONTALE LED GIMA 5W - NOIR
LED STIRNLEUCHTE GIMA 5W - SCHWARZ
CASCO LED GIMA 5W - NEGRO
CAPACÊTE LED GIMA 5W - PRÊTO
ΦΩΤΙΣΜΟΣ ΚΕΦΑΛΗΣ LED GIMA 5W - ΜΑΥΡΟ
خوذة بصمام ثنائي باعث للضوء (LED) من صناعة GIMA - 5W - لون أسود

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Warnings and Cautions

This product must not be used in the presence of flammable gases.

This product should not be immersed in fluids.

Do not shine the light directly into the patient's eyes.

Ensure that the headband power cord is located in the clips provided. The headlight body should not be left covered with a cloth or similar.

Federal law restricts this device to sale or order of a physician.

1. Purpose and Features of the product

1.1 Purpose

The headlight is used for examination and surgical operation in multifarious microsurgery.

1.2 Features

- Light Source: Our LED without the heat of traditional Halogen headlight. With a color temperature of 5000 Kelvin, this light supplies bright, white, shadow-free light allowing the doctor to see tissue characteristics without distortion.
- Enjoy the flexibility of operating the headlight using the rechargeable Lithium battery, can use without the direct power source. It is one of the most comfortable lightweight headlights and suitable for all kinds of surgery.
- The headlight could be used when it is charging.
- Has low-battery indicator light.
- Adapts to any Binocular Loupes of GIMA.
- Smart design, convenient for mobile diagnosis.

2. Operating Instructions

Attach the power cord by means of the clip provided. Place the instrument on your head and adjust for width and height to that it fits comfortably.

Connect the Headlight to the power adapter. Switch it on and use the headband right or left brightness control knob to select the required brightness.

Loosen the fixing screw to adjust the light head up and down and backwards or forwards. After adjustment, re-tighten the screw. The illuminated field of view can be selected by means of the slide control.

3. Maintenance

3.1 Battery charging and replacing

During charging the LEDs operate sequentially, in use they remain switched on. If the LEDs remain switched off, the unit is non-operational or the battery is switched off. The number of illuminated or sequentially-switched LEDs indicates the state of charge of the Li-ion battery. The state of charge is shown whether an instrument is connected or not.

All four LEDs illuminated indicated a fully-charged battery,

If only the yellow LED illuminated, the battery must be recharged. Switch off the battery and recharge it by means of the transformer.

Caution: Please use the appropriate battery. The battery should never be placed in municipal waste. Please check local regulations for disposal of batteries.

3.2 Battery use and care

The capacity of Lithium battery will be changed by it's using condition. For maximum battery capacity, use the battery in room temperature (10 °C ~ 40 °C). Do not expose the battery to extreme temperatures, if the battery is used in high or low temperatures, the battery capacity will be reduced.

The old batteries should never be placed in municipal waste.

3.3 Cleaning and sterilization

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.

Caution: Be careful to prevent the solvent into the instrument.

4. Technical Specification

Optical specifications


420mm illumination light spot	20-80mm
Illumination in 200mm working distance	≥ 60000 lx
Adjustable range:	
Fluctuation distance	12.5mm
Angle range Vertical	±25°
Weight: Headlight (excluding head band)	69g

Electrical specifications

Illumination Source	Extremely bright light, 5W LED, life 20000 hours continuous burn, do not need to change the LED lamp.
Battery	7.4 V (4400 mAh) Lithium battery
Battery Run-time	Not less than 7 hours continuous runtime on a fully charged battery.
Battery Charging-time	Not more than 5 hours
Rating Adapter Voltage	100VAC-240VAC 50Hz/60Hz
Maximum Adapter Input Power	≤ 24 VA
Power Adapter Output	12V DC 2A
Electrical Safety	EN60601-1, EN60601-1-2 / CLASS II / Internal electrical power source

Environment requirements

Operation:	
Environment Temperature	+10°C - +40°C
Relative Humidity	30% - 75%
Atmospheric Pressure	700 hPa - 1060 hPa
Storage:	
Environment Temperature	-40°C - +55°C
Relative Humidity	10% - 90%
Atmospheric Pressure	500 hPa - 1060 hPa

	Caution: read instructions (warnings) carefully		Follow instructions for use
	Medical Device compliant with Regulation (EU) 2017/745		Date of manufacture
	Keep away from sunlight		Keep in a cool, dry place
	Product code		Lot number
	WEEE disposal		Class II applied
	Internal electrical power source		Medical Device
	Type BF applied part		

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.