VNYUS ED610

VIDEONYSTAGMOSCOPE SYSTEM

OPERATING MANUAL

Rev 06 – 26/05/2021

CONTENTS

<u>MANUFACTURER IDENTIFICATION</u>	<u>4</u>
<u>PRODUCT IDENTIFICATION.....</u>	<u>5</u>
<u>SYMBOL MEANING</u>	<u>6</u>
<u>CLASSIFICATION</u>	<u>7</u>
<u>DISPOSAL.....</u>	<u>7</u>
<u>CAUTIONS</u>	<u>8</u>
<u>USE DESTINATION.....</u>	<u>9</u>
EQUIPMENT PLACEMENT	9
ENVIROMENTAL COMPLIANCE	9
OPERATOR COMPLIANCE	10
KIND OF FUNCTIONNING	10
FUNCTIONING FAULTS	10
COUNTER EFFECTS	10
WARNING, POSSIBLE UNDESIRED EFFECTS	10
<u>DESCRIPTION OF INSTALLATION AND/OR ASSEMBLING PROCEDURE.....</u>	<u>11</u>
ENVIRONMENT COMPLIANCE	11
ACTIONS TO BE EXECUTED BEFORE SWITCHING THE SYSTEM ON.....	11
FIRST ACTIONS TO BE EXECUTED.....	12
EMERGENCY SIGNAL AND WARNING.....	12
<u>PERIODICAL START UP</u>	<u>13</u>
<u>MANUFACTURER RESPONSIBILITIES</u>	<u>14</u>
<u>USE INDICATIONS</u>	<u>14</u>

GENERAL DESCRIPTION.....	14
<u>CONTROLS DESCRIPTION.....</u>	<u>15</u>
CONTROL BOX	15
<u>CLEANING ED610</u>	<u>16</u>
<u>TECHNICAL FEATURERS</u>	<u>17</u>
<u>FAULTS RESEARCH</u>	<u>17</u>
THE EQUIPMENT DOES NOT START	17
NO IMAGES ON THE MONITOR.....	17
<u>GUIDE CHARTS AND CONSTRUCTOR DECLARATIONS.....</u>	<u>18</u>
<u>TECHNICAL ASSISTANCE SERVICE.....</u>	<u>22</u>

MANUFACTURER IDENTIFICATION

MEDI-CARE SOLUTIONS S.r.l.

Manufacturer:

Via della Zecca nr. 1 40121 Bologna (BO) Italy

Facility:

Via Pietro Nenni nr. 3 40026 Imola (BO) Italy

P.IVA/C.F. 03390661209

Phone ++39 0542 642046

Fax ++39 0542 642355

Web site: MEDI-CARE SOLUTIONS.it

e-mail: service@MEDI-CARE SOLUTIONS.it

PRODUCT IDENTIFICATION

1. Control Box

2. B/W IR Camera

3. USB type C - USB type C cable

4. Cable

5. Stopper

6. Cable Rca/Bnc

SYMBOL MEANING

1.3 Symbols

Classification of the "BF" type machinery. "BF" type machinery (machinery of the type I,II, III or machinery housing an internal electric source with a protection degree against direct and indirect contact and in particular for:
permitted current dispersion
reliability of the ground contact)
Class II machinery BF.

This symbol indicates: TO CLARIFICATION CONSULT THE ATTACHED DOCUMENTATION

This symbol indicates: ATTENTION! CONSULT THE ATTACHED DOCUMENTATION BEFORE USE THE UNIT

Symbol that indicate "generic warning".

Earthling contact of the machinery

Alternated current

Symbol discharge with normative

Symbol that indicate hot surface

Medical Device

CLASSIFICATION

CLASS II type machinery B:

DISPOSAL

The symbol of crossed out garbage bin (on the back of the equipment) indicates that the product at the end of its working life must be collected separately from other garbage. This refers to the following law:

art. 13 of D.L. 25 /07/2005, n.151: "Implementation of the regulations 2002/95/CE, 2002/96/CE and 2003/108/CE, about decrease of dangerous substances on electric and electronic equipments" and moreover it refers to garbage disposal.

The separate collection of rubbish of this equipment (felt into disuse) must be arranged and coordinated by the manufacturer. When the user wants to throw away the equipment, he will contact the Manufacturer Company that will indicate all procedures to follow for the disused equipment separate collection of rubbish.

Therefore, the appropriate separate collection of rubbish, aimed to actuate afterwards disused equipment recycling, treatment and to ambient garbage disposal, contributes to avoid possible bad effects on the environment and on the health. In addition, these procedures raise the reemployment and/or the recycling of the materials of which the equipment is composed.

The unauthorized garbage disposal, on the part of the holder, requires administrative law sanctions application.

CAUTIONS

1. Do not install the machinery in environments which are exposed to explosions or frost
2. The temperature of the environment in which the equipment is installed and used must be between 10 ° C and 30 ° C (max 40 ° C for storage).
3. The relative humidity environment of use or storage must be understood between 30% and 70%; pressure 500hpa at 1060hpa
4. During installing, check that the electrical tension on the sticker on the machine corresponds to that in use in order to avoid damaging.
5. The building, repairing, modifying, calibration and all the other operations that require the opening of the unit should be carried out by personnel authorized by MEDI-CARE SOLUTIONS.
6. The electric circuit in which the unit is connected to, should comply to the IEC 64.8 Sec. 710 regulation (regulation regarding electrical circuits in rooms used for medical purposes).
7. The machinery should be used exclusively following the instructions illustrated in this manual.
8. MEDI-CARE SOLUTIONS refuses any responsibility for the following: damages or accidents caused by the carelessness of the user, the non compliance to the safety rules, a non correct use of the machine, a use of the machine which is different to that of its purpose.
9. MEDI-CARE SOLUTIONS denies all responsibility which has nothing to do with faults in the material and non correct building.
10. Before installation and functioning check that there are not loose parts, burnt electric cables, or various anomalies and if they are evident please contact MEDI-CARE SOLUTIONS or authorized technical staff.
11. Mobile phones and radios can influence the good functioning of medical device
12. If during the use of the MD happen strange changes of performance or various anomalies we suggest to interrupt all operations, switch off the machine and contact authorised staff MEDI-CARE SOLUTIONS in order to avoid any kind of damage.
13. Do not install this unit in places exposed to the risk of explosions, rich in oxygen or exposed to frost
14. The DM is not protected against splashing water and foreign bodies: DEGREE IP20
15. WARNING!: To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
16. WARNING!: There shall be no modification of this device
17. WARNING! The DM in question should not be used adjacent to or stacked with other equipment, if necessary use it close to or stacked with other equipment it is necessary to observe the DM to check whether, in the configuration in which it is used, its operation remains normal
18. WARNING! In the event of a serious accident occurring in relation to the device, it is necessary to report to the Manufacturer and the competent authority of your country
19. The MD can be transported and stored with the environment :
Temperature: da -10°C a 60°C;
Humidity: 10% - 90%;
Pressure = by 500hpa at 1060hpa

USE DESTINATION

VNYUS ED610 Videonystagmoscope system facilitates view and registration of Nystagmus reached by standard vestibular tests (both position and instrumental tests). This device allows to view patient eyes images on the monitor during the tests.

It 's not the intended use of the device in the surgery room

EQUIPMENT PLACEMENT

- Please, put the equipment on an horizontal surface that must be sufficiently large, solid and stable.
- Connect the power-supply cable to the plug;

Connect the USB-C cable only to a USB-C port of your PC/Laptop

ENVIROMENTAL COMPLIANCE

The environment in which the equipment is used must be hospital-type exempt from strong electric interferences and therefore far from high tension conductors, from large transformers, from x-ray implants, from diathermy and Marconi therapy generators.

Further inconvenience elements are elevator's contactors, automatic gear-cases, some lacks of the electric implant, such as conductors not perfectly fixed to the lamp-holders, to the current plugs, to the switches, etc.

In normal conditions the MD is capable of functioning without particular precautions. In critical conditions, that is when the inconveniences menace to compromise the performance of the good performance of the test, it is necessary to connect all equipment to a single efficient ground outlet.

The installation of the equipment must be made respecting the instructions of the following booklet of instructions.

The user should always make sure that the electric circuit complies to IEC values (in other words: 220/240 Vac, alternate monophase 50 Hz, the plant should comply with IEC 64.8 sec. 710 regulation).

It is advisable to insert the plug in a socket with a bipolar switch compliant to IEC and having the following technical characteristics:

230 Vac - 50 A contact opening of 3 mm.

It's not the intended use of the device in the surgery room

OPERATOR COMPLIANCE

The operator, in charge of use VNYUS ED610 Videonystagmoscope system, must have read and understood the instruction manual, therefore he must be able to recognize a potential malfunctioning and capable of following the procedure recommended in the manual.

The operator , in charge of use ED610 must be an ENT professional.

KIND OF FUNCTIONNING

VNYUS ED610 Videonystagmoscope system is an equipment with a continuous kind of working.

FUNCTIONING FAULTS

Not being an explosion-proof system, it is necessary the user not to use it near flammable anaesthetic.

Do not let VNYUS ED610 videonystagmoscope system let out to atmospheric factors. Do not let it unnecessarily turned on if nobody is using it.

Functioning faults are represented by specific situations when AD9240 mask is not preferably usable and, as a result, the device must be in case used with a defined limit, with regard to benefits/costs ratio that resting with the system in point.

The DM is not protected against splashing water and foreign bodies: DEGREE IP20

COUNTER EFFECTS

Counter effects are those specific situations in which VNYUS ED610 is not preferably usable and, as a result, the device must be in case used with a defined limit, with regard to benefits/costs ratio linked to system use.

WARNING, POSSIBLE UNDESIRE EFFECTS

Before starting to use VNYUS ED610 videonystagmoscope system, please, do not forget to make always a careful medical history to identify situations where it is advised against mask use for Nystagmus research.

MEDI-CARE SOLUTIONS S.R.L.**DESCRIPTION OF INSTALLATION AND/OR ASSEMBLING
PROCEDURE****ENVIRONMENT COMPLIANCE**

The environment in which the equipment is used must be hospital-type exempt from strong electric interferences and therefore far from high tension conductors, from large transformers, from x-ray implants, from diathermy and Marconi therapy generators.

Further inconvenience elements are elevator's contactors, automatic gear-cases, some lacks of the electric implant, such as conductors not perfectly fixed to the lamp-holders, to the current plugs, to the switches, etc.

In normal conditions the MD is capable of functioning without particular precautions. In critical conditions, that is when the inconveniences menace to compromise the performance of the good performance of the test, it is necessary to connect all equipment to a single efficient ground outlet.

The installation of the equipment must be made respecting the instructions of the following booklet of instructions.

The user should always make sure that the electric circuit complies to IEC values (in other words: 220/240 Vac, alternate monophasic 50 Hz, the plant should comply with IEC 64.8 sec. 710 regulation).

It is advisable to insert the plug in a socket with a bipolar switch compliant to IEC and having the following technical characteristics:

230 Vac - 50 A contact opening of 3 mm.

ACTIONS TO BE EXECUTED BEFORE SWITCHING THE SYSTEM ON

Push away all non-used cables and those near the system, they might cause encumbrance.

Furthermore, not being an explosion-proof system, it is necessary not to use it near flammable anaesthetic. Check all system device components integrity.

FIRST ACTIONS TO BE EXECUTED

1. After connecting USB-C plug into the USB-C socket, insert USB connector into the proper socket of PC/NOTEBOOK. Then, install DRIVER of camera using the CD about installation/video-acquisition.
2. Run VIEW/ACQUISITION PROGRAM; select the device; check green LED placed on the panel to be turned on: it means that software can interact with the camera.
3. If you are using EDM_Videocapture software, the camera will be automatically recognized.
4. Pressing the button I/O (placed on the system panel), check the lighting of red LED on the camera higher part.
5. Check that the USB cable is connected to the Power Supply if you use Analog Monitor

**EMERGENCY SIGNAL AND WARNING**

The green LED (placed on the fore panel) allows to view functioning condition about the two B/W Ir camera system. After connecting the camera, green LED lighting informs that camera drivers have been recognized by HARDWARE and by the used Software.

PERIODICAL START UP

The user must pay attention to the actions that are necessary to keep the device in good conditions, in order to make a "periodical start up" it in total efficiency:

- Check at sight the integrity of device and of its components (control box, mask, video camera, cables)
- place the device in an environment with the following features:

Temperature from	20° to 45°C
Relative humidity from	30 to 75%
Atmospheric pressure from	700 to 1060 hPa

- After checking disconnection of USB connector, take away the dust from the control box using an humid cloth.
- in case of prolonged inactivity, it is preferable to cover VNYUS ED610 videonystagmoscope together with its accessories (mask, camera, etc.) in a hiding place, after providing to a simple removal of the possible dust.
- Check to have the equipped power supplier of system electricity connection
- Require periodical checks (every years) c/o MEDI-CARE SOLUTIONS S.p.A..

MANUFACTURER RESPONSIBILITIES

VNYUS ED610 Videonystagmoscope system, manufactured by Medi-Care Solutions srl, is guaranteed against potential manufacture and/or material defects for one year from the purchasing date and according to the following conditions:

- The room interior wiring must comply with the current IEC regulations.
- The equipment must be used following the instructions and regulations as in the present manual.
- Installation, modifications and repairs must be executed by technicians authorized by Medi-Care Solutions srl Company.

The following voices are not under warranty:

- Defects caused by normal deterioration
- Parts subject to wearing out
- Malfunction caused by damaging (accidental causes), by maintenance or modifications executed by non-authorized people.
- Defects caused by the machine's use for different purposes from those it's been built for.

USE INDICATIONS

GENERAL DESCRIPTION

This manual describes VNYUS ED610 videonystagmoscope system functioning.

The equipment is composed of a small control box, a mask together with a camera, a USB-mini cable and a stopper.

After camera driver installation on PC/Notebook, using a visualization software, it will be possible to view and register the images reached by the camera.



VNYUS ED610 videonystagmoscope system allows to start the red led inside the mask during tests, in order to make Visual Suppression Test using the button called I/O, placed on the control box.

If you use the VIDEONYSTAGMOSCOPE ED610 connected to an Analog Monitor you must be connect the USB-mini-USB cable to Power Supply. All you need is to connect the device analog output to the monitor using the video cable.

MEDI-CARE SOLUTIONS S.R.L.

CONTROLS DESCRIPTION

CONTROL BOX

CLEANING ED610

The choice of materials used for building the unit has been carefully taken in order to correspond to the cleaning and disinfecting needs. Given the big differences between chemical products used in a medical study, it is possible to ruin the surfaces. The major damages are caused by the amount of time that these chemical products are left on the tops. It is important to immediately dry the interested parts with a damp cloth. To clean the painted or metal parts, use a specific product for these materials which doesn't damage their composition.

IMPORTANT! Do not use products with alcohol, ammonia or petrol

DO NOT USE SOLVENTS OR ALCOHOL FOR CLEANING AND DISINFECTING. THE USE OF THESE PRODUCTS CAN DAMAGE THE PLASTIC PIECES OF THE SYSTEM

Cleaning and disinfecting products chart	
Disinfectants	Pharmaceutic products and cosmetics
+ bactol 5%	+ blood plasma
+ clorammina	+ sun tan lotion
- DDT	+ hydroplex
+ delegol 5%	= iodine dye
= dimammina 5%	+ lanolyn
= iodine dye	= mantolo 90% in alcohol
- fenic acid	-nail varnish
+ lysoform 2%	- solvent
- TB lysoform	+ mouth wash
+ bactol	+ periston R
+ merfen 2%	+ vaselyn
+ octozone	+ vicks rub
+ perydrol	
+ resorcina 1%	
= sagrotane 5%	
+ alcohol to be purely burnt	
+ mercury chloride	
+ trosilina G	
+ oxiginated water	
= zefirol	
CHART SIGNS	
+ compatible	
- non compatible	
= compatible but with caution	

TECHNICAL FEATUTURES

CONTROL BOX

Electrical tension:	5V DC (USB port)
Power absorbed:	max 2,5W
Weight	115 gr
Dimensions Max(in mm):	115 x80 x 23

MASK

Weight	115 gr
Dimensions Max(in mm):	180 x90 x 45

CAMERA

Weight	200 gr (with cable)
Dimensions Max(in mm):	56 x55 x 45

FAULTS RESEARCH

It is possible the equipment cannot reach the desired function.

Before calling assistance technician, please, try to solve the problem alone according to the following checking list, that could help to identify and eliminate faults.

THE EQUIPMENT DOES NOT START

- Make sure that USB cable plug is inserted on the USB socket of your PC/NOTEBOOK and check the socket to be active.
- Make sure to have opportunely installed camera drivers
- Call MEDI-CARE SOLUTIONS technical assistance service

NO IMAGES ON THE MONITOR

- Make sure that USB cable plug is inserted on the USB socket of your PC/NOTEBOOK and check the socket to be active.
- Make sure to have opportunely installed camera drivers
- Verify to have selected the device on the choices of your software

MEDI-CARE SOLUTIONS S.R.L.

- Call MEDI-CARE SOLUTIONS technical assistance service

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TECHNICAL ASSISTANCE SERVICE**OF****MEDI-CARE SOLUTIONS SRL****TEL. / PHONE: +39 0542 642046****FAX: +39 0542 642355****E-MAIL: info@MEDI-CARE SOLUTIONS.it****!!!!!! VERY IMPORTANT ADVICE !!!!!!**

Before sending back the equipment, even if it is covered by its warranty, it is necessary to forewarn and require MEDI-CARE SOLUTIONS S.p.A. authorisation. Otherwise, the warranty could fall into decay.

In case, due to any reason, the equipment must be delivered to the technical assistance service, we suggest to use the original packing, in order to avoid eventual damages due to transport.

Therefore Medi-Care Solutions srl declines any responsibility in relation to these damages and will charge the new packing to be used for sending again the equipment to the customer.