**MONOFILAMENTO
MONOFILAMENT
MONOFILAMENT
EINZELFADEN
MONOFILAMENTO
MONO-FILAMENTO
MONOMHMATIO**

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where you registered office is located.

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
ENGLISH

MONOFILAMENT

Show the filament to the patient and touch it to his/her hand or arm so that he/she knows it does not hurt.

Use the 10 gram filament to test sensation at the indicated sites on each foot as shown. Apply the filament along the perimeter of, and NOT directly on, an ulcer, callus, scar or necrotic tissue.

The monofilament is an instrument for a preliminary diagnosis of sensory neuropathy which can be caused, for example, by diabetes. The foot is one of the parts of the body which are affected by the damage caused by neuropathy as the peripheral sensory perception is compromised; the monofilament can detect a patient's ability to perceive moderate pressures.

The MONOFILAMENT bends at a pressure of 10 grams for at least 50/60 times before bending to a lower force; it is recommended to change the wire after the use on 6/15 patients and allowed to rest for at least 24 hours before being reused. It can work optimally at a temperature of 20/26°C

with humidity 30/50% without deforming. With temperature higher than 33°C and with humidity at 80%, the filament can be deformed and make the diagnosis less accurate. The filament was tested to determine the accuracy of bending to 10 grams and the deformation force has been recognized at +/- 10%. Substitute definitely the filament after 150/200 bendings.

Hold the filament by exercising a moderate pressure perpendicular to the skin and use a smooth motion when testing.

Hold the filament perpendicular to the skin and use a smooth motion when testing. Use a 3-step sequence that includes (1) touch the skin, (2) bend the filament, and (3) lift from the skin (See Figures 1-3). Do not use rapid movement. The approach, skin contact, and departure of the filament should be approximately 1.5 seconds in duration.

Ask the patient to respond "yes" when the filament is felt. If the patient does not respond when you touch a given point on the foot, continue on to another site.

When you have completed the sequence, REPEAT on the area(s) where the patient did not indicate feeling the filament.

Use the filament in a random sequence.

LOSS OF PROTECTIVE SENSATION AT ANY ONE OF THE EIGHT SITES (shown on the diagram above) INDICATES A FOOT AT HIGH RISK.

It is recommended to disinfect the product after each use on a patient by using a non-alcoholic antibacterial disinfectant.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

REF	GB - Product code FR - Code produit IT - Codice prodotto DE - Erzeugniscode ES - Código producto PT - Código produto GR - Καθώς προϊόντος PL - Numer katalogowy SE - Produktkod RO - Cod produs HU - Termékkód SA - كود المنتج	GB Medical Device compliant with Regulation (EU) 2017/745 FR Dispositif médical conforme au règlement (UE) 2017/745 IT Dispositivo medico conforme al regolamento (UE) 2017/745 DE Medizinprodukt im Sinne der Verordnung (EU) 2017/745 ES Producto sanitario conforme con el reglamento (UE) 2017/745 PT Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 GR Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 PL - Wyrob medyczny zgodny z Rozporządzeniem (UE) 2017/745 SE - Den medicintekniska produkten överensstämmer med förordning (EU) 2017/745 RO - Dispozitiv medical realizat în conformitate cu prevederile regulamentului (UE) 2017/745 HU - A 2017/745/EU rendeletheek megfelelő orvostechnikai eszköz SA - جهاز طبي يتوافق مع التوجيه (UE) 2017/745
LOT	GB - Lot number FR - Numéro de lot IT - Numero di lotto DE - Chargennummer ES - Número de lote PT - Número de lote GR - Αριθμός παρτίδας PL - Kod partii SE - Satsnummer RO - Număr de lot HU - Tételszám SA - رقم الدفعة	
	GB - Keep in a cool, dry place FR - À conserver dans un endroit frais et sec IT - Conservare in luogo fresco ed asciutto DE - An einem kühlen und trockenen Ort lagern ES - Conservar en un lugar fresco y seco PT - Armazenar em local fresco e seco GR - Διατηρείται σε δροσερό και στεγνό περιβάλλον PL - Przechowywać w suchym miejscu SE - Förvara på svalt och torrt ställe RO - A se păstra într-un loc răcoos și uscat HU - Száraz, hűvös helyen tárolandó SA - يحفظ في مكان بارد وجاف	
	GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil IT Conservare al riparo dalla luce solare DE Vor Sonneneinstrahlung geschützt lagern ES Conservar al amparo de la luz solar PT Guardar ao abrigo da luz solar GR Κρατήστε το μακριά από ηλιακή ακτινοβολία PL Przechowywać z dala od światła słonecznego SE Skyddas från solljus RO A se păstra ferit de razele soarelui HU Napfénytől védve tárolandó SA - يحفظ بعيدا عن أشعة الشمس	
	GB - Manufacturer FR - Fabricant IT - Fabricante DE - Hersteller ES - Fabricante PT - Fabricante GR - Παράγωγος PL - Producent SE - Tillverkare RO - Producător HU - Gyártó SA - الشركة المصنعة	MD
		SA - جهاز طبي GB - Consult instructions for use FR - Consulter les instructions d'utilisation IT - Leggere le istruzioni per l'uso DE - Gebrauchsanweisung beachten ES - Consultar las instrucciones de uso PT - Consulte as instruções de uso GR - Διαβάστε προσεκτικά τις οδηγίες χρήσης PL - Przeczytaj instrukcję użytkowania SE - Läs bruksanvisningen RO - Citii instrucțiunile de utilizare HU - Olvassa el a használati utasításokat SA - اقرأ بصفة وحرص تعليمات الاستخدام

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Gima S.p.A.