PROFESSIONAL MEDICAL PRODUCTS

LUCE DI WORTH**WORTH LIGHTTEST****TEST DE WORTH****WORTH-TEST****TEST DE WORTH****TESTE DE WORTH****ΤΕΣΤ ΦΩΤΟΣ ΤΗΣ WORTH****WORTH** ثروء عضو فحص

Manuale utente - User manual - Notice d'utilisation - Gebrauchsanweisung

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Worth light Test is used to assess a patient's flat Fusional ability.

Flat Fusional Testing is Indicated any time stereopsis falls between (50) Sec of Arc, on those patients with suspected strabismus and on preschool Children. The Worth Test should also be used when evaluating cases of reduced monocular visual acuity that does not improve with the Pinhole Test.

Procedure

1. With the best refraction correction worn by the patient, place the Anaglyph glasses over the patient's correction, with red filter over the right eye.
2. In a slightly dimmed room, turn on the flashlight and hold the Worth Test with the red dot orientated up at approximately 16" from the patient and slightly below the LOS.
3. Conduct the following Monocular Check first:
 - a. Cover the right eye; ask how many dots the patient sees. They should report 3 green.
 - b. Then Cover the left eye and ask how many dots they see now. They should report 2 red.
4. Next, Conduct the binocular test:
 - a. With both eyes uncovered ask a third time how many dots they see. If normal flat fusion is present they should report 4.

Pediatric Note: This test can reliably be conducted on preschoolers as young as 2 if they are allowed to point to the dots.

5. Abnormal responses:
 - a. If the patient reports only 2 red dots under binocular conditions, this indicates that they are suppressing the left eye.
 - b. If the patient reports seeing 3 green dots under binocular conditions, then they are suppressing the right eye.
 - c. If they report 5 dots, they are diplopic. The type of diplopia can then be determined by asking which side are the green dots.
If the green dots are located on the right, the patient has a eso deviation; to the left, an exo deviation.
 - d. If the green dots are reported above or below the red dots then a vertical deviation exists. A report of the green dots above the red dots would be seen with aright hyper deviation.
 - e. With a report of 6 or more dots, one should question the patient's reliability.
6. The binocular view is repeated at 5 and 10 feet, then repeat all these distances under greatly reduced light as both suppression and ocular deviation can be different under varying lighting conditions.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

	<p>IT Fabbricante GB Manufacturer FR Fabricant DE Hersteller ES Fabricante PT Fabricante GR Παραγωγός</p> <p style="text-align: right;">SA الشركة المصنعة</p>
	<p>IT Data fabbricazione GB Date of manufacture FR Date de fabrication DE Herstellungsdatum ES Fecha de fabricación PT Data de fabrico GR Ημερομηνία παραγωγής</p> <p style="text-align: right;">SA تاريخ صنع المنتج</p>
	<p>IT Codice prodotto GB Product code FR Code produit DE Erzeugniscode ES Código producto PT Código produto GR Κωδικός προϊόντος</p> <p style="text-align: right;">SA كود المنتج</p>
	<p>IT Numero di lotto GB Lot number FR Numéro de lot DE Chargennummer ES Número de lote PT Número de lote GR Αριθμός παρτίδας</p> <p style="text-align: right;">SA رقم الدفعة</p>
	<p>IT Conservare in luogo fresco ed asciutto GB Keep in a cool, dry place FR À conserver dans un endroit frais et sec ES Conservar en un lugar fresco y seco DE An einem kühlen und trockenen Ort lagern PT Armazenar em local fresco e seco GR Διατηρείται σε δροσερό και στεγνό περιβάλλον</p> <p style="text-align: right;">SA يحفظ في مكان بارد وجاف</p>
	<p>IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil DE Vor Sonneneinstrahlung geschützt lagern ES Conservar al amparo de la luz solar PT Guardar ao abrigo da luz solar GR Κρατήστε το μακριά από ηλιακή ακτινοβολία</p> <p style="text-align: right;">SA يحفظ بعيدا عن أشعة الشمس</p>
	<p>IT Seguire le istruzioni per l'uso GB Follow instructions for use FR Suivez les instructions d'utilisation DE Folgen Sie den Anweisungen ES Siga las instrucciones de uso PT Siga as instruções de uso GR Ακολουθήστε τις οδηγίες χρήσης</p> <p style="text-align: right;">SA استخدم المنتج كما هو مكتوب</p>
	<p>IT Dispositivo medico conforme al regolamento (UE) 2017/745 GB Medical Device compliant with Regulation (EU) 2017/745 FR Dispositif médical conforme au règlement (UE) 2017/745 DE Medizinprodukt im Sinne der Verordnung (EU) 2017/745 ES Producto sanitario conforme con el reglamento (UE) 2017/745 PT Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 GR Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 (UE) 2017/745 مع التوافق مع توجيه (UE) 2017/745</p> <p style="text-align: right;">SA جهاز طبي يتوافق مع التوجيه (UE) 2017/745</p>
	<p>IT Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB Caution: read instructions (warnings) carefully FR Attention: lisez attentivement les instructions (avertissements) DE Achtung: Anweisungen (Warnungen) sorgfältig lesen ES Precaución: lea las instrucciones (advertencias) cuidadosamente PT Cuidado: leia as instruções (avisos) cuidadosamente GR Προσοχή: διαβάστε προσεκτικά τις οδηγίες (εισατάσεις)</p> <p style="text-align: right;">SA التحذير: قراءة التعليمات (التحذيرات) بعناية</p>
	<p>IT Smaltimento RAEE GB WEEE disposal FR Disposition DEEE DE Beseitigung WEEE ES Disposición WEEE PT Disposição REEE GR Διάθεση WEEE</p> <p style="text-align: right;">SA التخلص WEEE</p>
	<p>IT Dispositivo medico GB Medical Device FR Dispositif médical DE Medizinprodukt ES Producto sanitario PT Dispositivo médico GR Ιατροτεχνολογικό προϊόν</p> <p style="text-align: right;">SA جهاز طبي</p>
	<p>IT Importato da GB Imported by FR Importé par DE Importiert von ES Importado por PT Importado por GR Εισάγεται από</p> <p style="text-align: right;">SA داري تاسا</p>
	<p>IT Rappresentante autorizzato nella Comunità europea GB Authorized representative in the European community FR Représentant autorisé dans la Communauté européenne ES Representante autorizado en la Comunidad Europea PT Representante autorizado na União Europeia GR Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα</p> <p style="text-align: right;">SA مُمثل مفوض في المجتمع الأوروبي</p>

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**Gima S.p.A.**