# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**OTOSCOPIO F.O. SIGMA**  
**SIGMA F.O. OTOSCOPE**  
**OTOSCOPE F.O. SIGMA**  
**FASEROPTIK-OTOSKOP SIGMA**  
**OTOSCOPIO F.O. SIGMA**  
**OTOSCÓPIO F.O. SIGMA**  
**OTOSKOP SIGMA F.O.**  
**SIGMA F.O. OTOSCOPI**  
**ΩΤΟΣΚΟΠΙΟ ΟΠΤΙΚΗΣ ΙΝΑΣ SIGMA**

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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



Lens Locking Point

**Fig-6**

Insufflation Bulb connecting with Adaptor.

**Fig-7**

Insufflation Bulb.

## **Instructions For Use**

The operating and maintenance instructions found in this manual should be followed to ensure years of reliable service.

Please read these instructions carefully before using the Sigma and keep them in a safe place for future or use reference.

### **Intended Use**

Sigma Oscopes and associated accessories are intended for the examination of ear, the external auditory canal and the eardrum.

This instrument has an optical magnifying lens and battery-powered illumination handle, and recommended for use by medically qualified personnel to assess the health of ear through transmitting light into the ear canal.

This device is intended for transit treatment and examination.

## **Setting Up the Product**

### **Head & Handle Attachment**

To setup the Sigma Oscope, align the otoscope thread with handle and rotate the head clockwise until head and handle firmly joined.(Fig-2, A ).

### **Sigma Battery Handle**

The Sigma battery handle should only be used as a power source for Sigma instruments.

Grip the handle attached to Oscope head. To turn ON the light, slide the switch with thumb towards the bottom of the handle.

To turn OFF the light, slide the switch with thumb towards the top of the handle. (Fig-2, B).

### **Insert Batteries**

Arrange two AA/LR6 size good quality alkaline batteries for power. Rotate end cap counterclockwise to open the handle end cap. Insert two batteries into handle with positive (+) terminal towards the head and close the cap by rotating clockwise (Fig-2, C).

**Note:** Batteries are not included in sets.

### **Handle with rechargeable system:**

- The handle with rechargeable system has all the characteristics of a standard Sigma.
- A metal cap (cod. 31588) converts standard handles into rechargeable ones.
- The metal cap can be used both with “AA” batteries and the rechargeable battery.

For grips with rechargeable system, a 2.5 V lithium-ion battery is used (code 31587), which is included in the set.

### **Note:**

The metal cap can be purchased separately, thereby converting a standard handle into a rechargeable one.

The Oscope is intended for transit treatment and the handle may gets warm by longer use.

So let it be cool down before next application.

A harmless max temperature of 60° C will not be exceeded.

## **Sigma xenon otoscope model**

### **Bulb Replacement**

Rotate otoscope head counterclockwise to detach from handle.

Pull down the bulb and insert new bulb into socket and rotate the head clockwise until head and handle firmly joined (Fig-3).

Before replacing the bulb read below instruction carefully.

Let the device cool down before changing the bulb.

Verify the bulb voltage / specification before changing.

### Operating

Sigma Oscope should only be inserted into the auditory canal when a ear tip is fitted. Place the ear tip on the otoscope in such a way the projection inside the ear tip fits into the slit (Fig- 4,1) and tighten the ear tip by twisting it slightly to the right.

- Do not use otoscope without ear tip that may cause ear injury and / or cross contamination.
- Do not reuse disposable ear tip that may result in contamination from one patient to another patient.
- Do not point light in patient eyes that may cause eye injury.
- Do not use excess force while locking the ear tip on head.

### Magnifying Lens

Sigma Oscope has 3x magnification that produces distortion-free images of the examined object. The magnification viewing window is permanently attached and can be swivelled aside for instrumentation. To facilitate the instrumentation rotate the viewing window clockwise as much as required (Fig-5).

- Peel off the lens protective sheet before use (Fig-5).

### *Sigma LED otoscope model*

The LED fixed and sealed into the otoscope head and cannot be changed because of its long life.

LED has 50,000 hours of extended lifetime.

Ensure no replacement of LED ever in the lifetime of the instrument.

Excellent light output of 50K LUX.

Extended battery life of 50 HOURS due to low current consumption of LED.

With a colour temperature of 4000K and CRI>80 colours are visible as they are.

### *Sigma-C led otoscopes*

- Light weight alternative to the standard Desk top direct light Oscope for easy examination of ear and nose, and small procedures.
- Made from impact resistant re-enforced plastic and highly corrosion resistant stainless steel material for exterior finish.
- Swivel type window with 3X magnification lens for sharp image and minimal reflection.
- Insufflation port for pneumatic testing of tympanic mobility.
- 50,000 hours of extended lifetime ensures no replacement of LED ever in the lifetime of the instrument.
- Brighter than Xenon model. 7500 LUX as compared to 4000 on the Standard model Oscopes.
- Extended battery life of 50 hours due to low current consumption of the LED.
- With a color temperature of 4300K, CRI > 95 and R9 = 88 ensures that colors are visible as they are.

### Ear Tip

Ear tip is used in otoscope as viewing passage for direct examination of the external acoustic duct and tympanic membrane.

Choose the size of Ear tip suitable for the examination of the acoustic duct.

### Insufflation Port (optional)

Sigma Oscope have a connection port to facilitate the tympanic mobility test (Fig-7).

This can be carried out when the viewing window is closed by using optionally available insufflation bulb and adaptor.

- Dose the pressure carefully with the insufflation bulb.

### *Cleaning / Sterilization and Maintenance*

The Sigma F.O. Pocket otoscope is not suitable for treatment in the ultrasonic bath, for sterilization, spray disinfection or treatment by immersion in liquids.

### *Manual Cleaning Process*

Use a damp cloth for cleaning the Sigma pocket otoscope.

For cleaning and disinfection, a disinfectant agent suitable for plastic medical products must be used.

## Handle Cleaning

The handle can be cleaned with a cloth soaked with disinfectant, however it cannot be dipped into liquids.

Do not use the equipment in case it is damaged. Apply to your retailer.

Periodically check the battery conditions, making sure that no sign of corrosion of oxidation is present. In case of necessity replace them with new alkaline ones.

Carefully handle the batteries as the liquids they contain can irritate skin and eyes.

Before being used, thoroughly check the product. The same operation shall be carried out after cleaning.

Check that the connection between the head and the handle is perfect and that the On/Off button works correctly.

If the light is intermittent or in case it does not turn on, Check the bulb, the batteries and the electrical contacts.

### Precautions

Read the label before losing packing. Do not use if packaging has been open or damaged. Only qualified person should use it.

Handle the device with care and keep the instrument away from the range of insects and rodents.

Autoclaving is not permissible for this type of handle.

Ultrasonic cleaning is not recommended.

Store the device in clean environment at normal temperature.

Please ensure that the batteries, Handles and Ear Tips and other Sub-assemblies are discarded as clinical waste in accordance with local policies.

In accordance with local regulations this product should be disposed of as an electronic device separately.

## Transport and Storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks bumps or falls during the transport itself.

Damage to the appliance caused during transport and handling is not covered by the guarantee.

The device must be stored in a dry, Cool area away from direct sunlight.

It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

### Operation

Temperature 10° C - 40° C

Humidity 30% - 75%

Air Pressure 700 hPa - 1060 hPa

### Storage & Transport

Temperature -20° C - 70° C

Humidity 10% - 90% (without condensation)

Air Pressure 500 hPa - 1060 hPa

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

## GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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