GIMA

PROFESSIONAL MEDICAL PRODUCTS

BRACCIALE PER INFUSIONE
INFUSION CUFF
BRASSARD DE PERFUSION A PRESSION
INFUSIONSMANSCHETTE
BRAZAL PARA INFUSIÓN
MANGUEIRA DE INFUSÃO
ΕΓΧΥΜΑΤΙΣΜΕΝΟ ΜΑΝΣΑΚΙ
تسريب الكفة

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Infusion Cuff Set

This device is applicable to be used as the ancillary of infusion in clinic, and it is recommended to be used together with the infusion bag. Comparing with the traditional infusion method with pressure provided by gravity, this device can avoid the contamination of the ambient air as no ambient air will contact with the liquid medicine. In the meantime, it makes the paramedic more convenient as the paramedic don't need to hold the infusion set at given level accompanying the patient, especially when the patients go for eating, outdoor activity or toilet.

The device consists of infusion cuff, inflation bulb, manometer and knot. Please connect these parts before use.

The instruction of using this instrument

This device can provide necessary pressure for the infusion bag so as to ensure the liquid medicine flows from the infusion bag to the patients' body.

- Put the infusion bag into the infusion cuff, and the infusion bag will be just below the transparent net. Put the infusion cuff in safe area to avoid extrusion or falling down, fix the infusion bag by the pressure of infusion cuff;
- Pump the inflation bulb to inflate air into the inflation cuff. Observing the manometer carefully to make sure that the pressure reaches the average of the patient's systolic pressure and diastolic pressure;
- Connect the disposable infusion set with the infusion bag, and the liquid medicine will flow into the whole disposable infusion set. Make sure that there is not bubble in it, and close the infusion switch. Then it is ready for other ordinary infusion operation;
- During the infusion, the pressure will decrease gradually. Please pump the inflation bulb to increase the pressure. If the infusion speed is too low or too high, please modify the pressure in the infusion cuff by pumping the bulb to inflate or releasing deflation valve to deflate so as to get expected infusion speed;

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

| | | | |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>IT Codice prodotto GB Product code FR Code produit DE Erzeugniscode ES Código producto PT Código produto GR Κωδικός προϊόντος</p> <p>SA كود المنتج</p> |  | <p>IT Dispositivo medico GB Medical Device FR Dispositif médical DE Medizinprodukt ES Producto sanitario PT Dispositivo médico GR Ιατροτεχνολογικό προϊόν</p> <p>SA جهاز طبي</p> |
|  | <p>IT Numero di lotto GB Lot number FR Numéro de lot DE Chargennummer ES Número de lote PT Número de lote GR Αριθμός παρτίδας</p> <p>SA رقم الدفعة</p> |  | <p>IT Leggere le istruzioni per l'uso GB Consult instructions for use FR Consulter les instructions d'utilisation DE Gebrauchsanweisung beachten ES Consultar las instrucciones de uso PT Consulte as instruções de uso GR Διαβάστε προσεχτικά τις οδηγίες χρήσης</p> <p>SA اقرأ بدقة وحرص تعليمات الاستخدام</p> |
|  | <p>IT Conservare in luogo fresco ed asciutto GB Keep in a cool, dry place FR À conserver dans un endroit frais et sec DE An einem kühlen und trockenen Ort lagern ES Conservar en un lugar fresco y seco PT Armazenar em local fresco e seco GR Διατηρείται σε δροσερό και στεγνό περιβάλλον</p> <p>SA يحفظ في مكان بارد وجاف</p> |  | <p>IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil DE Vor Sonneneinstrahlung geschützt lagern ES Conservar al amparo de la luz solar PT Guardar ao abrigo da luz solar GR Κρατήστε το μακριά από ηλιακή ακτινοβολία</p> <p>SA يحفظ بعيداً عن ضوء الشمس</p> |
|  | <p>IT Fabricante GB Manufacturer FR Fabricant DE Hersteller ES Fabricante PT Fabricante GR Παραγωγός</p> <p>SA الشركة المصنعة</p> |  | <p>IT Data di fabbricazione GB Date of manufacture FR Date de fabrication DE Herstellungsdatum ES Fecha de fabricación PT Data de fabrico GR Ημερομηνία παραγωγής</p> <p>SA تاريخ التصنيع</p> |
|  | <p>IT Dispositivo medico conforme al regolamento (UE) 2017/745 GB Medical Device compliant with Regulation (EU) 2017/745 FR Dispositif médical conforme au règlement (UE) 2017/745 DE Medizinprodukt im Sinne der Verordnung (EU) 2017/745 ES Producto sanitario conforme con el reglamento (UE) 2017/745 PT Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 GR Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 SA جهاز طبي يتوافق مع التوجيه (UE) 2017/745</p> |  | <p>IT Contenuto o presenza di lattice naturale GB Contains or presence of natural rubber latex FR Contient ou présence de latex de caoutchouc naturel DE Enthält oder Anwesenheit von Naturlatex ES Contiene o presencia de látex de caucho natural PT Contém ou presença de látex de borracha natural GR Περιέχει ή περιέχει λατέξ φυσικού καουτσούκ</p> <p>SA يحتوي على أو وجود مطاط طبيعي من المطاط</p> |

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Gima S.p.A.