

ERKA.

THE ORIGINAL

Made in Germany since 1889

USER MANUAL CUFFS



360° SERVICE
globale Kompetenz
aus einer Hand

Sie haben Fragen zu einem unserer Produkte,
zu dessen Anwendung oder Instandhaltung?

Von A wie Auftragserteilung bis Z wie Zertifiziert
nach DIN EN ISO 13485: unser kompetentes
ERKA.-Serviceteam berät und unterstützt Sie gerne
persönlich unter der kostenfreien Rufnummer
+8000-CALL-ERKA (008000-2255-3752).

Get in touch with us! International free call
+8000-CALL-ERKA (008000-2255-3752)* oder
per E-mail an service@erka.org

*Unsere kostenlose Telefonberatung ist für Sie
werktags von 07:00 bis 19:00 Uhr erreichbar

360° SERVICE
Global expertise
from a single source

Do you have any questions about our products,
their use or maintenance?

For all your needs concerning ERKA. products and
services, our expert ERKA. Service Team is happy
to provide advice and support at the following free
phone number:
+8000-CALL-ERKA (008000-2255-3752).

Get in touch with us! International free call
+8000-CALL-ERKA (008000-2255-3752)* or
by E-mail to service@erka.org

*Our free telephone consultation service can
be reached on weekdays between 7:00 a.m. and
7:00 p.m.

EXPLANATION OF SYMBOLS | LÉGENDE |
ERKLÄRUNG DER SYMBOLE | EXPLICACIÓN DE SÍMBOLOS |
SPIEGAZIONE DEI SIMBOLI:



Please pay attention to the accompanying User Manual.
Veuillez lire le mode d'emploi ci-joint.
Bitte beiliegende Bedienungsanleitung beachten.
Por favor, tenga en cuenta el manual de instrucciones adjunto.
Si prega di osservare le allegate istruzioni per l'uso.



This product complies with the Council Directive 93/42/EEC
concerning medical devices

Ce produit est en conformité avec la directive 93/42/CEE
du Conseil pour les dispositifs médicaux

Dieses Produkt entspricht der Richtlinie 93/42/EWG
des Rates über Medizinprodukte

Este producto se corresponde con la directiva 93/42/CEE
del Consejo relativa a productos sanitarios.

Questo prodotto è conforme alla Direttiva 93/42/CEE
del Consiglio in materia di prodotti medicali.



Manufacturer
Fabricant
Hersteller
Fabricante
Produttore



Catalogue number
Numéro de catalogue
Katalognummer
Número de catálogo
Numero di catalogo



Serial number
Numéro de série
Seriennummer
Número de serie
Numero di serie



Lot number
Numéro de lot
Lot-Nummer
Número de lote
Numero di lotto



Pack of ten
Pack de dix
Zehnerpack
Paquete de 10
Pacco da dieci



Caution
Attention
Achtung
Atención
Attenzione

EXPLANATION OF SYMBOLS | LÉGENDE |
ERKLÄRUNG DER SYMBOLE | EXPLICACIÓN DE SÍMBOLOS |
SPIEGAZIONE DEI SIMBOLI:



Latex-free
Ne contient pas de latex
Latexfrei
Sin látex
Privo di lattice



Humidity threshold
Limite d'humidité de l'air
Luftfeuchte Begrenzung
Limitación de la humedad del aire
Limitazione umidità dell'aria



20,5-28 cm
Arm circumference and cuff size
Circonférence du bras et taille du brassard
Armumfang und Manschettengröße
Perímetro del brazo y tamaño del brazalete
Misura braccio e misura polsino



PVC-free
Ne contient pas de PVC
PVC-frei
Sin PVC
Privo PVC



Temperature threshold
Limite de température
Temperaturbegrenzung
Limitación de la temperatura
Limitazione della temperatura



Store in a dry place
À stocker dans un endroit sec
Trocken aufbewahren
Conservar en seco
Conservare in luogo asciutto



Fragile: Handle with care!
Fragile, à manipuler avec précaution
Zerbrechlich, mit Sorgfalt handhaben
Fragil, manipular con cuidado
Fragile, manipolare con cura

R_x only


Only sold to physicians or upon their request
La vente ne s'effectue qu'à des médecins,
resp. sur l'ordre de ses derniers
Der Verkauf erfolgt nur an Ärzte bzw. auf
deren Anordnung
La venta solo se realiza a médicos o si estos
así lo indican
La vendita avviene solo a medici o su loro
disposizione

USER MANUAL BLOOD PRESSURE CUFFS

TABLE OF CONTENTS:

1. In General
2. Product description
3. Putting on the cuff
4. Selection of the cuff size
5. Important information
6. Maintenance and care
7. Warranty

1. IN GENERAL

 Dear Customer,
With this blood pressure cuff, you have purchased a quality product from ERKA. We are delighted that you have expressed your confidence in ERKA. quality and we thank you for doing so. We have a wide range of blood pressure cuffs in various models. High quality materials that are suitable to the task guarantee optimal operation. If you have questions regarding the right choice or the operating of your blood pressure cuff, you can direct them to made@erka.org. We would be happy to help you.



Please read this user's manual before use carefully and thoroughly and follow the instructions for maintenance and keep the manual for further use. Important for home users: Consult your doctor at regular intervals. He will tell you your individual values for normal blood pressure, as well as the value at which the level of the blood pressure is considered hazardous.

2. PRODUCT DESCRIPTION

PURPOSE

ERKA. cuffs are used along with blood pressure monitors of the company ERKA. Kallmeyer Medizintechnik GmbH & Co.KG.

SIDE EFFECTS

No side effects are known so far.

BLOOD PRESSURE CUFFS

GREEN CUFF - SUPERB RAPID CUFF

- Waterproof surface for easy and thorough cleaning
- Very soft materials for maximum comfort against the skin
- Durable cover with removable bladder
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener

GREEN CUFF - SUPERB HOOK CUFF

- Waterproof surface for easy and thorough cleaning
- Very soft materials for maximum comfort against the skin
- Durable cover with removable bladder
- The hook cuff has metal pins sewn in and is closed by a clasp

GREEN CUFF - SUPERB D-RING

- Waterproof surface for easy and thorough cleaning
- Very soft materials for maximum comfort against the skin
- Durable cover with removable bladder
- Suitable particularly for application at clearly conical upper arms
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener

GREEN CUFF - SUPERB HOME CUFF

- Waterproof surface for easy and thorough cleaning
- Very soft materials for maximum comfort against the skin
- Durable cover with removable bladder
- Suitable particularly for application at clearly conical upper arms
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener

GREEN CUFF - SMART RAPID CUFF

- Waterproof surface for easy and thorough cleaning
- With the ERKA quick releasing fastener, the cuff can be changed in seconds
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener

GREEN CUFF - SMART RAPID CUFF ONE PIECE

- Waterproof surface for easy and thorough cleaning
- With the ERKA quick releasing fastener, the cuff can be changed in seconds
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener
- Cover with integrated bladder

GREEN CUFF - SMART RAPID CUFF SINGLE PATIENT SOFT

- With the ERKA quick releasing fastener, the cuff can be changed in seconds
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener
- Shelf life for the entire hospital stay of patients

GREEN CUFF - SMART RAPID CUFF SINGLE USE

- With the ERKA quick releasing fastener, the cuff can be changed in seconds
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener
- For single use

CLASSIC RAPID CUFF

- Cotton cover makes for a comfortable fit; especially suitable for sensitive patients with paper skin or aged skin
- With removable bubble
- The super soft adhesive tape and Velcro tape closes the cuff simply like a hook and loop fastener

3. PUTTING ON THE CUFF

Select appropriate cuff size for the circumference of the extremities. The cuff is put on so that the lower edge of the cuff is about 2 to 3 cm above the elbow or about 5 cm above the knee joint. Note the marking for the location of the artery when doing so. Do not wrap the cuff too loosely and not too tight around the limb and make sure that the edge of the marking is within the permissible range.

4. SELECTION OF THE CUFF SIZE

The following table provides information on the relevant area of application and meets the standards of the latest AAMI and AHA guidelines on the proper size of the blood pressure cuff:

Size 1	Infant	10.0 to 15.0 cm circumference	Free of PVC, DEHP and latex
Size 2	Toddler	14.0 to 21.5 cm circumference	Free of PVC, DEHP and latex
Size 3	Adults Small	20.5 to 28.0 cm circumference	Free of PVC, DEHP and latex
Size 4	Grown-ups	27.0 to 35.0 cm circumference	Free of PVC, DEHP and latex
Size 5	Adults Large	34.0 to 43.0 cm circumference	Free of PVC, DEHP and latex
Size 6	Thigh	42.0 to 54.0 cm circumference	Free of PVC, DEHP and latex

5. IMPORTANT INFORMATION



Check the stethoscope before each use as to its proper function. Do not use this product if it shows visible damage.

- The Medical Products Operator Ordinance is applicable for Germany
- The interpretation of the results of the measurement should only be done by qualified persons

- Do not pump up the manometer above 300 mmHg by any means
- Infants have a lower blood pressure than adults, do not pump too strong
- The measurement time should not exceed two minutes
- The cuff should be used with hygiene protective strips or cleaned before it is used for the next patient
- Only use original ERKA. spare parts
- Avoid damaging the rubber parts with sharp or pointed objects

6. MAINTENANCE AND CARE

CARE INSTRUCTIONS	Green Cuff Superb	Green Cuff Smart	Classic Cuff
	<ul style="list-style-type: none"> ● Rapid ● Hook ● D-Ring ● Home 	<ul style="list-style-type: none"> ● Rapid ● Rapid one piece 	<ul style="list-style-type: none"> ● Rapid
WASHING	After removing the bladder up to 30° C mild detergent	After removing the bladder up to 30° C mild detergent (not Green Cuff Smart Rapid one piece)	After removing the bladder up to 90° C mild detergent
Fold covers with Rapid stopper before washing, so that the Velcro tape and the hook tape are closed to each other.			

CLEANING BY SPRAY DISINFECTANT OR WIPE DISINFECTANT

We recommend the following cleaning agents and disinfectants:
 Mild Soap, Bacillol AF, Meliseptol, Oxivir Plus Spray, Terralin Protect
 70% isopropyl alcohol, Fugaten® Spray (Lysoform-Spray)

After cleaning/disinfecting, wipe with clean water to remove any residue and allow to air dry. Please only use soft cotton or microfiber cloths.

Never use autoclaving!

ENVIRONMENTAL CONDITIONS

- Optimal is a storage/use at temperatures between 0° C and +40° C and a relative humidity of less than 85%.
- Do not expose the cuff and the appliance directly to direct UV irradiation, this is particularly applicable for all plastic parts.

7. WARRANTY

The warranty is subject to the statutory provisions and is 2 years from the date of delivery.

Within this period, we will remove all the defects, which relate to the material or manufacturing defects, by repair or by replacing the product, free of cost. The warranty does not cover damage caused by wear and tear as well as damage due to improper treatment/handling, use of parts/replacement parts that are not originally from ERKA. or in case of non-compliance of the instructions.

If you encounter any questions, one of our friendly staff will be happy to help you at **8000-CALL ERKA (008000-2255 3752)** at any time.

Warranty - Made in Germany.

Opting for an **ERKA** product means you appreciate unrivalled quality, performance and durability. Take advantage of the free **ERKA. XXL warranty for an extra five years of coverage** - with no ifs, ands or buts.

Garantie - Made in Germany.

Wer sich für ein **ERKA** Produkt entscheidet sucht bewusst höchste Qualität und Langlebigkeit. Registrieren Sie Ihr Produkt unter www.erka.org/garantie und sichern Sie sich eine zusätzliche Garantiezeit von bis zu fünf Jahren - ohne Wenn und Aber.



ERKA.
XXL GARANTIE
www.erka.org/garantie

International Freecall +8000-CALL ERKA
008000-22553752



Made in
Germany

ERKA.

Kallmeyer Medizintechnik

Im Farchet 15

83646 Bad Tölz

Germany

www.erka.org

Date / Stand: 2015