BRACCIALI PER TOURNIQUET

TOURNIQUET CUFFS

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MANGUITOS PARA TORNQUETES

Manuale d'uso - User Manual - Notice d'utilisation - Manual del usuario

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto. - **ATTENTION:** The operators must carefully read and completely understand the present manual before using the product. - **AVIS:** Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit. - **ATENCIÓN:** Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.

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INSTRUCTIONS FOR USE (IFU) TOURNIQUET CUFFS

1. Product Description and Application

Tourniquet cuffs are medical accessories used in conjunction with automatic tourniquet devices. They are used to stop blood flow during surgical procedures, allowing for better visibility of the surgical site, reducing operation time, and minimizing blood loss.

2. Usage Conditions and Application

- Should only be used by trained medical personnel.
- The device must be compatible with single or dual-hose systems.
- Before use, the cuff should be tested up to 300 mmHg to check for any leaks.
- The cuff should be placed as proximally (closest to the body) as possible or at the widest diameter of the extremity.
- A gauze pad or protective sleeve should be used before applying the cuff.
- The cuff's tubes should be positioned either on the side or front.
- The cuff should be secured with a Velcro system and should not be too tight (there should be enough space for a finger to pass through).
- After securing with Velcro, a fastening strap (safety strap) should be tied over the cuff to prevent loosening.
- The cuff should overlap by at least 1/4 but should not exceed 1/2 of its width.
- To ensure venous drainage, the extremity should be elevated for at least 30 seconds before cuff inflation.
- Cuff pressure should be approximately 50-70 mmHg above the patient's systolic blood pressure, and about twice the systolic pressure for lower extremities.
- The maximum usage time is 120 minutes, with a recommended break of 20 minutes. In elderly and diabetic patients, this duration should not exceed 60 minutes.

3. Contraindications (Usage Restrictions)

- Peripheral artery disease, severe atherosclerosis.
- Sickle cell disease.
- Open fractures or areas with skin grafts.
- Severe traumatic injuries and circulatory disorders.
- Should be used with caution in elderly, diabetic, and alcohol-dependent patients.

4. Warnings and Precautions

- Improper positioning of the tourniquet cuff may result in nerve damage or paralysis.
- The tourniquet should not be twisted or folded on the skin.
- Excessively low or high pressure may cause intraoperative bleeding.
- Special caution is needed for patients at risk of vascular occlusion.
- After use, the cuff should be fully deflated before removal.

5. Sterilization and Maintenance

- Single-use cuffs should not be reused.
- Reusable cuffs should be inspected before each use, and any damaged cuffs should be replaced.

- Cuffs should not be autoclaved or sterilized at high temperatures.
- Mild disinfectant solutions should be used for cleaning, and the cuff should be thoroughly dried afterward.
- During storage, the cuff should be kept in a dry and cool environment.

6. Side Effects

- Pain, stiffness, and weakness in the extremity.
- Skin discoloration (reactive hyperemia).
- In rare cases, paralysis, vascular occlusion, or embolism.
- Severe adverse effects may occur in IVRA/Bier block procedures.

7. Model Variations

Model Code	Size (Inches)	Size (Cm) ~	Chamber Count
AVR-CUFF S008	8"L x 2"W	20 cm x 5 cm	Single Chamber
AVR-CUFF S012	12"L x 3"W	30 cm x 8 cm	Single Chamber
AVR-CUFF S015	15"L x 4"W	38 cm x 10 cm	Single Chamber
AVR-CUFF S018	18"L x 4"W	46 cm x 10 cm	Single Chamber
AVR-CUFF S024	24"L x 4"W	61 cm x 10 cm	Single Chamber
AVR-CUFF S030	30"L x 5"W	76 cm x 13 cm	Single Chamber
AVR-CUFF S034	34"L x 5"W	86 cm x 13 cm	Single Chamber
AVR-CUFF S038	38"L x 5"W	97 cm x 13 cm	Single Chamber
AVR-CUFF S044	44"L x 6"W	112 cm x 15 cm	Single Chamber
AVR-CUFF 2XD014	14"L x 4"W	36 cm x 10 cm	Dual Chamber
AVR-CUFF 2XD020	20"L x 6"W	51 cm x 15 cm	Dual Chamber
AVR-CUFF 2XD026	26"L x 6"W	66 cm x 15 cm	Dual Chamber

For more information, please contact your authorized distributor or manufacturer.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

**Indice dei simboli - Symbol index - Index des symboles
- Índice de símbolos**

	IT - Data di fabbricazione GB - Date of manufacture FR - Date de fabrication ES - Fecha de fabricación
	IT - Fabbricante GB - Manufacturer FR - Fabricant ES - Fabricante
	IT - Dispositivo medico GB - Medical Device FR - Dispositif médical ES - Producto sanitario
	IT - Numero di lotto GB - Lot number FR - Numéro de lot ES - Número de lote
	IT - Codice prodotto GB - Product code FR - Code produit ES - Código producto
	IT - Dispositivo medico conforme al regolamento (UE) 2017/745 GB - Medical Device compliant with Regulation (EU) 2017/745 FR - Dispositif médical conforme au règlement (UE) 2017/745 ES - Producto sanitario conforme con el reglamento (UE) 2017/745
	IT - Identificatore univoco del dispositivo GB - Unique device identifier FR - Identifiant unique de l'appareil ES - Identificador de dispositivo único PT - Identificador exclusivo do dispositivo
	IT - Leggere le istruzioni per l'uso GB - Consult instructions for use FR - Consulter les instructions d'utilisation ES - Consultar las instrucciones de uso
	IT - Attenzione - Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB - Caution - read instructions (warnings) carefully FR - Attention - lisez attentivement les instructions (avertissements) ES - Precaución - lea las instrucciones (advertencias) cuidadosamente
	IT - Distribuito da GB - Distributed by FR - Distribué par ES - Distribuido por