

IT - ELETTRODI ECG - ISTRUZIONI PER L'USO
EN - ECG ELECTRODES - INSTRUCTIONS FOR USE
FR - ELECTRODES ECG - MODE D'EMPLOI
DE - EKG-ELEKTRODEN - BEDIENUNGSANLEITUNG
ES - ELECTRODOS DE ECG - INSTRUCCIONES DE USO
PT - ELETRODOS DE ECG - INSTRUÇÕES DE UTILIZAÇÃO
RU - ЭЛЕКТРОДЫ ДЛЯ ЭКГ - ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
BG - ЕКГ ЕЛЕКТРОДИ - ИНСТРУКЦИИ ЗА УПОТРЕБА
CS - EKG ELEKTRODY - NAVOD K POUZITI
DA - EKG-ELEKTRODER - BRUGERVEJLEDNING
EL - ΗΛΕΚΤΡΟΔΙΑ ΗΚΓ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ET - EKG ELEKTROODID - KASUTUSJUHEND
FI - EKG-ELEKTRODIT - KÄYTTÖOHJEET
HR - EKG ELEKTRODE - UPUTE ZA UPOTREBU
HU - EKG ELEKTRODAK - HASZNÁLATI UTMUTATÓ
LV - EKG ELEKTRODI - LIETOSANAS PAMACĪBA
LT - EKG ELEKTRODAI - NAUDOJIMO INSTRUKCIJA
NL - ECG ELEKTRODEN - GEBRUIKSAANWIJZING
PL - ELEKTRODY EKG - INSTRUKCJA UZYCIA
RO - ELECTROZI ECG - INSTRUCȚIUNI DE UTILIZARE
SK - EKG ELEKTRODY - NAVOD NA POUZITIE
SL - EKG ELEKTRODE - NAVODILA ZA UPORABO
SV - EKG-ELEKTRODER - BRUKSANVISNING
TR - EKG ELEKTROTLARI - KULLANIM TALIMATLARI



Via Costoli, 4 – 50039 - Vicchio
(Florence, Italy) – www.fiab.it

52502408IU3B
2021-05



IT - ELETTRODI ECG - ISTRUZIONI PER L'USO
EN - ECG ELECTRODES - INSTRUCTIONS FOR USE
FR - ELECTRODES ECG - MODE D'EMPLOI
DE - EKG-ELEKTRODEN - BEDIENUNGSANLEITUNG
ES - ELECTRODOS DE ECG - INSTRUCCIONES DE USO
PT - ELETRODOS DE ECG - INSTRUÇÕES DE UTILIZAÇÃO
RU - ЭЛЕКТРОДЫ ДЛЯ ЭКГ - ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
BG - ЕКГ ЕЛЕКТРОДИ - ИНСТРУКЦИИ ЗА УПОТРЕБА
CS - EKG ELEKTRODY - NAVOD K POUZITI
DA - EKG-ELEKTRODER - BRUGERVEJLEDNING
EL - ΗΛΕΚΤΡΟΔΙΑ ΗΚΓ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ET - EKG ELEKTROODID - KASUTUSJUHEND
FI - EKG-ELEKTRODIT - KÄYTTÖOHJEET
HR - EKG ELEKTRODE - UPUTE ZA UPOTREBU
HU - EKG ELEKTRODAK - HASZNALATI UTMUTATÓ
LV - EKG ELEKTRODI - LIETOSANAS PAMĀCĪBA
LT - EKG ELEKTRODAI - NAUDOJIMO INSTRUKCIJA
NL - ECG ELEKTRODEN - GEBRUIKSAANWIJZING
PL - ELEKTRODY EKG - INSTRUKCJA UZYCIA
RO - ELECTROZI ECG - INSTRUC IUNI DE UTILIZARE
SK - EKG ELEKTRODY - NAVOD NA POUZITIE
SL - EKG ELEKTRODE - NAVODILA ZA UPORABO
SV - EKG-ELEKTRODER - BRUKSANVISNING
TR - EKG ELEKTROTLARI - KULLANIM TALIMATLARI

Rx ONLY - CAUTION:
U.S. FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR
ON THE ORDER OF A PHYSICIAN

- The equipment and connection cables to be used with these products must comply with the regulations in force.
- Do not use the products on patients with previous episodes of skin irritation or contact allergies.
- The electrodes comply with the biocompatibility standard ISO 19993-1

STORAGE

The product must be stored in the original package at the environmental conditions (temperature and relative humidity) specified on the pouch's label. Putting external heavy weights on the package, the product could be damaged.

GENERAL NOTES

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

WASTE DISPOSAL

Refuses deriving from health structures must be disposed in according to the regulation in force.



Via Costoli, 4 – 50039 - Vicchio
(Florence, Italy) – www.fiab.it

52502408IU3B
2021-05



EN INTENDED USE

The electrode is applied to the surface of the body to transmit the electrical signal at the body surface to a recording and monitoring device that produces an electrocardiogram or vectorcardiogram.

CONTRA INDICATIONS

None known.

DIRECTIONS FOR USE

- Check that the package, the electrodes and the gel are intact; in case of any damage do not use the product.
- If needed, clean application areas with a mild soap and water and dry thoroughly.
- For adults only: Clip off hair on excessively hairy areas if necessary. Only gently abrade if skin impedance is too high.
- The correct placement of electrodes is typically determined by the hospital/physician and/or local protocol/procedures.
- Peel the electrode away from the protective liner and place it on the skin. Dispose of the protective liner immediately after removal of electrode.
- Place the electrode on the skin by gently pressing around the adhesive edge with a circular smooth and firm motion.
- Connect the leadwire from the ECG monitoring device to the electrode. If the electrode is supplied with a pre-attached leadwire it may be connected directly to the ECG monitoring device.
- Electrode Removal: Gently lift the tab or the edge of the electrode, maintaining a low profile, and slowly pull off the electrode from the skin, holding and supporting down the newly exposed skin as the electrode is removed. Keep the electrode close to the skin surface when pulling back over itself. Avoid removing the electrode at 45-degree angle, to reduce the risk of mechanical trauma.
- Gel residue can be removed using paper towels or water.
- Close the bag carefully so that the remaining electrodes do not dry out.
- As a general precept, replacement of the electrodes is recommended after 48 hours and mandatory after 72 hours.

WARNINGS

- Surface electrodes should be applied only to intact, clean skin (e.g., not over open wounds, lesions, infected, or inflamed areas).
- For single use only. Reuse of the product may cause contamination, leading to infections.
- After removal of electrodes from the protective liner, the protective liner should be disposed of in assigned waste bin.
- During surgical procedures electrodes should be placed as far as possible from the electro-surgical area to minimize unwanted RF current flow which could lead to skin burns.
- Electrodes should only be used by or in consultation with a health care provider familiar with their proper placement and use.
- Do not soak, rinse, or sterilize the electrodes as these procedures may leave harmful residues or cause malfunction of the electrodes.
- To minimize the risk of skin irritation, when replacing electrodes find a new area of skin on which to place them.

