



Figure 13 Delete selection interface

(4) Denote value

Select "Denote Value" in Data Management interface to enter its sub-menu as shown in Figure 14, after selecting the parameter, it will automatically return to Data Management interface.



Figure 14 Denote value setting interface

Note: when GLI or SBPT is selected, there is no PEF option in denote value setting interface.

(5)Exit

In Data Management interface, select "Exit" or press RETURN to return to Menu interface.

c.Settings

Select "Settings" in Menu interface to enter the setting interface as shown in Figure 15, where the language, time, and calibration can be set, and device information can be viewed.




Figure 15 Settings interface

(1) Language

Select "Language" in Settings interface, then press UP or DOWN key to select "中文", "English", "Español", "Português", "Italiano", "Deutsch", "Français" or "pyck". (this operation is invalid if the device does not have the built-in language selection function.)

(2) Time setting

Select "Time" to enter its setting interface, select "Year" to display current year as shown in Figure 16, press UP or DOWN key to change the value, after selecting, press CONFIRM key to save.

The operation steps of "Month", "Day", "Hour", "Minute" and "Second" are the same to the "Year".




Figure 16 Time setting interface

(3) Calibration

Select "Calibration" in Settings interface to enter its sub-menu as shown in Figure 17, 2L and 3L are optional, after selecting, it will enter the calibration interface as shown in Figure 18.



Figure 17 Calibration selection interface




Figure 18 Calibration interface

Under Calibration interface, push the syringe once, the device will display "Please repeat", then push the syringe once again. After continuous three correct operations, the calibration is succeed, and the device will display "OK!". Finally the interface will jump to the former interface before calibration (The former interface: if calibrating after measuring, it will return to Settings interface; if calibrating before measuring, it will return to Testing interface.).

If the device displays "Error!", it indicates something wrong with the operation or the syringe selects improper volume, please confirm that the calibration volume is correct, then repeat calibrating until succeeding. If you need to stop calibrating, just press the CONFIRM key to exit to the interface before calibrating.

Select "Adjust" in Calibration interface to display the current calibration value as shown in Figure 19. Press UP or DOWN key to change the value, press CONFIRM key to save.

Note:

| Number | Model | Cable length (m) | Mask or no | Remark |
|--------|---------------------|------------------|------------|--------|
| 1 | Power adapter cable | 1.0 | YES | / |

The value determines the accuracy of measurement, please do NOT change it randomly.

After replacing the turbine, calibration shall be applied for inputting parameters of new turbine, which guarantees the accuracy of measurement after replacing.

When replacing the turbine, please use the one recommended by our company.

Improper calibration may affect the measurement accuracy, please be careful.



Figure 19 Calibration adjustment interface

In Calibration selection interface, select "Exit" or press RETURN to return to Settings interface.

(4) About

Select "About" in Settings interface to enter its sub-menu to check the device name and software version, then press CONFIRM or RETURN key to return to Settings interface.

(5) Exit

In Settings interface, select "Exit" or press RETURN to return to Menu interface.

d.Power off

Select "Power Off" in Menu interface to turn off the device.

Note: If there is no operation within two minutes, the device will shut down automatically.

e.Exit

In Menu interface, select "Exit" or press RETURN to return to Main interface, if the measurement is not completed before entering Main interface, it will return to Testing interface.

4.2.5 Repeated measure

The device has the function of repeated measurement, long press CONFIRM key for 2 seconds to enter Testing interface, when the memory is full, the information "The memory is full!" will display on the screen, shown as Figure 20, select "Yes" to enter data delete interface, select "No" to enter Menu interface.



Figure 20 Memory full interface

4.2.6 Charging

The device will automatically enter the charging interface when it is charging. Under this interface, all keys are unfunctional, and the device can't be used.

Two methods for charging:

- Charge the device by connecting to a computer via USB cable.
- Charge the device by connecting to the power adapter.

Do not use the device during charging.

During charging, the message "Charging..." displays on the interface, the battery icon is a lighting symbol, and the indicator light is orange. It is green after fully charged.

When charging, do not position the device so that it is difficult to operate the disconnection device. After charging, remove the power adapter and disconnect the device from the mains supply.

4.2.7 Data transmission

1) Install PC software into a computer, after that, connect the device with the computer by the equipped USB cable, open the software and turn on the device, then data transmission is available.

2) The device has Bluetooth transmission function. After power on, the Bluetooth is always ON, which can be searched and connected. After the connection is established, the device can communicate.

4.3 Attention

- Please check the device before using to confirm that it can work normally.
- Automatic power off when there is no operation in two minutes.
- It is power supplied by rechargeable lithium battery.
- It is recommended that the device should be measured in room.
- Excessive ambient light may affect measurement accuracy. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight, etc.
- Intense activity of the subject or electrosurgical interference may also affect the accuracy.
- Please clean and disinfect the device after using according to the User Manual (7.1).
- Please use the USB cable recommended by our company if it is necessary to replace the USB cable.

Chapter 5 Maintenance, Transportation and Storage

5.1 Cleaning and disinfection

Use medical alcohol to wipe the device enclosure, nature dry or clean it with a clean and soft cloth. It's necessary to clean the turbine periodically for accuracy, keep the diaphaneity of the lucency part, and keep it away from sundries(such as hair or lesser sediment). Immerse the turbine in disinfectant after use, after a few minutes, clean it with clean water and air dry (but don't make the turbine rinsed with water directly), this disinfection method will not bring pollution to environment. (Note: The disinfectant is 75% alcohol).

5.2 Maintenance

1) Please clean and disinfect the device before using according to the User Manual(5.1).

2) Please charge the device when the screen displays low voltage(the battery power is ).

3) Charge the battery in time after it is fully discharged. If the device is not used for a long time, it should be charged every 6 months, which could greatly extend the battery service life. Users are forbidden to replace the battery by themselves, if necessary, place contact the local service center or our company.

4) The device needs to be calibrated once a year(or according to the calibrating program of hospital). It can be performed at the state-appointed agent or just contact us for calibration.

5.3 Transportation and storage

1) The packed device can be transported by ordinary conveyance or according to transport contract. The device can not be transported mixed with toxic, harmful, corrosive materials.

2) The packed device should be stored in room with no corrosive gas and good ventilation. Temperature: -30°C~+55°C; Relative Humidity: ≤95%.

Chapter 5 Date of manufacture, service life and accessory list

6.1 Date of manufacture: see the label.

6.2 Service life: ten years from the date of manufacture.

6.3 List of accessories































| Accessories | Quantity | Replacement cycle | Size | Replacement method | Remark |
|--------------------------|----------|---------------------------------|------------------------|-----------------------|----------------------|
| User Manual | 1 pc | No need to replace. | --- | --- | --- |
| USB cable | 1 pc | Ten years or when it is damaged | --- | --- | Contact the supplier |
| Mouthpiece | 2 pcs | Single-use | 30 mm (outer diameter) | Refer to section 4.1. | Contact the supplier |
| Power adapter (optional) | 1 pc | Ten years or when it is | --- | --- | Contact the supplier |

| | | | | | |
|--|------|---------------------|------------------------|-----|----------------------|
| | | damaged | | | |
| PC software | --- | No need to replace. | --- | --- | --- |
| Nose clip (optional) | 1 pc | Single-use | --- | --- | Contact the supplier |
| Disposable filter for respiratory (optional) | 1 pc | Single-use | 30 mm (outer diameter) | --- | Contact the supplier |

Note: If other power adapters are used, the following requirements should be met: output voltage is DC 5 V, current is no less than 1A, and the power adapter should comply with IEC 60950 or IEC 60601-1.

Chapter 7 Symbols

7.1 Symbols

| Symbol | Meaning | Symbol | Meaning |
|---|---|---|---|
|  | Full battery |  | Covering Protection rate |
|  | Low battery |  | Non-ionizing radiation |
|  | Health status indicator bar |  | Serial number |
|  | Anticlockwise rotate to unlock the turbine |  | Manufacturer |
|  | Clockwise rotate to lock the turbine |  | Type BF applied part |
|  | Do not re-use |  | For indoor use only |
|  | Do not insert |  | Class II applied |
|  | Atmospheric pressure limit |  | WEEE disposal |
|  | Temperature limit |  | Follow instructions for use |
|  | Humidity limit |  | Lot number |
|  | Fragile, handle with care |  | Date of manufacture |
|  | This way up |  | Use-by date |
|  | Keep in a cool, dry place |  | Medical device |
|  | Medical Device compliant with Directive 93/42/EEC |  | Authorized representative in the European community |
|  | Product code |  | Imported By |

7.2 Measured parameters

| Parameter | Description | Unit |
|-----------|--|------|
| FVC | Forced vital capacity (total expiratory volume) | L |
| FEV1 | Forced Expiratory Volume in one second | L |
| FEV6 | Forced Expiratory Volume in six seconds | L |
| PEF | Peak expiratory flow | L/s |
| FEV1/FVC | Forced expiratory rate in one second, FEV1/FVC×100 | % |
| FEF25 | Forced expired flow at 25% of FVC | L/s |
| FEF50 | Forced expired flow at 50% of FVC | L/s |
| FEF2575 | Forced expiratory flow between 25% and 75% of FVC | L/s |
| FEF75 | Forced expired flow at 75% of FVC | L/s |

Remark:

time zero: at the PEF (peak expiratory flow) point on the volume-time chart, draw a tangent line with the same slope as the PEF, and the intersection point between the tangent line and the time axis is the time zero.

Chapter 8 Troubleshooting

| Trouble | Possible Reason | Solution |
|---|---|--|
| The device can't finish measurement for a long time, and the data can't be displayed. | The start speed is too low, the device does not measure. | Remeasure according to the User Manual. |
| | Device malfunction. | Remeasure or restart the device. |
| | Aging of the sensor. | Please contact the local service center. |
| Data error | Operate the device falsely. | Operate the device according to the User Manual. |
| | Device malfunction. | Please contact the local service center. |
| The device can not be powered on. | Low voltage or no voltage. | Please charge the device. |
| | Aging or damage of the battery electrodes. | Please contact the local service center. |
| | Device damaged. | Please contact the local service center. |
| The display disappears suddenly. | The device is set to automatic power off when there is no operation in 2 minutes. | Normal |
| | Low voltage | Please charge the device. |
| The use time is too short after charging. | The device is not fully charged. | Please charge the device. |
| | Device battery damaged. | Please contact the local service center. |
| The device can not be fully charged after charging more than 10 hours. | Device battery damaged. | Please contact the local service center. |

Appendix I

1. Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warnings: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warnings: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

2. Instructions for use

all necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.

Guidance and manufacturer's declaration -electromagnetic emissions and Immunity.

Table 1

| Guidance and manufacturer's declaration – electromagnetic emission | |
|--|------------|
| Emission test | Compliance |
| RF emissions CISPR 11 | Group 1 |
| RF emission CISPR 11 | Class B |
| Harmonic emissions IEC 61000-3-2 | Class A |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Comply |

Table 2

| Guidance and manufacturer's declaration – electromagnetic immunity | | |
|--|--|--|
| Immunity test | IEC 60601-1-2 test level | Compliance level |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air | ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency | ±2 kV for power supply lines Not applicable 100 kHz repetition frequency |
| Surge IEC 61000-4-5 | ±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode | ±0.5 kV, ±1 kV differential mode Not applicable |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT, 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT, 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle | 0 % UT, 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT, 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle |
| Power frequency magnetic field IEC 61000-4-8 | 30 A/m 50Hz/60Hz | 30 A/m 50Hz/60Hz |
| Conducted RF IEC61000-4-6 | 3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz | 3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz |
| Radiated RF IEC61000-4-3 | 10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz | 10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz |

NOTE U_i is the a.c. mains voltage prior to application of the test level.

Table 3

| Guidance and manufacturer's declaration – electromagnetic immunity | | | | | | |
|--|----------------------|------------|---|----------------------------------|--------------------------------|------------------------|
| Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment) | Test Frequency (MHz) | Band (MHz) | Service | Modulation | IEC 60601-1-2 Test Level (V/m) | Compliance level (V/m) |
| | 385 | 380–390 | TETRA 400 | Pulse modulation 18 Hz | 27 | 27 |
| | 450 | 430–470 | GMRS 460, FRS 460 | FM ±5kHz deviation 1 kHz sine | 28 | 28 |
| | 710 | 704–787 | LTE Band 13, 17 | Pulse modulation 217 Hz | 9 | 9 |
| | 745 | | | | | |
| | 780 | | | | | |
| | 810 | 800–960 | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 | Pulse modulation 18 Hz | 28 | 28 |
| | 870 | | | | | |
| | 930 | | | | | |
| | 1720 | 1700–1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation 217 Hz | 28 | 28 |
| | 1845 | | | | | |
| | 1970 | | | | | |
| Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields) | 2450 | 2400–2570 | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7 | Pulse modulation 217 Hz | 28 | 28 |
| | 5240 | 5100–5800 | WLAN 802.11 a/n | Pulse modulation 217 Hz | 9 | 9 |
| | 5500 | | | | | |
| | 5785 | | | | | |

Table 4

| Guidance and manufacturer's declaration - electromagnetic immunity | | | | |
|--|----------------|-----------------------------|--------------------------------|------------------------|
| Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields) | Test Frequency | Modulation | IEC 60601-1-2 Test Level (A/m) | Compliance level (A/m) |
| | 30 kHz | CW | 8 | 8 |
| | 134,2 kHz | Pulse modulation 2.1 kHz | 65 | 65 |
| | 13,56 kHz | Pulse modulation 50 kHz | 7,5 | 7,5 |

Attention: With the exception of energy exchange and cables sold by manufacturers of lung function devices as spare parts for internal components, the use of accessories and cables other than those specified will result in increased product emission or reduced anti-interference.

The following cable types must be used to ensure compliance with interference radiation and immunity standards.

Table:Cable overview



Smaltimento: Il prodotto non deve essere smaltito assieme agli altri rifiuti domestici. Gli utenti devono provvedere allo smaltimento delle apparecchiature da rottamare portandole al luogo di raccolta indicato per il riciclaggio delle apparecchiature elettriche ed elettroniche.

CONDIZIONI DI GARANZIA GIMA

Si applica la garanzia B2B standard Gima di 12 mesi