

ADULT PADS**Instructions**


REF : CUA05080

Step-1


Tear open the packaging of the pads

Step-2

Remove the pads from the packaging

Step-3

Peel off the protective cover of the conducting gel of the pads

Step-4

Attach the pads to the patient. See the illustration at the back of the pads for proper placement

See separate user instructions

Contains no Latex

CU Medical Systems, Inc.

**Multifunction Defibrillation ADULT PADS**

Defibrillation, Monitoring, Synchronized Cardioversion, Pacing

! CAUTIONS

- Do not open until ready for use.
- Do not use if gel is dry.
- Do not bend, fold, crush, or puncture.
- Do not use if electrode pad is damaged.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- For long term use, apply a new set of electrodes every 24 hours
- Do not apply any substance to the skin surface that will leave a residue e.g. lotion, oils, gels, water, etc. Make sure the skin surface is dry, i.e. wipe dry water, sweat alcohol, etc.

! WARNING

- Misuse or misapplication of any electrode pad may result in patient burns or ineffective therapy.
- Pacing requires separate leads and electrodes for monitoring.

Manufacturer

Cu Medical Systems, Inc.
Dongwha Medical Instrument Complex
1647-1 Dongwha-ri, Munmak-up, Wonju-si,
Gangwon-do, 220-801 Republic of Korea
Tel : +82 33 747 7657
Fax : +82 33 747 7659

Medical Device Safety Service, GmbH
Schiffgraben 41, D-30175 Hannover,
Germany

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GB 2 self-adhesive, polymer, disposable electrodes for monitoring, defibrillation synchronised cardioversion and pacing.

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USA 2 self-adhesive, polymer, disposable electrodes for monitoring, defibrillation synchronised cardioversion and pacing. Federal (USA) law restricts this device to sale by or on the order of a physician.

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