

Digital Color Doppler Ultrasound System

QBit 3/QBit 5

Operation Manual

V2.0

Mar. 1st, 2022

57-00850-01

CHISON Medical Technologies Co., Ltd.

We reserve the right to make changes to this manual without prior notice.

Regulatory Requirement

CE 0197 This product conforms to the essential requirements of the Medical Device Directive 93/42/EEC. Accessories without the CE mark are not guaranteed to meet the Essential Requirements of the Medical Device Directive.

This manual is a reference for the QBit 3/QBit 5. Please verify that you are using the latest revision of this document. If you need the latest revision, please contact your distributor.

Statement

1. No part of this manual may be reproduced, modified, copied or reprinted, in whole or in part, without written permission from CHISON.

2. The contents of this manual are subject to change without prior notice and without our legal obligation.

3. Before operating the system, please read and understand this manual. After reading, keep this manual in an easily accessible place. If you have any question or doubt, please contact CHISON's authorized service engineer.

4. CHISON's Warranty only cover material and parts costs for repair, but does not cover any labor cost or onsite service cost at end user's side.

Meaning of the signal words

In this manual, the warning words NOTE, CAUTION and WARNING are used for regarding safety and other important instructions. Please understand their meanings clearly before reading this manual. The words and their meanings are defined as follows:

Warning words	Meaning		
	Indicates information of interest to users of the equipment as to exceptional		
	conditions or operating procedures and make the system work in good		
	condition.		
	Indicates an imminently hazardous situation which can not be avoided will		
	result in death, serious injury to the user or damage to the system.		
	Indicates a potentially hazardous situation which can not be avoided, may		
	result in death, serious injury to the user or make the system misoperation.		

Important information

1. Federal law restricts the device to sale by or on the order of a licensed practitioner or therapist.

2. It is prohibited to use the device for fetal sex examination, except for necessary medical needs.

3. The device can only be sold to qualified medical institutions or doctors.

4. The users shall have got the qualification, and shall comply with the local laws and regulations, the local religion and customs, etc.

5. The users should read the operation manual carefully before operating the devices. Turning on the device means the users have read the operation manual and accept the listed cautions, warnings, and notes in the manuals. If the users disagree and cannot accept the cautions, the users can ask for returning the device.

<u> Моте:</u>

1. It is the customer's responsibility to maintain and manage the system after delivery.

2. The warranty does not cover the following items, even during the warranty period:

a. Damage or loss due to misuse or abuse with system and probes, for example, drop the probe, the liquid or the metal part fall into the system.

b. Damage or loss caused by Acts of God such as fires, earthquakes, floods, lightning, etc.

c. Damage or loss caused by failure to meet the specified conditions for this system, such as inadequate power supply, improper installation or environmental conditions.

d. Damage or loss caused by non-approved transportation by CHISON.

e. Damage or loss due to use the system outside the region where the system was originally sold.

f. Damage or loss involving the system purchased from a source other than CHISON or its authorized agents.

3. Do not make changes or modifications to this system and probes. The System modified or repaired by people other than CHISON's qualified service engineers, CHISON shall not be liable for the system.

4. The system is to provide physicians with data for clinical diagnosis. It is the physician's responsibility for diagnostic procedures. CHISON shall not be liable for the results.

5. This manual contains warnings regarding foreseeable potential dangers, but user shall always be alert to dangers other than those indicated as well. CHISON shall not be liable for damage or loss that results from negligence or from ignoring the precautions and operating instructions described in this operation manual.

6. Important data must be backed up on external memory media. CHISON shall not be liable for loss of data stored in the memory of this system caused by operator error or accidents.

7. Please put this manual with the system to ensure operator and manager can reach it at any time. Due to negligence not following operation manual, CHISON shall not be liable for the results.

8. LED display screen may have some dark or light dots, it is normal for the LED. It does not mean that LED screen is defective.

Service Responsibility

If users install, use and maintain the system fully according to CHISON's installation manual, operation manual and service manual, then the main unit has a life time of 5 years and probes have life time of 5 years after ex-work.

The warranty of the system and probes after ex-work is as the time in the warranty card.

The system is a precise electronic system. Only the CHISON's authorized service engineer could replace the defective parts. Any assembly, disassembly, handling, repair, or replacement by any other people may have adverse impact on the safety and effectiveness of the systems and probes, and thus will reduce the life time of the system and probes, and such systems and probes will not be covered by CHISON warranty after the above improper handling. Standard maintenance must be performed by CHISON's authorized service engineer during the life time of the product.

<u> Acaution:</u>

When the above life time is expired, the effectiveness and safety of system and probes maybe greatly affected, so it's NOT suggested to continue using the system and probes even the system and probes seem work properly. But if user still wants to continue using the system and probes, user should first contact CHISON service center at CHISON headquarter to arrange the necessary safety check and calibration by CHISON's authorized service engineer. If CHISON headquarter service center provides the calibration certificate for the related system or probe, then user could continue use the system or probes according to the calibration certificate. However, if CHISON headquarter service center concludes that the system or probe is no longer complied to the safety and effectiveness standard, then user should immediately stop using the system or probe. User understands that such check and calibration cost will be born by the user.

Systems and probes keep on using after the life time may also be difficult to repair and maintain, so it's suggested to renew the product after the life time.

Content

Digital Color Doppler Ultrasound System1
Chapter 1 Introduction 1
1.1 System Overview1
1.2 Contact Information1
Chapter 2 System Safety3
2.1 Safety Overview
2.2 Electrical Safety
2.3 Label7
2.4 Patient Environmental Devices
2.5 Biological Safety
2.6 Scanning Patients and Education13
Chapter 3 System Introduction21
3.1 Console View
3.2 Physical Specification
3.3 External Interface View22
3.4 Key System Features24
3.5 Installation Procedures
Chapter 4 Control Panel
4.1 Keyboard Appearance
4.2 Alphanumeric Keyboard
4.3 Function Keys/Knobs
4.4 Central Control
4.5 Information Area Indicating Machine Status 42
Chapter 5 Operation and Exam Mode 43
5.1 Preparing the System for Use
5.2 Choose Exam Mode43
5.3 Patient Data Entry 44
5.4 Image Interface Display45
5.5 Image Mode45

	5.6 B Image Menu & Parameters	. 52
	5.7 M Image Menu & Parameters	.55
	5.8 CFM/CPA/DPD/TDI Image Menu & Parameters	.56
	5.9 PW/CW Image Menu & Parameters	.57
	5.10 Image and Cine Disposition	.57
	5.11 Edit Comment	. 59
	5.12 Set Body Mark	61
	5.13 Set the Direction of Arrow	. 62
	5.14 Image Browse	. 62
	5.15 Archive Management	.64
	5.16 Report	.65
	5.17 DICOM	. 66
Cha	apter 6 Measurement and Calculation	69
	6.1 Keyboard for Measurement	. 69
	6.2 B Mode General Measurement Methods	.70
	6.3 B Fast Measurement	.74
	6.4 B General Measurement	. 75
	6.5 B Abdomen Measurement	. 76
	6.6 B OB Measurement	.77
	6.7 B Pediatrics Measurement	. 80
	6.8 B ORTH Measurement	. 81
	6.9 B GYN Measurement	.82
	6.10 B Small Parts Measurement	.83
	6.11 B Vascular Measurement	.83
	6.12 B Urology Measurement	.83
	6.13 B Cardiology Measurement	. 84
	6.13 Normal Measurement in M, B/M mode	. 87
	6.14 General Measurement in M mode	. 88
	6.15 M Cardiology Measurement	.88
	6.16 Other Urology Measurement in M Mode	. 90

6.17 PW mode measurement methods9	0
6.18 PW Fast Measurement9	2
6.19 PW General Measurement9	2
6.20 PW OB Measurement9	4
6.21 PW GYN Measurement9	4
6.22 PW Cardiology Measurement9	5
6.23 PW Vascular Measurement10	0
6.24 Other Urology Measurement in PW Mode10	0
Chapter 7 Preset 10	1
7.1 General setting10	1
7.2 Measurement	4
7.3 Comment11	3
7.4 Body marks11	6
7.5 Exam Mode11	7
7.6 Report12	1
7.7 DICOM	4
7.8 Network	5
7.9 System	6
Chapter 8 System Maintenance12	8
8.1 Machine Cleaning12	8
8.2 Safety Check	8
8.3 Malfunction Check	9
Chapter 9 Probes	0
9.1 General Description13	0
9.2 Care and Maintenance	0
9.3 Probe Operation Instructions14	3
Appendix A: Acoustic Output Report Table14	6
Appendix B: Guidance and Manufacturer's Declaration	2
Appendix C: Measurement Results Summary23	6
Appendix D: Display Accuracy and Acoustic Measurement Uncertainties	7

Appendix E: Transducer Maximum Surfa	ce Temperature	238
Appendix F: Procedures of setting netwo	ork sharing	

Chapter 1 Introduction

This manual contains necessary information for safe system operation.

Read and understand all instructions in this manual before operating the system. Always keeping this manual with the equipment, and periodically review the procedures for operation and safety precautions.

1.1 System Overview

Indications for Use

The device is a general-purpose ultrasonic imaging instrument intended for use by a qualified physician for evaluation of Fetal, Abdominal, Pediatric, Small Organ (breast, thyroid, testes), Adult Cephalic, Cardiac(adult, pediatric), Musculo-skeletal(Conventional, Superficial), Peripheral Vascular, Transrectal, Transvaginal, Urology.

Contraindication

The system is NOT intended for ophthalmic use or any use that causes the acoustic beam to pass through the eye.

1.2 Contact Information

For additional information or assistance, please contact your local distributor or the appropriate support resource shown below:

CHISON website	www.chison.com
Service Support	CHISON Medical Technologies Co., Ltd.
	Tel:0086-0510-85311707
	Fax: 0086-0510-85310726
	E-mail: service@chison.com.cn
Placing an Order	CHISON Medical Technologies Co., Ltd.
	Tel: 0086-0510-8531-0593/0937
	Fax: 0086-0510-85310726
	Email: export@chison.com.cn

Manufacturer	CHISON Medical Technologies Co., Ltd.		
	No.3 Changjiang South Road, Xinwu District, Wuxi, 214028 Jiangsu, P.R.		
	China		
US Agent	Mr. Marco Mu, 2219 Rimland Drive,Suite 301,Bellingham, Barkley Villiage		
	Bellingham,Washington,98226, UNITED STATES		
	Phone: 360-3257028, Fax: 360-9253199, Email: us.agent@mid-link.net		
	MID-LINK INTERNATIONAL CO., LTD		
EC Representative	Shanghai International Holding Corp.GmbH (Europe)		
	Address: Eiffestrasse 80, 20537 Hamburg, Germany		
	Tel: 0049-40-2513175		
	Fax: 0049-40-255726		
	E-mail: antonsissi@hotmail.com shholding@hotmail.com		

Chapter 2 System Safety

2.1 Safety Overview

This section discusses the measures to ensure the safety of both the operator and patient. To ensure the safety of both operator and patient, please read the relevant details in this chapter carefully before operating this system. **Disregarding the warnings or violation of relevant rules may result in personal injury for operator or patient, or even loss of life.**

Users should observe the following precautions:

> This system complies with Type BF general equipment, and the IEC standard. Please follow this operation manual to use this system properly.

> Please do not modify this system in any way. If modifications are necessary, please contact the manufacturer first to get more information and permission.

> This system has been fully adjusted at the factory. Do not adjust any fixed adjustable parts.

In the event of a malfunction, turn off the system immediately and inform the manufacturer or its designated agents.

> The power cord of the system should be connected to a grounded power socket. Do not remove the ground cable for any reason.

> Only connect this system, either electronically or mechanically, with devices that comply with the IEC/EN60601-1 standard. Recheck the leakage current and other safety performance indices of the entire system to avoid potential system damage caused by leakage from a current superposition.

> The system does not incorporate any specialized protective measures in the event it is configured with high-frequency operation devices. The operator should use caution in these types of applications.

> The system should be installed only by personnel authorized by the manufacturer. Do not attempt to install the system by yourself.

> Only a CHISON's authorized service engineer can perform maintenance.

> Only a qualified operator, or someone under qualified supervision, can use the system.

> Do not use this system in the presence of flammable substances, otherwise an explosion may

occur.

> Do not continuously scan the same part of a patient or expose the patient to prolonged scanning. Otherwise, it may harm the patient.

> When using the system for ultrasound testing, only use qualified ultrasound gel that complies with system standards.

> Do not unplug probe when the system is in active operation. Always go to transducer Selection screen when need to remove the probe.

> To prevent from arm or neck injury, the operator should not stay at the same position for too long during patient scanning without taking break.

> Do not put liquid on top of the main unit.

If there is any liquid or metal to enter to the system, please power off the system and stop using it immediately. Please first contact CHISON's authorized service engineer to make sure it's safe before restarting to use it.c

MNOTE:

* The system has built-in screen saver to avoid the tic mark on the display. It is not recommended to constantly turn on and off the unit.

* To dispose of this product properly, please contact the local CHISON's Authorized Service Representative.

2.2 Electrical Safety

Type of protection against electric shock

Class I Equipment

CLASS I EQUIPMENT in which protection against electric shock does not rely on basic insulation only, but which includes an additional safety precaution in that accessible conductive parts are connected to the protective earthing conductor in the electrical installation in such a way that accessible parts cannot become live in the event of a failure of the basic insulation.

Degree of protection against electric shock

• Type BF Applied part (for Probes marked with BF symbol)

TYPE BF APPLIED PART providing a specified degree of protection against electric shock, with

particular regard to allowable LEAKAGE CURRENT.

Level of protection against harmful ingress of water

• Parts of probe likely to come into contact with operator or patient meet the requirements of drip-proof equipment (IPX1)

Parts of probe intended to be immersed in normal use meet the requirements of watertight equipment (IPX7).

• The IP Classification of System is Ordinary Equipment (IPX0)

The Equipment is not suitable for use in the presence of a flammable anesthetic mixed with air (with oxygen or with oxide)

Mode of operation

• Continuous Operation

For maximum safety, always follow these guidelines:

> Proper grounding of the system is critical to avoid electric shock. For protection, ground the chassis with a three-wire cable and plug the system into a hospital-grade, three-hole outlet.

> Do not remove or circumvent the grounding wire.

> Do not remove the protective covers on the system. These covers protect users against hazardous voltages. Cabinet panels must remain in place while the system is in use. A qualified electronic technician must make all internal replacements.

> Do not operate this system in the presence of flammable gases or anesthetics.

All peripheral devices (unless certified as medical grade) that are connected to the system must be powered through the electrical outlet with an optional isolation transformer.

Notice upon Installation of Product

Separation distance and effect from fixed radio communications equipment: field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast transmitter cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength

in the location in which the ultrasound system is used exceeds the applicable RF compliance level as stated in the immunity declaration, the ultrasound system should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as re-orienting or relocating the ultrasound system or using an RF shielded examination room may be necessary.

• Use either power supply cords provided or designated by CHISON. Products equipped with a power source plug should be plugged into the fixed power socket which has the protective grounding conductor. Never use any adaptor or converter to connect with a power source plug (e.g. three-prong-to-two-prong converter).

- Locate the equipment as far away as possible from other electronic equipment.
- Be sure to only use the cables provided by or designated by CHISON. Connect these cables

following the installation procedures (e.g. wire power cord separately from signal cables).

• Lay out the main equipment and other peripherals following the installation procedures described in this manual.

Notice against User Modification

The user should never modify this product.

User modifications may cause degradation in Electrical Safety. Modification of the product includes changes in:

- Cables (length, material, wiring, etc.)
- System configuration/components

User modifications may cause degradation in EMC performance. Modification of the product includes changes in:

- Cables (length, material, wiring, etc.)
- System installation/layout
- System configuration/components
- Securing system parts (cover open/close, cover screwing)

2.3 Label



Rear panel label

2.3.1 Warning Symbols on Label

lcon	Meaning
	Do not use the following devices near this equipment: cellular phone, radio receiver, and mobile radio transmitter, radio controlled toy, etc. Use of these devices near this equipment could cause this equipment to perform outside the published specifications. Keep power to these devices turned off when near this equipment.
	This mark indicates that this product contains a limited amount of hazardous substances in the Chinese Standard GB/T 26572-2011 "Limited Requirements for Restricted Substances in Electrical and Electronic Products". The numbers in the logo are the environmental protection use period of the product, indicating that under the normal use conditions, the harmful substances will not leak or be abrupt. The use of the product will not cause serious pollution to the environment or cause personal or property serious damage, the term unit is year.
6	Refer to instruction manual/booklet.
⚠	Caution, consult accompanying documents. This symbol advises the reader to consult the accompanying documents for important safety related information such as warnings and pre-cautions that cannot be presented on the device itself.
X	WASTE OF ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE): This symbol is used for Environment Protection, it indicates that the waste of electrical

	and electronic equipment must not be disposed as unsorted waste and must be				
	collected separately. Please contact your local Authority or distributor of the				
	manufacturer for information concerning the decommissioning of your equipment.				
(F 0197	The CE mark of Conformity indicates this equipment conforms to the Council				
	Directive 93/42/EEC.				
Rx only	This symbol indicates that federal law restricts the device to sale by or on the				
	order of a licensed practitioner or therapist in the united states.				
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY: This				
EC REP	symbol is accompanied by the name and the address of the authorized				
	representative in the European Community.				
	MANUFACTURER: This symbol is accompanied by the name and the address of				
	the manufacturer.				
	Manufacturing date of the device in the form YYYY-MM.				
SN	Serial number of the device.				
~	The "Alternating current" symbol indicates that the equipment is suitable for				
	alternating current only.				

2.3.2 Symbols used in the system

lcon	Meaning	
人	Power On/off	
	CAUTION: This Power Switch cannot isolate Mains Supply completely.	
Ϊ	Type-BF applied part	
	Main power switch ON	
\bigcirc	Main power switch OFF	
∇	Potential equilibrium connection	
IPX0	No protection against ingress of water (system)	
IPX7	Protection against the effects of immersion (probes)	

2.4 Patient Environmental Devices

Right side:

- 2 Probe ports (standard)
- 2 USB ports on control panel

Rear panel:

- 4 USB ports
- 1 VIDEO OUT port
- 1 DVI port
- 1 VGA port
- 1 Remote port
- 1 ECG port
- 2 Footswitch ports
- 1 LAN port
- 1 S-VIDEO port
- 1 Power in port
- 1 Power out port

Acceptable Devices

The Patient Environmental devices shown above are specified to be suitable for use within the PATIENT ENVIRONMENT.

• DO NOT connect any probes or accessories without approval by CHISON within the PATIENT ENVIRONMENT.

• DO NOT touch patient and devices without IEC/EN 60601-1 approval to avoid the leakage current risk within the PATIENT ENVIRONMENT.

Unapproved Devices

- DO NOT use unapproved devices.
- If devices are connected without the approval of CHISON, the warranty will be INVALID.

 The system can't be used with HF surgical equipment; otherwise the burns to patient may occur.

Any device connected to this system must conform to one or more of the requirements listed below:

- IEC standard or equivalent standards appropriate to devices.
- The devices shall be connected to PROTECTIVE EARTH (GROUND).

ACAUTION:

Unsafe operation or malfunction may occur. Use only the accessories, options and supplies approved or recommended in these instructions for use.

Peripheral used in the patient environment

The system has been verified for overall safety, compatibility and compliance with the following on-board image recording devices:

B/W video printer: SONY UP-X898MD

The system may also be used safely while connected to devices other than those recommended above if the devices and their specifications, installation, and interconnection with the system conform to the requirements of IEC/EN 60601-1.

The connection of equipment or transmission networks other than as specified in the user instructions can result in an electric shock hazard or equipment malfunction. Substitute or alternate equipment and connections require verification of compatibility and conformity to IEC/EN 60601-1 by the installer. Equipment modifications possible resulting malfunctions and electromagnetic interference are the responsibilities of the owner.

General precautions for installing an alternate off-board, remote device or a network would include:

- The added device(s) must have appropriate safety standard conformance and CE Marking.
- There must be adequate mechanical mounting of the device and stability of the combination.
- Risk and leakage current of the combination must comply with IEC/EN 60601-1.
- Electromagnetic emissions and immunity of the combination must conform to IEC/EN

60601-1-2.

Peripheral used in the non-patient environment

The system has been verified for compatibility, and compliance for connection to a local area network (LAN) via a wire LAN. The provided LAN components are IEC/EN 60950 compliant.

General precautions for installing an alternate off-board, remote device or a network would include:

- The added device(s) must have appropriate safety standard conformance and CE Marking.
- The added device(s) must be used for their intended purpose having a compatible interface.

CAUTION: Make sure using ONLY the dedicated USB disk or removable media to save or back up data. Before connecting to the ultrasound system, make sure using the latest antivirus software on the USB disk or removable media to clean any virus. It is user's responsibility to ensure the USB disk or removable media is virus-free. Improper use of USB disk or removable media may cause the virus infections of system and eventually malfunction may occur. Such malfunction may impact the stability, effectiveness and safety of the system and probes, and users should immediately stop using the system and probes until CHISON authorized engineer has checked the system and confirm the effectiveness and safety of the system and probes.

CAUTION: Use only secure Local Area Network connection. Don't connect the ultrasound system to Internet. Make sure your hospital's firewall software is configured correctly, thus blocking incoming connection requests from Internet. Improper use of network connection may cause the virus infections of system and eventually malfunction may occur.

2.5 Biological Safety

This product, as with all diagnostic ultrasound equipment, should be used only for valid reasons and should be used both for the shortest period of time and at the lowest power settings necessary (ALARA - As Low As Reasonably Achievable) to produce diagnostically acceptable images. The AIUM offers the following guidelines:

Clinical Safety Quoted from AIUM

Approved March 26, 1997

Diagnostic ultrasound has been in use since the late 1950s. Given its known benefits and recognized efficacy for medical diagnosis, including use during human pregnancy, the American Institute of Ultrasound in Medicine herein addresses the clinical safety of such use:

There are no confirmed biological effects on patients or instrument operators caused by exposures from present diagnostic ultrasound instruments. Although the possibility exists that such biological effects may be identified in the future, current data indicate that the benefits to patients of the prudent use of diagnostic ultrasound outweigh the risks, if any that may be present.

<u>Heating</u>: Elevating tissue temperature during obstetrical examinations creates medical concerns. At the embryo development stage, the rise in temperature and the length of time exposed to heat combine to determine potential detrimental effects. Exercise caution particularly during Doppler/Color exams. The Thermal Index (TI) provides a statistical estimate of the potential temperature elevation (in centigrade) of tissue temperature. Three forms of TI are available: Soft Tissue Thermal Index (TIS), Bone Thermal Index (TIB) and Cranial Bone Thermal Index (TIC).

Soft Tissue Thermal Index (TIS). Used when imaging soft tissue only, it provides an estimate of potential temperature increase in soft tissue.

Bone Thermal Index (TIB). Used when bone is near the focus of the image as in the third trimester OB examination, it provides an estimate of potential temperature increase in the bone or adjacent soft tissue.

Cranial Bone Thermal Index (TIC). Used when bone is near the skin surface as in transcranial examination, it provides an estimate of potential temperature increase in the bone or adjacent soft tissue.

<u>Cavitations:</u> Cavitations may occur when sound passes through an area that contains a cavity, such as a gas bubble or air pocket (in the lung or intestine, for example). During the process of cavitations, the sound wave may cause the bubble to contract or resonate. This oscillation may cause the bubbles to explode and damage the tissue. The Mechanical Index (MI) has been created to help users accurately evaluate the likelihood of cavitations and the related adverse effects.

MI recognizes the importance of non-thermal processes, cavitations in particular, and the Index is an attempt to indicate the probability that they might occur within the tissue.

2.6 Scanning Patients and Education

The Track-3 or IEC60601-2-37 output display standard allows users to share the responsibility for the safe use of this ultrasound system. Follow these usage guidelines for safe operation:

> In order to maintain proper cleanliness of the probes, always clean them between patients.

> Always use a disinfected sheath on all EV/ER probes during every exam.

> Continuously move the probe, rather than staying in a single spot, to avoid elevated temperatures in one part of the patient's body.

> Move probe away from the patient when not actively scanning.

> Understand the meaning of the TI, TIS, TIB, TIC and MI output display, as well as the relationship between these parameters and the thermal/cavitation bio effect to the tissue.

> Expose the patient to only the very lowest practical transmit power levels for the shortest possible time to achieve a satisfactory diagnosis (ALARA - As Low As Reasonably Achievable).

2.6.1 Safe Scanning Guidelines

• Ultrasound should only be used for medical diagnosis and only by trained medical personnel.

• Diagnostic ultrasound procedures should be done only by personnel fully trained in the use of the equipment, in the interpretation of the results and images, and in the safe use of ultrasound (including education as to potential hazards to the patient and the operator).

• Operators should understand the likely influence of the machine controls, the operating mode (e.g. B-mode, color Doppler imaging or spectral Doppler) and probe frequency on thermal and cavitations hazards.

• Select a low setting for each new patient. Output should only be increased during the examination if penetration is still required to achieve a satisfactory result, and after the Gain control has been adjusted to its maximum value.

• Maintain the shortest examination time necessary to produce a useful diagnostic result.

• Do not hold the probe in a fixed position for any longer than is necessary. It should be removed from the patient whenever there is no need for real-time imaging or spectral Doppler acquisition. The frozen frame and Cine loop capabilities allow images to be reviewed and discussed without exposing the patient to continuous scanning.

• Do not use endo-cavitary probes if there is noticeable self heating of the probe when operating in the air. Although applicable to any probe, take particular care during trans- vaginal exams during the first eight weeks of gestation.

• Take particular care to reduce output and minimize exposure time of an embryo or fetus when the temperature of the mother is already elevated.

• Take particular care to reduce the risk of thermal hazard during diagnostic ultrasound when exposing: an embryo less than eight weeks after gestation; or the head, brain or spine of any fetus or neonate.

• Operators should continually monitor the on-screen thermal index (TI) and mechanical index (MI) values and use control settings that keep these settings as low as possible while still achieving diagnostically useful results. In obstetric examinations, TIS (soft tissue thermal index) should be monitored during scans carried out in the first eight weeks after gestation, and TIB (bone thermal index) thereafter. In applications where the probe is very close to bone (e.g. trans-cranial applications), TIC (cranial bone thermal index) should be monitored.

MI> 0.3 There is a possibility of minor damage to neonatal lung or intestine. If such exposure is necessary, reduce the exposure time as much as possible.

MI> 0.7 There is a risk of cavitations if an ultrasound contrast agent containing gas micro-spheres is being used. There is a theoretical risk of cavitations without the presence of ultrasound contrast agents. The risk increases with MI values above this threshold.

TI> 0.7 The overall exposure time of an embryo or fetus should be restricted in accordance with Table 2-2 below as a reference:

TI	Maximum exposure time (minutes)	
0.7	60	
1.0	30	
1.5	15	
2.0	4	
2.5	1	

Maximum recommended exposure times for an embryo or fetus

• Non-diagnostic use of ultrasound equipment is not generally recommended. Examples of non-diagnostic uses of ultrasound equipment include repeated scans for operator training, equipment demonstration using normal subjects, and the production of souvenir pictures or videos of a fetus. For equipment of which the safety indices are displayed over their full range of values, the TI should always be less than 0.5 and the MI should always be less than 0.3. Avoid frequent repeated exposure of any subject. Scans in the first trimester of pregnancy should not be carried out for the sole purpose of producing souvenir videos or photographs, nor should their production involve increasing the exposure levels or extending the scan times beyond those needed for clinical purposes.

• Diagnostic ultrasound has the potential for both false positive and false negative results. Misdiagnosis is far more dangerous than any effect that might result from the ultrasound exposure. Therefore, diagnostic ultrasound system should be performed only by those with sufficient training and education.

2.6.2 Understanding the MI/TI Display

Track-3 follows the Output Display Standard for systems that include fetal Doppler applications. The acoustic output will not be evaluated on an application-specific basis, but the global maximum de-rated lspta must be \leq 720 mW/cm² and either the global maximum MI must be \leq 1.9 or the global maximum de-rated lsppa must be \leq 190 W/cm². An exception is for ophthalmic use, in which case the TI=max (TIS_as, TIC) is not to exceed 1.0; lspta.3 \leq 50mW/cm², and MI \leq 0.23. Track-3 gives the user the freedom to increase the output acoustic power for a specific exam, and still limit output acoustic power within the global maximum de-rated lspta \leq 720 mW/cm² under an Output Display Standard.

For any diagnostic ultrasonic systems, Track-3 provides an Output Indices Display Standard. The diagnostic ultrasound systems and its operation manual contain the information regarding an ALARA (As Low As Reasonably Achievable) education program for the clinical end-user and the acoustic output indices, MI and TI. <u>The MI describes the likelihood of cavitations, and the TI offers</u> the predicted maximum temperature rise in tissue as a result of the diagnostic examination. In general, a temperature increase of 2.5°C must be present consistently at one spot for 2 hours to cause fetal abnormalities. Avoiding a local temperature rise above 1°C should ensure that no thermally induced biologic effect occurs. When referring to the TI for potential thermal effect, a TI equal to 1 does not mean the temperature will rise 1 degree C. It only means an increased potential for thermal effects can be expected as the TI increases. A high index does not mean that bioeffects are occurring, but only that the potential exists and there is no consideration in the TI for the scan duration, so minimizing the overall scan time will reduce the potential for effects. These operator control and display features shift the safety responsibility from the manufacturer to the user. So it is very important to have the Ultrasound systems display the acoustic output indices correctly and the education of the user to interpret the value appropriately.

RF: (De-rating factor)

In Situ intensity and pressure cannot currently be measured. Therefore, the acoustic power measurement is normally done in the water tank, and when soft tissue replaces water along the ultrasound path, a decrease in intensity is expected. The fractional reduction in intensity caused by attenuation is denoted by the de-rating factor (RF),

RF=10 (-0.1 a f z)

Where a is the attenuation coefficient in dB cm-1 MHz-1, f is the transducer center frequency, and z is the distance along the beam axis between the source and the point of interest.

De-rating factor RF for the various distances and frequencies with attenuation coefficient 0.3dB cm-1 MHz-1 in homogeneous soft tissue is listed in the following table. An example is if the user uses 7.5MHz frequency, the power will be attenuated by .0750 at 5cm, or 0.3x7.5x5=-11.25dB. The De- rated Intensity is also referred to as '.3' at the end (e.g. Ispta.3).

Distance	Frequency (MHz)			
(cm)	1	3	5	7.5
1	0.9332	0.8128	0.7080	0.5957
2	0.8710	0.6607	0.5012	0.3548
3	0.8128	0.5370	0.3548	0.2113
4	0.7586	0.4365	0.2512	0.1259
5	0.7080	0.3548	0.1778	0.0750
6	0.6607	0.2884	0.1259	0.0447
7	0.6166	0.2344	0.0891	0.0266
8	0.5754	0.1903	0.0631	0.0158

I'=I*RF Where I' is the intensity in soft tissue, I is the time-averaged intensity measured in water.

Tissue Model:

Tissue temperature elevation depends on power, tissue type, beam width, and scanning mode. Six models are developed to mimic possible clinical situations.

	Thermal Models	Composition	Mode	Specification	Application
1	TIS	Soft tissue	Unscanned	Large aperture (>1cm²)	Liver PW
2	TIS	Soft tissue	Unscanned	Small aperture (<1cm ²)	Pencil Probe
3	TIS	Soft tissue	Scanned	Evaluated at surface	Breast color
4	TIB	Soft tissue and bone	Scanned	Soft tissue at surface	Muscle color
5	TIB	Soft tissue and bone	Unscanned	Bone at focus	Fetus head PW
6	TIC	Soft tissue and bone	Unscanned/scann ed	Bone at surface	Transcranial

Soft tissue:

Describes low fat content tissue that does not contain calcifications or large gas-filled spaces.

Scanned: (auto-scan)

Refers to the steering of successive burst through the field of view, e.g. B and color mode.

Unscanned:

Emission of ultrasonic pulses occurs along a single line of sight and is unchanged until the transducer is moved to a new position. For instance, the PW, and M mode.

<u>TI:</u>

TI is defined as the ratio of the In Situ acoustic power (W.3) to the acoustic power required to raise tissue temperature by 1°C (Wdeg), TI=W.3/Wdeg.

Three TIs corresponding to soft tissue (TIS) for abdominal; bone (TIB) for fetal and neonatal cephalic; and cranial bone (TIC) for pediatric and adult cephalic, have been developed for applications in different exams.

An estimate of the acoustic power in milli-watts necessary to produce a 1°C temperature elevation in soft tissue is:

Wdeg=210/fc,for model 1 to 4, where fc is the center frequency in MHz.

Wdeg=40 K Dfor model 5 and 6, where K (beam shape factor) is 1.0, D is the aperture diameter in cm at the depth of interest.

<u>MI:</u>

Cavitation is more likely to occur at high pressures and low frequencies in pulse ultrasound wave in the tissue, which contains the bubble or air pocket (for instance, the lung, intestine, or scan with gas contrast agents). The threshold under optimum conditions of pulsed ultrasound is predicted by the ration of the peak pressure to the square root of the frequency.

MI=Pr'/sqrt(fc)

Pr' is the de-rated (0.3) peak rare-fractional pressure in Mpa at the point where PII is the maximum, and fc is the center frequency in MHz. PII is the Pulse Intensity Integral that the total energy per unit area carried by the wave during the time duration of the pulse. The peak rare- fractional pressure is measured in hydrophone maximum negative voltage normalized by the hydrophone calibration parameter.

Display Guideline:

For different operation modes, different indices must be displayed. However, only one index needs to be shown at a time. Display is not required if maximum MI is less than 1.0 for any setting of the operating mode, or if maximum TI is less than 1.0 for any setting of the operating mode. For TI, if the TIS and TIB are both greater than 1.0, the scanners need not be capable of displaying both indices simultaneously. If the index falls below 0.4, no display is needed. The display increments

are no greater than 0.2 for index value less than one and no greater than 1.0 for index values greater than one (e.g. 0.4, 0.6, 0.8, 1, 2, and 3).

Display and Report in Different Mode

Located on the upper middle section of the system display monitor, the acoustic output display provides the operator with real-time indication of acoustic levels being generated by the system.

For B-Scan Mode

Only display and report MI, and start from 0.4 if maximum MI > 1.0, display in increments of 0.2.

For Color Mode

Only display and report TIS or TIB and start from 0.4 if maximum TI > 1.0, display in increments of 0.2 for values of indices of 2.0 or less, and 0.5 for values of indices greater than 2.0.

For Doppler Mode

Only display and report TIS or TIB and start from 0.4 if maximum TI > 1.0, display in increments of 0.2 for values of indices of 2.0 or less, and 0.5 for values of indices greater than 2.0.

Below is a simple guideline for the user when TI exceeds one limit exposure time to 4(6-TI) minutes based on the 'National Council on Radiation Protection. Exposure Criteria for Medical Diagnostic Ultrasound: I. Criteria Based on Thermal Mechanisms. Report No.113 1992'.

Operator Control Features:

The user should be aware that certain operator controls may affect the acoustic output. It is recommended to use the default (or lowest) output power setting and compensate using Gain control to acquire an image. Other than the output power setting in the soft-menu, which has the most direct impact on the power; the PRF, image sector size, frame rate, depth, and focal position also slightly affect the output power. The default setting is normally around 70% of the allowable power depending on the exam application mode.

Controls Affecting Acoustic Output

The potential for producing mechanical bioeffects (MI) or thermal bioeffects (TI) can be influenced by certain controls.

Direct: The Acoustic Output control has the most significant effect on Acoustic Output.

Indirect: Indirect effects may occur when adjusting controls. Controls that can influence MI and TI are detailed under the bio effect portion of each control in the Optimizing the Image chapter. Always observe the Acoustic Output display for possible effects.

Best practices while scanning

HINTS: Raise the Acoustic Output only after attempting image optimization with controls that have no effect on Acoustic Output, such as Gain and STC.

WARNING: Be sure to have read and understood control explanations for each mode used before attempting to adjust the Acoustic Output control or any control that can affect Acoustic Output.

Use the minimum necessary acoustic output to get the best diagnostic image or measurement during an examination. Begin the exam with the probe that provides an optimum focal depth and penetration.

Acoustic Output Default Levels

In order to assure that an exam does not start at a high output level, the system initiates scanning at a reduced default output level. This reduced level is preset programmable and depends upon the exam icon and probe selected. It takes effect when the system is powered on or New Patient is selected. To modify acoustic output, adjust the Power Output level on the Soft Menu.

Chapter 3 System Introduction

3.1 Console View



Fig. 3-1 Console View

NOTE: The photograph above is adapted for 15 inch monitor. For device of 19 inch monitor please be subject to the actual product.

3.2 Physical Specification

Dimensions of main unit (approx.):

- 1. 15 inch monitor: 747mm (Length)*517mm (Width)*1283mm (Height)
- 2. 19 inch monitor (option): 747mm (Length) ×517mm (Width) ×1358.23mm (Height)

Net weight of main unit (approx.): 50 kg (no probe included)

3.3 External Interface View




Fig. 3-2 Front Side View 1. Monitor 2. Speaker 3. Keyboard 4. Probe holder




Fig. 3-4 Rear panel 1.USB 2.0 ports 2.Remote port 3.LAN port 4.DVI port 5.Foot Switch ports 6.Video Out port 7.S-Video port 8.VGA port




Fig. 3-5 Power panel

1. Power In 2. Power Switch 3. Potential equilibrium connection 4. Power Out

3.4 Key System Features

1. Zoom and depth adjustment.

2. Set the total gain, contrast, frequency band, 8 segments of STC, dynamic range, persistence.

3. 256 gray-scale image display technology, Q-image technology, stable performance, high resolution.

4. Image freezing and storage function; the stored images can be recalled for analysis.

5. Scanning direction can be changed and the image can be reversed in left/right, up/down direction.

6. Distance, area, circumference, volume, fetal weight, heart rate etc. measurements are available and automatic calculation of OB, cardiology are available. Direct display of gestation age and expected date of child delivery.

7. Elliptical method and tracing method are provided for area/circumference measurement.

8. Many kinds of body marks can be displayed together with corresponding probe position indication.

9. Comment function in image area of the screen, special comment terms for different exam mode can be added according to user's requirement.

10. Display of Patient ID, Time and Date display according to real-time clock.

11. Trackball available for operation and measurement. Characters can be input directly by keyboard.

12. When one function is under operation, the corresponding key on the control panel will be brightly lit. When exiting from the function, the corresponding key on the control panel will be slightly lit.

13. Measure the percentage of stenosis, blood flow velocity, velocity ratio, blood flow volume and pressure gradient. Automatically measure the values of maximum velocity, minimum velocity, time interval, pulsatility index and resistance index.

14. Possess multi-language interface display User interface change, shear plate, printing, DICOM

3.0, biopsy guided functions.

15. It has permanent storage for image and cine on HDD. It can also be connected to removable storage via USB 2.0 port. To realize mass storage ,can recall saved image for analysis.

16. Cine loop storage 256 frames real time image.

17. Output standard PAL or NTSC video signal and VGA/DVI signal.

18. Print or export graphic report.

3.4.1 Image Modes

Mode	QBit 3	QBit 5
B Mode	Standard	Standard
2 B Mode	Standard	Standard
B/M Mode	Standard	Standard
M Mode	Standard	Standard
4B Mode	Standard	Standard
2D Steer	Standard	Standard
CFM Mode	Standard	Standard
CPA(PD) Mode	Standard	Standard
DPD Mode	Standard	Standard
PW Mode	Standard	Standard
Triplex Mode	Standard	Standard
CW Mode	N/A	Option
TDI	N/A	Option
Color M Mode	N/A	Option
Trapezoidal Mode	Standard	Standard
ECG	N/A	Option
B/BC Mode	Standard	Standard
Super Needle	Option	Option
HD 3D	Option	Option

3.4.2 Accessories

Probes:

Probe	Frequency	QBit 3	QBit 5
D3C60L Convex probe	3.5MHz, 2.0-6.8MHz	Option	Option
D7L40L Linear probe	7.5MHz, 4.0-15.0MHz	Option	Option
Micro convex probe D6C12L	6.0MHz, 4.0-12.0MHz	Option	Option
Micro convex probe D7C10L	7.5MHz, 4.0-15.0MHz	Option	Option
Linear probe D7L40L-REC	7.5MHz, 4.0-15.0MHz	Option	Option
Micro convex probe D3C20L	3.0MHz, 2.0-6.8MHZ	Option	Option

Micro convex probe D6C15L	6.0MHz, 4.0-12.0MHz	Option	Option
Micro convex probe D5C20L	5.0MHz, 4.0-10.7MHz	Option	Option
Phased array probe D3P64L	2.5MHz, 1.5-5.3MHz	N/A	Option
Linear probe D12L40L	12.0MHz, 7.0-15.0MHz	Option	Option

Peripherals

S-VIDEO, VGA, DVI output for external monitor

VIDEO output for B&W video printer

LAN port output

LAN for DICOM and image review station

USB 2.0 for flash drive

Foot Switch

3.5 Installation Procedures

NOTE: Please do not turn on the power switch until finishing all the installation and necessary preparation.

3.5.1 Environment Condition

The system should be operated under the following environment.

3.5.1.1 Operation Environment Requirement

Ambient Temperature: 10 °C~40 °C

Relative Humidity: 30%~75%RH

Atmospheric Pressure: 700hPa~1060hPa

Strong radiation sources or powerful electromagnetic waves (e.g. electro-magnetic waves from radio broadcasting) may result in image ghosting or noise. The system should be isolated from such radiation sources or electromagnetic waves.

To prevent damage to the system, do not use in the following locations:

- Exposed to direct sunlight
- > Subject to sudden changes in temperature
- Dusty
- Subject to vibration
- Near heat generators
- ➢ High humidity

This equipment generates, uses and can radiate radio frequency energy. The equipment may cause radio frequency interference to other medical and non-medical devices and radio communications. To provide reasonable protection against such interference, this product complies with emissions limits for a Group 1, Class A Medical Devices Directive as stated in IEC/EN 60601-1-2. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment is found to cause interference (which may be determined by turning the equipment on and off), the user (or qualified service personnel) should attempt to correct the problem by one or more of the following measure(s):

- Reorient or relocate the affected device(s).
- Increase the separation between the equipment and the affected device.
- Power the equipment from a source different from that of the affected device.
- Consult the point of purchase or service representative for further suggestions.

3.5.1.2 Transport and Storage Environmental Requirement

The following environmental transport and storage conditions are within system tolerances:

Temperature: -5°C ~ 40°C

Relative Humidity: ≤ 80% non-condensing

Atmosphere Pressure: 700hPa ~ 1060hPa

3.5.1.3 Electrical Requirements

Power Requirements: AC 100-240V, 50-60Hz.

Fuse Requirements: Fuse specification is 250V, 5.0 A (time-lag), the model is 50T T5AL 250V.

CAUTION: Please use the fuse provided by manufactory if the fuse need replace. User can't buy and exchange the fuse by their self.

Power Consumption: 600 VA.

Voltage Fluctuation

WARNING:

Maintain a fluctuation range of less than $\pm 10\%$ of voltage labeling on rear panel of the system, otherwise the system may be damaged.

Grounding

Before connecting the power cable, connect the attached ground protection cable from Equipotentiality terminal on system rear panel to a specialized grounding device.

• Please follow the outlined power requirements. Only use power cables that meet the system guidelines-failure to follow these procedures may produce system damage.

• Line power may vary in different geographic locations. Refer to the detailed ratings on the rear panel of the system for detailed information.

3.5.1.4 Operation Space

Please leave enough free space for the system to ensure well ventilation.

CAUTION: Leave enough free space for the system; otherwise, with the increasing of the temperature inside the unit, malfunction may occur.

3.5.2 System Positioning & Transporting

When moving or transporting the system, take the precautions described below to ensure maximum safety for personnel, the system and other equipment.

Before Moving the System

- > Power off the system, refer to power off section.
- > Switch off the power and unplug the power cord (if the system is plugged into wall outlet).
- > Disconnect all cables from off-board peripheral devices (external printer, etc.) from the console.

- To prevent damage to the power cord, DO NOT pull excessively on the cord or sharply bend the cord while wrapping it.
- Store all probes in their original cases or wrap them in soft cloth or foam to prevent damage.
- Replace gel and other essential accessories in the appropriate storage case.
- Ensure that no loose items are left on the main unit.

When Moving the System

Carry the system with handle, or put the system on the cart to move it. Use extra care when crossing door or elevator thresholds.

Always use the handle to move the system. The system weighs approx. 50 kg. For avoid possible injury or equipment damage.

- Walk slowly and carefully when moving the system.
- Do not let the system strike walls or doorframe.

Transporting the System

> Use extra care when transporting the system in a vehicle. After preparing the system as described above, take the following additional precautions:

Powering the System

- > Before transporting, place the system in its original storage case.
- > Ensure that the system is firmly secured while inside the vehicle.
- > Load the unit abroad the vehicle carefully and over its center of gravity.

> Keep the storage case still and upright. Secure that the system firmly with straps or as directed within the vehicle to prevent movement during transport.

> Any movement, coupled with the weight of the system, could cause it to break loose. Drive carefully to prevent damage from vibration.

> Avoid unpaved roads, excessive speeds, and erratic stops or starts.

3.5.3 Powering the System

3.5.3.1 Acclimation Time

After being transported, the unit requires one hour for each 2.5 ° increment if its temperature is below 10 $^{\circ}$ C or above 40 $^{\circ}$ C.

Please keep at least 20 to 30 cm spare space away from the back of the system to ensure well ventilation. Otherwise, with the increasing of the temperature inside the unit, malfunction may occur.

3.5.3.2 Connecting the electric power

To connect the system to the electrical supply:

- > Check the power voltage input labeling at rear panel of the system.
- > Ensure that the wall outlet is the appropriate type and well grounded.
- > Ensure that the system powers off.
- > Unwrap the power cord, and allow sufficient slack in the cable so that the plug will not be pulled

out of the wall outlet if the system is moved slightly.

- > Attach the power plug to the system and secure it in place by using the retaining clamp.
- > Push the power plug securely into the wall outlet.

- Only use the power cord provided by Manufacturer.
- Use caution to ensure that the power cord does not disconnect during system use.
- If the system is accidentally unplugged, data may be lost.

WARNING:

• To avoid risk of fire, the system power must be supplied from a separate, properly rated outlet.

• Under no circumstances should the AC power plug be altered, changed, or adapted to a configuration rated less than specified. Never use an extension cord or adapter plug.

• To help assure grounding reliability, connect to a "hospital grade" or "hospital only" grounded power outlet.

3.5.3.3 Power on



Turn on the green power switch (main power circuit breaker switch) at the back of the main

unit, and then press the Power button 0 on the control panel to turn on the system.

Starting Sequence

The system is initialized and start-up status is reflected on the monitor:

- Control panel initiate
- Boot up the operation system
- > Load the software
- Entering into examination status

HINTS

The starting procedure takes about approx. 100 seconds. If a problem occurs, take a picture and record the error information for service reference.

NOTE:

- After powering off the system, please wait for 1 minute to power on again.
- When the system is powered on, for safety reason, please avoid moving the main unit.

3.5.3.4 Power off

To power off the system:

Press () on control panel. When the screen shows shut down dialog box, press "shutdown" to shut down the system.

Or press 0 for 4 seconds to shut down the system directly.

After powering off the system

- Disconnect the probes. Clean or disinfect all probes as necessary. Store them in their original cases to avoid any damage.
- To ensure the system is disconnected from the power source, disconnect power plug from the wall outlet.

3.5.4 Probes

Only use the probes approved by Manufacturer.

Selecting probes

- > Choose the probe according to the different examination.
- > Begin the scanning session by choosing the correct application and preset for the examination.

Connecting the Probe

When you connect the probes, please ensure that the probe ports are not active. Freeze the image by pressing [**Freeze**] key to deactivate the probe ports.

To connect a probe:

- > Place the probe's carrying case on a stable surface and open the case.
- > Carefully remove the probe and unwrap the probe cord.

> DO NOT allow head of the probe hang freely. Impact to head of the probe could result in irreparable damage.

Inspect the probe before and after each use for damage or degradation to the housing, strain relief, lens, seal and connector. DO NOT use a probe that appears damaged until its functional and safe performance is verified. A thorough inspection should be performed during the cleaning process.


> Align the connector with the probe port and carefully push into place with the cable facing the back of the system.

> Turn the probe connector locking lever to "lock" status.

- > Carefully position the probe cord so it is free to move and is not resting on the floor.
- > When the probe is connected, the system will be automatically recognized.

• Fault conditions can result in electric shock hazard. DO NOT touch the surface of probe connector that is exposed when the probe is removed. DO NOT touch the patient when connecting or disconnecting a probe.

• Take precautions with probe cables. DO NOT bend the cable acutely.






Fig. 3-6 Probe connector "Unlock" status

Fig. 3-7 Probe connector "Lock" status

Deactivating the Probe

When deactivating the probe, the probe is automatically placed in a standby mode.

To deactivate a probe:

- > Ensure the system is in freezing mode.
- > Gently wipe off the excess gel from the probe surface.
- > Carefully slide the probe toward the probe holder, and place the probe gently in the probe

holder.

Disconnecting the Probe

Probes can be disconnected when the system is in freezing mode.

To disconnect a probe:

- > Turn the connector locking lever to the "Unlock" position.
- > Pull the probe and connector straight out of the probe port.
- > Carefully slide the probe and connector away from the probe port.
- > Ensure that the head of the probe is clean before placing the probe in its storage box.

Transporting the Probe

When transporting a probe a long distance, store it in its original carrying case.

Storing the Probe

It is recommended that all probes should be stored in the original carrying case.

Place the probe connector into the carrying case.

Carefully wind the cable into the carrying case.

Carefully place the probe head into the carrying case. DO NOT use excessive force or impact on the probe head.

3.5.5 Accessories Installation

CAUTION: Please only use the optional parts provided or suggested by manufacturer! Using other types of optional devices may cause the damage to the system and the connected optional devices.

3.5.5.1 Video printer installation

1. Put video printer stably.

2. Connect cable of video printer to video port in the back of the device. And connect the other side to video signal output port in the rear side.

3. Connect the printer line to print control port in the printer rear side, and connect the other side to the print control port in the unit rear side.

4. Connect power cable of video printer to power system.

5. Adjust printer parameter preset according to the type of printing paper.

CAUTION: Do not use any other power cable to replace 3-wire power cable manufacturer provides, otherwise there is a danger of electric shock.

Video printer sign introduction

- Video signal input port



: Video signal output port



: Print control port

O: Video printer switch

3.5.5.2 Graphic & PC printer installation

Put the printer stably; connect printer cable to USB port of the system (right side of control panel or

back of the main unit).

Connect the power cable of the printer to power system.

CAUTION: Please see packing list for fundamental configuration!

Chapter 4 Control Panel

4.1 Keyboard Appearance



Fig. 4-1 keyboard appearance



4.2 Alphanumeric Keyboard

Fig. 4-2 Alphanumeric Keyboard

The alphanumeric keys are used for inputting patient number, name, character and figure etc.


4.3 Function Keys/Knobs

lcon	Name	Function description
C	Switch	Turn on or off the system.
ญิมิ	Patient	Set up a new patient data, input name and other information.
A	Probe	Press this key for selecting probe. It can only select the connected probe.
ARROW	Arrow	Press this key to enter into arrow status. Add arrows or cross icon to the image area.
e)	Comment	Press this key to enter into annotation status, and add annotations in the image area on the screen.
	Bodymark	Press this key to enter into body mark working status, select the body mark and confirm the probe scanning position on the screen.
	Report	Produce/Save/Recall an examination report.
PRINT1 PRINT2	Print 1 & Print 2	PRINT1: print the screen image by video printer connected to the system. PRINT2: print the report by printer connected to the system (Only report page works). Or print the image in the scanning page; Or print the image in the review page.

lcon	Name	Function description
ANGLE ZOOM	Angle_Zoom	Press the knob to switch the knob function between Angle and Zoom. Enter into image zoom or adjusting angle. Rotate the knob to adjust the activated function.
DEPTH FOCUS	Depth_Focus	Press the knob to switch the knob function between Depth and Focus Position. Rotate the knob to adjust the activated function.
B	B knob	Press the knob to enter into B mode, rotate the knob to adjust the gain of B mode.
C	C knob	Press the knob to enter into CFM, rotate the knob to adjust the gain of CFM. In 3D mode, rotate this knob to spin the 3D image around X axis.
	D knob	Press the knob twice to enter into PW mode, rotate the knob to adjust the gain of PW mode.In 3D mode, rotate this knob to spin the 3D image around Y axis.
M	M knob	Press the knob once to enter into B/M mode, and press it again to enter into M mode. Rotate the knob to adjust the gain of M mode.In 3D mode, rotate this knob to spin the 3D image around Z axis.
СРА	CPA	Press this key to enter into CPA (PD) mode.
CW	CW	Press this key or [UPDATE] key to enter into CW mode. (only for model QBit 5)
TDI	TDI	Press this key to enter into TDI mode. (only for model QBit 5)
O x-contrast	X-CONTRAST	Press this key to switch the X-CONTRAST parameter.
AIO	AIO	Press this key to activate the AIO function.
FHI	FHI	Press this key to turn on or off the FHI function.
	Full screen	Press this key to enter into or quit full screen display mode.
	2B	Press this key to enter into 2B mode.
INVERT	Invert	Press this key to invert the image from up and down. Or invert the color map of CFM. Or invert the

lcon	Name	Function description
		spectrum of PW/CW.
	Cine	Press this key to save the current cine loop.
	Save	Press this key to save the current image.
FREEZE	Freeze	Press this key to freeze or unfreeze the current image.
ESC	ESC	Press this key can exit measurement, dialog, and menu.
END	END	Press this key to finish the current exam.
CLEAR	Clear	Press this key to clear all the measurement lines, body mark, and annotations.
	Biopsy	Press this key to activate biopsy.
	Single B	This key's function is the same as key B. Display B mode.
	4B	Press this key to enter into 4B mode.
	Volume	Press the volume key to adjust PW/CW volume.
SYS	SYS	The reserved key.
HELP	Help	The reserved key.
	Review	File management of system, you can view and edit the patient data.
	Setup	Press this key to enter into or quit the system setting page.
Q-BEAM	Q-BEAM	Press this key to turn on or off the Q-beam function.
Q-FLOW	Q-FLOW	Press this key to turn on or off the Q-flow function.
Ê	Change	Press this key to change the menu.
Del	Del	Press this key to delete the measurement, comment step by step.

4.3.1 STC




STC can be used for adjusting gain compensation in different image depth.

4.3.2 Parameter Control Keys/Knobs



Increase/decrease the corresponding parameters of the screen, or activate/deactivate the function.

4.4 Central Control



1.UPDATE 2.EXIT 3.CURSOR 4.CLEAR 5.MENU

Name	Function description
	This multifunctional key is work with trackball. The function switches with the unit
ENTER	status. Such as, set the cursor position, body mark position, comment position,
	toggle trackball function, selected the menu, and confirm the input.
UPDATE	This multifunctional key is work with trackball. The function switches with the unit
OFDATE	status. Such as, call the annotation and back in measuring.
EXIT	Press this key can exit measurement, dialog, and menu.
CURSOR	Press this key to show or hide the cursor.
CLEAR	Press this key to clear all the measurement lines, body mark, and annotations.
	Press the MENU knob for second time to select the item and adjust the
MENU	parameters. Press this knob for third time to exit from current item.
	Rotate the MENU knob to select the item.
DIST	Press this key to enter into distance measurement.
TRACE	Press this key to enter into trace measurement.
CALC	Press this key to enter into measurement software package.
The slab sli	Trackball is the main operation tool on screen. Position calipers in measurement,
Trackball	the function of the trackball is different under diverse working status.

6.DIST 7.TRACE 8.CALC 9.ENTER 10.Trackball

4.5 Information Area Indicating Machine Status



Left-to-Right of the up row: hard disk, cable network, USB

• Hard disk: press this icon to show the capacity of disc has been used to save data or USB flash disk in current system.

• Cable network: show the present situation of cable network; press this icon to show the IP address of current system.

• USB: show whether this system connects USB flash disk or not, press this icon to show USB safely remove interface.

Left-to-Right of the down row: input method, task queue, power

- Input method: press this icon to switch between uppercase and lowercase characters.
- Task queue: press this icon to show task and its situation. To terminate the task, delete, and so

on.

• Power: show the connecting situation of the power.

Chapter 5 Operation and Exam Mode

This chapter mainly describes the process of the normal operation of the device, including the preparation before examination, how to get the image, optimize the image, add comments, body mark and so on.

5.1 Preparing the System for Use

5.1.1 The Device Inspection

- 1. The device is placed stability.
- 2. The grid voltage AC 100-240V, 50-60Hz.
- 3. Cable is properly connected, firm and ground, the adapter is properly connected to the device.
- 4. Probe is connected and fixed.

5.1.2 Start the system

Press O to start the system, wait for the system to enter into the user interface and activate the probe into the B mode.

5.2 Choose Exam Mode

5.2.1 The Probe Identification

The system default automatically identifies the current probe type, when the probe is inserted,

Press \checkmark to switch the probe.

CAUTION: Please connect or disconnect the probe only after the system is freezing, in order to ensure stability and extend the service life of the probe.

5.2.2 Mode Selection

In probe selection interface, probe and clinical application selection page is displayed, you can choose needed probe and inspection part, and press an application preset into the B mode, start scan detection.

igta NOTE: The system has been set clinical application preset before leaving factory, each

probe has its own presets.

The detailed operation steps of the clinical application preset for the probe please refer to the exam mode part of preset section.

5.3 Patient Data Entry

Press to display the patient screen.			
CHISON	Patient Information	$\overline{\mathbf{x}}$	
ID 2018-04-29-0001 Last Name	BirthDay DD-MM-YYYY Age 0 Y 0 M Sex	Search In Archive	
Middle Name	GYN CARD URO SMP PED VAS	WorkList	
Height 0.00	cm Weight 0.00 kg BSA 0.00	m²	
Study Description Primary Indications	CPT4 Code		
Secondary Indications	CPT4 Description		
Ref.Physician Comment	▼ Diagnostician		
New Patie	ent New Study End Study OK Cancel	≅ 5- ¥ En 🖵 Ϋ	

Function Buttons on Patient screen:

Archive: Operation on the patient information which has already existed.

Worklist: Recall patient information in worklist. And need to activate the DICOM function.

New Patient: Create a new patient information identity.

New Study: Choose exam applications (OB, GYN, CARD and so on) for the current patient.

End Study: End the current patient's exam item.

OK: Save patient information.

Cancel: Cancel the operation of new patient information.

Operation Methods:

1. Move the Trackball to the position of inputting character, and then input patient information by character keyboard.

2. Use the Trackball and the **[ENTER]** key to switch between different input options: ID, patient name, doctor's name, birthday (It can be automatically calculated when input age), age (It can be automatically calculated when input birthday), gender.

3. Select the exam items, and input the regular inspection information.

4. After inputting the required information, click on the **[OK]** button to save the patient information, the system will return to the B mode.

5. Recall information of the previous patient, you can use the Archive or Worklist to recall patient information to exam.

CAUTION: Creating a diagnostic record, you should check the accuracy of the patient information before saving measurement or image; otherwise, it will be stored in the wrong patient records. After checking the patient, press to save the patient information in the system.

5.4 Image Interface Display

Logo 2, Control menu 3, Image status prompt 4, Image parameter area
 System state prompt 6, Image parameter area 7, Gray-scale strip
 8, Start point of scanning 9, Image Region 10, Cine loop

5.5 Image Mode

5.5.1 B Mode

Press **[B]** key, and display the single B Mode image, B Mode is the basic mode for two-dimensional scanning and diagnosis.

5.5.2 2B Mode

Press or key to display double B mode images side by side. One image is in real-time status; the other is in frozen status. The real-time image has start scan marker and ruler

marker .Press key in 2B mode, the original active image is frozen while the original frozen image is activated.

In frozen status, Press to choose a B mode image to be activated when unfreezes the image.

5.5.3 4B Mode

Press key to enter into 4B mode, the screen will display four B mode images side by side, but only one image is in real-time status. Pressing it again can switch the real-time status among four images.

In frozen status, Press to choose a B mode image to be activated when unfreezes the image.

5.5.4 B/M Mode

Press **[M]** key, a real time B mode image and a real-time M-mode image will be displayed at the same time. And a sample line will appear in the B mode image area, which indicates the active sample position for M image on the B image area. Move the sampling line by trackball.

5.5.5 M Mode

Press **[M]** key again, B mode image will disappear; M mode image is still active on the whole screen. M mode image stands for the tissue movement status at the sampling line. The M mode image varies with time, so it is mainly used for cardiac applications.

5.5.6 CFM Mode

CFM is a Doppler mode intended to add color-coded qualitative information concerning the relative velocity and direction of fluid motion within the B mode image.

CFM is useful to see flow in a broad area. It allows visualization of flow in the CROI, whereas Doppler mode provides spectral information in a smaller area. CFM is also used a stepping stone to Doppler mode. You can use CFM to locate flow and vessels prior to activating Doppler.

In CFM mode, move the trackball to change the position of sampling box. Activate the **[STEER]** knob and rotate the knob to adjust the angle of color sampling box (if current probe is linear probe). Press **[ENTER]** key to set the position of color sampling box. At this time adjusts the size of color sampling box through moving trackball. Press **[ENTER]** key again and move trackball to change the color sampling position again.

Press **[C]** key to enter into CFM mode; after **[C]** key light is on, rotate the **[GAIN]** knob for adjusting the gain of CFM.

CFM mode Exam Procedure:

> Follow the same procedure as described under B mode to locate the anatomical area of interest.

> After optimizing the B mode image, add Color Flow.

> Move the color region of interest CROI as close to the center of the image as possible.

> Optimize the color flow parameters so that a high frame rate can be achieved and appropriate flow velocity can be visualized.

> Press [FREEZE] key to hold the image in cine memory.

Record color flow image as necessary.

CFM Scanning Hints:

PRF: increase/decrease the PRF on the color bar. Imaging of higher velocity flow requires increased velocity scale values to avoid aliasing

Wall Filter: affect low flow sensitivity versus motion artifact

Color Map: allow you to select a specific color map. It shows the direction of the flow and highlights the higher velocity flows.

Color Gain: amplify the overall strength of echoes processed in the CROI

Persistence: affect temporal smoothing and color Doppler 'robustness'.

5.5.7 B/BC Mode

In active color mode, turn the **[B/BC]** item on to display a real B mode image at the right side of the screen and active Color mode image at the left side of the screen.

5.5.8 CPA (PD) Mode

Power Doppler Imaging (PD) is a color flow mapping technique used to map the strength of the Doppler signal coming from the flow rather than the frequency shift of the signal. Using this technique, the ultrasound system plots color flow based on the number of reflectors that are moving, regardless of their velocity. PD does not map velocity; therefore it is not subject to aliasing.

Press [CPA] key to enter into the CPA (PD) mode.

Direction PD mode

In Power Doppler (CPA) mode, press **[MENU]** knob to pop up the PD Menu. Rotate **[MENU]** knob to select the DPD Mode and press **[MENU]** knob to enter into DPD mode.

If you need go back to PD mode from DPD mode, you could press **[CPA]** key or select the PD mode item in the DPD mode.

5.5.9 PW Mode

Doppler is intended to provide measurement data concerning the velocity of moving tissues and fluids. PW Doppler lets you examine blood flow data selectively from a small region called the Sample Volume.

The X axis represents time while the Y axis represents velocity in either a forward or reverse direction.

PW Doppler is typically used for displaying the speed, direction, and spectral content of blood flow at selected anatomical sites.

PW Doppler can be combined with B mode for quick selection of the anatomical site for PW Doppler examination. The site where PW Doppler data is derived appears graphically on the B mode image (Sample Volume Gate). The Sample Volume Gate can be moved anywhere within B

mode image.

PW mode Exam Procedure:

> Get a good B mode image. Press [C] key to help locate the vessel you wish to examine.

> Press **[D]** key to display the sample volume cursor and gate.

> Position the sample volume cursor by moving the Trackball left and right. Position or re-size the sample volume gate by moving the Trackball up and down, then press **[ENTER]** key.

Press [UPDATE] key to display PW Doppler spectrum and the system will run in combined B+Doppler mode. The Doppler signal can be heard through the speakers.

> Optimize the PW Doppler spectrum as necessary.

Ensure that the sample line is parallel to the blood flow.

- Press reaction bold the trace in cine memory and stop imaging.
- Perform measurements and calculations, as necessary.
- Record results with your recording devices.

> Press (B FREEZE) to resume imaging.

- Repeat the above procedure until all relevant flow sites have been examined.
- Replace the probe in its respective holder.

When entering Duplex mode for the first time, the Doppler spectrum is not activated. The Doppler Sample Volume appears in the default position, and the B mode image or 2D (either B or Color) mode are active. Moving the Trackball will change the Sample Volume position. Press the **[ENTER]** key to toggle the Trackball function between Sample Volume Gate position and size. Press the **[UPDATE]** key after the Sample Volume Gate is defined to activate the Spectral Doppler mode. Press the **[UPDATE]** key for second time to toggle back to 2D (B or Color) update and deactivate the Spectral Doppler.

Doppler Sample Volume Gate Position and Size (Trackball and ENTER)

Move the sample volume on the B mode's Doppler cursor. The gate is positioned over a specific position within the vessel.

- > To move Doppler cursor position, turn the trackball left or right until positioned over the vessel.
- > To move sample volume gate position, move the trackball up or down until positioned inside the

vessel.

> To size sample volume gate, press **[ENTER]** key to toggle trackball function from sample volume gate positioning to sizing, then move the trackball to change sample volume gate size.

5.5.10 CW Mode

Continuous Wave Doppler allows examination of blood flow data all along the Doppler cursor rather than from any specific depth. Gather samples along the entire Doppler beam for rapid scanning of the heart. Range gated CW allows information to be gathered at higher velocities. It works with a phased array or pediatric probe.

If the velocity of the blood flow is even too high for the HPRF mode to detect, you have to try CW mode. Press **[CW]** key to enter CW mode when the probe supports CW mode.

5.5.11 TDI Mode

TDI mode is tissue Doppler mode, which is intended to provide information of low-velocity tissue motion, specifically for cardiac movement. Only phased array probe is available for TDI function.

5.5.12 Color M Mode

Color M mode is used for fetal cardiac applications. Color flow overlays color on the M mode image using velocity and variance color maps. The color flow wedge overlays the B mode image and M mode timeline. The color flow maps available in M mode are the same as in CFM mode.

Color M mode is a Doppler mode intended to add color coded qualitative information concerning the relative velocity and direction of fluid motion within the M mode image.

If the system is in color mode and the probe supports Color M mode (e.g. phased array probe), press **[M]** key to enter Color M mode.

5.5.13 2D Steer

2D Steer is available for linear probes. It can steer the beam to obtain the left or right image and enlarge the area without to rotate the probes.

In B mode, press menu knob to pop up B Menu, adjust 2D Steer item to change 2D steer angle.

5.5.14 Trapezoidal Mode

Trapezoidal image is available for linear probes. In B mode, press menu knob to pop up B Menu, turn Trapezoidal Mode menu on to enter to Trapezoidal Mode.

5.5.15 ECG

The ECG module is a device that provides the 3 lead ECG signal acquisition for cardiac application. It is not intent for the ECG diagnostic purpose as in the 12-lead module. In the cardiac application, the ECG trace is displayed on the bottom of the screen. For echo-stress, the R-wave triggering is used to gate or synchronize the image acquisition. The ECG has 3 leads: LL (left leg, RED), LA (left arm, BLACK), RA (right arm, WHITE). LA is for reference, which usually provides a bias voltage from the ECG module, and the LL, LA are the two signals from the body and going to the differential input of the ECG isolation amplifier.

5.5.16 HD 3D

Before using the HD 3D function, connect the 3D module and the ultrasonic device with a network cable.

Operation Procedure:

a. In B mode, press the **[Freeze]** key to enter the frozen status. Rotate the corresponding parameter knob of **[3D ROI]** to turn on the 3D ROI. System will display the ROI box.

b. Press the **[ENTER]** key, the ROI box changes to the dotted line and move the trackball to adjust the size of the ROI box. Then press **[ENTER]** key again to confirm the size.

c. Move the trackball to drag the ROI box to the target area, and press the **[UPDATE]** key to confirm the location.

d. Then move the trackball and press the **[ENTER]** key to choose the start pos and end pos of the designated image.

e. Press the parameter knob of **[HD 3D]** to enter the HD 3D mode. The system will display the rebuilt 3D image. Rotate the **[C]**, **[D]** or **[M]** key on the control panel, the 3D image will spin around X axis, Y axis and Z axis. You can press the **[ANGLE_ZOOM]** knob to reset.

f. Rotate the parameter knob of **[Vol Spacing]** to set the physical distance of each frame of image movement.

g. After finishing these operation, you can press **[Freeze]** key or **[B]** knob to go back to the real B mode. Or press **[Exit]** key to go back to the frozen B mode.

Tips: User also can operate 3D mode directly according to the step d to g, without adjusting ROI box.

5.5.17 Biopsy and Super Needle

1. How to enter into Biopsy

Activate the [Biopsy] and press the [MENU] knob to show or hide biopsy line.

2. How to adjust the biopsy

After the biopsy line shows, press the **[ENTER]** key to activate the adjustment function of biopsy line, horizontal rolling the trackball can translate the biopsy line, vertical rolling the trackball can adjust the line angle, press the **[UPDATE]** key to set the default biopsy line position.

3. Super Needle

Super needle is used for enhance the needle image in the B mode image. After turning on the super needle, super needle and needle angle function will be active and user can adjust the needle angle to optimize the image for needle only (The angle is 5 degree per step.).

5.6 B Image Menu & Parameters

Menu & Parameters	Function description
Keyboard Keys	
Gain	In real status, rotate [B] knob to adjust the Gain.
STC	STC curves can be used for adjusting gain compensation in different image depth.Drag the slide of STC to adjust the value.STC curve will disappear automatically 1 second later after stopping adjustment.
Depth	Press [DEPTH_FOCUS] knob until the indicator of [DEPTH] is lit, and then rotate the knob to change the depth of image.
Focus	Press [DEPTH_FOCUS] knob until the indicator of [FOCUS] is lit, and then rotate the knob to change the focus position.
Zoom	Press [ANGLE_ZOOM] knob to appear the zoom box, press [UPDATE] key and move trackball or rotate the [ANGLE_ZOOM] knob to change the zoom ratio.

Menu &		
Parameters	Function description	
	B mode image and B/M mode image can be reversed horizontally and	
	vertically.	
	Press, the displayed image is reversed in the up-down direction.	
	The horizontal flip status indicators of the upper-left corner of the image	
Invert	window have the following meanings:	
	The meaning of the symbol " \bigcirc " indicating the probe initiative scanning	
	position " \bigcirc " situated in the left indicates that the first scanning line in the left of	
	the screen is corresponding to the initiative scanning position of the probe, " \bigcirc "	
	situated in the right indicates that the first scanning line in the right of the	
	screen is corresponding to the initiative scanning position of the probe.	
FHI	Press To turn on or off FHI function.	
AIO	AIO	
	Press V to activate the AIO function.	
	Full screen the image area. Press to activate the function;	
Full Screen Show	Press, [EXIT] or again to exit full screen show.	
	When full screen shows, press [MENU] knob to show the menu of current	
	mode (except bodymark, annotation, measurement).	
Bottom Menu		
Freq	In real status, rotate the corresponding knob of [Freq] to adjust the frequency.	
	Dynamic range is used for adjusting the contrast resolution of B mode image	
Dynamic	and mode image, compressing or enlarging the display range of gray scale.	
	At the real-time status, rotate the corresponding knob of [Dynamic] to adjust	
	dynamic.	
	In B mode, 9 focus points can be selected simultaneously, and the number	
Focus Num	controlled by the depth, SRA and Compound.	
	Rotate the corresponding knob of [Focus Num] to adjust.	
Focus Pos.	Rotate the corresponding knob to change Focus Position.	
	In real status, rotate the corresponding knob of [Compound] to turn it on or	
Compound	off.	
	The SRA can't be edited after opening the compound.	
SRA	In real status, rotate the corresponding knob of [SRA] to turn it on or off.	
Q-Image	In real status, rotate the corresponding knob of [Q-Image] to adjust.	
Persistence	In real status, adjust the contrast and resolution.	
	In real status, rotate the corresponding knob of [Persistence] to adjust.	

Menu &	Function description
Parameters	
	Scan Line Density function is only valid for the image in B mode, 2B mode,
	B/M mode or 4B mode image. The line density has two types: high density and
Density	low density. High density means better image quality while low density image
	has higher frame rate.
	Rotate the corresponding knob of [Density] adjust the line density.
	Adjust image gray scale inhibition parameters.
B Rejection	Activate [B Rejection], and rotate [MENU] knob to adjust B Rejection
	parameters.
Left Menu	
Scan Width	Rotate the corresponding knob of [Scan Width] to adjust the scan width of the
	corresponding size.
	Smoothness function is used for restraining the image noise and performing
Smooth	axial smooth processing to make the image smoother.
	Rotate the corresponding knob of [Smooth] to adjust.
	Edge enhancement is used for enhancing the image outline. In this way the
Edge Enhance	user can view the tissue structure more clearly.
	Rotate the corresponding knob of [Edge Enhance] to adjust.
A Power	Acoustic power means the acoustic power transmitting from the probe.
	Rotate the corresponding knob of [A power] to adjust.
L/R Flip	Activate the [L/R Flip] and press the [MENU] knob to invert the image from left
_/p	and right.
	Menus of Advance
Zoom Coef	Adjust the size of ruler & image.
	Rotate the corresponding knob of [Zoom Coef] to adjust.
Trapezoidal Mode	Activate the [Trapezoidal Mode] and press the [MENU] knob to turn on or off
	the trapezoidal function.
2D Steer	Activate the [2D Steer] and press the [MENU] knob to turn on or off the 2D
	Steer function.
ECG	Activate the [ECG] and press the [MENU] knob to turn on or off the ECG
	function.
UD Invert	Activate the [UD Invert] and press the [MENU] knob to turn on/off the ECG
	invert.
ECG GAIN	Activate the [ECG GAIN] and rotate the [MENU] knob to increase or decrease
	the ECG gain.
ECG POS	Activate the [ECG POS] and rotate the [MENU] knob to set the ECG trace
	position.

Menu &	Eurotian depariation	
Parameters	Function description	
ECG Velocity	Activate the [ECG Velocity] and rotate the [MENU] knob to set the ECG	
	Velocity.	
ECG Dyn	Activate the [ECG Dyn] and rotate the [MENU] knob to set the ECG dynamic.	
	Show or hide biopsy.	
	Activate the [Biopsy] and press the [MENU] knob to show or hide biopsy line.	
Biopsy	After the biopsy line shows, press the [ENTER] key to activate the adjustment	
ыорзу	function of biopsy line, horizontal rolling the trackball can translate the biopsy	
	line, vertical rolling the trackball can adjust the line angle, press the [UPDATE]	
	key to set the default biopsy line position.	
Super Needle	Activate the [Super Needle] and press the [MENU] knob to turn on or off the	
Super Needle	Super Needle function.	
Needle Angle	Activate the [Needle Angle] and rotate the [MENU] knob to set the needle	
	angle.	
Center Line	Show or hide Center Line.	
	Press the [MENU] knob to show or hide center line.	

5.6.1 Utility Menu

This function includes post processing, slide show and other items. Press **[MENU]** knob, then choose Utility selection, it will appear Utility options.

Menu & Parameters	Function description	
	Chroma:Adjust the type of the chroma.	
De et Des es es	2D Map: Select the type of the scale curve.	
Post Process	Gamma:Adjust image gray value parameters.	
	B Rejection: Adjust image gray scale inhibition parameters	
Slide Show	Activate the [Slide Show] to enter into slide show mode.	

5.7 M Image Menu & Parameters

Menu & Parameters	Function description
Keyboard Keys	
	Refer to B image adjustment for keyboard keys.
Bottom Menu	
Speed	Rotate the corresponding knob of [Speed] to adjust M speed.
M Chroma	Rotate the corresponding knob of [M Chroma] to select the Chroma type.
Layout	Rotate the corresponding knob of [Layout] to select the layout of B/M image.
Left Menu	

Menu & Parameters	Function description
M 2D Map	Rotate the corresponding knob of [M 2D Map] knob to select the scale curve
	type.

5.8 CFM/CPA/DPD/TDI Image Menu & Parameters

Menu &	Function description					
Parameters						
Keyboard Keys						
Refer to B image adjustment for keyboard keys.						
Bottom Menu						
Freq	In real status, rotate the corresponding knob of [Freq] to adjust the frequency.					
Steer	Press [ANGLE_ZOOM] selection knob until the indicator of [ANGLE] is lit, and					
	then rotate the knob to change the steer of linear image ROI. Or rotate the					
	corresponding knob of [Steer] to change the steer of linear image ROI.					
Wall Filter	Rotate the corresponding knob of [Wall Filter] to adjust wall filter.					
Color Map	Rotate the corresponding knob of [Color Map] knob to select the color type.					
	(CFM mode only)					
O flow						
Q-flow	Press to turn on or off the Q-flow function.					
PRF	Rotate the corresponding knob of [PRF] to Adjust PRF, adjustment range					
	depends on probe.					
Color Invert	Press the corresponding knob of [Color Invert] to invert the color of flow.					
Persistence	Rotate the corresponding knob of [Persistence] to adjust.					
Q-beam	Press to turn on or off the Q-beam function.					
Base Line	Adjust the position of baseline.					
	Rotate the corresponding knob of [Baseline] to change position of baseline.					
Left Menu						
CF Mode	Press [MENU] knob on [CF Mode] item to adjust CF mode, velocity or					
	variance.					
Wall Thre.	Activate [Wall Thre.], and rotate [MENU] knob to adjust wall thre.					
Blood Effection	Activate [Blood Effection], and rotate [MENU] knob to adjust blood effection,					
	smooth or resolution.					
Density	Press the corresponding knob of [Density] to adjust the line density, high or					
	low.					
B/BC	Press [MENU] knob on [B/BC] item to enter into or quit the B/BC mode.					

Menu &					
Parameters	Function description				
Keyboard Keys					
Refer to B image adjustment for keyboard keys.					
Steer	When the indicator of [ANGLE] is lit, rotate [ANGLE_ZOOM] knob to change				
	the steer of linear image ROI. Or rotate the corresponding knob of [Steer] to				
	change the steer of linear image ROI.				
Bottom Menu					
Freq	In real status, rotate the corresponding knob of [Freq] to adjust the frequency.				
PW/CWD Chroma	Rotate the corresponding knob of [Chroma] to select the Chroma type.				
Speed	Rotate the corresponding knob of [Speed] to adjust PW speed.				
Audio	Adjust audio volume. Rotate the corresponding knob of [Audio] or press & & & to adjust audio volume, adjustment range from 0~100%.				
Triplex Mode	Press the corresponding knob of [Triplex Mode] to turn on or off triplex mode.				
PRF	Rotate the corresponding knob of [PRF] to Adjust PRF, adjustment range depends on probe.				
Base Line	Adjust the position of baseline. Rotate the corresponding knob of [Baseline] to change position of baseline.				
Left Menu					
D/CW 2D Map	Activate [2D Map] and rotate [MENU] knob to select the 2D map.				
Spectrum Enhance	Activate [Spectrum Enhance] and rotate [MENU] knob to adjust spectrum enhance.				
Dynamic Range	Activate [Dynamic Range] and rotate [MENU] knob to adjust dynamic range.				


5.9 PW/CW Image Menu & Parameters

5.10 Image and Cine Disposition

5.10.1 The Principle of Cine Storage

In real image status, the image can be stored in the movie memory in chronological order, maximum frames can be set. The maximum number of frames of the film storage can be set, please refer to preset chapter.

If the movie memory is full, the recent frame saved into memory, the previous frame removed from memory.



Cine loop indicate diagram

5.10.2 Manual Loop

Press to freeze image, pop cine playback bar, at this time, move trackball to play by hand.

5.10.3 Automatic Loop

After freezing image, press [Play/Pause] to play, press it again to stop.

Press [ENTER] key to choose the area of automatic playback area needed.

5.10.4 Loop Range Set

Rotate the corresponding knob of **[Start Pos]** and **[End Pos]** to set the loop start position and end position, the loop range settled.

Rotate the corresponding knob of [Reset Ranges] to reset the loop range to maximum.

5.10.5 Save and Recall Image

Press to save current image, the image will be displayed below the screen. If you need to recall images that have been stored, move cursor to needed image, press **[ENTER]** key to recall it; or you can recall archived patient's information to recall image, please refer to archive chapter.

5.10.6 Save and Recall Cine

In freeze status, press (ENTER) to save cine, then it will be displayed below the screen, move cursor to needed cine, press [ENTER] key to recall cine.

5.10.7 Send image

After recalling the images, press the corresponding knob of [Send] to send images to USB flash

disk, net storage, DICOM storage and print.

Hint1: Activate the DICOM before DICOM storage and print.

Hint2: Click the "Cine transform to AVI" to save the cine to AVI.


ExportInfo						
UDisk	Drive	U disk C				
DicomPrint) JPG 🛛 DCM 📝 Cine transform t			
DicomStorage	Image Type	ВМР О ЈРС	🔘 DCM 🔛	Cine transform to	D AVI	
Network Storage	Folder Name	2016-12-13-0002_Kim				
		Export		Close		

Sending Interface

5.11 Edit Comment

5.11.1 Overview

The comment is to enter text or symbols on the image.



- 1. Press or , then system will go into the comment process; or input characters with alphanumeric keyboard, the system will enter into comments status too.
- 2. Move the cursor to the position where need to comments.
- 3. Input characters at cursor position by keyboard then press [ENTER] key to confirm.
- 4. Press or again or [EXIT], to exit

5.11.3 Input Comment Library Characters

- 1. In comment status, move trackball to image area to edit.
- 2. Rotate the corresponding knob of **[Font size]** to adjust front size of comments, the range is 10~20.
- 3. Rotate [MENU] knob to select needed comments, then press [MENU] knob to exit.

5.11.4 Edit Quick Comments

- 1. Rotate the corresponding knob of [Edit] to pop quick comment edit box.
- 2. Input customized comments.
- 3. Press [Done] to finish edit, press [X] to cancel edit.

5.11.5 Input Quick Comments

- 1. Press [Text] to select needed quick comments.
- 2. Adjust front size of comments.
- 3. Press the corresponding knob of [Input] to place comments in the image area.

5.11.6 Move Comments

- 1. In comment status, move trackball to the comment, press [ENTER] key to activate it.
- 2. Move trackball to place the comment to target area.
- 3. Press [ENTER] key again to confirm the comment.

5.11.7 Edit Comments

- 1. In comment status, move trackball to the comment, press [ENTER] key to activate it.
- 2. Press [BackSpace] key to delete unnecessary characters.

3. Press [ENTER] key to confirm.

5.11.8 Delete Comments

5.11.8.1 Delete Characters

In comment status, activate the comment that need to be deleted, then press the **[ENTER]** key, it will display "|" on the screen, press **[BackSpace]** key to delete the character.

5.11.8.2 Delete Single Comment

Activate the comment that needs to be deleted, press to delete comment.

5.11.8.3 Delete All Contents of the Comment

Don't activate the single comment, and then press $\underbrace{\square}_{\square EAR}$ to delete all characters that has input. **CAUTION:** Press $\underbrace{\square}_{\square EAR}$ will delete the measurement and body mark at the same time.

5.11.9 Set the Position of Default Comment

Operation:

1. Press the corresponding knob of **[Save Home Pos.]** to move cursor to the initial position.

2.Press the corresponding knob of [Load Home Pos.] to set the initial position .

5.12 Set Body Mark

5.12.1 General Description

The body mark indicates patient's examination position and the direction of probe scan on the image.

Body marks are divided into: obstetric, abdomen, gynecology, heart and small parts, each has different body mark. Each type of body mark automatically is corresponding to current examination mode.

5.12.2 Body Mark Operation

Operation:

ESC

1. Press or to enter into the body status.

2. Select the body mark that you need.

3. Move trackball after adding the body mark image and then adjust the position of probe. Rotate

[MENU] or **[ANGLE]** knob can adjust the probe direction. Press **[ENTER]** key to confirm when adjustment was finished.

4. Move the Trackball to change the position of the body mark.

5. If you want to exit from the body mark function, press \checkmark

🔟 or 🕅 again or [EXIT].

to exit; and the body mark is fastened to the screen.

6. Press (or (bel to delete the body mark.

5.13 Set the Direction of Arrow

Operation:

- 1. Press or know to display arrow.
- 2. Press to select cursor type: arrow or cross.
- 3. Rotate the corresponding knob of [Size] to change the size of the cursor.
- 4. Rotate [MENU] or [ANGLE] knob can adjust the arrow direction.
- 5. Press **[ENTER]** key to confirm when adjustment was finished.
- 6. Press [**EXIT**] key or to exit the arrow settings.
- 7. Press to delete the arrows one by one.
- 8. Press to clear all the arrows inputted.

5.14 Image Browse

Press **[ENTER]** key on an item to activate the item's function.



- ID: ID of Current patient.
- Name: Current patient's name.
- Information: Enter into current patient's information interface.
- Report: Enter into current patient's report interface.
- Send images: Send image to USB disk , DICOM storage and print.
- Print images: print the image which be chosen, it will be printed as the arrangement set.
- Delete images: Delete selected image.
- Row*Column: Select image's format.
- Pre page: Page up.
- Next page: Page back.
- Select All: Select all the images of this patient.
- **Deselect All:** Deselect all the images of this patient.
- New Exam: Exit current examination and open a new dialog box.
- Continue Exam: Exit image browsing interface and go on checking current patient.
- Archive: Open up archive management interface.
- Exit: Turn off image browsing interface.

5.15 Archive Management

Archive management can search for patient's information which has been stored in system. Press

C	HIS	01	N			Arc	nive					×
Item	Name	V	Keyword	kim		Period	Six Mont	1s 🔻 [Multiple Choice	Data Source	Local	•
ID				Name		Birth	Date	Sex	StoreTime			
= 20	16-09-12	-0591		Kim		198	4-08-15	Male	2016-09-12T03:3	6:36		nfo
												eport
												xam
												xam
											Send Ex	
		Patie	nt View								Delete I	
												▲ 0/0 ▼
			Ne	w Exam	Continue Ex	am	Review		Ex	it.	in 📮	ې مې

[Archive] to archive management interface, all process can be opened up by moving cursor.

- Item: Type selection, select Patient's ID or Name.
- Keyword: Search for key words.
- **Period:** Time filter, select today, one week, one month, three months, six months, recent one year and all.
- Multiple Choice: Multiple choice.
- Data Source: Path choice, select hard disk or U disk.
- Patient info: Enter into patient's information interface.
- **Review Report:** Enter into report interface.
- **Backup Exam:** Select examination information to USB hard disk.
- **Restore Exam:** Recover examination information from USB hard disk.
- Send Exam: Send selected examination information remotely to USB hard disk or DICOM

Storage/Print (Need to activate the DICOM).

- **Delete Exam:** Delete selected examination information.
- Patient View: Change display mode of information.

- Expand All: Select Patient View, it will display sub-directory.
- Collapse All: Exit sub-directory.
- Select All: Select all examination information.
- New Exam: Exit current patient's examination.
- Continue Exam: Exit archive management interface and go on checking current patient.
- Easy View: Exit archive management interface and open up image browsing interface.
- Exit: Exit archive management interface and go on checking current patient.

5.16 Report

Move cursor to the images and press **[ENTER]** key to add the image into the report page. The report can be saved and printed. It is convenient for the doctor to view and edit the patient's information.

Reports contain normal report, abdominal report, cardiac report, small part report etc. Move the cursor to the required report page and press **[ENTER]** key to choose.

Press and the system pop up the report page of the current exam mode. Change the other exam mode report by the drop-down box.

CHISON		Report		8
Report Title Cardiol	ogy Report	•		
	hospita	1		
	Cardiology I	Report		
Patient Name:	Age:0	Sex:		2
Patient ID:2018	-04-29-0001			
Diagnostic				
				3
				A AND AND A
				4
Ultrasound Ir	nage:			
7000-75.45-110 - 5	Toppo 9000 6311 Johnson Ganta 2000 564 65213 A 8	7002154-000 10000 10000 0000 00000 0000 Grade 20054-000 00000		
	A A A A A A A A A A A A A A A A A A A			
ľ				
				1/1
				t All Deselect All
-11		-14	▲ Select	CAIL Deselect All
Description			▼ Ing	out Template
] ‰ ¥
	Clear All Print	Export Save	Exit	<u> </u>
				n 📙 😗

• Report Title: Report options, different kinds of report can choose, such as Normal, OB/GYN

etc.

- Hosp: Display the hospital name.
- Abdomen Report: Display the kind of report.
- Patient Name: Display the patient name.
- Age: Display the patient age.
- Sex: Display the patient sex.
- Patient ID: Display the patient ID
- **Diagnostic:** Input the diagnostic instructions.
- **Description:** Input the description of symptom.
- Tips: Input note information.
- Clear All: Clear all the data including the selected images the measurement result and patient's information etc.
- **Print:** Print the report with image.
- **Export:** Export the PDF report to the U disk.
- **Save:** Save the report in system.
- The image on the right side: Press [ENTER] key on the image to add the image into the report.
- Select All: Select all the images of the image on the right side.
- **Deselect All:** Cancel to select all the images of the image on the right side.

5.17 DICOM

5.17.1 DICOM Worklist

Press [Worklist] button in Patient Interface, pop up the following dialog box

			Work	List				×
Name		Se	arch By	Name		-		
ID	/ Name	BirthDate	Sex	Docto	r	Acc#		
							Update	

- Name: input ID or some characters , fuzzy query needs server.
- Search By: select term , ID or name.

ID: display the ID of patients.

Name: display the names of patients.

- **BirthDate:** display the birth date of patients.
- Sex: display sex of patients.
- **Doctor:** display names of doctors.
- Acc#: display the accession number of patients.
- Update: press this button to do search operation.

• Apply: select the searched patient and press this button , input all patient information into the new patient interface.

• Clear: clear all searched content.

5.17.2 DICOM Storage

Check the "Send while saving" in setting, then DICOM storage when saving cine and images.

Press send button in archive or freeze interface. DICOM Send interface as follows:

	Detail >>
ort Close	

Select DICOM Storage in left, choose DICOM server and press Export button to DICOM storage. Enter Task Queue and watch or edit DICOM process.

5.17.3 DICOM Print

DICOM Print operation is the same as DICOM storage.

5.17.4 **DICOM SR**

Press [Send DICOM SR] button in report interface, this task is added into Task Queue.

Chapter 6 Measurement and Calculation

Main content of this chapter:

Normal calculation and measurement on B mode image and M mode image, OB calculation and Urology measurement etc., system can enter into corresponding measurement mode depend on current exam mode, and enter into the corresponding report depend on the measurement mode. System has built-in the default measurement according to the exam mode, the change of measurement please refers to the chapter of preset settings.

CAUTION: Please select the most appropriate ultrasound images, measurement tools and measurement methods for measurements according to your diagnosis needs. The final measurement results must be determined and verified by a physician. Measurement accuracies are affected by many non-technical factors, for example operator's experience, patient's status. Please do not only use the ultrasound measurement results as the sole basis for diagnosis, please always use other clinical information to do integrated diagnostics.

6.1 Keyboard for Measurement

6.1.1 Trackball

Trackball is used to move the cursor, main functions are as follows:

1. Before starting a measurement, use the trackball to choose the menu options.

2. After starting a measurement, move the trackball to move the cursor, during the measurement, the cursor should not be moved out image area.

3. During the Ellipse method measurement, use trackball to change the length of short axis.

4. Update the moving of the measurement result, move the trackball to change the position of the measurement result.

6.1.2 [ENTER] key

During the measurement, the functions of **[ENTER]** key are as follows:

- 1. When cursor is on the menu, press the key to choose the options and start the measurement.
- 2. During the measurement, press the key to anchor the start point and end point.

6.1.3 [UPDATE] key

1. Before the measurement, press **[UPDATE]** key to change the measurement method, such as ellipse, trace. The changeable measurement item has "< >".

2. During the measurement, **[UPDATE]** key is used to switch the start point and end point, long axis and short axis when the measurement is not finished.

3. During the distance measurement, press **[ENTER]** key to fix the start point, when the end point is not fixed, press **[UPDATE]** key to switch the start point and end point.

4. During the Ellipse measurement, when fix the long axis, but the short axis is not fixed, press

[UPDATE] key to switch the long axis and short one.

6.1.4 [CLEAR] key

6.1.5 [Del] key

6.1.6 [CHANGE] key

Press **[EXIT**] key or **ESC** to exit.

6.1.7 [Exit] key

Press [**EXIT**] key or to exit the measurement menu.

6.1.8 Parameters control key

Press the corresponding key to update the function and use the function.

6.2 B Mode General Measurement Methods

The system B mode contains Distance, Ellipse, and Trace.

6.2.1 Meas. Distance

Measurement steps:

1. Press the **[CALC]** key to enter into measurement. Update the **[distance]** item in the menu or press the quick measure key **[DIST]**, it will display a segment "+"icon.

2. Move the "+"icon by trackball to fit the one point of the line. Press **[ENTER]** key to fix the start point and the cursor can be moved to the next position.

3. Press [UPDATE] key can change the activated point, and fit the other point of the line.

4. Move the cursor to the end-point, press **[ENTER]** key again to complete the measurement.

5. After the measurement, the result will display in the measurement results area.

6. Repeat the steps from 1 to 4 to start next "distance" measurement.

Each group of measurement is limited, if the measurement results beyond, it will begin a new group of measurement automatically.

6.2.2 Ellipse

Measurement steps:

1. Press the **[CALC]** key to enter into measurement. Update the **[Ellipse]** item on the menu; it will display a segment "+"icon.

2. Move the "+"icon by cursor, Press **[ENTER]** key to fix the point and the cursor can be moved to form a round.

3. Press [UPDATE] key can exchange the activated point and the fixed point.

4. Move the cursor to the end-point of the ellipse, press [ENTER] key to fix the axis, at the same

time, the next axis be updated, and can change the size of the axis by the cursor.

5. Now press [UPDATE] key can exit to the step 4.

6. After fixing the next axis, can press [ENTER] key to complete the measurement.

7. After the measurement, the result will display in the measurement results area.

8. Repeat the steps from 1 to 6 to start next "ellipse" measurement.

Each group of measurement is limited, if the measurement results beyond, it will begin a new group of measurement automatically.

6.2.3 Trace

Measurement steps:

1. Press the **[CALC]** key to enter into measurement. Update the **[Trace]** item in the menu or press the quick measure key **[TRACE]**, it will display a segment "+"icon.

2. Move the "+"icon by cursor, Press **[ENTER]** key to fix the point and the cursor can be moved to the next position.

3. Make the cursor tracing along the edge of required area, the traced line can be not closed.

4. Now press [UPDATE] key to cancel the tracing.

5. Press **[ENTER]** key again in the endpoint, the start point and end point of trace line will be closed by a straight line.

6. After the measurement, the result will display in the measurement results area.

7. Repeat the steps from 1 to 5 to start next "trace" measurement.

<u> MNOTE:</u>

Each group of measurement is limited, if the measurement results beyond, it will begin a new group of measurement automatically.

6.2.4 Histogram

Histogram is used to calculate the gray distribution of the ultrasound echo signals within a specified area. Use the rectangle, ellipse or trace method to draw along the desired measurement area. The result is shown in the form of histogram.

Histogram can be measured only on the frozen image.

Measurement steps by rectangular method:

1. Press to freeze the image.

- 2. Press [ENTER] key on [Histogram] menu to enter into measurement status.
- 3. Press [ENTER] key to fix one apex of the rectangle.
- 4. Move the trackball to change the cursor position and fix the diagonal point of the rectangle

5. Move the trackball to change the cursor position, fix the diagonal point of the rectangle, and press **[ENTER]** key again to confirm the measurement area. The result will display on the measurement result area.

◆ Measure the histogram by ellipse or trace method: The method is the same as that to measure ellipse or trace method, press **[UPDATE]** key to change the measurement between ellipse and trace.

The horizontal axis represents the gray scale of the image ranging from 0 to 255.

The vertical axis represents the distribution ratio of each gray scale. The value shown on the top of vertical axis represents the percentage of the maximally distributed gray in the whole gray distribution.

6.2.5 Profile

Profile is used to measure the gray distribution of the ultrasound signals in the vertical or horizontal direction on a certain profile (section).

This measurement is only available in the frozen mode.

Measurement steps:

- 1. Press (FREEZE) to freeze the image.
- 2. Press [CALC] key, and choose [B General Meas.].
- 3. Draw a straight line at the measuring position. The method is the same as that to measure distance.
- 4. The calculated result of the profile will be displayed at the center of the screen.

1-The horizontal (or vertical) axis represents the projection of the profile line on the horizontal direction.

2-The vertical (or horizontal) axis represents the gray distribution of the corresponding points on the profile line. The range is 0 to 255.

6.2.6 Volume of A4CTrace method

A4CTrace tool is used to measure the cardiology item, Single Plane, Simpson BP etc.. These measurements are used to calculate the volumes.

Measurement steps:

1. Press

to freeze the image.

2. Press **[CALC]** key, and choose **[Cardiology]**. Activate an A4CTrace measurement such as Single Plane-EDV.

3. Move the cursor to a point and press [ENTER] key to set the start point, then move the cursor

to next point and press [ENTER] key until the last point, and then set the straight line across the whole area.

4. The Volume of the chosen area will display on the measurement result area.

6.3 B Fast Measurement

Press **[DIST]** key to enter B Fast measurement in B mode. Press the corresponding parameter control key to switch the fast measurement item.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
	Distance	cm	Refer to Distance Meas.	
Distance	Ratio(Distance)		Refer to Distance Meas. Formula: R=D1/D2	D1: First Distance D2: Second Distance
	Angle	deg	Refer to Distance Meas.	Angle Range: 0°~ 180°
Area	Perimeter/Area	Area:cm ² Perimeter: cm	Refer to Ellipse and trace meas.	Ellipse and trace.
Alea	Ratio(Area)		Refer to ellipse Meas. Formula: R=A1/A2	A1: First Area A2: Second Area
	Volume(1Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D ³	D means: Depth
Volume	Volume(1Ellipse)	ml	Refer to ellipse meas. Formula: V=(π/6)*A*B ²	A: Long Axis B: Short Axis
volume	Volume(2Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D1*D2 ²	D1: the longer distance D2: the shorter distance
	Volume(3Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D1*D2*D3	D1, D2, D3: Distance

Меа	as. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
		Volume(1Dis1Ellip)	ml	Refer to distance and ellipse Meas. Formula: V=(π/6)*A*B*M	A: Long Axis B: Short Axis M: Distance

6.4 B General Measurement

Press [B], [2B] or [4B] to enter into B, 2B or 4B mode, then click [CALC] key to enter into

measurement status. Or press to choose the General measurement.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
	Distance	cm	Refer to Distance Meas.	
	Perimeter/Area	Area:cm ² Perimeter: cm	Refer to Ellipse and trace meas.	Ellipse and trace. Press [UPDATE] key to change.
	Volume(1Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D ³	D means: Depth
General	Volume(1Ellipse)	ml	Refer to ellipse meas. Formula: V=(π/6)*A*B ²	A: Long Axis B: Short Axis
Measurement	Volume(2Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D1*D2 ²	D1: the longer distance D2: the shorter distance
	Volume(3Distance)	ml	Refer to Distance Meas. Formula: V=(π/6)*D1*D2*D3	D1, D2, D3: Distance
	Volume(1Dis1Ellip)	ml	Refer to distance and ellipse Meas. Formula: V=(π/6)*A*B*M	A: Long Axis B: Short Axis M: Distance
Ratio	Ratio(distance)		Refer to Distance Meas. Formula: R=D1/D2	D1: First Distance D2: Second Distance
	Ratio(area)		Refer to ellipse Meas. Formula: R=A1/A2	A1: First Area A2: Second Area

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Angle		deg	Refer to Distance Meas.	Angle Range: 0°~ 180°
Histogram			Refer to histogram	
Profile			Refer to Profile	

6.5 B Abdomen Measurement

Choose Abdomen exam mode. Freeze the required image, then Press [CALC] key to enter into

ABD measurement status. Or press to choose the Abdomen measurement.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Distance		cm	Refer to Distance Meas.	
CBD		cm	Refer to Distance Meas.	
GB wall		cm	Refer to Distance Meas.	
Liver Length		cm	Refer to Distance Meas.	
	Height	cm	Refer to Distance Meas.	
	Width	cm	Refer to Distance Meas.	
	StD%	%	Refer to Distance Meas. Formula: ((D1-D2)/D1)*100 %	D1: Length of Normal D2: Length of Stenosis
Prox Aorta	StA%	%	Refer to ellipse Meas. Formula: ((A1-A2)/A1)*100 %	A1: Area of Normal A2: Area of Stenosis
	Vessel Area	cm ²	Refer to Ellipse and Trace Meas.	Ellipse and trace. Press [UPDATE] key to change.
	Vessel Distance	cm	Refer to Distance Meas.	
Mid Aorta	The same as above	The same as above	The same as above	The same as above
Distal Aorta	The same as above	The same as above	The same as above	The same as above

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Spleen		cm cm cm	Refer to Distance Meas. Formula: V=(π/6)*L*H*W	L: Length H: Height W: Width
Renal Volume (Right/Left)		cm	Refer to Distance Meas. Formula: V=(π/6)*L*H*W	L: Length H: Height W: Width
	Height	cm	Refer to Distance Meas.	
	Width	cm	Refer to Distance Meas.	
	StD%	%	Refer to Distance Meas. Formula: ((D1-D2)/D1)*100 %	D1: Length of Normal D2: Length of Stenosis
Lliac (Right/Left)	StA%	%	Refer to ellipse Meas. Formula: ((A1-A2)/A1)*100 %	A1: Area of Normal A2: Area of Stenosis
	Vessel Area	cm ²	Refer to Ellipse and Trace Meas.	Ellipse and trace. Press [UPDATE] key to change.
	Vessel Distance	cm	Refer to Distance Meas.	

6.6 B OB Measurement

Choose OB exam mode. Freeze the required image, then Press [CALC] key to enter into OB

measurement status. Or press to choose the OB measurement.								
Meas. Menu	Submenu	Unit	Meas. Method	Comment				
Distance		cm	Refer to Distance					
Distance		cm	Meas.					
				Formula to choose: CFEF,				
GS			Refer to Distance	Campbell, Hadlock,				
63		cm	Meas.	Hansmann, Korean, Merz,				
				Shinozuka				

to choose the OB measurement.

Meas. Menu	Submenu	Unit	Meas. Method	Comment
CRL		cm	Refer to Distance Meas.	Formula to choose: Hadlock, Hansmann, Korean, Nelson, Osaka, Rempen, Robinson, Shinozuka
BPD		cm	Refer to Distance Meas.	Formula to choose: Bessis, CFEF, Campbell, Chitty, Hadlock, Hansmann, Jeanty, Johnsen, Korean, Kurtz, Merz, Osaka, Rempen, Sabbagha, Shinozuka
нс		cm	Refer to Ellipse and Trace Meas.	Formula to choose: CFEF, Campbell, Chitty, Hadlock, Hansmann, Johnsen, Korean, Merz
AC		cm	Refer to Ellipse and Trace Meas.	Formula to choose: CFEF, Campbell, Hadlock, Hansmann, Korean, Merz, Shinozuka
FL		cm	Refer to Distance Meas.	Formula to choose: Bessis, CFEF, Campbell, Chitty, Doubilet, Hadlock, Hansmann, Hohler, Jeanty, Johnsen, Korean, Merz, Osaka, Shinozuka
	YS	cm	Refer to Distance Meas.	
	OFD	cm	Refer to Distance Meas.	Formula to choose: Hansmann, Korean
	APD	cm	Refer to Distance Meas.	Formula: Bessis
	TAD	cm	Refer to Distance Meas.	Formula: CFEF
Fetal Biometry	FTA	cm	Refer to Distance Meas.	Formula: Osaka
	SL	cm	Refer to Distance Meas.	
	APTD	cm	Refer to Distance Meas.	Formula: Hansmann
	TTD	cm	Refer to Distance Meas.	Formula: Hansmann
	ThC	cm	Refer to Distance Meas.	

Meas. Menu	Submenu	Unit	Meas. Method	Comment
			Refer to Distance	Formula to choose: Jeanty,
	Humerus	cm	Meas.	Korean, Merz, Osaka
	Ulna Long	cm	Refer to Distance Meas.	Formula: Jeanty
Fetal Long	Tibia Long	cm	Refer to Distance Meas.	Formula to choose: Jeanty, Merz
Bones	Radius Long	cm	Refer to Distance Meas.	
	Fibula Long	cm	Refer to Distance Meas.	
	Clavicle Long	cm	Refer to Distance Meas.	Formula: Yarkoni
	Cerebellum	cm	Refer to Distance Meas.	Formula to choose: Chitty, Hill
	СМ	cm	Refer to Distance Meas.	
	NF	cm	Refer to Distance Meas.	
	NT	cm	Refer to Distance Meas.	
Fetal Cranium	OOD	cm	Refer to Distance Meas.	Formula: OOD
	IOD	cm	Refer to Distance Meas.	
	NB	cm	Refer to Distance Meas.	
	LVent	cm	Refer to Distance Meas.	Formula: Tokyo
	HC Width	cm	Refer to Distance Meas.	
	LtKid	cm	Refer to Distance Meas.	
	RtKid	cm	Refer to Distance Meas.	
	LtRenalAP	cm	Refer to Distance Meas.	
OB Others	RtRenalAP	cm	Refer to Distance Meas.	
	LVWrHEM	cm	Refer to Distance Meas.	
	MAD	cm	Refer to Distance Meas.	
AFI		cm	Refer to Distance Meas.	AFI=AFI1+AFI2+AFI3+AFI 4

Meas. Menu	Submenu	Unit	Meas. Method	Comment
CX L		cm	Refer to Distance	
		cm	Meas.	

6.6.1 Twins Measurement

1. In the new patient OB page, choose the number of gestation from one to four.

2. In the measurement menu, press baby A, and then press **[ENTER]** key to switch babies, which could measure the babies separately.

6.6.2 EDD (estimated date of delivery) Estimation

6.6.2.1 Calculating EDD by LMP (Last menstrual period)

1. In the new patient OB page, update the LMP input box.

2. Choose the LMP from the date dialog box or input the LMP date directly.

3. The calculated EDD value will appear in the result measurement area of OB page.

6.6.2.2 Calculating EDD by BBT (Basal body temperature)

1. In the new patient OB page, update the Ovul.Date input box and input the bbt date.

2. The method is the same with the LMP method.

6.6.3 Growth curves

Function: Growth curves comparison is used to compare the measured data of the fetus with the normal growth curve in order to judge whether the fetus grows normally.

Measurement steps:

- 1. Finish the measurement of the OB item and get into the report page.
- 2. Choose the growth curve at the right list and press **[ENTER]** key to display the growth curve.
- 3. Choose the growth curve need to display, and check it to show the growth curve on report.
- 4. Click [*] icon on the dialog box to exit.

Tips: The abscissa of Growth curves is the gestational weeks calculated according to the LMP in patient information.

6.7 B Pediatrics Measurement

Select OB mode, press [CALC] to enter OB mode, then enter into pediatrics mode. Or press

to switch to pediatrics measurement menu.

6.7.1 HIP Angle

٢Ì

HIP function is used for evaluating the fetal hip growth. In order to make calculation, three lines need to be added on the image, which is to conform to the fetal anatomic structure. The system will calculate and display two angles for doctor's reference.

Measurement steps:

1. Choose [HIP Angle] menu item, and click it to enter into measurement.

2. Click on the line image region, and appear one line with"+" .Move the line to the target measurement region.

3. Rotate [MENU] knob to adjust the line angle, press [ENTER] key to fix the line.

4. Then appear the second line, adjust the line according to the step 3, and fix the line.

5. Fix the 3 lines; the measurement result of the angle appears in the district.


ACAUTION:

D 3 shows bias line between protruding of conjunction and acetabular bone

D 2 shows direct line between osileum and acetabular bone

D 1 shows base line between cotyle, joint purse, gristle periosteum and ilium.

 β is the angle between D1 and D 2 (acute angle); α is the angle between D 1 and D 3(acute angle).



HIP angle

6.8 B ORTH Measurement

Refer to pediatrics measurement.

6.9 B GYN Measurement

GYN measurement includes measurement of UT-D (uterus diameter), ENDO (endometrium), CX-L (Uterine cervix length), LEFT OV and RIGHT OV (volume of left and right ovary) and LEFT FO and RIGHT FO (left and right follicle). The result will be calculated and displayed automatically on the screen by measuring relevant parameters.

Freeze the required image under GYN examination, then press **[CALC]** key to enter into GYN measurement status.

Meas. Menu	Submenu	Unit	Meas. Method	Comment
Distance		cm	Refer to Distance	
Distance		CIII	Meas.	
	UT_L	cm	Refer to Distance	
		om	Meas.	
	CX_L	cm	Refer to Distance	
UT			Meas.	
	UT_W	cm	Refer to Distance	
			Meas.	
	UT_H	cm	Refer to Distance	
	01_11		Meas.	
			Refer to Distance	L: Cervix Length
Cervix Vol.		ml	Meas.	H: Cervix_Length
			Formula:	W: Cervix Width
			V=(π/6)*L*H*W	
ENDO		cm	Refer to Distance	
			Meas.	
Ovary			Refer to Distance	L: Ovary_Length
Volume		ml	Meas.	H: Ovary Height
(Left/Right)			Formula:	W: Ovary Width
(V=(π/6)*L*H*W	···· • • • • • • • • • • • • • • • • •
	Follicle width	cm	Refer to Distance	
			Meas.	
	Follicle height	cm	Refer to Distance	
			Meas.	
Follicle (Left/Right)				2distance
			Formula of 2distance	A: the longer distance
	Follicle Volume r		V=(π/6)*A*B ²	B: the shorter distance
		ml	Formula of 3distance	3distance
			V=(π/6)*L*H*W	L: follicle length
				H: follicle height
				W: follicle width

6.10 B Small Parts Measurement

Freeze the required image under Small parts examination, then press **[CALC]** key to enter into small parts measurement status.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Distance		cm	Refer to Distance Meas.	
Thyroid (Left/Right)		ml	Refer to Distance Meas. Formula: V=(π/6)*L*H*W	L: Thyroid_ Length H: Thyroid_Heigh W: Thyroid_Width
Angle		deg	Refer to Distance Meas.	Angle Range: 0°∼180°
Ratio			Refer to Distance Meas. Formula: R=D1/D2	D1: First Distance D2: Second Distance

6.11 B Vascular Measurement

Refer to normal measurement in B mode.

6.12 B Urology Measurement

Normally urology measurements are performed in B and 2B mode.

Freeze the required image under Urology examination, then press **[CALC]** key to enter into Urology measurement status.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Kidney Volume (Left/Right)		ml	Refer to Distance Meas. Formula: V=(π/6)*L*H*W	L: Kidney_L H: Kidney_H W: Kidney_w
Bladder Volume	V(L*W*H)	ml	Refer to Distance Meas. Formula: V=0.497*L*H*W	L: Bladder _L H: Bladder _H W: Bladder _w
Prostate	Volume	ml	Refer to Distance Meas. Formula : V=(π/6)*L*H*W	L: Prostate _L H: Prostate _H W: Prostate _w
	PPSA	ng/ml	Formula: PPSA=0.12*V	

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
	PSAD	ng/ml	Formula: PSAD=SPSA/V	SPSA: input the SPSA when create a new patient
Residual	Volume	ml	Refer to Distance Meas. Formula: V=0.7*L*H*W	L: RVU_L H: RVU_H W: RVU_w

6.13 B Cardiology Measurement

Normally urology measurements are performed in B and 2B mode.

Freeze the required image under Cardiac examination, then press **[CALC]** key to enter into Cardiac status.

Hint: input the BSA when creating a new patient.
--

Meas. Menu	Submenu	Meas. Method/	Comment
weas. wenu	Submenu	Meas. Formula	Comment
		Refer to Distance Meas.	
		Formula:	
		EDV=(7*LVIDd ³)/(2.4+LVIDd)	
		ESV=(7*LVIDs ³)/(2.4+LVIDs)	
Teichholz	LVIDd	SV= EDV-ESV	
Telchinoiz	LVIDs	CO=SV*HR/1000	
		EF=SV/EDV*100	
		FS=(LVIDd-LVIDs)/LVIDd*100	
		SI=SV/BSA	
		CI=CO/BSA	
		Refer to Volume of A4CTrace method	
		Formula:	
	EDV(A2C/A4C	SV= EDV-ESV	
Single Plane)	CO=SV*HR/1000	
	ESV(A2C/A4C)	EF=SV/EDV*100	
		SI=SV/BSA	
		CI=CO/BSA	
		Refer to Volume of A4CTrace method	
	EDV(A2C)	Formula:	
	ESV(A2C)	SV= EDV-ESV	
Simpson BP	EDV(A4C)	CO=SV*HR/1000	
	ESV(A4C)	EF=SV/EDV*100	
		SI=SV/BSA	
		CI=CO/BSA	

Meas. Method/				
Meas. Menu	Submenu	Meas. Formula	Comment	
Modi Simpson	LVLd LVLs LVAMd LVAMs LVAPd LVAPs	Refer to Distance Meas. & Ellipse Mea. Formula: EDV= $\frac{LVLd}{9} \times (4 \times LVAMd + 2 \times LVAPd + \sqrt{LVAMd \times LVAPd})$ ESV= $\frac{LVLs}{9} \times (4 \times LVAMs + 2 \times LVAPs + \sqrt{LVAMs \times LVAPs})$ SV= EDV-ESV CO=SV*HR/1000 EF=SV/EDV*100 SI=SV/BSA CI=CO/BSA		
Cube	LVSd LVIDd LVPWd IVSs LVIDs LVPWs	Refer to Distance Meas. Formula: EDV=LVIDd ³ ESV=LVIDs ³ SV= EDV-ESV CO=SV*HR/1000 EF=SV/EDV*100 FS=(LVIDd-LVIDs)/LVIDd*100 SI=SV/BSA CI=CO/BSA		
Bullet Volume	LVLd LVLs LVAMd LVAMs	Refer to Distance Meas. Formula: EDV=(5/6.0)*LVLd*LVAMd ESV=(5/6.0)*LVLs*LVAMs SV= EDV-ESV CO=SV*HR/1000 EF=SV/EDV*100 SI=SV/BSA CI=CO/BSA		
Gibson	LVIDd LVIDs	Refer to Distance Meas. Formula: EDV=π/6*(0.98*LVIDd+5.9)*LVIDd*LVIDd ESV=π/6*(0.98*LVIDs+5.9)*LVIDs*LVIDs SV= EDV-ESV CO=SV*HR/1000 EF=SV/EDV*100 FS=(LVIDd-LVIDs)/LVIDd*100 SI=SV/BSA CI=CO/BSA		
Mitral Valve	MV Diam MV Area LA/AO	Refer to Distance Meas. & Ellipse Mea. Formula: LA/AO=LAD/AOD		

Meas. Menu	Submenu	Meas. Method/	Comment
Aortic Valve	AV Diam	Meas. Formula Refer to Distance Meas. & Ellipse Mea.	
Pulmonary	AV Area PV Diam		
Valve	PV Area	Refer to Distance Meas. & Ellipse Mea.	
Tricuspid Valve	TV Diam TV Area	Refer to Distance Meas. & Ellipse Mea.	
LVOT	LVOT Diam LVOT Area	Refer to Distance Meas. & Ellipse Mea.	
RVOT	RVOT Diam RVOT Area	Refer to Distance Meas. & Ellipse Mea.	
RV/LV	RVIDd LVIDd	Refer to Distance Meas. Formula: RV/LV=RVIDd/LVIDd*100	
PISA	PISA MR(Rad, Als Vel) PISA AR(Rad, Als Vel) PISA TR(Rad, Als Vel) PISA PR(Rad, Als Vel)	Refer to Distance Meas. Formula: Flow Rate=2π*Rad*Rad*Als Vel	
	Cube(LVSd, LVIDd, LVPWd)	Refer to Distance Meas. Formula: LV Mass=1.04*((LVSd+ LVIDd+ LVPWd) ³ - LVIDd ³)-13.6 LV Mass Index=LV Mass/BSA	
LV Mass	A-L(LVAd sax Epi, LVAd sax Endo, LVLd apical)	Refer to Distance Meas. & Ellipse Mea. Formula: $1.05 \times \frac{5}{6} \left\{ \left[LVAd \ sax \ Epi \times \left(LVLd \ apical + \sqrt{\frac{LVAd \ sax \ Epi}{\pi}} - \sqrt{\frac{LVAd \ sax \ Endo}{\pi}} \right) \right] - (LVAd \ sax \ Endo \times LVLd \ apical) \right\}$	
	T-E(LVAd sax Epi, LVAd sax Endo, a, d)	Refer to Distance Meas. & Ellipse Mea. Formula: $1.05\pi \times \left\{ (b+t)^2 \left[\frac{2(a+t)}{3} + d - \frac{d^3}{3(a+t)^2} \right] - b^2 \left(\frac{2a}{3} + d - \frac{d^3}{3a^2} \right) \right\}$	
Qp/Qs	AV Diam PV Diam	Refer to Distance Meas.	

6.13 Normal Measurement in M, B/M mode

At real-time status, press **[M]** key twice to enter M mode, press **[CALC]** key to enter into M mode measurement status.

OR

At real-time status, click [M] key to enter B/M mode, press [CALC] key to enter into M mode measurement status.

6.13.1 Distance

Measurement steps:

1. Select menu item [Distance]to enter into measurement.

2. Click on the M image area, it will display an orange dotted line with two horizontal short line. The orange dotted line represents the position need to be measured. The distance between the two short lines is the distance you want to measure. The yellow short line represents it's in active status. Click it and drag the short line to anywhere you want to put.

3. Press **[UPDATE]** key to activate the two short lines in turns and drag them to change the distance between them. The measurement result will be displayed on the result area.

6.13.2 Time

Measurement steps:

1. Select menu item [Time] to enter into measurement.

2. Click on the M image area, it will display two orange straight dotted line. The orange dotted line with one yellow short line on it represents it is in active status. The distance between the two straight lines stands for time you want to measure. You can drag the active straight line to anywhere you want to change the measured time.

3. Press **[UPDATE]** key to activate the two straight lines in turns and drag them to change the distance between them. The measurement result will be displayed on the result area.

6.13.3 Heart rate

Heart rate is used to calculate the number of heart beats per minute from cardiac image.

Measurement steps:

- 1. Choose [Heart Rate] menu item to enter into measurement.
- 2. The method is the same as Time.
- 3. After the above measurement, the calculated heart rate result is displayed in the measurement result area.
- 4. Repeat the steps form1 to 3 to start next measurement.

6.13.4 Velocity

Measurement steps:

- 1. Click [Velocity] menu item to enter measurement condition.
- 2. Select the starting point and press **[ENTER]** key, the starting point and the removable cursor display, drag cursor to the end point.
- 3. Press **[ENTER]** key again, measurement completes, the result displays in the region of measurement
- 4. Repeat 1~3, and enter the next measurement of velocity.

NOTE: The maximum number of the measurement result on the image area is one. The second measurement result will cover the first one. The measurement result area will list all the measurement values.

6.14 General Measurement in M mode

At real time status, press [M] key twice to enter M mode, click [CALC] key to enter into M mode cardiology measurement status.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Distance		cm	Refer to M Distance Meas.	
Time		S	Refer to M Time Meas.	
Velocity		cm/s	Refer to M Velocity Meas.	
Lleart Date	One Cycle	bpm	Refer to M Time Meas.	
Heart Rate	Double Cycles	bpm	Refer to M Time Meas.	

6.15 M Cardiology Measurement

Freeze the required image under Urology examination, then press [CALC] key to enter into

Meas.	Submenu	Meas. Method/	Comment
Menu	Submenu	Meas. Formula	Comment
Distance		Refer to M Distance Meas.	
Time		Refer to M Time Meas.	
Slope		Refer to M Velocity Meas.	
		Refer to M Distance Meas.	
		Formula:	
		EDV=LVIDd ³	
	Cube(LVSd,	ESV=LVIDS ³	
	LVIDd, LVPWd,	SV= EDV-ESV	
	LVSs, LVIDs,	CO=SV*HR/1000	
	LVPWs)	EF=SV/EDV*100	
		FS=(LVIDd-LVIDs)/LVIDd*100	
		SI=SV/BSA	
		CI=CO/BSA	
		Refer to M Distance Meas.	
		Formula:	
		EDV=7*LVIDd ³ /(2.4+ LVIDd)	
		ESV=7*LVIDs ³ /(2.4+ LVIDs)	
Left	Teichholz(LVIDd,	SV= EDV-ESV	
Ventricle	LVIDs)	CO=SV*HR/1000	
ventricie		EF=SV/EDV*100	
		FS=(LVIDd-LVIDs)/LVIDd*100	
		SI=SV/BSA	
		CI=CO/BSA	
		Refer to M Distance Meas.	
		Formula:	
		EDV=π/6*(.98*LVIDd+0.59)*LVIDd ²	
		ESV=π/6*(1.14*LVIDs+4.18)*LVIDs ²	
	Gibson(LVIDd,	SV=EDV-ESV	
	LVIDs)	CO=SV*HR/1000	
		EF=SV/EDV*100	
		FS=(LVIDd-LVIDs)/LVIDd*100	
		SI=SV/BSA	
		CI=CO/BSA	
	IVSd/LVPWd	Refer to M Distance Meas.	
	EPSS	Refer to M Distance Meas.	
	MV E Amp	Refer to M Distance Meas.	
Mitral	MV A Amp	Refer to M Distance Meas.	
Valve	MV D-E Exc Dist	Refer to M Velocity Meas.	
	MV E-F Slope	Refer to M Velocity Meas.	
	MV D-E Slope	Refer to M Velocity Meas.	
	A-C Int Slope	Refer to M Velocity Meas.	

Urology measurement status. Or press to choose the cardiac measurement.

Meas.	Submenu	Meas. Method/ Meas. Formula	Comment
Menu	E Dumatian		
	E Duration	Refer to M Time Meas.	
	A Duration	Refer to M Time Meas.	
	AOD	Refer to M Distance Meas.	
	Ao Sinus Diam	Refer to M Distance Meas.	
	Ao Asc Diam	Refer to M Distance Meas.	
Aortic	Ao Arch Diam	Refer to M Distance Meas.	
Valve	Ao Desc Diam	Refer to M Distance Meas.	
	LVOT Diam	Refer to M Distance Meas.	
	LAD	Refer to M Distance Meas.	
	LVPEP/LVET	Refer to M Time Meas.	
	AA	Refer to M Distance Meas.	
	RVOT Diam	Refer to M Distance Meas.	
Tricuspid	RA Diam	Refer to M Distance Meas.	
Valve	D-E Exc Dist	Refer to M Velocity Meas.	
Valve	E-F Slope	Refer to M Velocity Meas.	
	A-C Int Time	Refer to M Velocity Meas.	
	RVPEP	Refer to M Time Meas.	
Pulmonar	RVET	Refer to M Time Meas.	
y Valve	A wave Amp	Refer to M Distance Meas.	
	B-C Slope	Refer to M Velocity Meas.	
RV/LV	RVIDd	Refer to M Distance Meas.	
	LVIDd		

6.16 Other Urology Measurement in M Mode

M Abdomen, OB, GYN, Urology, Small Parts, Pediatric and ORTH measurement refer to general measurement in M mode.

6.17 PW mode measurement methods

Press [D] key to enter PW mode, and then press [CALC] key to enter PW mode measurement.

- In order to get accurate result, the PW image must be clear and high quality.
- Insure that you have fixed the cursor at the exact place of cardiac systole and diastole.

6.17.1 Velocity

Measurement steps:

1. Select menu item [Velocity] to enter into measurement.

2. Click on the PW image area, the measurement result will be displayed on the result area.

6.17.2 Time

Measurement steps:

1. Select menu item [Time] to enter into measurement.

2. Click on the PW image area, it will display two orange straight dotted line. The orange dotted line with one yellow short line on it represents it is in active status. The distance between the two straight lines stands for time you want to measure. You can drag the active straight line to anywhere you want to change the measured time.

3. Press **[UPDATE]** key to activate the two straight lines in turns and drag them to change the distance between them. The measurement result will be displayed on the result area.

6.17.3 Slope

Measurement steps:

1. Click [Slope] menu item to enter measurement condition.

2. Select the starting point and press **[ENTER]** key, the starting point and the removable cursor display, drag cursor to the end point.

3. Press **[ENTER]** key again, measurement completes, the result displays in the region of measurement

4. Repeat 1~3, and enter the next measurement of slope.

6.17.4 Auto Trace

Measurement steps:

1. Move the trackball to select the start point of the one cycle and press [ENTER] key to fix it.

2. A second cursor "^" will appear, and move the trackball to the end point of the cycle, press **[ENTER]** key to fix it.

3. The measurement results will be displayed on the monitor and calculate other values of parameters.

6.17.5 Manual Trace

Measurement steps:

1. Move the trackball to select the start point of the one cycle and press [ENTER] key to fix it.

2. Move the trackball along the spectrum and press **[ENTER]** key to complete.

3. The measurement results will be displayed on the monitor and calculate other values of parameters.

6.18 PW Fast Measurement

Press [**DIST**] key to enter PW Fast measurement in PW mode. Press the corresponding parameter control key to switch the fast measurement item.

Meas.	Culture and	Unit	Meas. Method/	Commont
Menu	Submenu	Unit	Meas. Formula	Comment
	Vs	cm/s	Refer to PW Slope Meas.	
	Pressure(s)	ma ma l l ai	Formula:	
	Flessule(s)	mmHg	Pressure=4*Vs*Vs/10000	
	Vd	cm/s	Refer to PW Slope Meas.	
Peak	Pressure(d)	mmHg	Formula:	
	Flessule(u)	i i i i i i i g	Pressure=4*Vd*Vd/10000	
	SD		Formula: SD=Vs/Vd	
	RI		Formula: SD=(Vs-Vd)/Vs	
	Time	s	Refer to PW Time Meas.	
HR	Single wave	bpm	Refer to PW Time Meas.	
	Vs	cm/s	Refer to PW Slope Meas.	
	Pressure(s)	mmHg	Formula:	
			Pressure=4*Vs*Vs/10000	
	Vd	cm/s	Refer to PW Slope Meas.	
	Pressure(d)	mmHg	Formula:	
		i i i i i i i g	Pressure=4*Vd*Vd/10000	
Auto	VMean	cm/s	Refer to PW Slope Meas.	
Trace/Manu	Pressure(VMea		Formula:	
al Trace	n)	mmHg	Pressure=4*VMean*VMean/100	
	,		00	
	TVI	cm		
	SD		Formula: SD=Vs/Vd	
	RI		Formula: SD=(Vs-Vd)/Vs	
	PI		Formula: SD=(Vs-Vd)/VMean	
	HR(Single	bpm		
	wave)	~~~···		

6.19 PW General Measurement

Press [CALC] to enter PW measurement in PW mode. Press the corresponding parameter control

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Velocity		cm/s	Refer to PW Slope Meas.	
Distance		cm	Refer to B Distance Meas.	
	Vs	cm/s	Refer to PW Slope Meas.	
	Pressure(s)	mmHg	Formula: Pressure=4*Vs*Vs/10000	
	Vd	cm/s	Refer to PW Slope Meas.	
Peak	Pressure(d)	mmHg	Formula: Pressure=4*Vd*Vd/10000	
	SD		Formula: SD=Vs/Vd	
	RI		Formula: SD=(Vs-Vd)/Vs	
	Time	s	Refer to PW Time Meas.	
	Vs	cm/s	Refer to PW Slope Meas.	
	Pressure(s)	mmHg	Formula: Pressure=4*Vs*Vs/10000	
	Vd	cm/s	Refer to PW Slope Meas.	
	Pressure(d)	mmHg	Formula: Pressure=4*Vd*Vd/10000	
	VMean	cm/s	Refer to PW Slope Meas.	
Trace <autotrac e></autotrac 	Pressure(VMea n)	mmHg	Formula: Pressure=4*VMean*VMean/100 00	
	TVI	cm		
	SD		Formula: SD=Vs/Vd	
	RI		Formula: SD=(Vs-Vd)/Vs	
	PI		Formula: SD=(Vs-Vd)/VMean	
	HR(Single wave)	bpm		
	Distance1	cm	Refer to B Distance Meas.	
StD%	Distance2	cm	Refer to B Distance Meas.	
	StD%	%	Formula: StD%=((D1-D2)/D1)*100%	D1:Distance1, D2:Distance2
	Area1	cm ²	Refer to B Ellipse Meas.	
StA%	Area2	cm ²	Refer to B Ellipse Meas.	
SIA 70	StA%	%	Formula: StA%=((A1-A2)/A1)*100%	A1:Area1,A2:Area 2
	ICA	cm/s	Refer to PW Velocity Meas.	
10.1/02.5	Pressure(ICA)	mmHg	Formula: Pressure=4*ICA*ICA/10000	
ICA/CCA	CCA	cm/s	Refer to PW Velocity Meas.	
	Pressure(CCA)	mmHg	Formula: Pressure=4*CCA*CCA/10000	

key to switch the general measurement item.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
	ICA/CCA		Formula: ICA/CCA	
	Diam	cm	Refer to B Distance Meas.	
	TVI	cm		
	Time	s	Refer to PW Time Meas.	
Flow Volume	HR(Single wave)	bpm		
	SV	ml	Formula: SV=0.785*Diam*Diam* TVI	
	СО	l/min	Formula: CO=SV*HR(Single wave)/1000	

6.20 PW OB Measurement

Press [CALC] key to enter PW measurement in PW OB mode. Or press the corresponding parameter control key to switch the OB measurement item.

Meas.	Submenu	Unit	Meas. Method/	Comment
Menu	Gubinenu	Onic	Meas. Formula	Comment
	Vs	cm/s	Refer to PW Trace Meas.	
	Pressure(s)	mmHg	Formula:	
Umb A	Flessule(s)	minig	Pressure=4*Vs*Vs/10000	
Aorta,	Vd	cm/s	Refer to PW Trace Meas.	
Descendi	Pressure(d)	mmHg	Formula:	
ng Uterine	Flessule(u)	i i i i i i g	Pressure=4*Vd*Vd/10000	
Artery(Lt),	VMean	cm/s	Refer to PW Trace Meas.	
Uterine	Pressure(VMea		Formula:	
Artery(Rt),	n)	mmHg	Pressure=4*VMean*VMean/100	
Pulmonar	11)		00	
y Artery,	TVI	cm		
MCA	SD		Formula: SD=Vs/Vd	
Duct Veno	RI		Formula: SD=(Vs-Vd)/Vs	
FHR	PI		Formula: SD=(Vs-Vd)/VMean	
	HR(Single	hnm		
	wave)	bpm		

6.21 PW GYN Measurement

Press [CALC] to enter PW measurement in PW GYN mode. Or press the corresponding

Meas.	Submenu	Unit	Meas. Method/	Comment
Menu		Onit	Meas. Formula	
Umb A,	Vs	cm/s	Refer to PW Trace Meas.	

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
MCA,	Pressure(s)	mmHg	Formula:	
Uterin A,		minig	Pressure=4*Vs*Vs/10000	
Fetal AO	Vd	cm/s	Refer to PW Trace Meas.	
	Proceuro(d)	mmHa	Formula:	
	Pressure(d)	mmHg	Pressure=4*Vd*Vd/10000	
	VMean	cm/s	Refer to PW Trace Meas.	
	Pressure(VMea n)	mmHg	Formula:	
			Pressure=4*VMean*VMean/100	
			00	
	TVI	cm		
	SD		Formula: SD=Vs/Vd	
	RI		Formula: SD=(Vs-Vd)/Vs	
	PI		Formula: SD=(Vs-Vd)/VMean	
	HR(Single wave)	bpm		

6.22 PW Cardiology Measurement

Press [CALC] to enter PW measurement in PW cardiology mode. Or press the corresponding

parameter control key to switch the cardiology measurement item.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
Velocity		cm/s	Refer to PW Velocity Meas.	
Acceleratio n		cm/s ²	Refer to PW Slope Meas.	
Time		S	Refer to PW Time Meas.	
Slope		cm/s	Refer to PW Slope Meas.	
HR(Single Wave)		bpm	Refer to PW Time Meas.	
ED/PS	ED PS	cm/s	Refer to PW Velocity Meas. Formular: RI=(ED-PS)/ED	
Mitral Valve	MV E Vel	cm/s	Refer to PW Velocity Meas. Formular: MV E PG=4*MV E Vel*MV E Vel/10000 E/Ea=MV E Vel/Ea	

			Refer to PW Velocity Meas.	
			Formular:	
	MV A Vel	cm/s	MV A PG=4*MV A Vel*MV A	
			Vel/10000	
			E/A=MV E Vel/MV A Vel	
			A/E=MV A Vel/MV E Vel	
			Refer to PW Trace Meas.	
			Formular:	
			MV PGmax=4*MV Vmax*MV	
			Vmax/10000	
	MV VTI		MV PGmean=4*MV	
			Vmean*MV Vmean/10000	
			MV SV=0.785*MV Diam*MV	
			Diam* MV VTI	
			MR Fraction=MR Flow/MV	
			SV*100	
			Formular:	
	MVA(PHT)		MV PHT=(-1)*0.3*VPeak/Slope	
			MVA PHT=220/MV PHT/1000	
			Formular:	
	MVA(VTI)		MV(VTI)=π* LVOT VTI*LVOT	
			Diam*LVOT Diam /MV VTI/4	
	MV E Dur	S	Refer to PW Time Meas.	
	MV A Dur	S	Refer to PW Time Meas.	
	MV DecT	S	Refer to PW Time Meas.	
			Refer to PW Velocity Meas.	
	MR Vmax	cm/s	Formula:	
			MR PGmax=4*MR Vmax*MR	
			Vmax/10000	
			Refer to PW Trace Meas.	
			Formula:	
	MR VTI		MR PGmax=4*MR Vmax*MR	
			Vmax/10000	
			MR PGmean=4*MR	
			PGmean*MR PGmean/10000	
	dp/dt		Formula:	
	-		dp/dt=32/Time	
	IVRT	S	Refer to PW Time Meas.	
	IVCT	S	Refer to PW Time Meas.	
			Refer to PW Velocity Meas.	
Aortic	AV Vmax	cm/s	Formular:	
			AV PGmax=4*AV Vmax*AV	
			Vmax/10000	

			Refer to PW Trace Meas.	
			Formular:	
			AV PGmax=4*AV Vmax*AV	
			Vmax/10000	
			AV PGmax=4*AV Vmean*AV	
	AV VTI		Vmean/10000	
			AV SV=0.785*AV Diam*AV	
			Diam* AV VTI	
			AR Fraction=AR Flow/AV	
			SV*100	
			Formular:	
	AVA(VTI)		AVA(VTI)=π* LVOT VTI*LVOT	
			Diam*LVOT Diam//AV VTI/4	
			Refer to PW Velocity Meas.	
			Formular:	
	LVOT Vmax	cm/s	LVOT PGmax=4*LVOT	
			Vmax*LVOT Vmax/10000 Refer to PW Trace Meas.	
	LVOT VTI		Formular:	
			LVOT PGmax=4*LVOT	
			Vmax*LVOT Vmax/10000	
			LVOT PGmax=4*LVOT	
			Vmean*LVOT Vmean/10000	
	LVPEP	S	Refer to PW Time Meas.	
	LVET	S	Refer to PW Time Meas.	
			Refer to PW Velocity Meas.	
	AR Vmax	cm/s	Formular:	
			AR PGmax=4*AR Vmax*AR	
			Vmax/10000	
			Refer to PW Trace Meas.	
			Formular:	
	AR VTI		AR PGmax=4*AR Vmax*AR	
			Vmax/10000	
			AR PGmean=4*AR Vmean*AR	
			Vmean/10000	
	AR DecT	S	Refer to PW Time Meas.	
	AR PHT			
			Refer to PW Velocity Meas.	
Tricuspid	T) () (Formular:	
	TV Vmax	cm/s	TV PGmax=4*TV Vmax*TV	
			Vmax/10000	
Valve			Refer to PW Velocity Meas.	
			Formular:	
	TV E Vel	cm/s	TV E PG=4*TV E Vel*TV E	
			Vel/10000	
L				

			Refer to PW Velocity Meas. Formular:	
	TV A Vel	cm/s	TV A PG=4*TV A Vel*TV A	
		011/5	Vel/10000	
			E/A=TV E Vel/TV A Vel	
			A/E=TV A Vel/TV E Vel	
			Formular:	
	TVA(VTI)		TVA(VTI)=π* RVOT VTI*LVOT	
			Diam*LVOT Diam /TV VTI/4	
			Formular:	
			TR PGmax=4*TR Vmax*TR	
	RVSP		Vmax/10000	
			RVSP=RAP+4*TR Vmax*TR	
			Vmax/10000	
			Refer to PW Velocity Meas. Formular:	
	RVOT Vmax		RVOT PGmax=4*RVOT	
			Vmax*RVOT Vmax/10000	
			Formular:	
	RVOT VTI		RVOT PGmax=4*RVOT	
			Vmax*RVOT Vmax/10000	
			RVOT PGmean=4*RVOT	
			Vmean*RVOT Vmean/10000	
			Refer to PW Velocity Meas.	
			Formular:	
	PV Vmax	cm/s	PV PGmax=4*PV Vmax*PV	
			Vmax/10000	
Dulmonom			Formular:	
Pulmonary Valve			PV PGmax=4*PV Vmax*PV	
Valve	PV VTI		Vmax/10000	
			PV PGmean=4*PV Vmean*PV	
			Vmean/10000	
			Formular:	
	PVA(VTI)		PVA(VTI)=π* RVOT VTI*LVOT	
			Diam*LVOT Diam /PV VTI/4	
			Refer to PW Velocity Meas.	
	MPA Vmax	cm/s	Formular:	
			MPA PGmax=4*MPA	
			Vmax*MPA Vmax/10000	
			Refer to PW Velocity Meas.	
	RPA Vmax	cm/s	Formular: RPA PGmax=4*RPA	
			Vmax*RPA Vmax/10000	

			Refer to PW Velocity Meas.
	LPA Vmax	cm/s	Formular:
			LPA PGmax=4*LPA Vmax*LPA
			Vmax/10000
	RVPEP	S	Refer to PW Time Meas.
	RVET	S	Refer to PW Time Meas.
			Formular:
	RAEDP		PR PGmax=4*PR Vmax*PR
			Vmax/10000
			RAEDP=RAP+4*PR Vmax*PR
			Vmax/10000
	RVET		Refer to PW Velocity Meas.
			Formular:
			PR PGmax=4*PR Vmax*PR
			Vmax/10000
			Formular:
			PR PGmax=4*PR Vmax*PR
	PR VTI		Vmax/10000
			PV PGmean=4*PR Vmean*PR
			Vmean/10000
			Refer to PW Velocity Meas.
Pulmonary	PVein D Vel		Refer to PW Velocity Meas.
	PVein A Vel		Refer to PW Velocity Meas.
Vein	PVein A Dur		Refer to PW Time Meas.
-	PVein S VTI		Refer to PW Trace Meas.
	PVein D VTI		Refer to PW Trace Meas.
PISA	MR,AR,TR,PR		Refer to PW Trace Meas.
			Formular:
			EROA= Flow Rate/Vmax
			Flow=EROA/VTI
			Fraction=Flow/SV*100
Qp/Qs	AV VTI		Refer to PW Trace Meas.
			Formular:
			AV SV=0.785*AV Diam*AV
			Diam/ AV VTI
			CO=AV SV*AV HR/1000
			Refer to PW Trace Meas.
			Formular:
			PV SV=0.785*PV Diam*PV
	PV VTI		Diam/ PV VTI
			CO=PV SV*PV HR/1000
			Qp/Qs
			Qp-Qs
Tei Index	IRT		Refer to PW Time Meas.

	ICT	Refer to PW Time Meas.
	ET	Refer to PW Time Meas. Formular:
		Tei Index=(ICT+IRT)/ET
TDI	MV medial (Sa,Ea,Aa,ARa,DR a)	Refer to PW Velocity Meas. & Slope Meas.
	MV lateral	Refer to PW Velocity Meas. &
	(Sa,Ea,Aa,ARa)	Slope Meas.

6.23 PW Vascular Measurement

Press [CALC] to enter PW measurement in PW vascular mode. Or press the corresponding parameter control key to switch the vascular measurement item.

Meas. Menu	Submenu	Unit	Meas. Method/ Meas. Formula	Comment
	Vs	cm/s	Refer to PW Slope Meas.	
	Pressure(s)	mmHg	Formula:	
	Flessule(s)	minig	Pressure=4*Vs*Vs/10000	
Subclavian A,	Vd	cm/s	Refer to PW Slope Meas.	
Prox CCA,	Pressure(d)	mmHg	Formula:	
Mid CCA,	Flessule(u)	minig	Pressure=4*Vd*Vd/10000	
Distal CCA,	VMean	cm/s	Refer to PW Slope Meas.	
Bulb,	Pressure(V		Formula:	
Prox ICA,	Mean)	mmHg	Pressure=4*VMean*VMean/100	
Mid ICA,			00	
Distal ICA,	TVI	cm		
ECA,	SD		Formula: SD=Vs/Vd	
Vertebral A	RI		Formula: SD=(Vs-Vd)/Vs	
	PI		Formula: SD=(Vs-Vd)/VMean	
	HR(Single	bpm		
	wave)			

6.24 Other Urology Measurement in PW Mode

PW Abdomen, Urology, Small Parts and Pediatric measurement refer to PW general measurement.

Chapter 7 Preset

This chapter introduces the operation to make settings of the system through preset menu at preset mode.

Preset function is used to set up working environment and status, parameters of each examination mode. The setting will be stored in the memory of system and not be lost even after the system is switched off. When the system is switched on, it will work automatically with the status which is required by the operator.

In preset interface, all operation relies on moving trackball to required Function key position. Press **[ENTER]** key to start operation.

7.1 General setting

7.1.1 General Interface

CHISON		Setup			×
General Measuren	nent Comment BodyMark	ExamMode Repo	rt DICOM	Net System	
General Image	&Cine Keyboard				
Hospital			Screen Controlle		
Hospital	hospital		Enable ScreenSave	r Wait	
Dept.	neike				
Date and ti			ScreenSaver		
Date	29-04-2018	Undo	Default		
Date format	DD-MM-YYYY			Import Previe	≌W
Date format			Enable Screen Off	Wait	
Time	02:59:27	Undo			
Time Zone	UTC+8:00 Chinese Standard Time	2	Auto Export		
Language			Export files to U disk	k after storage	
	-				
Language	English	▼	Export files to Netwo	ork after storage	
				🗟 🏹	5 ¥
		Default	Ok Can		

Press the [SETUP] to enter system setting interface. User can do user-defined setting.

Click [X] in the title bar or the exit key on the button can exit the system setting interface.

Function name	Setting method	Function description			
Hospital, Department	Input freely	Set up the hospital name which is shown at top left corner of "General Setting" dialog box, 20 characters Max can be input.			
Date and Time	Input freely	Set up the system date (calendar format), select current date directly. Date format can be changed by format setting.			
Date Format	Set up freely	Set up date format: Year/Month/Date, Month/Date/Year, Date/Month/Year.			
Time zone	Set up freely	Set up the working clock of the system.			
		Enable screen saver, user can custom screen			
Screen	Select the function	saver images, a custom picture named			
Controller	and set its start time	"screensaver", JPG, PNG, BMP format, size does			
		not exceed 512 * 384 pixels.			
Language	Select language	Select the language of operation interface(Simplified Chinese English and so on)			
Auto Export	Select the need option	Select to export files to U disk or Network after storage.			
Default	Press button	Recover all preset to factory setting.			

7.1.2 Image & Cine setting

User can select Image&Cine option to set the screen type, parameters of image, cine or image area. See the following sheet 2 for detailed description.

CHISON		Setup					×
General Measurement Comment E	BodyMark	ExamMode	Report	DICOM	Net	System	
General Image&Cine Keyboard							
Screen Type			Options	for Image Area			n II
Only image			STC curve	Always hide			
Image and patient information							
O Full screen			Options	for Transmitted	Images		
Cine Frames					Opt		
Real Time saving 4 s	▼ fr	rames	PC Print	Image Option			
After freeze saving 256	▼ fr	rames	Only image) In	age and info		
			Images columr	n for print in revi	ew 2	-	
							-
		Defaul	t C	0k	Cancel	🗟 妃 En 🔲	∳ ÿ

Function name	Setting method	Function description
	Coloct the need type	Set the content which screen picture contain: only
Screen Type	Select the need type	picture, image and patient information, full screen.
Cine Frames	Select the parameter	Set the default frame or time when saving the film.
Options for	Click the button to	Set the STC curve, including always shows, always
Image Area	open the setting box	hide and hide for 1 to 8 seconds.
Options for	Click the button to	Adjust the parameters of transmitted images:
Transmitted	open the setting box	brightness, contrast, and gamma and select the need
Images	open the setting box	image type.
PC Print Image	Select the need	Including print area and arrangement condition in the
Option	option	image foresee interface.

7.1.3 Keyboard setting

Users can set the different functions for number keys and other optional buttons according to the using habits.

CHISO	Ν			Setup					×
General Mea	surement	Comment	BodyMark	ExamMode	Report	DICOM	Net	System	
General	mage&Cine	Keyboard							
				Select Fur					
Key		unction		Output I	mage Mode	Print M	leasurement	Other	
Foot SW1	Video Print	t		· · · · · · · · · · · · · · · · · · ·					
Foot SW2	Freeze			store Image					
	GS			Store Cine					
	CRL			Store Image To					
	BPD			Store Cine To U	Disk				
	AC		٩	lone					
	FL								
	NT								
	AFI								
	M Fetal HR								
	HIP Angle								
	Trapezoida	alMode							
Save Image	Store Imag	je							
Save Cine	Store Cine								
								8 2	ψ
				Defaul		Ok	Cancel		
								En 💻	Ŷ

- Key: display the settable number keys and buttons and selected function.
- **Output:** including store images, store films, store images to U disk, store film to U disk, etc. function options.
- Image mode: including image full screen, biopsy, chroma, etc.
- Print: including video print, pc print etc.
- Measurement: including GS, CR, BPD, HC, AC, etc.
- Other: including arrow, archive.

7.2 Measurement

Measurement includes general measurement setting and measurement formula setting.

7.2.1 General measurement setting

General settings can only change the display of measurement unit.

CHISON		Setup					×
General Measurement Co	mment BodyMark	ExamMode	e Report	DICOM	Net	System	
General Configuration	OB Table						
Unit			Measure	Option			
Distance	cm		🛃 Clear results	upon unfreeze			
			🛃 Freeze upon	measure			
Area	cm²		Measure	Result			
Volume	ml		Font Color	Yellow			
Time	s		Follicle Method	3 Distances			
Speed	cm/s		Heart Ra	te 1 Cycle			
Acceleration			HIP GRAF				
Acceleration	cm/s²		Display HIP C				
Angle	deg						
Wataba							
Weight	g						
						5	₽ ₩
		Defau	ilt O	k C	Cancel	En 🛛	₽ ₽

- Distance: cm, mm
- Area: cm², mm²
- Volume: ml, l
- Time: s, ms
- Speed: mm/s, cm/s, m/s
- Acceleration: mm/s², cm/s², m/s²
- Angle: deg, rad
- Weight: g, kg
- Measure Option: whether clear results upon unfreeze, and the automatic froze image.
- Measure Result: the color of the result font is alternative, including yellow/white/orange/green.
- Follicle Method: the ways to measure follicle, you can choose two distances and three distances.
- Heart Rate: cycle options for calculating the number of heart beats per minute of cardiac image.
- HIP GRAF: select to display HIP GRAF in result or not.

7.2.2 Measurement formula setting

CHIS			Setu	P				
General	Measurement Con	iment BodyM	ark ExamM	Node Report	DICOM	Net	System	n
General	Configuration (0B Table						
Region	General 🔻	Manage						
B	○ M ○ D							
	ral Measurement stance		Folde Caption	r Edit General Measurement	Measure Orde	r None		
	cumference/Area	Down	Var Name	GeneralMeasurement	Default Item	None		
	lume(3Distances) lume(2Distances)	Add Folder	Unique Name	GeneralMeasurement	🛃 Display In	Menu 🔲 Co	ontinue	
	lume(1Distance)	Add Measurement	📝 Distance	e	🥑 Volume(1	1Ellipse)		
	lume(1Ellipse) lume(1Ellip1Dis)	Add Calc Measurement	🛃 Circumf	ference/Area	🛃 Volume(:	1Ellip1Dis)		
+ Ratio		Remove	📝 Volume	(3Distances)				
+ Angle		Cut	🛃 Volume	(2Distances)				
Histor			🛃 Volume	(1Distance)				
							Save	
							5	
			D	efault Ok	Car		En 🖸	⊒ ¢

7.2.2.1 Interface Description-Measurement Menu

• Region: pull down and select needed measurement menu.

• Manage: pop up measurement software edition interface, add modify delete change marshalling sequence in measurement menu.

- **B**, **M**, **D**: display measurement of each Exam mode.
- Up: press this button to move selected measurement term up.
- **Down:** press this button to move selected measurement term down.
- Add Folder: add a measurement item. In the left column when the term is fold there is "+" otherwise "-".

• Add Measurement: add a measurement item for a term in the right column there is selected

item and detailed parameter.

- Add Calc Measurement: add a calc item for a measurement term.
- **Remove:** remove selected measurement term or item.
- **Default:** restore all measurement term as factory setting.
- Save: save measurement item modification users did.

Folder Edit content description

Caption	Display	names	of	all	items	that	is	displayed	name	in
Caption	measure	ement m	ienu	J.						

	The name of built-in selected measurement menu, user
Var Name	don't need modify while display order according to the
	names.
Unique Name	Built-in code, user doesn't need modify.
Measurement	None: Disable rule, Repeat: Repeat this item , Next:
Order	measure by sequence.
	After choosing the Repeat and Sequential, choose one
Default item	measurement or calculation to activate the measurement
	rule.

7.2.2.2 Interface Description-Measurement Manipulation

CHISON		Setup					×
General Measurement Com	ment BodyMark	ExamMode	Report	DICOM	Net	System	
General Configuration O Region General ▼ ● B M D	B Table Manage						
 General Measurement Distance Circumference/Area Volume(3Distances) Volume(2Distances) Volume(1Ellipse) Volume(1Ellip1Dis) Ratio Angle Histogram Profile 	Down Va Add Folder Un Add Measurement Add Calc	ir Name Circu	umference/Area umferenceArea1 cumferenceArea1 Jnit Display In Re n Yes	Tool I 💌 Display	Circumferenc Ellipse In Menu 😨 I In Report srage Avg Avg	▼ On Line Select	
		Default	Ok		ancel	En 🗌	لا ق

Measurement Edit content description

Caption	Display the name of selected term that is display names in the measurement menu.						
Var Name	The name of built-in selected measurement menu, user don't need modify while display order according to the names.						
Unique Name	Built-in code, users don't need modify.						
Display In Menu	Check the required item and it will display on the measurement menu. The item without checking will not display on the measurement menu.						
Choose	Check the methods in measurement menu, press						
measurement	[update] key to switch the method; Otherwise, the						
methods	measurement is not available.						

Display in Result	Check and display the item in measurement menu,						
	otherwise, the item is not displayed.						
Name	Measurement operation of specific measurement display						
Name	in results.						
Unit	Data unit which measurement operation produces.						
Display In Report	Whether display in the report or not.						
Average Method	The average rule of data.						
Modify	Press this button to pop up interface to modify						
Modify measurement operation.							
Bomovo	Press this button to delete selected measurement						
Remove	operation.						

7.2.2.3 Interface Description-Measurement Calculation

CHISON		Setup					×	:
General Measurement Com	ment BodyMa	ark ExamMode		Report DI	COM Net		ystem	
General Configuration O	B Table							
Region General v	Manage							
● B ● M ● D								
+ General Measurement - Ratio	Up	Calc Mea Caption Ra	sureme atio(Dis		Var Name Ratio	Distance2		
Distance1	Down	Unique Name B	_Ratio1	_RatioDistance2	📝 Display In Me			
Distance2	Add Folder	Name	Unit	Display In Resul	t Display In Repor	terage Meth		
Ratio(Distance)	Add	Distance1		Yes	No	Avg		
Areal	Measurement	Distance2		Yes	No	Avg	Modify	
Area2 Ratio(Area)	Add Calc Measurement	Ratio(Distance)		Yes	Yes	Calc Af		
+ Angle	Remove							
Histogram Profile	Cut						Remove	
							Save	
		Defau	ilt	Ok	Cancel	En	₽ ¥	

Calc Measurement Edit content description

Caption	Display the name of selected term that is display names in				
Caption	the measurement menu.				
	The name of built-in selected measurement menu, user				
Var Name	don't need modify while display order according to the				
	names.				
Unique Name	Built-in code, users don't need modify.				
	Select available measurement tool type.				
Tool Type/Tool	B distance (B line), area/circumference (ellipse, trace), M				
	distance (M vertical line), time (M horizontal line), M slope				
	(M slant).				

Name	Needed measurement operation of specific measurement						
Name	and calculation.						
Unit	Data unit which measurement operation produces.						
Display In Result	Whether display in the report or not.						
Up	Press this button to move measurement operation up.						
Down	Press this button to move measurement operation down.						
Average Method	The average rule of data.						
Add	Press this button to pop up interface to add measurement						
Add	operation.						
Modify	Press this button to pop up interface to modify						
wouny	measurement operation.						
Remove	Press this button to delete selected measurement						
Kelliove	operation.						

7.2.2.4 Create Measurement Operation

Press [Add] button in measurement operation interface, pop up the following dialog box.

lte	m Configuration						
Caption	Volume(2Distance)						
Var Name	Volume2DistanceItem1						
Unique Name	istance21_Volume2DistanceItem1						
Ruler Result	Calculation Edit						
Min	0						
Max	9999						
Unit	ml 🔹						
Method	Calc After Mean						
📝 Display In Result							
🛃 Display In Report							
	OK Cancel						

Operation content description of creating new measurement

Caption	Display the name of selected term that is display names in		
	the measurement menu.		
	The name of built-in selected measurement menu, user		
Var Name	don't need modify while display order according to the		
	names.		
Unique Name	Built-in code, user doesn't need modify.		
Dular Decult	Needed measurement operation of specific measurement		
Ruler Result	and calculation.		
	Enter into interface to edit formula when selecting		
Edit	calculation item.		
Мах	The maximum value displays in result zone and report.		
Min	The minimum value displays in result zone and report.		

Unit Data unit which measurement operation produces.					
Display In Result	Whether display in the result or not.				
Display In Report	Whether display in the report or not.				
Method	The average rule of data.				

7.2.2.5 Formula Edit-Normal

It is necessary to enter into the following interface when creating measurement operation except OB.

				Formula Config	uration				
	Formul								
Formula	{ SV }	/{ EDV(A2)	C/A4C) }*100						Validate
									Clear
Unit	%								
	Measur	ement Iter	ns			Functi	on		
Na	ame	Unit	Measure	Folder		Standard	Sys	tem	
LVIDd		cm	LVIDd	Teichholz					
EDV		ml	LVIDd	Teichholz		Sin	Cos	Tan	Atan
LVIDs		cm	LVIDs	Teichholz		Min	Max	Sqrt	Avg
ESV		ml	LVIDs	Teichholz					
sv		ml	LVIDs	Teichholz		ехр	Ln	Log2	Log10
со		l/min	LVIDs	Teichholz		()	abs	PI
EF		%	LVIDs	Teichholz					
FS		%	LVIDs	Teichholz			8		
SI			LVIDs	Teichholz		4		6	
сі			LVIDs	Teichholz					
LVALd		cm^2	LVALd	Single Plane					
LVLd		cm	LVLd	Single Plane					
LVALs		cm^2	LVALs	Single Plane	T				
								ancel	

- Formula: edit formula in input box via keyboard and built-in formula.
- Validate: press this button to check whether the formula is right or not after editing formula.
- Clear: clear the content in the input box.
- Unit: select the unit of calculation consequence.
- **Measurement Items:** display all available measurement operation in the measurement menu.
- Function: built-in formula, number input and some parameters that system needs such as BSA, SPSA etc.
- **Cancel:** cancel editing formula and close the interface.
- **OK:** save edited operation and close the interface.

7.2.2.6 Formula Edit-OB

It is necessary to call built-in OB formula sheet when creating OB measurement operation the

following function interface is required.

CAUTION: the results of GA and EDD don't require unit, the unit of this class has been built-in.

B Tables	OB Ta			ns	ment Iter	Measure
	AC	Folder		Measure	Unit	Name
	AFI		В	Distance	cm	Distance
	APD		В	GS	cm	GS
	APTD		В	GS		GS_GA
D	APTDxTTD		В	GS		GS_EDD
	BPD		В	CRL	cm	CRL
	CEREB		В	CRL		CRL_GA
	CLAV		В	CRL		CRL_EDD
	СМ		В	BPD	cm	BPD
	CRL		В	BPD		BPD_GA
	FIB		В	BPD		BPD_EDD
	FL.		В	нс	cm	нс
	FLdivHC		В	НС		HC_GA
			В	нс		HC_EDD
	FTA		в	AC	cm	AC
	GS		В	AC		AC_GA
	HC		в	AC		AC_EDD

- Measurement Items: display created measurement term by now.
- **OB Tables:** built-in OB formula table.

OB List

CHISON	Setup		8
General Measurement Commer	nt BodyMark ExamMode	Report DICOM Net	System
General Configuration OB Ta	ble		
Measure Item			
AC AFI APD APTD	GA Tables Growth Table CFEF	Fetal Weight EFW Equation EFW Growth Campbell	
APTDXTTD BPD CEREB	Campbell Hadlock Hansmann Korean	Ferrero Hadlock1 Hadlock2	
CLAV CM CRL	Merz Shinozuka	Hadlock3 Hadlock4 Hadlock	
FIB FL FL divHC		Hansmann Higginbottom Merz	Cancel
Info	Cancel Save	Osaka	Save
	Default	Ok Cancel	Image: Image

- Measure Item: OB Measurement Item.
- GA Tables: Gestational list for the current measuring project.
- Growth Table: Growth table for the current measuring.
- Fetal Weight: Fetal Weight calculation formula.
- **EFW Equation:** Fetal weight calculation for the current measuring.
- EFW Growth: Fetal weight growth curve for the current measuring.
- Info: Display the gestational age and fetal weight for the current measuring.
- **Cancel:** Cancel the operation of choosing the formula.
- Save: Save the users' choice of formulas.

7.3 Comment

CHISON	Setup			×
General Measurement Comment Bo	dyMark ExamMode R	eport DICOM	Net Syste	em
Comment Lib Abdomen 🔻	Comment Type	Abdomen 🔻	Edit Comment Type	
Comment Lib		nt Selected	Edit Comment	
Adrenal gland	Adrenal gland		Abbreviation	
Bile duct	Bile duct		Full Name	
Bladder	Bladder			
Body of pancreas	Body of pancre			
Common bile duct	Common bile o		Comment&Arro	
Gall bladder	< Gall bladder		Clear comments an upon unfreeze	d arrows
Head of pancreas	Head of pancro	Head of pancreas		ents and
Hepatic artery	Hepatic artery		arrows	
Hepatic vein	<< Hepatic vein			
Inferior vena cava	Inferior vena c	ava		
Kidney	Kidney			
Kidney vein	Kidney vein			
Left hepatic lobe	e From Lib Left hepatic loi	e 🔺		
Liver	Liver			
				₽. ¥
	Default	Ok	Cancel	
			En	므

7.3.1 Comment Library

The annotation database of the system is classified as: abdomen, OB, GYN, Cardiac, small parts, Pathological change Annotation can be made by inputting characters from the soft keyboard or recalling the terms saved in annotation database.

Press **[CommentLib]** pull down button, pop up annotation name within system, through trackball and **[ENTER]** key to show required annotation status.

7.3.1.1 Edit Comment Library

Operation:

1. At the annotation status, Move the cursor to the **[Edit Comment Type]** button then press **[ENTER]** key, the annotation will be updated, and can be edited.

2. Input name into the new created annotation status box, move the cursor to the **[Create]** button then press **[ENTER]** key, then create new annotation status and appear in selected annotation status list.

3. Move the cursor to the [Delete] button, press [ENTER] key, then delete current annotation

status in the selected annotation list.

4. Alter name of current annotation status list in **[Current Type Name]** input box, press **[ENTER]** key on the **[Rename]** button, then rename the selected annotation status name.

MNOTE:

Factory Settings comments cannot be deleted or renamed.

7.3.2 Edit Comment

Operator uses only current annotation instead of all content annotation status provides inlay common annotation. If necessary, user can import annotation or self-compiled annotation into common annotation.

7.3.2.1 Add comment from comment library

Operation:

1. Select needed source annotation status via Trackball and [ENTER] key.

2. Select needed annotation at **[CommentLib]** column then press the **[ENTER]** key to activate this annotation.

3. Press [ENTER] key on [>] button to import selected annotation into user-selected annotation status; press [ENTER] on [>] button to move selected annotation in [Comment Selected] column into source annotation.

4. Press **[ENTER]** on **[>>]** button to import all annotation in source into user-selected annotation status; Press **[ENTER]** on **[>>]** button to move all annotation in **[Comment Selected]** column into source annotation.

7.3.2.2 Add comment manually

Operation:

1. Activate **[Edit Comment]** input box via Trackball and **[ENTER]** key, and then input needed abbreviation and full name of annotation.

2. Press **[ENTER]** key on **[Add]** button, meanwhile this handout will be added into source and user-selected annotation status.

7.3.2.3 Alter comment

Operation:

1. Alter annotation in user-selected status, the abbreviation and full name of annotation will be displayed in **[Edit Comment]** box.

2. Activate needed abbreviation and full name via [ENTER] key and alter via keyboard.

3. Press **[ENTER]** key on **[Modify]** button, modify the annotation in both source and user-selected status.

7.3.2.4 Delete comment in library

Operation:


Select needed annotation in source status, press **[ENTER]** key on **[Delete from Lib]** button, then the annotation deleted.

7.3.3 Comment and Arrow Option

Optional: whether clear comments and arrows upon unfreeze.

Optional: whether freeze upon comments and arrows.

7.4 Body marks



7.4.1 Body Marks Library

Built-in body marks: Abdomen OB, GYN, Cardiac, Small part, Urology, Vessel.

Press [BodyMarkLib] pull down button, pop up name of built-in body marks, show needed body marks via Trackball and [ENTER] key.

7.4.1.1 Edit body marks library

Operation:

1. Press [ENTER] key on [Edit BodyMark Type], pop edit box.

2. Input name into the new created body marks box, Move the cursor to the **[Create]** button then press **[ENTER]** key, then the new body marks will be created and appear in selected body marks list

3. Move the cursor to the **[Delete]** button then press **[ENTER]** key, then delete current body marks in the selected list.

4. Alter current annotation status list name in **[Current Type Name]** input box, press **[ENTER]** key on the **[Rename]** button, then rename the selected body marks.

7.4.2 Body mark edition

Operation:

1. Select needed source body marks via Trackball and [ENTER] key.

2. Select needed body marks at **[BodyMarkLib]** column then press the **[ENTER]** key to activate it.

3. Press [ENTER] key on [>] button to import selected body marks into user-selected status; press [ENTER] key on [>] button to move selected body marks in [BodyMark Selected] column into source body marks.

4. Press [ENTER] key on [>>] button to import all body marks in source into user-selected status; press [ENTER] key on [>>] button to move all body marks in [BodyMarks Selected] column into source.

7.4.3 Bodymark Option

Optional: whether clear bodymark upon unfreeze.

Optional: whether freeze upon bodymark.

7.5 Exam Mode

7.5.1 Exam Mode Edit

Press Utility in the submenu of **[MENU]**, the parameters control area shows as following. Press the corresponding control key to turn on or off the function.

CHISON	2016-09-12-0591 Kim	HOPE Carotid	MI:0.2 TIS:0.1 D7L40L/H9.0M	2016-09-13 10:08:23
Utility Menu PostProcess	- 0			B FPS: 41 GN: 110 PWR: 100 % FRQ: H9.0M
Slide Show Cancel	C			DN: 102 D: 3.9 cm FHI: On X-con: Nor
	- 1			
	▶2			
	- 3			
	*		1/1	
15	16			3/3
Utility Menu No	.1 No.2 Rename	No.3 No. Load Preset	4 No.5 Save Save As	🖬 ^r o 🐇
TB:Cine Play	eset Car otid			En 📃 Ϋ

- **Preset:** Display the current preset.
- **Rename:** Rename the current preset.
- Load Preset: Load the preset displayed.
- Save: Save the current preset.
- Save As: Save the current preset as others.

7.5.2 Exam Mode Selection

Choose Utility and press **[MENU]**, the following interface pops up. Open the relate function by corresponding control area.

CHISON		Setup			×
General Measurement	Comment BodyMarl	c ExamMode P	eport DICOM	Net Syste	m
ExamMode Selection	ExamMode Config	Image Freeze Co	nfig		
Probe D3P64					
ExamMode Lib					
MSK		>	ExamMode Selected		
Prostate		Card	iac [Default] iacDifficult		
Thyroid		ABD	lacDimcult		
Chest		>>	ostal		
UT & Ovary		RUQ	ostai		
UpperExtVein					
Cardiac		<			
Bladder		Peivi	C		
Vertebral Vessel		<<			
Musculoskeletal					
Breast					
ABDKidney	•				
					₽ \$
		Default	Ok	Cancel En	🗖 Å

- Probe: select needed probe and ExamMode fit for it.
- ExamModeLib: show all existed Exam modes.
- ExamMode Selected: show Exam modes in probe column.
- >: import selected Exam modes from ExamModeLib column to ExamMode Selected column.
- >>: import all Exam modes from ExamModeLib column into ExamMode Selected column.
- <: delete selected Exam modes in ExamMode Selected column.
- <<: delete all Exam modes in ExamMode Selected column.
- Set as default: set selected Exam modes in ExamMode Selected column as default
- Move up: move selected Exam modes in ExamMode Selected column up.
- Move down: move selected Exam modes ExamMode Selected column down.
- Copy to: copy the exam mode selected in ExamMode to a specified preset.

7.5.3 Exam mode Config

Operator can define needed Exam mode in detail including annotation, body marks, measurement menu import and export etc.

HIS	SON			Setup					
eneral	Measurement	Comment	BodyMark	ExamMode	Report	DICOM	Net	System	
ExamM	lode Selection	ExamMod	ie Config	Image Freez	e Config				
E	ExamMode	Comm	ient	BodyMark		Meas	surement		
MSK		Small Parts	м	SK	Gen	eral			
Prostate		Small Parts		RO	Uro	ogy			
Thyroid		Small Parts		mall Parts	Sma	allPart			
Chest		Abdomen	A	bdomen	Abd	omen			
UT & Ov	vary	GYN	G	GYN G		Gynecology			
UpperE	xtVein	Small Parts		Vessel		Vascular			
Cardiac		Cardiac	c	ardiac Cardiology					
Bladder		Small Parts		RO	Uro	ogy			
Vertebr	al Vessel	Small Parts		essel	Vas	cular			
Musculo	oskeletal	Small Parts	м		Gen	eral			
Breast		Small Parts		reast	Sma	allPart			
ABDKid	ney	Abdomen		BDKidney	Uro	ogy			
Phantor		ОВ	0	B2/3	Obs	tetric			_
KidnevL	Jreter	Abdomen		RO	Uro	oqv			
				Export		Import		Restore	
								5	, ,
				Default		Ok	Cancel	En 🖵) ¢

- ExamMode: show all existed Exam modes in system.
- Comment: double press [ENTER] key to activate widget box, can choose existed annotation

status name. After setting, the default of the Exam mode is user-selected.

- **BodyMark:** same as Comment, select user-needed default body marks.
- Measurement: same as Comment, select user-needed default measurement menu.
- **Rename:** rename selected Exam mode.
- **Delete:** delete selected Exam mode.
- Export: export all built-in Exam modes into USB flash disk.
- Import: import all built-in Exam modes into USB flash disk.
- **Restore:** restore all Exam mode as factory setting.

7.5.4 Image Freeze Config

Operator can define status after freezing image. The statuses list as below.

CHISON		Setup				×
General Measurement	Comment BodyMark	ExamMode Re	eport DICOM	Net	System	
ExamMode Selection	ExamMode Config	Image Freeze Con	fig			
Status After Freez Cine Body Mark Comment						
 Trace Distance Calc 						
Status After Real	Time Saving Cine Frames					
		Default	Ok	Cancel	🗟 妃 En 🔲	\$ ©

- Status After Freeze: options for status after freezing.
- Status After Real Time saving Cine Frames: whether unfreezing upon save or not.

7.6 Report

Users can set the layout and diagnostic template for report according to the requirement.

7.6.1 Layout Setting

CHISON	Setup					×
General Measurement Comment	BodyMark ExamMod	e Report	DICOM	Net	System	
Layout Setting Diagnosic Templ	ate					
Report Name	Report Margir					
Report Name General Report	▼ Left 0 F	light 0	Тор О	Botton	n 0	
	Selected S			Fields	Dreparties	
Title	💌 Title			Fielus	Properties	
Patient Information	🖉 Patient Inform	ation				
Diagnostic	🗹 Diagnostic					
Measurement	Measurement					
Ultrasound Image	>> 🛛 🗹 Ultrasound Im	age				
Description	Z Description					
Tips	Tips					
Bottom Information	🛛 Bottom Inform	ation				
Growth Curve						
					Save	
					1	Ŷ
	Defa		Dk	Cancel	En 📃	Ÿ.

- Report Name: display the report type , select the report which you need to set .
- **Report Margin:** set the left, right, top and bottom space for the report margin.
- Section: all the optional contents you can select to display in the report.
- Selected Sections: show the selected contents displayed in the report.
- >: import the selected option from sections to the selected sections.
- >>: import all the content from sections to the selected sections.
- <: delete the selected content in the selected sections you does not want to display in the report.
- <<: delete all the content in the selected sections.
- UP/Down (right side of the selected sections): select [Up] or [Down] to set the displayed order of the selected contents in the report.
- Fields: display the corresponding content of options in the selected sections. For example, when you select title, it will display hospital name and title, you can select the content you want to display in the report or not.
- **Properties:** adjust the column or row layout in the vertical or horizontal direction for the selected content in selected sections.

• UP/Down (on the bottom of fields): adjust the displayed order in report for selected contents in fields.

7.6.2 Diagnostic Template setting

CHISON		Setup					×
General Measurement Comment	BodyMark	ExamMode	Report	DICOM	Net	System	
Layout Setting Diagnosic Templat	te						
Template			Name				ή.
– Template				escription			
IMAGE QUALITY							
+ STUDY TYPE							
+ STUDY COMDITION							
+ ECG RHYTHM							
+ LEFT VENTRICLE							
+ LEFT ATRIUM							
+ RIGHT ATRIUM							
+ RIGHT VENTRICLE							
+ AORTIC VALVE							
+ MITRAL VALVE							
+ TRICUSPID VALVE							
+ PULMONIC VALVE							
+ PERICARDIUM							
+ AORTA							
						0 0	
						8 🖓	Ŷ
		Default	c	ж	Cancel	En 📃	Ŷ

- **Template:** show the content of all the template.
- Name: show the name of the template you have selected.
- Description: the description display the detailed content of the submenu you have selected.
 when you select the option marked with "+", it will display the submenus and select an option you want to edit.
- Add Category: add a submenu for the option in the template you have selected. For example, when you select image quality, it will display a submenu named folder. Then you can select [Name] to modify the name for the folder.
- Add Item: add items for the submenus according to the user's requirement. After you add items, you can select [Name] to modify the items' name.
- Add Comment: add comment for the items under the submenu. After you add comment, it will display comment in the description area. Then you can select the bottom area in the description to modify the name for the comment.
- Move Up/Move Down: adjust the display order for the content in the template or in the description you have selected.
- **Delete:** delete the content in the template or in the description you have selected.

• Save: after you complete the setting, select this option to save all the modification.

7.7 DICOM

DICOM includes DICOM Storage, DICOM Worklist, DICOM Print and DICOM SR. If DICOM is to be applied, please make sure DICOM has been activated. In the system page of setting interface, you can check whether DICOM is open or not. If you want to activate DICOM, please contact with CHISON.

There must be DICOM SCP server which has been installed with PACS or other relative DICOM server software.

CHISON			Setup						
General Measurem	ent Comment	BodyMark	ExamMode	Repo	ort	ЫСОМ	Net	System	
Local AE Title									
Service	AE Title			Port	Timeo	ut(s) Default	Verify		
DICOM Worklist	worklist	192.168.	30.175	105		yes			
DICOM Storage	storage	192.168.	30.175	104	10	yes			
DICOM SR	SERVER	192.168.	30.175	104		yes			
DICOM Print	PRINT_SCP	192.168.	30.62	108	10	yes			
Add									
Add DICOM Storage				DIC	Se OM SR				

- Local AE Title: input local DICOM Title to separate the DICOM equipment in local network.
- Service: display the local DICOM function worklist,
- AE Title: display the name of local DICOM AE title,
- IP: display the IP of DICOM server,
- Port: display the port of DICOM server
- Timeout(s): display the delay time
- Available: display whether DICOM is default or not,
- Verify: Press verify button and display whether DICOM setting is correct or not
- Add: add DICOM function and pop up setting dialog,

- Delete: delete the existed DICOM function
- Set as default: set one DICOM service as default,
- Verify: affirm the DICON function.
- DICOM Storage: check this item and enable DICOM storage while saving image or cine, send

clip or image according to activated function.

• **DICOM SR:** select this option the system will send DICOM SR after the operator ends the exam.

7.7.1 Add/Edit DICOM Function

Done	DICOM Setting					
Net Setting Print	Setting Print Density Print Job					
	Default					
Service	DICOM Storage	•				
Server AE Title						
Server IP Address						
SCP Port						
Timeout	10					
Compression Mode	JPEG					
Compression Ratio	Low	•				
Max Framerate	Full	•				

- Default: check this option , set DICOM function as default.
- Service: choose DICOM Storage , DICOM Worklist , DICOM Print or DICOM SR.
- Server AE Title: input DICOM server AE name.
- Server IP address: input DICOM IP address.
- SCP port: input DICOM server SCP port.
- Timeout: set the delay time of DICOM.

 \bigtriangleup **Tips:** choosing the DICOM Print type must be before relative print setting.

7.8 Network

Set the unit's and target unit's IP and do the connection testing. And network storage settings, details see the appendix: Procedures of setting network sharing .

7.9 System

	CHIS	ON				Setup					×
ſ	General	Measurer	ment	Comment	BodyMark	ExamMode	Report	DICOM	Net	System	
2		tem Inforn									
	Software v	ersion	V1.0.	7.b687702e			Update	DICOM	▼ Opene	d Key Input	
	USB versio	n	7.0.6				Update				
	DSP versio		T7.2.5	52			Update			All Functions	5
				52							
	Keyboard version V1.1.2.3				Update	Installment					
	System ver	sion	3.0.12	2 XA2.00 UA2.	00 KA2.00		Update	There are 50 days maturity			
	SN.		QBit5								
										Rent Key Input	
		dware Fun				m Maintainar					
		xport Hard	lware S	SN				NTSC) pa	L-D	
						Maintain		😿 Video Opene	d 📝 Vo	GA Opened	
	Im	port Hard	ware K	(ey				USB Video Prin	nter Options		
L										~ n	
						Defau	It	ok c	ancel		Ŷ
										En 🛄	Ŷ

7.9.1 System information

Display the software version, hardware version, system version, etc.

7.9.2 Update

Software and Hardware can be updated by USB disk.

Software upgrades File Path: X:\update_machine code\AAA.

Hardware upgrades File Path: X:\fpga_machine code\AAA.

Keyboard upgrades File Path: X:\keyboard_update_machine code\AAA.

X means USB disk. AAA means upgrades content .It should restart manually after hardware update, and after software update, machine can be restarted automatically.

7.9.3 Function Setting

Show the status of current functions, and click the related button to turn on/off this function (need the key provided by CHISON).

7.9.4 Installment setting

Input relevant key to open trial function and the detail please contact CHISON Company.

7.9.5 Video VGA

Choose the video data: NTSC and PAL-D.

Video opened: Choose the item to open this function.

VGA opened: Choose the item to open this function.

7.9.6 Hardware Function

Export hardware SN and import hardware key, it only for engineer use.

7.9.7 System Maintenance

Only an authorized service engineer may perform the maintenance.

7.9.8 USB Video Printer Option

Adjust the parameters of Video Printer Option: Dark, Light, Sharpness, and Gamma.

Select the parameters needed to adjust, press **[Done]** button on the slider of the parameter, and move the trackball to change the parameter.

Chapter 8 System Maintenance

8.1 Machine Cleaning

CAUTION: Do turn off the power before cleaning and pulling out the cable from socket. There is possibility of electric shock if the device is on.

Cleaning methods:

Use the soft dry cloth to wrap the machine. If the device is quite dirty, use wet soft cloth. After wiping the blot, use soft dry cloth to wipe.

1. Don't use organic solvent such as alcohol, otherwise surface may be ruined.

2. When cleaning the machine, don't let the liquid inflow the machine, otherwise it may malfunction and there is a danger of electric shock.

3. When it is necessary to clean the probe connector and peripheral instrument, please contact Sales office contact customer service or agent of CHISON. Any self-cleaning may result in malfunction or degrading the function of device.

8.2 Safety Check

To make sure this device works normally, a piece of maintenance plan is suggested to make to check the device regularly If any abnormal phenomena, immediately stop operation and contact with Sale Office, Customer Service department or Agents of manufacturer.


If no image or menu but image, please check following malfunction table. If the malfunction can't be solved, please contact with Sale Office, Customer Service department or Agents of manufacturer.

8.3 Malfunction Check

Serial number	Malfunction	Reason	Measures
2	Power LED lights but LED no image	The interval time is too short to restart	Restart after 1 minute
3	LED display character menu but no scan image	Launch power, gain or STC control errors Not connect to probe or the probe connection is not correct Device is in freezing condition	Control launch power , gain or STC control Make sure of right connection Exit from freezing condition
4	Abnormal image	Exam mode errors Image processing setting errors	Whether Exam mode is proper or not Adjust image processing setting or set it as default
5	Probe works improperly	1.The plug plugs loosely 2.Internal circuit protects	1.Extract the probe and reinsert 2.Restart
6	No OB calculation package menu	Unselect the OB application before scanning.	Select the OB application
7	PRINT-button doesn't work	1.The connected printer is approved2.Printer power is not on3.Printer is not connect well	1.Change the approved printer2.Turn on the printer3.Connect the printer again

Chapter 9 Probes

9.1 General Description



Convex Probe Overview

The probes provide high spatial and contrast ultrasound imaging of frequencies from 2.0MHz to 11.0MHz. These probes operate by pulsing sound waves into the body and listening to the returning echoes to produce high-resolution brightness mode, and a real time display.

9.2 Care and Maintenance

The probes that come with the system are designed to be durable and dependable. These precision instruments should be inspected daily and handled with care. Please observe the following precautions:

1. Only person Received professional training can use the probes.

2. Probes can't receive pressure sterilizer, when operation in sterile area, disposable sterile probe hood should be applied.

3. Make sure not to drop the transducer on hard surface. This can damage the transducer elements and compromise the electrical safety of the transducer.

4. Be careful when doing operation, make sure not to scratch the probe surface.

5. Avoid kinking or pinching the transducer cable.

6. Make sure not to connect the probe to plug or put adjacent cable into any kind of liquid.

7. Keep the probe clean and dry. Power off or freeze when fixing or dismantling the probe.

8. Make sure not to use or deposit the probe in the environment above 50 degree.

9. If any abnormal phenomena of probe is found, immediately stop operation and contact with Sale Office, Customer Service department or Agents of manufacturer.

10. Use only approved ultrasonic coupling gels.

11. Follow the instructions for cleaning and disinfecting that come with each probe.

12. Each time before and after ultrasound examination, please check the probe surface, probe cable and sheath whether they are abnormal, such as cracking, peeling and deformation. Also check whether the lens is strongly fixed. Abnormal probes may cause electric shock and injure the patient. Once any abnormal, user must stop using and contact CHISON's authorized service engineer.

13. Please don't use solvents (such as paint thinner, benzine, or alcohol) or abrasive cleansers for cleaning the system (including monitor and probes, etc.). It may corrode the system and probes.

14. If the probe is dropped or scratched by hard part, please stop using the probe immediately. And contact CHISON's authorized service engineer to make sure the safety and effectiveness is in good condition before use.

9.2.1 Inspecting Probes

Before and after each use, inspect carefully the probe's lens, cable, casing, and connector. Look for any damage that would allow liquid to enter the probe. If any damage is suspected, do not use the probe until it has been inspected and repaired/replaced by an authorized Service Representative.

Keep a log of all probe maintenance, along with a picture of any probe malfunction.

The probes are designed to be used only with this ultrasound system. Use of these probes on any other system or a non-qualified probe may cause electrical shock or damage on the system/transducer.

9.2.2 Cleaning and Disinfecting

Cleaning

The cleaning procedure is fit for all probes. After operation every probe should do cleaning according to stated procedure of this passage. Inspection should be done for intracavity probe depends on condition of use.

• Cleaning procedure:

- 1. Disconnect the probe from the ultrasound console and remove all coupling gel from the probe by wiping with a soft cloth and rinsing with flowing water.
- 2. Wash the probe with mild soap in lukewarm water. Scrub the probe as needed using a soft sponge, gauze, or cloth to remove all visible residue from the probe surface. Scrub the probe as needed using a soft sponge, gauze, or cloth to remove all visible residue from the probe surface.
- 3. Rinse the probe with enough clean, potable water to remove all disinfectant residues.
- 4. Use a soft cloth to clean the cable and the user section of the probe with the cleaning disinfectant liquid. Make sure that the surface of the probe and cable is wetted thoroughly with the cleaning-disinfectant.
- 5. Allow probe to air dry completely.
- 6. Reconnect the probe to the ultrasound console and place the probe into its holder.

WARNING:

To avoid electrical shock, always turn off the system and disconnect the probe before cleaning the probe.

ACAUTION:

- Take extra care when handling the lens face of the Ultrasound transducer. The lens face is especially sensitive and can easily be damaged by rough handling. NEVER use excessive force when cleaning the lens face.
- Rinse the probe with enough clean potable water to remove all visible soap residue.
- Air dry or dry with a soft cloth.

To minimize the risk of infection from blood-borne pathogens, you must handle the probe and all disposables that have contacted blood, other potentially infectious materials, mucous membranes,

and non-intact skin in accordance with infection control procedures. You must wear protective gloves when handling potentially infectious material. Use a face shield and gown if there is a risk of splashing or splatter.

Disinfection

Disinfection procedure is fit for intracavity probe.

If it is necessary to use in surgery, please abide with instruction of professional disinfection person.

ACAUTION:

- After each use, please disinfect the probes. Ultrasound probes can be disinfected using liquid chemical germicides. The level of disinfection is directly related to the duration of contact with the germicide. Increased contact time produces a higher level of disinfection.
- In order for liquid chemical germicides to be effective, all visible residue must be removed during the cleaning process. Thoroughly clean the probe, as described earlier before attempting disinfection.
- You MUST disconnect the probe from the system prior to cleaning/disinfecting the probe.
 Failure to do so could damage the system.
- DO NOT soak probes in liquid chemical germicide for longer than is stated by the germicide instructions for use. Extended soaking may cause probe damage and early failure of the enclosure, resulting in possible electric shock hazard.

• Disinfecting procedure:

1. Obey the cleaning procedure to clean the probe totally.

2. Prepare and retreat the disinfection solution according to the instruction of manufacturer.

3. Put the insertion region of probe into disinfection solution, the inserting depth can't exceed insertion region.

4. Don't let the probe connector touch any liquid.

5. Soak the probe for 3 h.

6. Pull out the probe, immediately rinse totally with sterile water and saline. Make sure no any solution remaining. Please obey the rule to do right rinsing procedure including enough rinsing water and times.

7. When probe is used in sterile area, make sure to use disposable sterile probe hood.

Ultrasound transducers can easily be damaged by improper handling and by contacting with certain chemicals. Failure to follow these precautions can result in serious injury and equipment damage.


• Do not immerse the probe into any liquid beyond the level specified for that probe. Never immerse the transducer connector or probe adapters into any liquid.

• Avoid mechanical shock or impact to the transducer and do not apply excessive bending or pulling force to the cable.

- Transducer damage can result from contacting with inappropriate coupling or cleaning agents:
- Do not soak or saturate transducers with solutions containing alcohol, bleach, ammonium chloride compounds or hydrogen peroxide.
- Avoid contacting with solutions or coupling gels containing mineral oil or lanolin.

• Avoid temperatures above 60°C. Under no circumstances should the transducer be subjected to heat sterilization method. Exposure to temperatures above 60° C will cause permanent damage to the transducer.

Inspect the probe prior to use for damage or degeneration to the housing, strain relief, lens and seal. Do not use a damaged or defective probe.



Probe Immersion Levels

These transducers are not designed to withstand heat sterilization methods. Exposure to temperatures in excess of 60 ° C will cause permanent damage. The transducers are not designed to be totally submerged in fluid, as permanent damage will result if the entire transducer is submerged.

◆ <u>Recommended materials</u>

1. Cleaning Compatibility

Solution	Origin	Qualified Use	Active Ingredients
CIDEZYME	All	Soaking	Enzyme
Epizyme-Rapid(Rapid Multi-Enzyme cleaner)	All	Soaking	Enzyme
ANIOSYME DD1	FRA	Soaking	Enzyme
Neodisher mediclean	DEU	Soaking	Enzyme
MetriZyme	USA	Soaking	Enzyme
Endozime Xtreme Power	USA	Soaking	Enzyme
Sani-Cloth AF3	USA	Wiping	Quaternary ammonium chlorides

2. Disinfecting Compatibility

Solution	Origin	Qualified Use	Active Ingredients		
Aihujia Environmental Surface Wipes	All	Wiping	Quaternary ammonium chlorides		
WIP'ANIO	FRA	Wiping	Quaternary ammonium chlorides		
Sani-Cloth Bleach Germicidal Wipes	USA	Wiping	Sodium hypochlorite		
MetriCide OPA Plus	USA	Soaking	Orthophthalal-aldehyde		
CIDEX OPA	All	Soaking	Orthophthalal-aldehyde		

ACAUTION:

- Repeated disinfection will eventually damage the transducer, please check the transducer's performance periodically.
- Other agents such as quaternary ammonium compounds are not considered high level

disinfectants and should not be used. Isopropanol is not a high level disinfectant when used as a wipe and probe manufacturers generally do not recommend soaking probes in the liquid.

- The FDA has published a list of approved sterilants and high level disinfectants for use in processing reusable medical and dental devices. That list can be consulted to find agents that may be useful for probe disinfection.
- Practitioners should consult the labels of proprietary products for specific instructions. They should also consult instrument manufacturers regarding compatibility of these agents with probes. Many of the chemical disinfectants are potentially toxic and many require adequate precautions such as proper ventilation, personal protective devices (gloves, face/eye protection, etc.) and thorough rinsing before reuse of the probe.

9.2.3 Deposit

Please replace the probe in clean and dry environment, avoid direct sunshine.

Keep the environment to deposit the probe during -10~50℃, Do not put it in high pressure and vacuum environment.

When accessing probe be careful and avoid ruin.

During transportation or leisure, the probe should be deposited in probe box.

9.2.4 Probe Safety

9.2.4.1 Handling precautions

Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use carefully when handling and protecting from damage when not in use. DO NOT use a damaged or defective probe. Failure to follow these precautions can result in serious injury and equipment damage.

• Electrical shock hazard:

The probe is driven with electrical energy that can injure the patient or user if live internal parts are contacted by conductive solution:

• DO NOT immerse the probe into any liquid beyond the level indicated by the immersion level diagram. Never immerse the probe connector into any liquid.

• Prior to each use, visually inspect the probe lens and case area for cracks, cuts, tears, and other signs of physical damage. DO NOT use a probe that appears to be damaged until you verify functional and safe performance. You need to perform a more thorough inspection, including the cable, strain relief, and connector, each time you clean the probe.

• Before inserting the connector into the probe port, inspect the probe connector pins. If a pin is bent, DO NOT use the probe until it has been inspected and repaired/replaced by a CHISON Service Representative.

• Electrical leakage checks should be performed on a routine basis by CHISON Service or qualified hospital personnel.

Mechanical hazard:

A defective probe or excess force can cause patient injury or probe damage:

• Observe depth markings and do not apply excessive force when inserting or manipulating endo-cavitary probe.

• Inspect probes for sharp edges or rough surfaces that may injure sensitive tissue.

• DO NOT apply excessive force to the probe connector when inserting into the probe port. The pin of a probe connector may bend.

Special handling instructions

Using protective sheaths:

The use of market cleared probe sheaths is recommended for clinical applications. Reference FDA March 29, 1991 "Medical Alert on Latex Products".

Protective sheaths may be required to minimize disease transmission. Probe sheaths are available for use with all clinical situations where infection is a concern. Use of legally marketed, sterile probe sheaths is strongly recommended for endo-cavitary procedures.

DO NOT use pre-lubricated condoms as a sheath. In some cases, they can damage the probe. Lubricants in these condoms may not be compatible with probe construction.

Devices containing latex may cause severe allergic reaction in latex sensitive individuals. Refer to

FDA's March 29, 1991 Medical Alert on latex products.

DO NOT use an expired probe sheath. Before using a sheath, verify if it has expired.

Endo-cavitary Probe Handling Precautions:

If the sterilization solution comes out of the endo-cavitary probe, please follow the cautions below: **Sterilant Exposure to Patient (e.g., Cidex):** Contact with a sterilant to the patient's skin for mucous membrane may cause an inflammation. If this happens, refer to instruction manual of the sterilant.

Sterilant Exposure from Probe handle to Patient (e.g. Cidex): DO NOT allow the sterilant to contact the patient. Only immerse the probe to its specified level. Ensure that no solution has entered the probe's handle before scanning the patient. If sterilant comes into contact with the patient, refer to the sterilant's instruction manual.

Sterilant Exposure from Probe connector to Patient (e.g. Cidex): DO NOT allow the sterilant to contact the patient. Only immerse the probe to its specified level. Ensure that no solution has entered the probe's connector before scanning the patient. If sterilant comes into contact with the patient, refer to the sterilant's instruction manual.

Endo-cavitary Probe Point of Contact: Refer to the sterilant's instruction manual.

9.2.4.2 Probe Handing and Infection Control

This information is intended to increase user awareness of the risks of disease transmission associated with using this equipment and provide guidance in making decisions directly affecting the safety of the patient as well as the equipment user.

Diagnostic ultrasound systems utilize ultrasound energy that must be coupled to the patient by direct physical contact.

Depending on the type of examination, this contact occurs with a variety of tissues ranging from intact skin in a routine exam to recirculating blood in a surgical procedure. The level of risk of infection varies greatly with the type of contact.

One of the most effective ways to prevent transmission between patients is with single use or disposable devices. However, ultrasound transducers are complex and expensive devices that must be reused between patients. It is very important, therefore, to minimize the risk of disease transmission by using barriers and through proper processing between patients.

Risk of Infection:

ALWAYS clean and disinfect the probe between patients to the level appropriate for the type of

examination and use FDA-cleared probe sheaths where appropriate.

Adequate cleaning and disinfection are necessary to prevent disease transmission. It is the responsibility of the equipment user to verify and maintain the effectiveness of the infection control procedures in use. Always use sterile, legally marketed probe sheaths for intra-cavitary procedures.

Coupling gels:

DO NOT use unrecommended gels (lubricants). They may damage the probe and void the warranty. AQUASONIC Gel made by R. P. Kincheloe Company in USA is recommended.

In order to assure optimal transmission of energy between the patient and probe, a conductive gel must be applied liberally to the patient where scanning will be performed.

- ✓ Please do not use any gel or other materials which are not provided by CHISON. Un-authorized gel, lubricants and other materials may corrode probes and other parts of the device, for example the keyboard. This may reduce the safety and effectiveness of the system and probes, and may also reduce the life time of the systems and probes. Damages caused by such reason will not be covered by the warranty.
- ✓ DO NOT apply gel to the eyes. If there is gel contact to the eye, flush eye thoroughly with water.

Coupling gels should not contain the following ingredients as they are known to cause probe damage:

- Methanol, ethanol, isopropanol, or any other alcohol-based product.
- Mineral oil
- Iodine
- Lotions
- Lanolin
- Aloe Vera
- Olive Oil
- Methyl or Ethyl Parabens (para hydroxybenzoic acid)
- Dimethyl silicone

9.2.4.3 Planned and Special maintenance

The following maintenance plan is suggested for the system and probes to ensure optimum operation and safety.

- Daily: inspect the probes.
- After each use: clean the probes and disinfect the probes.
- As necessary: inspect the probes, clean the probes, and disinfect the probes.

Probe Covers

The transducer should be covered with a barrier. If the barriers used are condoms, these should be nonlubricated and nonmedicated. Practitioners should be aware that condoms have been shown to be less prone to leakage than commercial probe covers, and have a six-fold enhanced AQL (acceptable quality level) when compared to standard examination gloves. They have an AQL equal to that of surgical gloves. Users should be aware of latex-sensitivity issues and have available nonlatex-containing barriers.

Returning/Shipping Probes and Repair Parts

Transportation dept. and our policy require that equipment returned for service MUST be clean and free of blood and other infectious substances.

When you returning a probe or part for service, you need to clean and disinfect the probe or part prior to packing and shipping the equipment.

Ensure that you follow probe cleaning and disinfection instructions provided in this manual.

This ensures that employees in the transportation industry as well as the people who received the package are protected from any risk.

9.2.4.4 AIUM Outlines Cleaning the Endocavotary Transducer

Guidelines for Cleaning and Preparing Endocavitary Ultrasound Transducers between Patients from AIUM *Approved June 4, 2003*

The purpose of this document is to provide guidance regarding the cleaning and disinfection of transvaginal and transrectal ultrasound probes.

All sterilization/disinfection represents a statistical reduction in the number of microbes present on a surface. Meticulous cleaning of the instrument is the essential icon to an initial reduction of the microbial/organic load by at least 99%. This cleaning is followed by a disinfecting procedure to ensure a high degree of protection from infectious disease transmission, even if a disposable barrier covers the instrument during use.

Medical instruments fall into different categories with respect to potential for infection transmission. The most critical level of instruments are those that are intended to penetrate skin or mucous membranes. These require sterilization. Less critical instruments (often called "semi-critical" instruments) that simply come into contact with mucous membranes such as fiber optic endoscopes require high-level disinfection rather than sterilization.

Although endocavitary ultrasound probes might be considered even less critical instruments because they are routinely protected by single use disposable probe covers, leakage rates of 0.9% - 2% for condoms and 8%-81% for commercial probe covers have been observed in recent studies. For maximum safety, one should therefore perform high-level disinfection of the probe between each use and use a probe cover or condom as an aid in keeping the probe clean.

There are four generally recognized categories of disinfection and sterilization. Sterilization is the complete elimination of all forms or microbial life including spores and viruses.

Disinfection, the selective removal of microbial life, is divided into three classes:

High-Level Disinfection - Destruction/removal of all microorganisms except bacterial spores.

Mid-Level Disinfection - Inactivation of Mycobacterium Tuberculosis, bacteria, most viruses, fungi, and some bacterial spores.

Low-Level Disinfection - Destruction of most bacteria, some viruses and some fungi. Low-level disinfection will not necessarily inactivate Mycobacterium Tuberculosis or bacterial spores.

The following specific recommendations are made for the use of Endocavitary ultrasound transducers. Users should also review the Centers for Disease Control and Prevention document on sterilization and disinfection of medical devices to be certain that their procedures conform to the CDC principles for disinfection of patient care equipment.

9.2.4.5 Aspetic Tchnique

For the protection of the patient and the health care worker, all endocavitary examinations should be performed with the operator properly gloved throughout the procedure. Gloves should be used to remove the condom or other barrier from the transducer and to wash the transducer as outlined above. As the barrier (condom) is removed, care should be taken not to contaminate the probe with secretions from the patient. At the completion of the procedure, hands should be thoroughly washed with soap and water.

NOTE: Obvious disruption in condom integrity does NOT require modification of this protocol. These guidelines take into account possible probe contamination due to a disruption in the barrier sheath.

In summary, routine high-level disinfection of the endocavitary probe between patients, plus the use of a probe cover or condom during each examination is required to properly protect patients from infection during endocavitary examinations. For all chemical disinfectants, precautions must be taken to protect workers and patients from the toxicity of the disinfectant.

Amis S, Ruddy M, Kibbler CC, Economides DL, MacLean AB. Assessment of condoms as probe covers for transvaginal sonography. J Clin Ultrasound 2000; 28:295-8.

Rooks VJ, Yancey MK, Elg SA, Brueske L. Comparison of probe sheaths for endovaginal sonography. Obstet. Gynecol 1996;87:27-9.

Milki AA, Fisch JD. Vaginal ultrasound probe cover leakage: implications for patient care. Fertil Steril 1998;69:409-11.

Hignett M, Claman P. High rates of perforation are found in endovaginal ultrasound probe covers before and after oocyte retrieval for in vitro fertilization-embryo transfer. J Assist Reprod Genet 1995; 12:606-9.

Sterilization and Disinfection of Medical Devices: General Principles. Centers for Disease Control, Division of Healthcare Quality Promotion. http://www.cdc.gov/ncidod/hip/sterile/sterilgp.htm (5-2003).

ODE Device Evaluation Information--FDA Cleared Sterilants and High Level Disinfectants with General Claims for Processing Reusable Medical and Dental Devices, March 2003. http://www.fda.gov/cdrh/ode/germlab.html (5-2003).

9.3 Probe Operation Instructions

For details on connecting, activating, deactivating, disconnecting, transporting and storing the probes.

9.3.1 Scanning the Patient

In order to assure optimal transmission of energy between the patient and probe, a conductive gel must be applied liberally to the patient where scanning will be performed.

After the examination is complete, follow the cleaning and disinfecting, or sterilizing procedures as appropriate.

9.3.2 Operating Transvaginal probe

The transvaginal probe is an endo-cavity probe, for the operation safety, please refer to "Care and Maintenance" for cleaning and disinfection.

The temperature at the tip of the probe displays on the screen for monitoring. No temperature above 43° C is allowed. It also depends on the patient's body temperature. When the temperature of probe tip exceeds 43° C, the probe will stop working immediately to protect the patient.

Transvaginal probe should be used with FDA approved condom or probe cover. See the following instructions to put the probe into the condom:

ACAUTION:

• Some patients may be allergic to natural rubber or medical device with rubber contains. FDA suggests that the user identify these patients and be prepared to treat allergic reactions promptly before scanning.

• Only water-solvable solutions or gel can be used. Petroleum or mineral oil-based materials may harm the cover.

• When the transvaginal probe is activated outside patient's body, its acoustic output level should be decreased to avoid any harmful interference with other equipment.

Operation Procedure:

- Put on medical sterile glove
- Get the condom for the package.

➤ Unfold the condom.

> Load some ultrasound gel into condom.

- > Take the condom with one hand, and put the probe head into the condom.
- > Fasten the condom on the end of the probe handle.

Confirm the integrity of the condom, and repeat the above steps to the condom if any damage to the condom is found.

9.3.3 Cleaning and Disinfecting TV and TR Probes

We strongly recommend wearing gloves when cleaning and disinfecting any endocavitary probe.

- Every time before and after each exam, please clean the probe handle and disinfect the transvaginal and transrectal probes probe using liquid chemical germicides
- If the probe is contaminated with body fluids, you should disinfect the probe after cleaning.
- Regard any exam waste as potentially infectious and dispose of it accordingly.

\triangle CAUTION:

- Since the probe is not waterproof, you should disconnect it from the system before cleaning or disinfecting.
- Before and after each exam, please clean the probe handle and disinfect the transvaginal and transrectal probes using liquid chemical germicides.

Cleaning

You can clean the transvaginal and transrectal probes to remove all coupling gel by wiping with a soft cloth and rinsing with flowing water. Then wash the probe with mild soap in lukewarm water. Scrub the probe as needed and use a soft cloth to remove all visible residues from the transvaginal probe surface. Rinse the probe with enough clean potable water to remove all visible soap residues, and let the probe air dry.

- ✓ Please remove the cover (if any) before cleaning the probe. (The cover like condom is one time usable).
- ✓ When cleaning the TV and TR probes, it is important to be sure that all surfaces are

thoroughly cleaned.

Disinfecting

Cidex is the only germicide that has been evaluated for compatibility with the material used to construct the probes.

To keep the effectiveness of the disinfection solutions, a thoroughly cleaning must be done to the probe before the disinfecting, make sure no residues remain on the probe.

Disinfecting Procedure:

1. Following all precautions for storage, use and disposal, prepare the germicide solution according to the manufacturer's instructions.

2. Place the cleaned and dried probe to contact with the germicide, being careful not to let the probe drop to the bottom of the container and thus damage the probe.

3. After placing/immersing, rotate and shake the probe while it is below the surface of the germicide to eliminate air pockets. Allow the germicide to remain in contact with the fully immersed probe. For high level disinfection, follow the manufacturer's recommended time.

4. Following all precautions for storage, use and disposal, prepare the germicide solution according to the manufacturer's instructions.

5. After removing from the germicide, rinse the probe according to the germicide manufacturer's rinsing instructions.

6. Flush all visible germicide residues from the probe and allow to air dry.

Appendix A: Acoustic Output Report Table

Transducer Model: <u>D3C60L</u> Operating Mode: <u>B</u>

					TIS		TIB	T 1
	Index Label		МІ	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	5
Global Maxin	num Index Value		0.8	0.97				
	P _{r.3}	(MPa)	1.32					
	Wo	(mW)		511.48				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
i diametei	z@PII.3max	(cm)	65.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	2.75	2.75				
	Dim of A _{aprt}	X (cm)		4.61				
	DITT OF Aaprt	Y (cm)		1.1				
	PD	(µsec)	0.55					
	PRF	(Hz)	9677.3					
Other	pr@PII _{max}	(MPa)	2.35					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.26				
	Focal Length	FLy (cm)		0.32				
	I _{PA.3} @ MI _{max}	(W/cm ²)	75.45					
	Mode	NA	В	В				
One should be a	Focus	(cm)	7	7				
Operating Control	Depth	(cm)	20.45	20.45				
Conditions	Scanning width	%	100	100				
	Freq	MHz	3	3				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C60L</u> Operating Mode: <u>B+C</u>

					TIS		TIB	
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.56	0.68				
	P _{r.3}	(MPa)	0.95					
	Wo	(mW)		345.6				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Faranielei	z@PII.3max	(cm)	61.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	2.9	2.92				
		X (cm)		4.61				
	Dim of A _{aprt}	Y (cm)		1.1				
	PD	(µsec)	1.32					
	PRF	(Hz)	6012.3 2					
Other	pr@PII _{max}	(MPa)	1.65					
Information	d _{eq} @PII _{max}	(cm)						
	Food Longth	FL _x (cm)		0.22				
	Focal Length	FLy (cm)		0.31				
	I _{PA.3} @ MI _{max}	(W/cm ²)	35.25					
	Mode	NA	С	С				
One station of	Focus	(cm)	6	8				
Operating Control	Depth	(cm)	11.82	9.36				
Conditions	Scanning width	%	100	100				
COnditions	Freq	MHz	3.0/3.0	3.0/3.0				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer <u>D3C60L</u> Operating Mode: <u>PW</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	ті С
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.5			1.77	1.02	
	P _{r.3}	(MPa)	0.84					
	Wo	(mW)					588.55	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				276.05		
	Z1	(cm)				4.1		
Associated Acoustic	Z _{bp}	(cm)				3.8		
Parameter	Z _{sp}	(cm)					6.45	
Falametei	z@PII.3max	(cm)	44					
	d _{eq} (z _{sp})	(cm)					0.67	
	f _c	(MHz)	2.88			2.87	2.88	
		X (cm)				6.91	6.91	
	Dim of A _{aprt}	Y (cm)				1.1	1.1	
	PD	(µsec)	1.29					
	PRF	(Hz)	9658.7 1					
Other	pr@PII _{max}	(MPa)	1.45					
Information	d _{eq} @PII _{max}	(cm)					0.65	
	Focal Length	FL _x (cm)				0.32		
	Focal Length	FLy (cm)				0.41		
	I _{PA.3} @ MI _{max}	(W/cm ²)	189.78					
	Mode	NA	PW			PW	PW	
One sections	Focus	(cm)	6			7	10	
Operating Control	Depth	(cm)	9.36			8.13	11.82	
	Scanning width	%	100			100	100	
Conditions	Freq	MHz	3.0/3.0			3.0/3.0	3.0/3.0	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C60L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	num Index Value		0.95			1.02	1.1	
	P _{r.3}	(MPa)	1.57					
	Wo	(mW)					511.23	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				227.56		
	Z1	(cm)				4.2		
Associated Acoustic	Z _{bp}	(cm)				3.8		
Parameter	Z _{sp}	(cm)					5.1	
Falameter	z@PII.3max	(cm)	56					
	d _{eq} (z _{sp})	(cm)					0.81	
	f _c	(MHz)	2.74			2.78	2.78	
	Dim of A _{aprt}	X (cm)				6.91	6.91	
		Y (cm)				1.1	1.1	
	PD	(µsec)	0.55					
	PRF	(Hz)	4759.6 2					
Other	pr@PII _{max}	(MPa)	2.64					
Information	d _{eq} @PII _{max}	(cm)					0.81	
	Easel Longth	FL _x (cm)				0.23		
	Focal Length	FLy (cm)				0.18		
	I _{PA.3} @ MI _{max}	(W/cm ²)	112.05					
	Mode	NA	М			М	М	
Onenting	Focus	(cm)	7			6	6	
Operating Control	Depth	(cm)	8.13			6.89	6.89	
	Scanning width	%	100			100	100	
Conditions -	Freq	MHz				227.56		
	Power	(%)				4.2		

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L</u> Operating Mode: <u>B</u>

	<u> </u>				TIS		TIB	ті
	Index Label		МІ	Scan	non-	scan	non-sca	C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	num Index Value		0.97	0.53				
	P _{r.3}	(MPa)	2.41					
	Wo	(mW)		73.73				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falametei	z@PII.3max	(cm)	19.95					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	6.23	6.23				
	Dim of A _{aprt}	X (cm)		2.05				
	Dim Of Aaprt	Y (cm)		0.45				
	PD	(µsec)	0.23					
	PRF	(Hz)	2873.6					
Other	pr@PII _{max}	(MPa)	3.37					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.26				
	T OCAI Lengti	FLy (cm)		0.36				
	I _{PA.3} @ MI _{max}	(W/cm ²)	43.96					
	Mode	NA	В	В				
Operating	Focus	(cm)	2	2				
Operating	Depth	(cm)	4.93	4.93				
Control Conditions	Scanning width	%	100	100				
	Freq	MHz	5	5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L</u> Operating Mode: <u>B +C</u>

	<u> </u>				TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	U
Global Maxim	num Index Value		0.55	0.79				
	P _{r.3}	(MPa)	1.38					
	Wo	(mW)		106.5				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Parameter	z@PII.3max	(cm)	16.5					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	6.36	6.37				
	Dim of A _{aprt}	X (cm)		2.05				
	Dim Of Aaprt	Y (cm)		0.45				
	PD	(µsec)	0.77					
	PRF	(Hz)	6944.4					
Other	pr@PII _{max}	(MPa)	1.94					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.23				
		FLy (cm)		0.35				
	I _{PA.3} @ MI _{max}	(W/cm ²)	91.4					
	Mode	NA	С	С				
Operating Control Conditions	Focus	(cm)	2.5	1				
	Depth	(cm)	4.9	3.7				
	Scanning width	%	100	100				
	Freq	MHz	5.0/6.5	5.0/6.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L</u> Operating Mode: <u>PW</u>

					TIS		TIB	ті
	Index Label		мі	Scan	non-	scan	non-sca	C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	um Index Value		0.63		0.56		0.88	
	P _{r.3}	(MPa)	1.59					
	Wo	(mW)			221.18		73.73	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)					1.85	
	z@PII. _{3max}	(cm)	17					
	d _{eq} (z _{sp})	(cm)					0.44	
	fc	(MHz)	6.34		6.23		6.34	
	Dim of A _{aprt}	X (cm)			4.1		4.1	
		Y (cm)			0.45		0.45	
	PD	(µsec)	0.59					
	PRF	(Hz)	2873.6					
Other	pr@PII _{max}	(MPa)	2.14					
Information	d _{eq} @PII _{max}	(cm)					0.44	
Information	Focal Length	FL _x (cm)			0.24			
	T Ocar Length	FLy (cm)			0.33			
	I _{PA.3} @ MI _{max}	(W/cm ²)	118.84					
	Mode	NA	PW		PW		PW	
Operating	Focus	(cm)	4		2.5		2.5	
Operating	Depth	(cm)	7.39		3.7		7.39	
Control Conditions	Scanning width	%	100		100		100	
	Freq	MHz	5.0/6.5		5.0/6.5		5.0/6.5	
	Power	(%)	100		100		100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L</u> Operating Mode: <u>B +M</u>

	node. <u>D i m</u>				TIS		TIB	
	Index Label		мі	0	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.77		1.1		1	
	P _{r.3}	(MPa)	1.91					
	Wo	(mW)			139.3		139.3	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)					1.8	
i diameter	z@PII. _{3max}	(cm)	16.5					
	d _{eq} (z _{sp})	(cm)					0.89	
	f _c	(MHz)	6.19		6.22		6.19	
	Dim of A _{aprt}	X (cm)			4.1		4.1	
		Y (cm)			0.45		0.45	
	PD	(µsec)	0.22					
	PRF	(Hz)	2873.6					
Other	pr@PII _{max}	(MPa)	2.6					
Information	d _{eq} @PII _{max}	(cm)					0.89	
mornation	Focal Length	FL _x (cm)			0.22			
	T Ocal Length	FLy (cm)			0.36			
	I _{PA.3} @ MI _{max}	(W/cm ²)	158.13					
	Mode	NA	М		М		М	
Operating	Focus	(cm)	2		1.5		2	
Operating Control	Depth	(cm)	3.7		3.7		3.7	
Conditions	Scanning width	%	100		100		100	
Conditions	Freq	MHz	5.0/6.5		5.0/6.5		5.0/6.5	
	Power	(%)	100		100		100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L-REC</u> Operating Mode: <u>B</u>

					TIS		TIB	- 1
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	J
Global Maxim	um Index Value		0.84	0.52				
	P _{r.3}	(MPa)	2.08					
	Wo	(mW)		73.28				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falameter	z@PII.3max	(cm)	19.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	6.14	6.14				
	Dim of A _{aprt}	X (cm)		2.05				
	Dim Of Aaprt	Y (cm)		0.45				
	PD	(µsec)	0.22					
	PRF	(Hz)	2865.3					
Other	pr@PII _{max}	(MPa)	3.07					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.24				
		FLy (cm)		0.33				
	I _{PA.3} @ MI _{max}	(W/cm ²)	78.21					
	Mode	NA	В	В				
Operating	Focus	(cm)	3	3				
Operating Control Conditions	Depth	(cm)	4.93	4.93				
	Scanning width	%	100	100				
	Freq	MHz	5	5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L-REC</u> Operating Mode: <u>B +C</u>

					TIS		TIB	- ,
	Index Label		МІ	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.49	0.78				
	P _{r.3}	(MPa)	1.25					
	Wo	(mW)		105.6				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Faranielei	z@PII.3max	(cm)	16.5					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	6.36	6.35				
	Dim of A _{aprt}	X (cm)		2.05				
	Dim Of Aaprt	Y (cm)		0.45				
	PD	(µsec)	0.75					
	PRF	(Hz)	6855.2					
Other	pr@PII _{max}	(MPa)	1.85					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.19				
		FLy (cm)		0.13				
	I _{PA.3} @ MI _{max}	(W/cm ²)	88.5					
	Mode	NA	С	С				
Operating	Focus	(cm)	2	2.5				
Operating Control	Depth	(cm)	4.93	3.7				
Conditions	Scanning width	%	100	100				
	Freq	MHz	5.0/6.5	5.0/6.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L-REC</u> Operating Mode: <u>PW</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.56		1.02		0.8	
	P _{r.3}	(MPa)	1.39					
	Wo	(mW)			73.73		73.73	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)					1.65	
i diameter	z@PII. _{3max}	(cm)	16.5					
	d _{eq} (z _{sp})	(cm)					0.45	
	fc	(MHz)	6.32		6.32		6.3	
	Dim of A _{aprt}	X (cm)			4.1		4.1	
	DIIII OI Aaprt	Y (cm)			0.45		0.45	
	PD	(µsec)	0.35					
	PRF	(Hz)	5877					
Other	pr@PII _{max}	(MPa)	1.9					
Information	d _{eq} @PII _{max}	(cm)					0.45	
Information	Focal Length	FL _x (cm)			0.25			
	Focal Length	FLy (cm)			0.18			
	IPA.3@ MImax	(W/cm ²)	120.41					
	Mode	NA	PW		PW		PW	
Onenetin	Focus	(cm)	2		2		3	
Operating Control	Depth	(cm)	7.39		7.39		7.39	
	Scanning width	%	100		100		100	
Conditions	Freq	MHz	5.0/6.5		5.0/6.5		5.0/6.5	
	Power	(%)	100		100		100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7L40L-REC</u> Operating Mode: <u>B +M</u>

					TIS		TIB	
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.69		1.1		0.49	
	P _{r.3}	(MPa)	1.71					
	Wo	(mW)			136.5		139.3	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)					1.86	
Falanielei	z@PII.3max	(cm)	24.52					
	d _{eq} (z _{sp})	(cm)					0.95	
	f _c	(MHz)	6.05		6.18		6.05	
	Dim of A _{aprt}	X (cm)			4.1		4.1	
	Dim Of Aaprt	Y (cm)			0.45		0.45	
	PD	(µsec)	0.22					
	PRF	(Hz)	2655.3					
Other	pr@PII _{max}	(MPa)	2.52					
Information	d _{eq} @PII _{max}	(cm)					0.89	
mornation	Focal Length	FL _x (cm)			0.23			
	T Ocal Length	FLy (cm)			0.34			
	I _{PA.3} @ MI _{max}	(W/cm ²)	134.51					
	Mode	NA	М		М		М	
Operating	Focus	(cm)	2		2		2	
Operating Control	Depth	(cm)	7.39		3.7		7.39	
Control Conditions	Scanning width	%	100		100		100	
	Freq	MHz	5.0/6.5		5.0/6.5		5.0/6.5	
	Power	(%)	100		100		100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D5C20L</u> Operating Mode: <u>B</u>

					TIS		TIB	T 1
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.89	1.06				
	P _{r.3}	(MPa)	1.96					
	Wo	(mW)		350.21				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falametei	z@PII.3max	(cm)	24.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	4.84	5				
	Dim of A _{aprt}	X (cm)		1.46				
	Diff Of Aaprt	Y (cm)		1				
	PD	(µsec)	0.43					
	PRF	(Hz)	14706					
Other	pr@PII _{max}	(MPa)	2.69					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.28				
		FLy (cm)		0.21				
	I _{PA.3} @ MI _{max}	(W/cm ²)	80.63					
	Mode	NA	В	В				
Onerations	Focus	(cm)	4	1.5				
Operating Control	Depth	(cm)	11.77	5.61				
	Scanning width	%	100	100				
Conditions	Freq	MHz	3.5	3.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D5C20L</u> Operating Mode: <u>B+C</u>

					TIS		TIB	ті
	Index Label		МІ	Scan	non-	scan	non-sca	C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	5
Global Maxin	num Index Value		0.7	1.1				
	P _{r.3}	(MPa)	1.39					
	Wo	(mW)		321.02				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Farameter	z@PII.3max	(cm)	28.5					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	4.01	4.01				
	Dim of A _{aprt}	X (cm)		1.46				
		Y (cm)		1				
	PD	(µsec)	0.96					
	PRF	(Hz)	5952.4					
Other	pr@PII _{max}	(MPa)	1.7					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.18				
	T Ocar Length	FLy (cm)		0.26				
	I _{PA.3} @ MI _{max}	(W/cm ²)	58.37					
	Mode	NA	С	С				
Operating	Focus	(cm)	3.5	3.5				
Operating Control Conditions	Depth	(cm)	9.3	9.3				
	Scanning width	%	100	100				
	Freq	MHz	3.5/4.0	3.5/4.0				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D5C20L</u> Operating Mode: <u>PW</u>

	<u></u>				TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	1um Index Value		0.8			1.13	0.45	
	P _{r.3}	(MPa)	1.6					
	Wo	(mW)					349.43	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				256.69		
	Z1	(cm)				2.05		
Associated Acoustic	Z _{bp}	(cm)				2.04		
Parameter	Z _{sp}	(cm)					2.2	
Farameter	z@PII.3max	(cm)	23.5					
	d _{eq} (z _{sp})	(cm)					0.72	
	fc	(MHz)	3.93			3.96	3.96	
	Dim of A _{aprt}	X (cm)				2.92	2.92	
		Y (cm)				1	1	
	PD	(µsec)	0.95					
	PRF	(Hz)	6024.1					
Other	pr@PII _{max}	(MPa)	2.1					
Information	d _{eq} @PII _{max}	(cm)					0.7	
mornation	Focal Length	FL _x (cm)				0.16		
	T Ocar Length	FLy (cm)				0.23		
	I _{PA.3} @ MI _{max}	(W/cm ²)	101.33					
	Mode	NA	PW			PW	PW	
Operating	Focus	(cm)	2			3	3	
Operating Control Conditions	Depth	(cm)	4.37			4.37	4.37	
	Scanning width	%	100			100	100	
	Freq	MHz	3.5/4.0			3.5/4.0	3.5/4.0	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D5C20L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.93			1.15	1.01	
	P _{r.3}	(MPa)	2.02					
	Wo	(mW)					356.04	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				177.51		
Associated	Z1	(cm)				2.05		
Associated Acoustic	Z _{bp}	(cm)				2.04		
Parameter	Z _{sp}	(cm)					2.05	
Falameter	z@PII.3max	(cm)	20.5					
	d _{eq} (z _{sp})	(cm)					0.74	
	f _c	(MHz)	4.75			4.92	4.83	
	Dim of A _{aprt}	X (cm)				2.91	2.92	
	Diff Of Aaprt	Y (cm)				1	1	
	PD	(µsec)	0.42					
	PRF	(Hz)	14706					
Other	pr@PII _{max}	(MPa)	2.83					
Information	d _{eq} @PII _{max}	(cm)					0.74	
mornation	Focal Length	FL _x (cm)				0.33		
		FLy (cm)				0.26		
	I _{PA.3} @ MI _{max}	(W/cm ²)	111.39					
	Mode	NA	М			М	М	
Operating	Focus	(cm)	3.5			2.5	2	
Operating Control Conditions	Depth	(cm)	4.37			3.14	5.61	
	Scanning width	%	100			100	100	
	Freq	MHz	3.5			3.5	3.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C20L</u> Operating Mode: <u>B</u>

					TIS		TIB	T 1
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.72	1.27				
	P _{r.3}	(MPa)	1.3					
	Wo	(mW)		365.57				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falanetei	z@PII.3max	(cm)	49.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	3.27	3.44				
	Dim of A _{aprt}	X (cm)		2.05				
	Dim Of Aaprt	Y (cm)		1.3				
	PD	(µsec)	0.6					
	PRF	(Hz)	3846.2					
Other	pr@PII _{max}	(MPa)	2.22					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.26				
		FLy (cm)		0.25				
	I _{PA.3} @ MI _{max}	(W/cm ²)	55.73					
	Mode	NA	В	В				
Operating	Focus	(cm)	5	7				
Operating Control Conditions	Depth	(cm)	8.84	10.7				
	Scanning width	%	100	100				
	Freq	MHz	3.5	3.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C20L</u> Operating Mode: <u>B+C</u>

	<u> </u>				TIS		TIB	ті
	Index Label		МІ	Scan	non-	scan	non-sca	C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	num Index Value		0.3	0.5				
	P _{r.3}	(MPa)	0.52					
	Wo	(mW)		136.19				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falametei	z@PII.3max	(cm)	39					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	2.84	2.84				
	Dim of A _{aprt}	X (cm)		2.05				
		Y (cm)		1.3				
	PD	(µsec)	0.91					
	PRF	(Hz)	5952.4					
Other	pr@PII _{max}	(MPa)	0.8					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.23				
		FLy (cm)		0.21				
	I _{PA.3} @ MI _{max}	(W/cm ²)	98.48					
	Mode	NA	С	С				
Operating	Focus	(cm)	5	5				
Operating Control	Depth	(cm)	7.61	7.61				
Conditions	Scanning width	%	100	100				
	Freq	MHz	3.5/3.5	3.5/3.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C20L</u> Operating Mode: <u>PW</u>

	<u></u>				TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	5
Global Maxim	num Index Value		0.32			0.5	0.8	
	P _{r.3}	(MPa)	0.54					
	Wo	(mW)					78.8	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				36.09		
	Z1	(cm)				3.9		
Associated Acoustic	Z _{bp}	(cm)				2.76		
Parameter	Z _{sp}	(cm)					3.5	
Falametei	z@PII.3max	(cm)	56					
	d _{eq} (z _{sp})	(cm)					1.13	
	f _c	(MHz)	2.9			2.9	2.89	
	Dim of A _{aprt}	X (cm)				3.58	3.58	
	DITT OF Aaprt	Y (cm)				1.3	1.3	
	PD	(µsec)	0.9					
	PRF	(Hz)	3846.2					
Other	pr@PII _{max}	(MPa)	0.75					
Information	d _{eq} @PII _{max}	(cm)					1.07	
mornation	Focal Length	FL _x (cm)				0.28		
	Focal Length	FLy (cm)				0.32		
	I _{PA.3} @ MI _{max}	(W/cm ²)	110.43					
	Mode	NA	PW			PW	PW	
Operation	Focus	(cm)	5			5	5	
Operating Control Conditions	Depth	(cm)	8.84			8.84	17.46	
	Scanning width	%	100			100	100	
	Freq	MHz	3.5/3.5			3.5/3.5	3.5/3.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3C20L</u> Operating Mode: <u>B+M</u>

	<u></u>				TIS		TIB	ті
	Index Label		МІ	Scan	non-	scan	non-sca	C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	5
Global Maxim	num Index Value	_	0.7			1.4	0.5	
	P _{r.3}	(MPa)	1.27					
	Wo	(mW)					487.42	
Associated	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				219.13		
	Z1	(cm)				3.55		
Acoustic	Z _{bp}	(cm)				2.76		
Parameter	Z _{sp}	(cm)					4.3	
Parameter	z@PII.3max	(cm)	46.5					
	d _{eq} (z _{sp})	(cm)					0.93	
	f _c	(MHz)	3.26			3.26	3.26	
	Dim of A _{aprt}	X (cm)				3.58	3.58	
		Y (cm)				1.3	1.3	
	PD	(µsec)	0.5					
	PRF	(Hz)	7575.8					
Other	pr@PII _{max}	(MPa)	1.92					
Information	d _{eq} @PII _{max}	(cm)					0.91	
mornation	Focal Length	FL _x (cm)				0.36		
		FLy (cm)				0.26		
	I _{PA.3} @ MI _{max}	(W/cm ²)	162.39					
	Mode	NA	М			М	М	
Operating	Focus	(cm)	5			5	5	
Operating Control Conditions	Depth	(cm)	6.37			6.37	6.37	
	Scanning width	%	100			100	100	
	Freq	MHz	3.5			3.5	3.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C15L</u> Operating Mode: <u>B</u>

	<u></u>				TIS		TIB	T 1
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	J
Global Maxin	num Index Value		0.81	1.24				
	P _{r.3}	(MPa)	1.93					
	Wo	(mW)		409.6				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falametei	z@PII.3max	(cm)	27					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	5.66	5.88				
	Dim of A _{aprt}	X (cm)		1.22				
	Dim Of Aaprt	Y (cm)		0.9				
	PD	(µsec)	0.26					
	PRF	(Hz)	5924.2					
Other	pr@PII _{max}	(MPa)	2.77					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.32				
	Focal Length	FLy (cm)		0.35				
	IPA.3@ MImax	(W/cm ²)	83.82					
	Mode	NA	В	В				
Onenating	Focus	(cm)	2.5	2				
Operating Control Conditions	Depth	(cm)	5.69	6.92				
	Scanning width	%	100	100				
	Freq	MHz	4	4				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C15L</u> Operating Mode: <u>B+C</u>

					TIS		TIB	T 1
	Index Label		мі	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.75	1.11				
	P _{r.3}	(MPa)	1.72					
	Wo	(mW)		296.7				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Assesiated	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Farameter	z@PII.3max	(cm)	18					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	5.26	5.36				
	Dim of A _{aprt}	X (cm)		1.22				
	Diffi Of Aaprt	Y (cm)		0.9				
	PD	(µsec)	0.72					
	PRF	(Hz)	6024.1					
Other	pr@PII _{max}	(MPa)	2.39					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.26				
		FLy (cm)		0.34				
	IPA.3@ MImax	(W/cm ²)	94.87					
	Mode	NA	С	С				
Operating	Focus	(cm)	2.5	1				
Operating Control Conditions	Depth	(cm)	11.85	3.22				
	Scanning width	%	100	100				
	Freq	MHz	4.0/5.3	4.0/5.3				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C15L</u> Operating Mode: <u>PW</u>

	<u></u>				TIS		TIB	
	Index Label		МІ	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.69			1.68	0.66	
	P _{r.3}	(MPa)	1.58					
	Wo	(mW)					204.29	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				104.65		
	Z1	(cm)				1.8		
Associated Acoustic	Z _{bp}	(cm)				1.77		
Parameter	Z _{sp}	(cm)					1.8	
Parameter	z@PII.3max	(cm)	28.5					
	d _{eq} (z _{sp})	(cm)					0.7	
	f _c	(MHz)	5.21			5.38	5.38	
	Dim of A _{aprt}	X (cm)				2.43	2.43	
	Diff Of Aaprt	Y (cm)				0.9	0.9	
	PD	(µsec)	0.5					
	PRF	(Hz)	5924.2					
Other	pr@PII _{max}	(MPa)	2.4					
Information	d _{eq} @PII _{max}	(cm)					0.7	
mornation	Focal Length	FL _x (cm)				0.31		
		FLy (cm)				0.29		
	I _{PA.3} @ MI _{max}	(W/cm ²)	138.48					
	Mode	NA	PW			PW	PW	
Onenations	Focus	(cm)	2.5			3.5	3.5	
Operating Control Conditions	Depth	(cm)	3.22			4.46	4.46	
	Scanning width	%	100			100	100	
	Freq	MHz	4.0/5.3			4.0/5.3	4.0/5.3	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C15L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.8			1.33	0.9	
	P _{r.3}	(MPa)	1.79					
	Wo	(mW)					272.38	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				128.38		
Associated	Z1	(cm)				2		
Associated Acoustic	Z _{bp}	(cm)				1.77		
Parameter	Z _{sp}	(cm)					2	
Parameter	z@PII.3max	(cm)	35.5					
	d _{eq} (z _{sp})	(cm)					0.78	
	fc	(MHz)	4.92			5.44	4.99	
	Dim of A _{aprt}	X (cm)				2.43	2.43	
	DITI OF Aaprt	Y (cm)				0.9	0.9	
	PD	(µsec)	0.28					
	PRF	(Hz)	11236					
Other	pr@PII _{max}	(MPa)	2.44					
Information	d _{eq} @PII _{max}	(cm)					0.65	
mornation	Focal Length	FL _x (cm)				0.22		
	Focal Length	FLy (cm)				0.31		
	I _{PA.3} @ MI _{max}	(W/cm ²)	180.8					
	Mode	NA	М			М	М	
Operating Control Conditions	Focus	(cm)	4			2	3.5	
	Depth	(cm)	10.62			4.46	8.15	
	Scanning width	%	100			100	100	
	Freq	MHz	4			4	4	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7C10L</u> Operating Mode: <u>B</u>

Index Label Global Maximum Index Value			МІ	TIS			TIB	
				Scan	non-scan		non-sca	TI C
					A _{aprt} ≤1	A _{aprt} >1	n	
			0.48	0.2				
Associated	P _{r.3}	(MPa)	0.99					
	Wo	(mW)		29.7				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
	Z _{bp}	(cm)						
Acoustic Parameter	Z _{sp}	(cm)						
Parameter	z@PII.3max	(cm)	22.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	4.15	4.18				
	Dim of A _{aprt}	X (cm)		1.48				
		Y (cm)		0.7				
Other Information	PD	(µsec)	0.67					
	PRF	(Hz)	4347.8					
	pr@PII _{max}	(MPa)	1.3					
	d _{eq} @PII _{max}	(cm)						
	Focal Length	FL _x (cm)		0.22				
		FLy (cm)		0.36				
	IPA.3@ MImax	(W/cm ²)	41.39					
Operating Control Conditions	Mode	NA	В	В				
	Focus	(cm)	2.5	1.5				
	Depth	(cm)	3.93	2.7				
	Scanning width	%	100	100				
	Freq	MHz	4	4				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7C10L</u> Operating Mode: <u>B+C</u>

Index Label Global Maximum Index Value			МІ	TIS			TIB	
				Scan	non-scan		non-sca	TI C
					A _{aprt} ≤1	A _{aprt} >1	n	
			0.41	0.65				
Associated	P _{r.3}	(MPa)	0.95					
	Wo	(mW)		77.2				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
	Z _{bp}	(cm)						
Acoustic Parameter	Z _{sp}	(cm)						
Parameter	z@PII.3max	(cm)	24.5					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	5.27	5.28				
	Dim of A _{aprt}	X (cm)		1.48				
		Y (cm)		0.7				
	PD	(µsec)	0.7					
Other Information	PRF	(Hz)	6024.1					
	pr@PII _{max}	(MPa)	1.38					
	d _{eq} @PII _{max}	(cm)						
	Focal Length	FL _x (cm)		0.21				
		FLy (cm)		0.29				
	I _{PA.3} @ MI _{max}	(W/cm ²)	30.2					
Operating Control Conditions	Mode	NA	С	С				
	Focus	(cm)	2.5	2				
	Depth	(cm)	10.09	3.93				
	Scanning width	%	100	100				
	Freq	MHz	4	4				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7C10L</u> Operating Mode: <u>PW</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.47			1.34	1.18	
	P _{r.3}	(MPa)	1.07					
	Wo	(mW)					100.97	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				53.19		
Assesiated	Z1	(cm)				1.75		
Associated Acoustic	Z _{bp}	(cm)				1.72		
Parameter	Z _{sp}	(cm)					2	
Falanielei	z@PII.3max	(cm)	26					
	d _{eq} (z _{sp})	(cm)					0.51	
	fc	(MHz)	5.3			5.3	5.28	
	Dim of A _{aprt}	X (cm)				2.97	2.97	
		Y (cm)				0.7	0.7	
	PD	(µsec)	0.87					
	PRF	(Hz)	5555.6					
Other	pr@PII _{max}	(MPa)	1.52					
Information	d _{eq} @PII _{max}	(cm)					0.5	
Information	Focal Length	FL _x (cm)				0.26		
	Focal Length	FLy (cm)				0.24		
	IPA.3@ MImax	(W/cm ²)	143.33					
	Mode	NA	PW			PW	PW	
Operating Control Conditions	Focus	(cm)	3			2.5	2	
	Depth	(cm)	8.86			6.39	2.7	
	Scanning width	%	100			100	100	
	Freq	MHz	4.0/5.3			4.0/5.3	4.0/5.3	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D7C10L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		МІ	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		0.48			0.57	0.24	
	P _{r.3}	(MPa)	0.98					
	Wo	(mW)					47.51	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				28.7		
A	Z1	(cm)				1.75		
Associated Acoustic	Z _{bp}	(cm)				1.72		
Parameter	Z _{sp}	(cm)					2.3	
Farameter	z@PII.3max	(cm)	24					
	d _{eq} (z _{sp})	(cm)					0.46	
	f _c	(MHz)	4.13			4.17	4.13	
	Dim of A _{aprt}	X (cm)				2.97	2.97	
	DITI OF Aaprt	Y (cm)				0.7	0.7	
	PD	(µsec)	0.67					
	PRF	(Hz)	5555.6					
Other	pr@PII _{max}	(MPa)	1.3					
Information	d _{eq} @PII _{max}	(cm)					0.45	
mornation	Focal Length	FL _x (cm)				0.26		
	Focal Length	FLy (cm)				0.19		
	I _{PA.3} @ MI _{max}	(W/cm ²)	140.31					
	Mode	NA	М			М	М	
Operating Control Conditions	Focus	(cm)	3			2	3	
	Depth	(cm)	10.09			5.16	10.09	
	Scanning width	%	100			100	100	
	Freq	MHz	4			4	4	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C12L</u> Operating Mode: <u>B</u>

					TIS		TIB	T 1
	Index Label		МІ	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.66	0.42				
	P _{r.3}	(MPa)	1.46					
	Wo	(mW)		43.67				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Falametei	z@PII.3max	(cm)	25					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	4.9	5.13				
	Dim of A _{aprt}	X (cm)		1.21				
	DITT OF Aaprt	Y (cm)		0.9				
	PD	(µsec)	0.31					
	PRF	(Hz)	7286.5					
Other	pr@PII _{max}	(MPa)	1.93					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.26				
		FLy (cm)		0.24				
	IPA.3@ MImax	(W/cm ²)	62.56					
	Mode	NA	В	В				
Operating Control	Focus	(cm)	3	2				
	Depth	(cm)	4.29	3.06				
Conditions	Scanning width	%	100	100				
	Freq	MHz	5.3	5.3				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C12L</u> Operating Mode: <u>B+C</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.45	0.83				
	P _{r.3}	(MPa)	0.91					
	Wo	(mW)		76.87				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
Faranielei	z@PII.3max	(cm)	33.5					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	4.09	5.24				
		X (cm)		1.21				
	Dim of A _{aprt}	Y (cm)		0.9				
	PD	(µsec)	1.16					
	PRF	(Hz)	5000.2 4					
Other	pr@PII _{max}	(MPa)	1.26					
Information	d _{eq} @PII _{max}	(cm)						
	Focal Length	FL _x (cm)		0.24				
	Focal Length	FLy (cm)		0.33				
	I _{PA.3} @ MI _{max}	(W/cm ²)	34.36					
	Mode	NA	С	С				
Operating Control Conditions	Focus	(cm)	4	2.5				
	Depth	(cm)	6.75	3.06				
	Scanning width	%	100	100				
Conditions	Freq	MHz	5.3/4.0	5.3/4.0				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C12L</u> Operating Mode: <u>PW</u>

	<u></u>				TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	5
Global Maxim	num Index Value		0.49			0.93	0.87	
	P _{r.3}	(MPa)	1.11					
	Wo	(mW)					71.86	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				38.11		
Associated	Z1	(cm)				1.77		
Associated	Z _{bp}	(cm)				1.76		
Parameter	Z _{sp}	(cm)					1.93	
Falanielei	z@PII.3max	(cm)	24					
	d _{eq} (z _{sp})	(cm)					0.42	
	f _c	(MHz)	5.16			5.21	5.18	
	Dim of A _{aprt}	X (cm)				2.43	2.43	
	DIIII OI Aaprt	Y (cm)				0.9	0.9	
	PD	(µsec)	0.69					
	PRF	(Hz)	7282.6					
Other	pr@PII _{max}	(MPa)	1.56					
Information	d _{eq} @PII _{max}	(cm)					0.43	
mornation	Focal Length	FL _x (cm)				0.24		
	Focal Length	FLy (cm)				0.23		
	I _{PA.3} @ MI _{max}	(W/cm ²)	142.87					
	Mode	NA	PW			PW	PW	
Onersting	Focus	(cm)	3			3.5	3.5	
Operating Control Conditions	Depth	(cm)	3.06			9.22	5.52	
	Scanning width	%	100			100	100	
	Freq	MHz	5.3/4.0			5.3/4.0	5.3/4.0	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D6C12L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		мі	Coor	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.78			0.3	1.12	
	P _{r.3}	(MPa)	1.72					
	Wo	(mW)					24.26	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				13.21		
	Z1	(cm)				1.78		
Associated Acoustic	Z _{bp}	(cm)				1.77		
Parameter	Z _{sp}	(cm)					2.07	
Falametei	z@PII.3max	(cm)	23					
	d _{eq} (z _{sp})	(cm)					0.21	
	f _c	(MHz)	4.83			4.86	4.83	
	Dim of A _{aprt}	X (cm)				2.43	2.43	
	Diffi Of Aaprt	Y (cm)				0.9	0.9	
	PD	(µsec)	0.29					
	PRF	(Hz)	7287.6					
Other	pr@PII _{max}	(MPa)	2.27					
Information	d _{eq} @PII _{max}	(cm)					0.2	
mornation	Focal Length	FL _x (cm)				0.32		
		FLy (cm)				0.24		
	IPA.3@ MImax	(W/cm ²)	120.53					
	Mode	NA	М			М	М	
Operating Control Conditions	Focus	(cm)	3			3	3	
	Depth	(cm)	4.29			10.45	5.22	
	Scanning width	%	100			100	100	
	Freq	MHz	5.3			5.3	5.3	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display $\left(\left(a \right)^{2} \right)^{2}$

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>B</u>

					TIS		TIB	ті
	Index Label		МІ	Scan	non-	scan	non-sca	C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	num Index Value		1.08	1.24				
	P _{r.3}	(MPa)	1.63					
	Wo	(mW)		269.95				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
i arameter	z@PII. _{3max}	(cm)	45					
	$d_{eq}(z_{sp})$	(cm)						
	fc	(MHz)	2.28	2.48				
	Dim of A _{aprt}	X (cm)		1.92				
	Diffi Of Aaprt	Y (cm)		1.58				
	PD	(µsec)	0.97					
	PRF	(Hz)	3333.3					
Other	pr@PII _{max}	(MPa)	13.83					
Information	d _{eq} @PII _{max}	(cm)						
mornation	Focal Length	FL _x (cm)		0.33				
	T Ocal Length	FLy (cm)		0.24				
	I _{PA.3} @ MI _{max}	(W/cm ²)	98.54					
	Mode	NA	В	В				
Operating Control Conditions	Focus	(cm)	4	1				
	Depth	(cm)	25.87	29.57				
	Scanning width	%	100	100				
	Freq	MHz	2.5	2.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>B+C</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	num Index Value		0.9	0.53				
	P _{r.3}	(MPa)	1.35					
	Wo	(mW)		92.16				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
	Z1	(cm)						
Associated	Z _{bp}	(cm)						
Acoustic	Z _{sp}	(cm)						
Parameter	z@PII. _{3max}	(cm)	33.67					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	2.25	2.32				
	Dim of A _{aprt}	X (cm)		1.92				
	Dim Of Aaprt	Y (cm)		1.58				
	PD	(µsec)	1.69					
	PRF	(Hz)	5952.4					
Other	pr@PII _{max}	(MPa)	1.75					
Information	d _{eq} @PII _{max}	(cm)						
information	Focal Length	FL _x (cm)		0.27				
		FLy (cm)		0.32				
	I _{PA.3} @ MI _{max}	(W/cm ²)	57.39					
	Mode	NA	С	С				
Oneration	Focus	(cm)	4	1				
Operating Control Conditions	Depth	(cm)	20.94	7.39				
	Scanning width	%	100	100				
	Freq	MHz	2.5/2.5	2.5/2.5				
	Power	(%)	100	100				

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>PW</u>

					TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-sca	ті С
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxim	num Index Value		0.89			0.77	0.37	
	P _{r.3}	(MPa)	1.25					
	Wo	(mW)					122.88	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				80.74		
Accesieted	Z1	(cm)				3.05		
Associated	Z _{bp}	(cm)				2.94		
Acoustic	Z _{sp}	(cm)					5.35	
Parameter	z@PII. _{3max}	(cm)	59.16					
	d _{eq} (z _{sp})	(cm)					0.4	
	fc	(MHz)	1.99			1.99	1.99	
	Dim of A _{aprt}	X (cm)				1.92	1.92	
		Y (cm)				1.58	1.58	
	PD	(µsec)	1.83					
	PRF	(Hz)	5050.5					
Other	pr@PII _{max}	(MPa)	111.81					
Information	d _{eq} @PII _{max}	(cm)					0.4	
mornation	Focal Length	FL _x (cm)				0.32		
	Focal Length	FLy (cm)				0.39		
	I _{PA.3} @ MI _{max}	(W/cm ²)	150.4					
	Mode	NA	PW			PW	PW	
Onersting	Focus	(cm)	1			1	1	
Operating Control Conditions	Depth	(cm)	7.39			7.39	7.39	
	Scanning width	%	100			100	100	
	Freq	MHz	2.5/2.5			2.5/2.5	2.5/2.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>B+M</u>

					TIS		TIB	
	Index Label		МІ	Seen	non-	scan	non-sca	TI C
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	C
Global Maxin	1um Index Value		1.2			0.75	0.62	
	P _{r.3}	(MPa)	1.64					
	Wo	(mW)					89.2	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				30.72		
	Z1	(cm)				2.85		
Associated	Z _{bp}	(cm)				2.94		
Acoustic	Z _{sp}	(cm)					2.92	
Parameter	z@PII. _{3max}	(cm)	25.5					
	d _{eq} (z _{sp})	(cm)					0.44	
	fc	(MHz)	1.89			1.89	1.88	
	Dim of A _{aprt}	X (cm)				1.92	1.92	
	Diff Of Aaprt	Y (cm)				1.58	1.58	
	PD	(µsec)	0.45					
	PRF	(Hz)	5532.6					
Other	pr@PII _{max}	(MPa)	1.59					
Information	d _{eq} @PII _{max}	(cm)					0.33	
mornation	Focal Length	FL _x (cm)				0.28		
		FLy (cm)				0.34		
	I _{PA.3} @ MI _{max}	(W/cm ²)	129.2					
	Mode	NA	М			М	М	
Onenations	Focus	(cm)	5			5	6	
Operating Control Conditions	Depth	(cm)	8.62			8.62	7.39	
	Scanning width	%	100			100	100	
	Freq	MHz	2.5			2.5	2.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>CW</u>

	<u></u>				TIS		TIB	
	Index Label		МІ	Scan	non-	non-scan non-sca	TIC	
				Scall	A _{aprt} ≤1	A _{aprt} >1	n	
Global Maxin	num Index Value		0.12			0.8	1.98	0.9 5
	P _{r.3}	(MPa)	0.2					
	Wo	(mW)					178	
	min of [W.3(z1),	(mW)				97		
	I _{TA.3} (Z ₁)]	(11100)				07		
	Z1	(cm)				5.4		
Associated	Z _{bp}	(cm)				1.86		
Acoustic	Z _{sp}	(cm)					6.75	
Parameter	z@PII.3max	(cm)	6.8					
	d _{eq} (z _{sp})	(cm)					0.8	
	fc	(MHz)	2.8			2.75	2.8	2.7 6
		X (cm)				1.92	1.92	
	Dim of A _{aprt}	Y (cm)				1.4	1.4	
	PD	(µsec)	20					
	PRF	(Hz)	0					
	pr@PII _{max}	(MPa)	0.12					
Other	d _{eq} @PII _{max}	(cm)					0.79	
Information	Focal Length	FL _x (cm)				0.46		0.4 5
		FLy (cm)				0.35		0.3 3
	I _{PA.3} @ MI _{max}	(W/cm ²)	43.21					
	Mode	NA	CW			CW	CW	CW
	Focus	(cm)	5			4	6	6
Operating Control	Depth	(cm)	6.16			8.62	7.39	7.3 9
Conditions	Scanning width	%	100			100	100	100
	Freq	MHz	2.5			2.5	2.5	2.5
	Power	(%)	100			100	100	100

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display

Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D3P64L</u> Operating Mode: <u>CFM-M</u>

	<u></u>				TIS		TIB	
	Index Label		мі	Seen	non-	scan	non-sca	тіс
				Scan	A _{aprt} ≤1	A _{aprt} >1	n	
Global Maxin	num Index Value		1.02			0.12	0.56	
	P _{r.3}	(MPa)	1.57					
	Wo	(mW)					254	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				12.49		
	Z1	(cm)				6.15		
Associated	Z _{bp}	(cm)				2.75		
Acoustic	Z _{sp}	(cm)					6.14	
Parameter	z@PII.3max	(cm)	6.15					
	d _{eq} (z _{sp})	(cm)					3.08	
	f _c	(MHz)	2.36			2.38	2.36	
	Dim of A	X (cm)				1.92	1.92	
	Dim of A _{aprt}	Y (cm)				1.4	1.4	
	PD	(µsec)	0.71					
	PRF	(Hz)	252					
	pr@PII _{max}	(MPa)	2.63					
Other	d _{eq} @PII _{max}	(cm)					3.01	
Information		FL _x (cm)				0.52		
	Focal Length	FLy (cm)				0.5		
	I _{PA.3} @ MI _{max}	(W/cm ²)	79.5					
	Mode	NA	CFM-M			CFM- M	CFM-M	
Operating	Focus	(cm)	5			3	6	
Control	Depth	(cm)	6.16			8.62	8.62	
Conditions	Scanning width	%	100			100	100	
	Freq	MHz	2.5			2.5	2.5	
	Power	(%)	100			100	100	

Notes:

((a) This index is not required for this operating mode. See section 4.1.3.1. of the Output Display Standard (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic uses.

(c) This formulation for TIS is less than that for an alternate formulation in this mode.

Transducer Model: <u>D12L40L</u> Operating Mode: <u>B</u>

	<u>,</u>				TIS		TIB	
	Index Label		м	Scan	non-	scan	non-scan	TIC
				Scall	A _{aprt} ≤1	A _{aprt} >1	non-scan	
Global N	laximum Inde	ex Value	0.6	0.1				
	p _{r.3}	(MPa)	1.39					
	Wo	(mW)		8.19				
	min of [W.3(z1), ITA.3(Z1)]	(mW)						
Associated	Z1	(cm)						
Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
	z@PII.3max	(cm)	1.8					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	5.91	5.91				
	Dim of	X (cm)		2.05				
	A _{aprt}	Y (cm)		0.45				
	PD	(µsec)	0.23					
	PRF	(Hz)	5319					
	pr@PII _{max}	(MPa)	2.03					
Other	d _{eq} @PII _{max}	(cm)						
Information	Focal	FL _x (cm)		0.33				
	Length	FLy (cm)		0.24				
	I _{PA.3} @ MI _{max}	(W/cm²)	84.56					
Operating	Mode	NA	В	В				
Control	Focus	(cm)	2.5	2.5				
Conditions	Power	(%)	100	100				

- (b) This probe is not intended for transcranial or neonatal cephalic uses.
- (c) This formulation for TIS is less than that for an alternate formulation in this mode.
- # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: <u>D12L40L</u> Operating Mode: <u>THI-B</u>

	<u>,</u>				TIS		TIB	
	Index Label		МІ	Scan	non-	scan	non-scan	TIC
				Scall	A _{aprt} ≤1	A _{aprt} >1	non-scan	
Global N	laximum Inde	ex Value	0.6	0.1				
	p _{r.3}	(MPa)	1.39					
	Wo	(mW)		8.19				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Acoustic	Z _{bp}	(cm)						
Parameter	Z _{sp}	(cm)						
	z@PII.3max	(cm)	1.8					
	d _{eq} (z _{sp})	(cm)						
	f _c	(MHz)	5.89	5.90				
	Dim of	X (cm)		2.05				
	A _{aprt}	Y (cm)		0.45				
	PD	(µsec)	0.23					
	PRF	(Hz)	5319					
	pr@PII _{max}	(MPa)	2.03					
Other	d _{eq} @PII _{max}	(cm)						
Information	Focal	FL _x (cm)		0.33				
	Length	FLy (cm)		0.22				
	I _{PA.3} @ MI _{max}	(W/cm²)	84.49					
Operating	Mode	NA	THI-B	THI-B				
Control	Focus	(cm)	2.5	2.5				
Conditions	Power	(%)	100	100				

- (b) This probe is not intended for transcranial or neonatal cephalic uses.
- (c) This formulation for TIS is less than that for an alternate formulation in this mode.
- # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: <u>D12L40L</u> Operating Mode: <u>C</u>

	<u> </u>				TIS		TIB	
	Index Label		мі	Scan	non	scan	non-scan	TIC
				Scall	A _{aprt} ≤1	A _{aprt} >1	non-scan	
Global N	laximum Inde	ex Value	0.5	0.4				
	p _{r.3}	(MPa)	1.18					
	Wo	(mW)		57.34				
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)						
Associated	Z1	(cm)						
Acoustic	Zbp	(cm)						
Parameter	Z _{sp}	(cm)						
	z@PII.3max	(cm)	1.8					
	d _{eq} (z _{sp})	(cm)						
	fc	(MHz)	5.12	5.40				
	Dim of	X (cm)		2.05				
	A _{aprt}	Y (cm)		0.45				
	PD	(µsec)	0.86					
	PRF	(Hz)	7887					
	pr@PII _{max}	(MPa)	1.51					
Other	d _{eq} @PII _{max}	(cm)						
Information	Focal	FL _x (cm)		0.27				
	Length	FLy (cm)		0.32				
	I _{PA.3} @ MI _{max}	(W/cm²)	72.65					
Operating	Mode	NA	С	С				
Control	Focus	(cm)	2.5	2.5				
Conditions	Power	(%)	100	100				

- (b) This probe is not intended for transcranial or neonatal cephalic uses.
- (c) This formulation for TIS is less than that for an alternate formulation in this mode.
- # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: D12L40L Operating Mode: PW

	<u> </u>				TIS		TIB	
	Index Label		МІ	Scan	non	-scan	non-scan	TIC
				Scan	A _{aprt} ≤1	A _{aprt} >1	non-scan	
Global N	laximum Inde	ex Value	0.2			1.4	1.9	
	p _{r.3}	(MPa)	0.54					
	Wo	(mW)					106.50	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				54.21		
Associated	Z1	(cm)				1.8		
Acoustic	Z _{bp}	(cm)				1.8		
Parameter	Z _{sp}	(cm)					1.8	
	z@PII.3max	(cm)	3.0					
	d _{eq} (z _{sp})	(cm)					0.86	
	f _c	(MHz)	5.43			5.43	5.43	
	Dim of	X (cm)				4.10	4.10	
	A _{aprt}	Y (cm)				0.45	0.45	
	PD	(µsec)	0.83					
	PRF	(Hz)	8051					
	pr@PII _{max}	(MPa)	0.75					
Other	d _{eq} @PII _{max}	(cm)					0.86	
Information	Focal	FL _x (cm)				0,26		
	Length	FLy (cm)				0.21		
	I _{PA.3} @ MI _{max}	(W/cm²)	103.77					
Operating	Mode	NA	PW			PW	PW	
Control	Focus	(cm)	2.5			2.5	2.5	
Conditions	Power	(%)	100			100	100	

- (b) This probe is not intended for transcranial or neonatal cephalic uses.
- (c) This formulation for TIS is less than that for an alternate formulation in this mode.
- # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: <u>D12L40L</u> Operating Mode: <u>M</u>

	<u>,</u>	_			TIS		TIB	
	Index Label		мі	Scan	non-	-scan	non-scan	TIC
				Scall	A _{aprt} ≤1	A _{aprt} >1	non-scan	
Global N	laximum Inde	ex Value	0.5			0.6	1.0	
	p _{r.3}	(MPa)	1.21					
	Wo	(mW)					40.98	
	min of [W _{.3} (z ₁), I _{TA.3} (z ₁)]	(mW)				19.75		
Associated	Z1	(cm)				1.8		
Acoustic	Z _{bp}	(cm)				1.8		
Parameter	Z _{sp}	(cm)					1.8	
	z@PII.3max	(cm)	1.8					
	d _{eq} (z _{sp})	(cm)					0.67	
	fc	(MHz)	5.85			5.87	5.82	
	Dim of	X (cm)				4.10	4.10	
	A _{aprt}	Y (cm)				0.45	0.45	
	PD	(µsec)	0.23					
	PRF	(Hz)	5319					
	pr@PII _{max}	(MPa)	1.74					
Other	d _{eq} @PII _{max}	(cm)					0.67	
Information	Focal	FL _x (cm)				0.25		
	Length	FLy (cm)				0.18		
	I _{PA.3} @ MI _{max}	(W/cm²)	170.00					
Operating	Mode	NA	М			М	М	
Control	Focus	(cm)	2.5			2.5	2	
Conditions	Power	(%)	100			100	100	

- (b) This probe is not intended for transcranial or neonatal cephalic uses.
- (c) This formulation for TIS is less than that for an alternate formulation in this mode.
- # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Operating Mode: B

operating			МІ	T	IS	TI	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	dex value		0.60	0.	10	0.	10	-
Index compo	nent value			0.10	0.10	0.08	0.10	
	$p_{r,\alpha}$ at Z_{MI}	Mpa	1.39					
	Р	mW		8.	19	8.	19	_
	P _{1x1}	mW		2.	00	2.	00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	1.68					
	Z _{pii,α}	cm	-					
	f _{awf}	MHz	5.91	5.	91	5.9	5.91	
	prr	Hz	-					
	srr	Hz	41.55					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	82.21					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	109.55					
	p _r at z _{pii}	Мра	-					
Operating	Focus	cm	2.50	2.	50	2.	50	_
control conditions	Power	%	100	1(00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB. NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Operating Mode: THI-B

			МІ	Т	IS	TI	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	dex value		0.60	0.	10	0.	10	-
Index compo	nent value			0.10	0.10	0.07	0.10	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.39					
	Р	mW		8.	19	8.	19	_
	P _{1x1}	mW		2.	00	2.	00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					-	
	Z _{MI}	cm	1.69					
	Z _{pii,α}	cm	-					
	f _{awf}	MHz	5.89	5.90		5.90		_
	prr	Hz	_					
	srr	Hz	41.55					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	82.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	110.32					
	p _r at z _{pii}	Мра	_					
Operating	Focus	cm	2.50	2.	50	2.	50	_
control conditions	Power	%	100	1(00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB. NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Operating Mode: C

opolating			МІ	T	IS	TI	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	dex value		0.50		40	0.4	40	-
Index compo	nent value			0.40	0.40	0.33	.33 0.40	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.18					
	Р	mW		57	.34	57.	.34	_
	P _{1x1}	mW		14	.00	14.	.00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					-	
	Z _{MI}	cm	1.55					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	5.12	5.40		5.40		_
	prr	Hz	_					
	srr	Hz	61.62					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha} \text{ at } z_{pii,\alpha} \text{ or } z_{sii,\alpha}$	mW/cm ²	70.45					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	94.68					
	p _r at z _{pii}	Мра	_					
Operating	Focus	cm	2.50	2.	50	2.	50	_
control conditions	Power	%	100	10	00	10	00	-

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB. NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Operating Mo	ode: PW
---------------------	---------

			МІ	T	IS	TI	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	dex value			90	-			
Index compo	nent value			1.23	1.40	1.78	1.90	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.54					
	Р	mW		106	6.50	106	6.50	-
	P _{1x1}	mW		26	.00	26.	.00	
Acoustic	Zs	cm			1.80			
Parameters	Z _b	cm					1.80	
	Z _{MI}	cm	2.64					
	Z _{pii,α}	cm	3.00					
	f _{awf}	MHz 5.43 5.43 5.43		43	-			
	prr	Hz	8051.00					
	srr	Hz	_					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	103.77					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	100.56					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	126.54					
	p _r at z _{pii}	Мра	0.75					
Operating	Focus	cm	2.50	2.	50	2.	50	-
control conditions	Power	%	100	10	00	10	00	-

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB. NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Operating Mode: M

			МІ	Т	IS	TI	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	dex value		0.50	0.	60	1.00		-
Index compo	nent value			0.55	0.60	0.88	1.00	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.21					
	Р	mW		40	.98	40.	.98	_
	P _{1x1}	mW		10	.00	10.	.00	
Acoustic	Zs	cm			1.80			
Parameters	Z _b	cm					1.80	
	Z _{MI}	cm	1.65					
	Z _{pii,α}	cm	1.80					
	f _{awf}	MHz	5.85	5.87		5.82		_
	prr	Hz	5319.00					
	srr	Hz	41.55					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	170.00					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	166.32					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	194.51					
	p _r at z _{pii}	Mpa	1.74					
Operating	ng Focus cm 2.50 2.50 2.00		00	_				
control conditions	Power	%	100	10	00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB. NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Transducer Model: <u>D3C60L</u> Operating Mode: <u>B</u>

	<u>,</u>		МІ	Т	IS	ТІ	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ir	ndex value		0.80	0.9	97	0.9	97	-
Index comp	onent value			0.97	0.97	0.88	0.97	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.32					
	Р	mW		511	.48	511	.48	_
	P _{1x1}	mW		74.	.00	74.	.00	
Acoustic	Zs	cm			-			
Paramete rs	Z _b	cm					-	
	Z _{MI}	cm	6.33					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	2.75	2.	75	2.	75	_
	prr	Hz						
	srr	Hz	75.60					
	n _{pps}		1					
Other	I _{pa,α} at z _{pii,α}	W/cm ²	_					
Informatio n	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $Z_{sii,\alpha}$	mW/cm ²	72.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	98.54					
	p _r at z _{pii}	Мра	-					
	Focus	cm	7.00	7.	00	7.00		-
Operating Depth		cm	20.45	20.	.45	20	.45	_
control	Scanning width	%	100	1(00	1(00	_
	Freq	MHz	3.0	3.	.0	3.	.0	_
	Power	%	100	10	00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Transducer Model: <u>D3C60L</u> Operating Mode: B+C

			МІ	Т	IS	Т	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ir			0.56	0.		0.		-
Index comp	onent value			0.68	0.68	0.66	0.68	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.95					
	Р	mW		345	5.60	345	5.60	-
A	P _{1x1}	mW		70.	.36	70	.36	
Acoustic Paramete	Zs	cm		-				
raramete	Zb	cm					—	
	Z _{MI}	cm	6.11					
	Z _{pii,α}	cm	-					
	f _{awf}		2.90	2.9	92	2.9	92	-
	prr	Hz	_					
	srr	Hz	46.97					
Other	N _{pps}		1					
Other Informatio	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
n	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	33.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	56.51					
	p _r at z _{pii}	Мра	_					
	Focus	cm	6.00	8.	00	8.	00	_
Operation	Depth	cm	11.82	9.3	36	9.3	36	_
Operating control conditions	Scanning width	%	100	1(00	1(00	_
	Freq	MHz	3.0/3.0	3.0/	/3.0	3.0/	/3.0	-
	Power	%	100	10	00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Transducer<u>D3C60L</u> Operating Mode: <u>PW</u>

	$ \begin{array}{c c} I & Z_{sii,\alpha} \\ \hline I_{spta} \text{ at } Z_{pii} \text{ or} \\ Z_{sii} \\ \hline p_r \text{ at } Z_{pii} \\ \hline \hline Focus \\ \hline Depth \\ \end{array} $		МІ	Т	IS	TI	в	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
			0.50					-
Index comp	onent value	í		1.55	1.77	0.99	1.02	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.84					
	Р	mW		588	8.55	588	.55	_
A	P _{1x1}	mW		85	.15	85.	15	
	Zs	cm			4.10			
	Z _b	cm					6.45	
	Z _{MI}	cm	4.32					
	Z _{pii,α}	cm	4.40					
	f _{awf}	MHz	2.88	2.	87	2.8	38	-
	srr	Hz	-					
	n _{pps}		1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	189.78					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or	mW/cm ²	180 22					
n	Z _{sii,α}							
	I _{spta} at z _{pii} or	mW/cm ²	206.25					
	p _r at z _{pii}	Мра	1.45					
		cm						_
Operating	Depth	cm	9.36	1.55 1.77 0.99 1.02 0.84 58 58 58 58 58 58 58 58 58 85 85 85 85 15 4.10 4.10 6.45 4.32 4.32 4.10 $1-10$ 6.45 4.32 2.88 2.87 2.88 2.88 2.87 2.88 $1-10$ -1 $1-10$ $1-10$ $1-10$ 189.78 2.12 2.12 1.12 180.22 2.12 2.12 2.12 $2.66.25$ 2.12 2.12 2.12 1.45 2.12 2.12 2.12 1.45 2.12 $1.1.82$ 100 1.00 1.00	_			
control	Scanning width	%	100	1(00	10	00	-
	Freq	MHz	3.0/3.0)/3.0 3.0/3.0		3.0/3.0		-
	Power	%	100	1(00	10	00	_

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Transducer Model: <u>D3C60L</u> Operating Mode: <u>B+M</u>

-	<u>g modol <u>D - m</u></u>		МІ	ТІ	S	т	IB	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum in	dex value		0.95	1.(02	1.	10	-
Index compo	onent value			0.98	1.02	1.04	1.10	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.57					
	Р	mW		511	.23	511.23		-
	P _{1x1}	mW		73.	96	73	.96	
Acoustic Parameter	Zs	cm			4.20			
S	Z _b	cm					5.10	
	Z _{MI}	cm	5.48					
	Z _{pii,α}	cm	5.60					
	f _{awf}		2.74	2.78		2.	78	_
	prr	Hz	4759.62					
	srr	Hz	37.18					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	112.05					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	113.28					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	133.22					
	p _r at z _{pii}	Мра	2.64					
	Focus	cm	7.00	6.0	00	6.	00	-
Operating	Depth	cm	8.13	6.89		6.	89	_
control	Scanning width	%	100	10	00	1	00	_
	Freq	MHz	3.0	3.	0	3	.0	_
	Power	%	100	10	00	1	00	-

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D7L40L

Operating			МІ	Т	IS	Т	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.97	0.:	53	0.	53	-
Index compo	nent value	-		0.53	0.53	0.44	0.53	
	$p_{r,\alpha}$ at Z_{MI}	Мра	2.41					
	Р	mW		73	.73	73	.73	—
	P _{1x1}	mW		18	.00	18	.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	1.88					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	6.23	6.:	25	6.	25	_
	prr	Hz	_					
	srr	Hz	22.45					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	40.22					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	61.54					
	p _r at z _{pii}	Mpa	-					
	Focus	cm	2.00	2.	00	2.	00	
Operating	Depth	cm	4.93		93		93	_
control	Scanning width	%	100	100			00	_
conditions	conditions Freq MHz		5.0	5.0		5.0		_
Power		%	100	10	00	100		—

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

			МІ	Т	IS	T	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.55	0.	79	0.	79	-
Index compo	nent value	-		0.79	0.79	0.66	0.79	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.38					
	Р	mW		106	6.50	106	5.50	-
	P _{1x1}	mW		26	.00	26	.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	1.54					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	6.36	6.37		6.37		_
	prr	Hz	_					
	srr	Hz	54.25					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	88.59					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	112.25					
	p _r at z _{pii}	Мра	-					
	Focus	cm	2.50		00		00	_
Operating	Depth	cm	4.90		70	3.		_
control	Scanning width	%	100	100		100		_
conditions	Freq	MHz	5.0/6.5	6.5 5.0/6.5 5.0/6.5			_	
	Power	%	100	1(00	10	100	

Operating Mode: <u>B +C</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

Operating			МІ	Т	IS	Т	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.63	0.	56	0.	88	-
Index compo	nent value			0.45 0.56		0.77	0.88	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.59					
	Р	mW		221	.18	221	.18	_
	P _{1x1}	mW		54	.00	54	.00	
Acoustic	Zs	cm			1.47			
Parameters	Z _b	cm					1.85	
	Z _{MI}	cm	1.64					
	Z _{pii,α}	cm	1.70					
	f awf	MHz	6.34	6.23		6.34		-
	prr	Hz	2873.60					
	srr	Hz	_					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	118.84					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	115.22					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	138.56					
	p _r at z _{pii}	Мра	2.14					
	Focus	cm	4.00	2.	50	2.	50	-
Operating	Depth	cm	7.39	3.	70	7.	39	_
control	Scanning width	%	100	100 100			_	
conditions	Freq	MHz %	5.0/6.5	5.0/6.5		5.0/6.5		_
	Power		100	10	00	10	00	-

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

			МІ	T	IS	T	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	lex value		0.77	1.	10	1.	00	-
Index compo	nent value	-		0.98 1.10		0.88	1.00	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.91					
	Р	mW		139	9.30	139.30		-
	P _{1x1}	mW		34.	.00	34	.00	
Acoustic	Zs	cm			1.68			
Parameters	Z _b	cm					1.80	
	Z _{MI}	cm	1.55					
	Zpii,α	cm	1.65					
	f _{awf}	MHz	6.19	6.22		6.19		_
	prr	Hz	2873.60					
	srr	Hz	22.45					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	158.13					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	152.25					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	178.55					
	p _r at z _{pii}	Мра	2.60					
	Focus	cm	2.00	1.	50	2.	00	-
Operating	Depth	cm	3.70	3.	70	3.	70	_
control	Scanning width	%	100	100 100			-	
conditions Freq MHz		MHz	5.0/6.5	5.0/6.5		5.0/6.5		_
Power		%	100	10	00	10	00	-

Operating Mode: <u>B +M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

			МІ	Т	IS	Т	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.84	0.	52	0.	52	-
Index compo	nent value	-		0.52	0.52	0.44	0.52	
	$p_{r,\alpha}$ at Z_{MI}	Мра	2.08					
	Р	mW		73	.28	73	.28	_
	P _{1x1}	mW		17	.89	17	.89	
Acoustic	Zs	cm			_			
Parameters	Zb	cm					_	
	Z _{MI}	cm	1.88					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	6.14	6.14		6.14		-
	prr	Hz	-					
	srr	Hz	22.39					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	74.44					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	99.56					
	p _r at z _{pii}	Мра	_					
	Focus	cm	3.00		00		00	_
Operating	Depth	cm	4.93	4.	93	4.	93	_
control	Scanning width	%	100	1(00	1(00	-
conditions	Freq	MHz	5.0	5.0		5.0		_
	Power	%	100	1(00	1(00	-

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D7L40L-REC

			МІ	Т	IS	Т	B	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.49	0.	78	0.	78	-
Index compo	nent value			0.78	0.78	0.65	0.78	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.25					
	Р	mW		105	5.60	105	5.60	_
	P _{1x1}	mW		25	.78	25	.78	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	1.55					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	6.36	6.35		6.35		_
	prr	Hz	_				35	
	srr	Hz	53.56					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	88.59					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	112.25					
	p _r at z _{pii}	Мра	_					
	Focus	cm	2.00		50	2.50		_
Operating	Depth	cm	4.93		70		70	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	5.0/6.5	/6.5 5.0/6.5 5.0/6.5		/6.5	_	
	Power	%	100	10	00	10	00	-

Operating Mode: B +C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D7L40L-REC

			МІ	Т	S	Т	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.56	1.	02	0.	80	-
Index compo	nent value			0.88	1.02	0.77	0.80	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.39					
	Р	mW		73	.73	73	.73	_
	P _{1x1}	mW		18	.00	18	.00	
Acoustic	Zs	cm			1.47			
Parameters	Z _b	cm					1.65	
	Z _{MI}	cm	1.56					
	Zpii,α	cm	1.65					
	f _{awf}	MHz	6.32	6.32		6.30		-
	prr	Hz	5877.00				30	
	srr	Hz	_					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	120.41					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	115.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	138.56					
	p _r at z _{pii}	Мра	1.90					
	Focus	cm	2.00	2.	2.00 3.00		00	_
Operating	Depth	cm	7.39		39		39	_
control	Scanning width	%	100	100 100		00	_	
conditions	Freq	MHz	5.0/6.5	5 5.0/6.5 5.0/6		/6.5	_	
	Power	%	100	10	00	10	00	-

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D7L40L-REC

			МІ	Т	S	TI	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	dex value		0.69	1.	10	0.4	49	-
Index compo	nent value	1		0.98	1.10	0.38	0.49	
	$p_{r,\alpha}$ at Z_{MI}	Mpa	1.71					
	Р	mW		136	6.50	139	.30	_
	P _{1x1}	mW		33	.33	34	.00	
Acoustic	Zs	cm			1.68			
Parameters	Zb	cm					1.86	
	Z _{MI}	cm	2.33					
	Z _{pii,α}	cm	2.45					
	f _{awf}	MHz	6.05	6.18		6.05		-
	prr	Hz	2655.30)5	
	srr	Hz	20.74					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	134.51					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	130.25					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	156.22					
	p _r at z _{pii}	Мра	2.52					
	Focus	cm	2.00	2.00 2.00		00	_	
Operating	Depth	cm	7.39		70		39	_
control	Scanning width	%	100	100 100			_	
conditions	Freq	MHz	5.0/6.5	5 5.0/6.5 5.0/6.5				
	Power	%	100	10	00	10	00	_

Operating Mode: **B** +M

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D5C20L

			МІ	т	IS	Т	B	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.89	1.	06	1.	06	-
Index compo	nent value			1.06	1.06	0.88	1.06	
	$p_{r,\alpha}$ at Z_{MI}	Mpa	1.96					
	Р	mW		350).21	350).21	-
	P _{1x1}	mW		120	0.00	120	0.00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	2.33					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	4.84	5.00		5.00		_
	prr	Hz	-				00	
	srr	Hz	114.89					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	77.54					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	102.22					
	p _r at z _{pii}	Мра	-					
	Focus	cm	4.00	1.	50	1.	50	_
Operating	Depth	cm	11.77	5.	61	5.	61	_
control	Scanning width	%	100	100 100		00	_	
conditions Freq MHz		MHz	3.5	3.5		3.5		_
	Power		100	1(00	10	00	_

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D5C20L

Index label			МІ	TIS		TIB		TIC
				At	Below	At	Below	
				surface	surface	surface	surface	
Maximum index value			0.70	1.	1.10		1.10	
Index component value			1.10	1.10	0.99	1.10		
Acoustic Parameters	$p_{r,\alpha}$ at Z_{MI}	Мра	1.39					
	Р	mW		321.02		321.02		-
	P _{1x1}	mW		110.00		110.00		
	Zs	cm			-			
	Z _b	cm					_	
	Z _{MI}	cm	2.55					
	Zpii,α	cm	_					
	f _{awf}	MHz	4.01	4.01		4.01		_
Other Information	prr	Hz	_					
	srr	Hz	46.50					
	N _{pps}		1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	55.22					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	79.45					
	p _r at z _{pii}	Мра	-					
Operating control conditions	Focus	cm	35	3.50		3.50		_
	Depth	cm	9.30	9.30		9.30		_
	Scanning width	%	100	100		100		_
	Freq	MHz	3.5/4.0	3.5/4.0		3.5/4.0		_
	Power	%	100	100		100		-

Operating Mode: B+C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Transducer Model: D5C20L

Operating			МІ	Т	IS	т	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.80	1.	13	0.4	45	-
Index compo	nent value			1.04	1.13	0.33	0.45	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.60					
	Р	mW		349	9.43	349	.43	-
	P _{1x1}	mW		119	9.73	119	0.73	
Acoustic	Zs	cm			2.05			
Parameters	Z _b	cm					2.20	
	Z _{MI}	cm	2.22					
	Z _{pii,α}	cm	2.35			2.06		
	f awf	MHz	3.93	3.	96	3.96		-
	prr	Hz	6024.10					
	srr	Hz	_					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	101.33					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	99.98					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	124.10					
	p _r at z _{pii}	Мра	2.10					
	Focus	cm	2.00	3.	00	3.	00	-
Operating	Depth	cm	4.37	4.	37	4.37		_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	3.5/4.0	3.5	/4.0	3.5	/4.0	_
	Power	%	100	10	00	10	00	_

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	dex value		0.93	1.	15	1.0	01	-
Index compo	nent value			1.08	1.15	0.88	1.01	
	$p_{r,\alpha}$ at Z_{MI}	Мра	2.02					
	Р	mW		356	õ.04	356	6.04	-
	P _{1x1}	mW		122	2.00	122	2.00	
Acoustic	Zs	cm			2.05			
Parameters	Zb	cm				2.05		
	Z _{MI}	cm	1.89			2.05 2.05 4.83		
	Zpii,α	cm	2.05				4.83	
	f _{awf}	MHz	4.75	4.	92	4.8	83	_
	prr	Hz	14706.00					
	srr	Hz	114.89					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	111.39					
Information	I _{spta,α} at z _{pii,α} or z _{sii,α}	mW/cm ²	106.25					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	128.55					
	p _r at z _{pii}	Мра	2.83					
	Focus	cm	3.50	2.	50	2.0	00	_
Operating	Depth	cm	4.37	3.	14	5.0	61	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	3.5		.5		.5	_
	Power	%	100	10	00	10	00	_

Operating Mode: <u>B+M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Operating			МІ	Т	IS	Т	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.72	1.:	27	1.:	27	-
Index compo	nent value			1.27	1.27	1.22	1.27	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.30					
	Р	mW		365	5.57	365	5.57	_
	P _{1x1}	mW		102	2.00	102	2.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					-	
	Z _{MI}	cm	4.44					
	Z _{pii,α}	cm	_				3.44	
	f _{awf}	MHz	3.27	3.4	44	3.4	44	_
	prr	Hz	_					
	srr	Hz	30.05					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	49.51					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	72.21					
	p _r at z _{pii}	Мра	-					
	Focus	cm	5.00		00		00	_
Operating	Depth	cm	8.84		.70		.70	_
control	Scanning width	%	100		00		00	_
conditions	Freq	MHz	3.5		.5		.5	_
	Power	%	100	1(00	10	00	_

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	T	IB	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.30	0.	50	0.	50	-
Index compo	nent value	-		0.50	0.50	0.41	0.50	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.52					
	Р	mW		136	6.19	136	136.19	
	P _{1x1}	mW		38	.00	38	.00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm				-		
	Z _{MI}	cm	3.88					
	Zpii,α	cm	-					
	f _{awf}	MHz	2.84	2.	84	2.84		_
	prr	Hz	-					
	srr	Hz	46.50					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	95.65					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	121.54					
	p _r at z _{pii}	Мра	-					
	Focus	cm	5.00		00		00	_
Operating	Depth	cm	7.61		61		61	_
control	Scanning width	%	100		00	10	00	_
conditions	Freq	MHz	3.5/3.5	3.5	/3.5	3.5/	/3.5	_
	Power	%	100	10	00	10	00	-

Operating Mode: B+C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	lex value		0.32	0.	50	0.8	80	-
Index compo	nent value	-		0.46	0.50	0.74	0.80	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.54					
	Р	mW		78.	.80	78	.80	-
	P _{1x1}	mW		21.	.99	21	.99	
Acoustic	Zs	cm			3.90			
Parameters	Z _b	cm					3.50	
	Z _{MI}	cm	5.45					
	Zpii,α	cm	5.60			2.89		
	f _{awf}	MHz	2.90	2.9	90	2.	89	_
	prr	Hz	3846.20					
	srr	Hz	_					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	110.43					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	105.52					
	l _{spta} at z _{pii} or z _{sii}	mW/cm ²	132.25					
	p _r at z _{pii}	Мра	0.75					
	Focus		5.00	5.	00	5.	00	_
Operating	Depth			-				
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	3.5/3.5		/3.5	3.5/3.5		_
	Power	%	100	10	00	10	00	_

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	TI	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	lex value		0.70	1.4	40	0.:	55	-
Index compo	nent value			1.23	1.40	0.41	0.50	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.27					
	Р	mW		487	' .42	487	'.42	-
	P _{1x1}	mW		136	6.00	136	6.00	
Acoustic	Zs	cm			3.55			
Parameters	Z _b	cm					4.30	
	Z _{MI}	cm	4.44					
	Z _{pii,α}	cm	4.65					
	f awf	MHz	3.26	3.:	26	3.:	3.26	
	prr	Hz	7575.80					
	srr	Hz	59.19					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	162.39					
Information	$I_{spta,\alpha} at z_{pii,\alpha} or z_{sii,\alpha}$	mW/cm ²	152.25					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	178.55					
	p _r at z _{pii}	Мра	1.92					
	Focus cm 5.00 5.00 5.00		00	-				
Operating	Depth	cm	6.37	6.	37	6.	37	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	3.5	3	.5	3.	.5	_
	Power	%	100	10	00	10	00	_

Operating Mode: <u>B+M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	T	IS	T	B	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.81	1.	24	1.:	24	I
Index compo	nent value			1.24	1.24	1.15	1.24	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.93					
	Р	mW		409	9.60	409	9.60	_
	P _{1x1}	mW		168	3.42	168	3.42	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	2.55					
	Z _{pii,α}	cm - Image: Constraint of the second						
	f _{awf}	MHz	5.66	5.	88	5.	88	-
	prr	Hz	_					
	srr	Hz	46.28					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	81.25					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	105.56					
	p _r at z _{pii}	Мра	_					
	Focus	cm	2.50	2.	00	2.	00	_
Operating	Depth	cm	5.69	6.	92	6.9	92	_
control	Scanning width	%	100	10	00	10	00	-
conditions	Freq	MHz	4.0	4	.0	4.	.0	_
	Power	%	100	10	00	10	00	_

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.75	1.	11	1.	11	-
Index compo	nent value	-		1.11	1.11	1.05	1.11	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.72					
	Р	mW		296	6.70	296	6.70	_
	P _{1x1}	mW		122	2.00	122	2.00	
Acoustic	Zs	cm			-			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	1.77					
	Zpii,α	cm	-					
	f _{awf}	MHz	5.26	5.	36	5.	36	-
	prr	Hz	-					
	srr	Hz	47.06					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	91.15					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	106.22					
	p _r at z _{pii}	Мра	-					
	Focus	cm	2.50	1.	00	1.	00	_
Operating	Depth	cm	11.85	3.:	22	3.	22	
control	Scanning width	%	100	10	00	1(00	_
conditions	Freq	MHz	4.0/5.3	4.0	/5.3	4.0	/5.3	_
	Power	%	100	10	00	10	00	_

Operating Mode: B+C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	T	S	TI	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.69	1.0	68	0.	66	_
Index compo	nent value			1.55	1.68	0.55	0.66	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.58					
	Р	mW		204	.29	204	.29	_
	P _{1x1}	mW		84.	.00	84	.00	
Acoustic	Zs	cm			1.80			
Parameters	Z _b	cm					1.80	
	Z _{MI}	cm	2.77					
	Zpii,α	cm	2.85					
	f _{awf}	MHz	5.21	5.3	38	5.3	5.38	
	prr	Hz 5924.20						
	srr	Hz	-					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	138.48					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	134.25					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	159.54					
	p _r at z _{pii}	Mpa	2.40					
	Focus	cm	2.50	3.	50	3.	50	_
Operating	Depth	cm	3.22	4.4	46	4.4	46	-
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	4.0/5.3	4.0/	/5.3	4.0	/5.3	_
	Power	%	100	10	00	10	00	_

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	dex value		0.80	1.3	33	0.9	90	-
Index compo	nent value	-		1.12	1.33	0.84	0.90	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.79					
	Р	mW		272	2.38	272	2.38	-
	P _{1x1}	mW		112	2.00	112	2.00	
Acoustic	Zs	cm			2.00			
Parameters	Zb	cm				2.00		
	Z _{MI}	cm	3.33			112.00 2.00 4.99		
	Zpii,α	cm	3.55			4.99		
	f _{awf}	MHz	4.92	5.4	44	4.9	99	_
	prr	Hz	11236.00					
	srr	Hz	87.78					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	180.80					
Information	I _{spta,α} at z _{pii,α} or Z _{sii,α}	mW/cm ²	175.55					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	201.23					
	p _r at z _{pii}	Mpa	2.44					
	Focus	cm	4.00		00	3.		_
Operating	Depth	cm	4.00		46	8.		_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	4	4	.0	4.	.0	_
	Power	%	100	10	00	10	00	-

Operating Mode: <u>B+M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	T	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.48	0.:	20	0.:	20	-
Index compo	nent value			0.20	0.20	0.11	0.20	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.99					
	Р	mW		29	.70	29	.70	_
	P _{1x1}	mW		10	.00	10	.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	2.11					
	Z _{pii,α}	cm	_					
	f _{awf}	MHz	4.15	4.	18	4.18		-
	prr	Hz	_					
	srr	Hz	33.97					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	40.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	61.54					
	p _r at z _{pii}	Мра	_					
	Focus	cm	2.50	1.	50	1.	50	—
Operating	Depth	cm	3.93	2.	70	2.	70	_
control	Scanning width	%	100	1(00	1(00	_
conditions	Freq	MHz	4.0	4	.0	4	.0	-
	Power	%	100	10	00	10	00	—

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	T	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.41	0.	65	0.	65	-
Index compo	nent value			0.65	0.65	0.55	0.65	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.95					
	Р	mW		77	.20	77.	.20	_
	P _{1x1}	mW		26	.00	26	.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	2.33					
	Zpii,α	cm	_					
	f _{awf}	MHz	5.27	5.:	28	5.2	5.28	
	prr	Hz	_					
	srr	Hz	47.06					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	28.89					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	50.22					
	p _r at z _{pii}	Mpa	_					
	Focus	cm	2.50	2.	00	2.	00	_
Operating	Depth	cm	10.09	3.	93	3.	93	_
control	Scanning width	%	100	1(00	10	00	_
conditions	Freq	MHz	4	4	1	4	1	_
	Power	%	100	10	00	10	00	-

Operating Mode: <u>B+C</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ind	dex value		0.47	1.3	34	1.	18	-
Index compo	nent value	-		1.22	1.34	1.05	1.18	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.07					
	Р	mW		100).97	100.97		_
	P _{1x1}	mW		34	.00	34.	.00	
Acoustic	Zs	cm			1.75			
Parameters	Z _b	cm					2.00	
	Z _{MI}	cm	2.55					
	Zpii,α	cm	2.60					
	f _{awf}	MHz	5.30	5.	30	5.2	28	_
	prr	Hz	5555.60					
	srr	Hz	_					
0.1	n _{pps}		1					
Other Information	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	143.33					
mornation	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	140.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	165.55					
	p _r at z _{pii}	Мра	1.52					
	Focus	cm	3.00	2.	50	2.0	00	_
Operating	Depth	cm	8.86	6.	39	2.1	70	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	4.0/5.3	4.0	/5.3	4.0/5.3		_
	Power	s cm 1.75 s cm 1.75 b cm 2.55 MI cm 2.55 Dii, α cm 2.60 wf MHz 5.30 5.30 rrr Hz 5555.60 5.2 rrr Hz 5555.60 1 pps 1 1 1 a, a t zpii, a W/cm² 143.33 1 pta, a t zpii, a or zsii, a mW/cm² 140.22 1 pta at zpii or zsii mW/cm² 165.55 1 occus cm 3.00 2.50 2.00 epth cm 8.86 6.39 2.77 canning width % 100 100 100	00	-				

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for

transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	T	S	T	В	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	lex value		0.48	0.	57	0.:	24	-
Index compo	nent value			0.48	0.57	0.19	0.24	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.98					
	Ρ	mW		47.	.51	47	.51	-
	P _{1x1}	mW		16	.00	16	.00	
Acoustic	Zs	cm			1.75		surface	
Parameters Z _b cm						2.30		
	Z _{MI}	cm	2.22					
	Z _{pii,α}	cm	2.40					
	f _{awf}	MHz	4.13	4.	17	4.13		_
	prr	Hz	5555.60					
	srr	Hz	43.40					
	n _{pps}	•	1					
Other	I _{pa,α} at z _{pii,α}	W/cm ²	140.31					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	133.22					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	159.54					
	p _r at z _{pii}	Mpa	1.30					
	Focus	cm	3.00	2.	00	3.	00	-
Operating	Depth	cm	10.09	5.	16	10	.09	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	4	4.	.0	4.0		_
	Power	%	100	10	00	10	00	_

Operating Mode: <u>B+M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Operating			МІ	T	IS	TI	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		1.08	1.:	21	1.:	21	-
Index compo	nent value	-		1.21	1.21	1.11	1.21	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.63					
	Р	mW		269	9.95	269	9.95	-
	P _{1x1}	mW		140	0.60	140	.60	
Acoustic	tic Z _s cm				-			
Parameters	Parameters Z _b cm						_	
	Z _{MI}	cm	1.33					
	Z _{pii,α}	cm	-					
	f _{awf}	MHz	2.28	2.4	48	2.4	48	_
	prr	Hz	-					
	srr	Hz	26.04					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	94.45					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	118.55					
	p _r at z _{pii}	Мра	-					
	Focus	cm	4.00	1.	00	1.	00	-
Operating	Depth	cm	25.87	29	.57	29	.57	_
control	Scanning width	%	100	1(00	10	00	-
conditions	Freq	MHz	2.5	2	.5	2	.5	-
	Power	%	100	1(00	10	100	

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ind	lex value		0.90	0.	53	0.	53	-
Index compo	nent value	-		0.53	0.53	0.48	0.53	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.35					
	Р	mW		92	.16	92	.16	-
	P _{1x1}	mW		48	.00	48	.00	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	3.22					
	Zpii,α	cm	_					
	f _{awf}	MHz	2.25	2.	32	2.	2.32	-
	prr	Hz	_					
	srr	Hz	46.50					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_				surface 53 0.53 0.53 0.0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	55.54					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	79.84					
	p _r at z _{pii}	Мра	_					
	Focus	cm	4.00	1.	00	1.	00	_
Operating	Depth	cm	20.94	7.	39	7.	39	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	2.5/2.5	2.5	/2.5	2.5/2.5		-
	Power	%	- 2.25 2.32 2.32 2.32 $-$ 2.32 2.32 2.32 $ 2.32$ 2.32 2.32 $ 2.32$ 2.32 2.32 $ 2.32$ 2.32 2.32 46.50 $ 2.55.54$ $ 4.00$ 1.00 1.00 1.00 1.00 $ 100$ 100 100 100 $ -$	-				

Operating Mode: B+C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	T	IS	T	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ind	dex value		1.20	0.	75	0.	62	-
Index compo	nent value	-		0.64	0.75	0.55	0.62	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.64					
	Р	mW		89	.20	89	.20	_
	P _{1x1}	mW		46	.46	46	.46	
Acoustic	Zs	cm			2.85			
Parameters	Z _b	cm					2.92	
	Z _{MI}	cm	2.45					
	Z _{pii,α}	cm	2.55					
	f _{awf}	MHz	1.89	1.	.89 1.88	_		
	prr	Hz	5532.60					
	srr	Hz	43.22					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	129.20					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	125.45					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	149.85					
	p _r at z _{pii}	Мра	1.59					
	Focus	cm	5.00	5.	00	6.	00	_
Operating	Depth	cm	8.62	8.	62	7.3	39	_
control	Scanning width	%	100	10	00	10	00	_
Parameters Other Information	Freq	MHz	2.5	2	.5	2.5		_
	Power	%	100	10	00	10	00	-

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	T	S	T	В	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ind	lex value		1.20	0.	75	0.	62	-
Index compo	nent value	_		0.64	0.75	0.55	0.62	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.64					
	Р	mW		89	.20	89	.20	_
	P _{1x1}	mW		46	.46	46	.46	
Acoustic	Zs	cm			2.85			
Parameters	Z _b	cm					2.92	
	Z _{MI}	cm	2.45					
	Z _{pii,α}	cm	2.55					
	f _{awf}	MHz	1.89	1.3	89	1.8	88	_
	prr	Hz	5532.60					
	srr	Hz	43.22					
	n _{pps}		1				surface 62 0.62 0.20 0.20	
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	129.20					
Information	$I_{spta,\alpha} \text{ at } z_{pii,\alpha} \text{ or } z_{sii,\alpha}$	mW/cm ²	125.45					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	149.85					
	p _r at z _{pii}	Мра	1.59					
	Focus	cm	5.00	5.	00	6.	00	_
Operating	Depth	cm	8.62	8.	62	7.3	39	_
control	Scanning width	%	100	1(00	10	00	_
conditions	Freq	MHz	2.5	2	.5	2	.5	_
	Power	%	100	1(00	10	00	_

Operating Mode: B+M

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			MI	Т	IS	TI	B	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum ind	dex value		0.12	0.	80	1.9	98	0.95
Index compo	nent value			0.74	0.80	1.88	1.98	
	$p_{r,\alpha}$ at Z_{MI}	Мра 0.20						
	Р	mW		178	8.00	178	8.00	178.00
	P _{1x1}	mW		92	.71	92.	.71	
Acoustic	Zs	cm			5.40			
Parameters	Z _b	cm					6.75	
	Z _{MI}	cm	6.66					
	Z _{pii,α}	cm	6.80					
	f _{awf}	MHz	2.80	2.	75	2.8	80	2.76
	prr	Hz	0.00					
	srr	Hz	-					
	n _{pps}	•	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	43.21					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	45.52					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	78.55					
	p _r at z _{pii}	Мра	0.12					
	Focus	cm	5.00	4.	00	6.	00	6.00
Operating	Depth	cm	6.16	8.	62	7.3	39	7.39
control	Scanning width	%	100	1(00	10	00	100
conditions	Freq	MHz	2.5	2	.5	2.5		2.5
	Power	%	100	1(00	10	00	100

Operating Mode: CW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		1.02	0.	12	0.	56	-
Index compo	nent value	-		0.08	0.12	0.44	0.56	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.57					
	Р	mW		254	1.00	254	l.00	_
	P _{1x1}	mW		132	2.29	132	2.29	
Acoustic	Zs	cm			6.15			
Parameters	Z _b	cm					6.14	
	Z _{MI}	cm	6.11					
	Zpii,α	cm	6.15					
	f _{awf}	MHz	2.36	2.3	2.38 2.36		36	_
	prr	Hz	252.00					
	srr	Hz	1.97					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	79.50					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	80.23					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	105.22					
	p _r at z _{pii}	Мра	2.63					
	Focus	cm	5.00	3.	00	6.	00	_
Operating	Depth	cm	6.16	9.	62	9.	62	_
control	Scanning width	%	100	10	00	1(00	_
conditions	Freq	MHz	2.5	2.5		2.5		_
	Power	%	100	10	00	10	00	-

Operating Mode: CFM-M

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	IB	TIC
	Index label			At	Below	At	Below	
				surface	surface	surface	surface	
Maximum inc	lex value		0.66	0.4	42	0.4	42	-
Index compo	nent value	-		0.42	0.42	0.38	0.42	
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.46					
	Р	mW		43	.67	43	.67	-
	P _{1x1}	mW		17	.96	17	.96	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	2.34					
	Z _{pii,α}	cm	-					
	f _{awf}	MHz	4.90	5.	13	5.13		-
	prr	Hz	_					
	srr	Hz	56.93					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	_					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	60.23					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	84.44					
	p _r at z _{pii}	Mpa	_					
	Focus	cm	3.00	2.	00	2.	00	_
Operating	Depth	cm	4.29	3.	06	3.	06	-
control	Scanning width	%	100	1(00	1(00	_
conditions	Freq	MHz	5.3	5.3		5.3		_
	Power	%	100	10	00	10	00	

Operating Mode: B

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	T	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum inc	lex value		0.45	0.	83	0.8	83	-
Index compo	nent value	-		0.83	0.83	0.77	0.83	
	$p_{r,\alpha}$ at Z_{MI}	Мра	0.91					
	Р	mW		76	.87	76	.87	_
	P _{1x1}	mW		31	.61	31	.61	
Acoustic	Zs	cm			_			
Parameters	Z _b	cm					_	
	Z _{MI}	cm	3.22					
	Z _{pii,α}	cm	-					
	f _{awf}	MHz	4.09	5.24		5.24		-
	prr	Hz	_					
	srr	Hz	39.06					
	n _{pps}		1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	-					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	33.22					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	55.45					
	p _r at z _{pii}	Мра	-					
	Focus	cm	4.00	2.	50	2.	50	_
Operating	Depth	cm	6.75	3.	06	3.	06	_
control	Scanning width	%	100		00		00	_
conditions	Freq	MHz	5.3/4.0		/4.0	5.3/		_
	Power	%	100	10	00	10	00	-

Operating Mode: B+C

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			МІ	Т	IS	Т	B	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	dex value		0.49	0.9	93	0.8	87	-
Index compo	nent value			0.88	0.93	93 0.77 0.87		
	$p_{r,\alpha}$ at Z_{MI}	Мра	1.11					
	Р	mW		71.	.86	71	.86	_
	P _{1x1}	mW		29.	.55	29	.55	
Acoustic	Zs	cm			1.77			
Parameters	Z _b	cm					1.93	
	Z _{MI}	cm	2.22				Below surface 0.87 0.87 1.86 0.55 1.1.93 1.1.93 1.18 1.193 1.18 1.193 1.193 1.18 1.18 1.18 1.18 1.18 1.18 1.18 1.18 1.18 1.18 1.18 1.193 1.18 1	
	Z _{pii,α}	cm	2.40					
	f _{awf}	MHz	5.16	5.21		5.18		-
	prr	Hz	7282.60					
	srr	Hz	_					
	N _{pps}	-	1					
Other	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	142.87					
Information	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$	mW/cm ²	140.22					
	I_{spta} at z_{pii} or z_{sii}	mW/cm ²	168.66					
	p _r at z _{pii}	Мра	1.56					
	Focus	cm	3.00	3.	50	3.	50	_
Operating	Depth	cm	3.06		22			-
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	5.3/4.0	5.3/4.0		5.3/4.0		-
	Power	%	100	10	00	100		-

Operating Mode: PW

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

			MI	Т	IS	Т	IB	TIC
	Index label			At surface	Below surface	At surface	Below surface	
Maximum ind	dex value		0.78	0.	30	1.	12	-
Index compo	nent value	-		0.20	0.30	1.05	1.12	
	$p_{r,\alpha}$ at Z_{MI}	Mpa	1.72					
	Р	mW		24	.26	24	.26	-
	P _{1x1}	mW		9.	98	9.	98	
Acoustic	Zs	cm			1.78			
Parameters	Z _b	cm					2.07	
	Z _{MI}	cm	2.22					
	Z _{pii,α}	cm	2.30					
	f _{awf}	MHz	4.83	4.	4.86 4.8		83	-
	prr	Hz	7287.60					
	srr	Hz	56.93				.98 2.07 .83	
	n _{pps}	_	1					
Other Information	$I_{pa,\alpha}$ at $z_{pii,\alpha}$	W/cm ²	120.53					
mormation	$I_{spta,\alpha} \text{ at } z_{pii,\alpha} \text{ or } z_{sii,\alpha}$	mW/cm ²	118.55					
	I _{spta} at z _{pii} or z _{sii}	mW/cm ²	142.22					
	p _r at z _{pii}	Мра	2.27					
	Focus	cm	3.00	3.	00	3.	00	_
Operating	Depth	cm	4.29	10	.45	5.	22	_
control	Scanning width	%	100	10	00	10	00	_
conditions	Freq	MHz	5.3	5.3		5.3		_
	Power	%	100	10	00	10	00	_

Operating Mode: <u>B+M</u>

NOTE 1 Only one operating condition per index.

NOTE 2 Data should be entered for "at surface" and "below surface" both in the columns related to TIS or TIB.

NOTE 3 Information need not be provided regarding TIC for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

NOTE 4 If the requirements of 201.12.4.2a) are met, it is not required to enter any data in the columns related to TIS, TIB or TIC.

NOTE 5 If the requirements of 201.12.4.2b) are met, it is not required to enter any data in the column related to MI.

NOTE 6 The depths z_{pii} and $z_{pii,\alpha}$ apply to NON-SCANNING MODES, while the depths z_{sii} and $z_{sii,\alpha}$ apply to SCANNING MODES.

Appendix B: Guidance and Manufacturer's Declaration

1. Guidance and manufacture	er's declaration – electromagnetio	c emissions
	-	c environment specified below. The is used in such an environment.
Emissions test	Compliance	Electromagnetic environment –
		guidance
RF emissions	Group 1	The QBit 3/QBit 5 uses RF energy
CISPR 11		only for its internal function.
		Therefore, its RF emissions are
		very low and are not likely to
		cause any interference in nearby
		electronic equipment.
RF emissions	Class A	The QBit 3/QBit 5 is suitable for
CISPR 11		use in medical establishments
		including domestic
Harmonic emissions	Class A	establishments and those directly
IEC 61000-3-2		connected to the public
Voltage fluctuations/	Complies	high-voltage power supply
flicker emissions		network
IEC 61000-3-3		that used for non-domestic
		purposes.

2. Guidance and manufacturer's declaration - electromagnetic immunity

The QBit 3/QBit 5 is intended for use in the electromagnetic environment the QBit 3/QBit 5 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Electrical fast	±2 kV for power	±2 kV for power	Mains power quality
transient/burst	supply lines	supply lines	should be that of a typical
IEC 61000-4-4	±1 kV for input/output	±1 kV for input/output	commercial or hospital
	lines	lines	environment.
Surge	±1 kV line(s) to line(s)	±1 kV line(s) to line(s)	Mains power quality
IEC 61000-4-5	±2 kV line(s) to earth	±2 kV line(s) to earth	should be that of a typical
			commercial or hospital
			environment.
interruptions	<5 % UT	<5 % UT	Mains power quality
and	(>95 % dip in UT)	(>95 % dip in UT)	should be that of a typical
voltage	for 0,5 cycle	for 0,5 cycle	commercial or hospital
variations	40 % UT	40 % UT	environment. If the user
on power	(60 % dip in UT)	(60 % dip in UT)	of the QBit 3/QBit 5
supply	for 5 cycles	for 5 cycles	requires continued
input lines	70 % UT	70 % UT	operation during power
IEC	(30 % dip in UT)	(30 % dip in UT)	mains interruptions, it is
61000-4-11	for 25 cycles	for 25 cycles	recommended that the
	<5 % UT	<5 % UT	QBit 3/QBit 5 be powered
	(>95 % dip in UT)	(>95 % dip in UT)	from an uninterruptible
	for 5 sec	for 5 sec	power supply or a battery.
Power	3 A/m	3 A/m	Power frequency
frequency			magnetic fields should be
frequency			at levels characteristic of
(50-60 Hz)			a typical location in a
magnetic field			typical commercial or
IEC 61000-4-8			hospital environment.
NOTE UT is the a.c. m	ains voltage prior to app	lication of the test level.	

3 Guidance and ma	anufacturer's decl	aration – electrom	nagnetic immunity
The QBit 3/QBit 5	is intended for u	se in the electror	nagnetic environment specified below. The
customer or the us	er of the QBit 3/Q	Bit 5 should assur	e that it is used in such an environment.
3.1. Immunity	IEC 60601 test	IEC 60601 test	Electromagnetic environment – guidance
test	level	level	

		1	
Conducted RF	3 Vrms	3 Vrms	Portable and mobile RF communications
IEC 61000-4-6	150 kHz to 80	3 V/m	equipment should be used no closer to
Radiated RF	MHz		any part of the QBit 3/QBit 5, including
IEC 61000-4-3	3 V/m		cables, than the recommended separation
	80 MHz to 2,5		distance calculated from the equation
	GHz		applicable to the frequency of the
			transmitter.
			Recommended separation distance
			$d = 1,2 \sqrt{P}$
			$d = 1, 2 \sqrt{P}$ 80 MHz to 800 MHz
			$d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			where P is the maximum output power rating o
			the transmitter in watts (W) according to the
			transmitter manufacturer and <i>d</i> is the
			recommended separation distance in metres (
			Field strengths from fixed RF transmitters, as
			determined by an electromagnetic site survey,
			should be less than the compliance level in ea
			frequency range. _b
			Interference may occur in the vicinity of
			equipment marked with the following symbol:
			$((\cdot,\cdot))$
		a latada a u fua au cana	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the QBit 3/QBit 5 is used exceeds the applicable RF compliance level above, the QBit 3/QBit 5 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the QBit 3/QBit 5.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between

portable and mobile RF communications equipment and the QBit 3/QBit 5

The QBit 3/QBit 5 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the QBit 3/QBit 5 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the QBit 3/QBit 5 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum	Separation distance	e according to frequency of t	ransmitter
output	m		
power of transmitter W	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Appendix C: Measurement Results Summary

Measurement	Useful Range	Accuracy
Distance	Full Screen	<±5%
Circumference:	Full Screen	<±5%
trace method, ellipse method		
Area:	Full Screen	<±10%
trace method, ellipse method		
Volume	Full screen	<±10%
Angle	Full screen	<±5%
Time	Full Screen	<±5%
Heart rate	Full Screen	<±5%
Velocity	Full Screen	<±10%

Appendix D: Display Accuracy and Acoustic Measurement Uncertainties

According to IEC60601-2-37 and NEMA UD-3 2004, the display accuracy and acoustic measurement uncertainties are summarized in the table below.

ltem	Measurement Uncertainty
nem	(Percentage, 95% Confidence Value
Center Frequency	±15%
Acoustic Power	±30%
Acoustic Intensity	±30%
Peak Rarefactional	+4 = 9/
Pressure	±15%

Display accuracy of MI is ±20%, and TI is ±40% or <0.1, if MI, TI below 0.5.

Appendix E: Transducer Maximum Surface Temperature

the According section 201.11 standard to the requirements of in the IEC60601-2-37:2007/AMD1:2015, the transducer surface temperature has been tested in two kinds of conditions: the transducer suspended in still air or transducer contacting human-tissue mimicking material. The calculation of the expanded uncertainty is based on the ISO Guide tout ye Expression of uncertainty in measurement. Three transducer samples have been tested and the confidence coefficient is at 95%, the value of t.975 is 4.30.

The measurement data were obtained under the test conditions employed at CHISON.

	Maximum surface temperature(℃)	Maximum surface temperature(℃)
Transducer model	Contacting human-tissue mimicking material	Suspending in air
D3C60L	41±1	48±1
D7L40L	40±1	47±1
D6C12L	41±1	43±1
D3P64L	41±2	48±1
D7L40L-REC	38±2	47±2
D7C10L	41±1	43±2
D5C20L	41±1	46±1
D3C20L	38±2	47±2
D6C15L	41±1	46±1
D12L40L	39±2	45±1

NOTE: Values following the "±"mark indicate the expanded uncertainty with a confidence lever of 95%, t.975=4.30.

Appendix F: Procedures of setting network sharing

For Windows set up, set up a shared document

1. Choose the file you want to share, as the "test" file

1 test	2015/5/20 19:28	File folder	
System Volume Information	2014/12/1 20:41	File folder	
4388_1485765_MVM_7.tmp	2015/5/18 4:27	TMP File	0 KB
4388_1485765_MVM_10.tmp	2015/5/18 4:27	TMP File	0 KB
7208_1754125_MVM_7.tmp	2015/5/13 18:36	TMP File	0 KB
7208_1754125_MVM_10.tmp	2015/5/13 18:36	TMP File	0 KB
Name	Date modified	Туре	Size

2. Right click this file, choose "properties", and click "share".

	Sharing	Security	Customi				
Netwo	ork File and	d Folder S	haring				
	test						
	Not S	Shared					
	ork Path:						
Not S	Shared						
5	Share						
Advar	nced Shari	ng					
	ustom per			ltiple sh	nares, ar	d set ot	ner
adva	nced shar	ing option	S.				
4	Advanc	ed Sharin	g				
Pace	vord Prote	ction					
500500		5					
	ole must ha ss shared		account	and pas	sword to	ir this co	mputer to
	ando this	setting, us	e the <u>Net</u>	work an	d Sharin	g Center	
acce	lange uns						

3. Then you can see the sharing setting interface, as you can see in the picture, choose "everyone", and click "add".

Choose people to share with			
Type a name and then click Add, or click the arrow to	find someone.		
Everyone	~	Add	
Name	Permission L	evel	
L junye zhang (zhanngjunye@hotmail.com)	Owner		

4. Choose "read/write" in the permission level in everyone, then click "share", after that, confirm.


Choose people to share with			
Type a name and then click Add, or click the arrow to	find someone.		
	~	Ad	d
Name	Permission Lev	el	
Severyone	Read/Write 🔻		Read
👗 junye zhang (zhanngjunye@hotmail.com)	Owner	~	Read/Write
			Remove

5. If the windows have not set the code, then turn off the password protected sharing is necessary. Instructions as the image shows.

a) Click "network and sharing center" in the sharing properties.

	test Properties
General	Sharing Security Customize
Netwo	ork File and Folder Sharing
	test Shared
1000	vork Path: MES-ZHANG\test
	Share
adva	anced sharing options.
adva	anced sharing options.
Pass	Advanced Sharing
Passy Peop acce	Advanced Sharing word Protection ble must have a user account and password for this computer to
Passy Peop acce	Advanced Sharing word Protection ble must have a user account and password for this computer to ess shared folders.

b) In the network and sharing center interface, choose "public", in the password protected sharing, choose turn off password protected sharing.



Set up in system

IP set up

1. First confirm the service address of shared files, you can get the IP address in windows

interface.

In windows "start"-"run" type in "cmd" and enter, then type in " ipconfig" and enter, you can see the IP address of local service.

Connection-specific DN	s	S	ufi	Fi>	<		
Link-local IPv6 Addres	s						fe80::d91e:e058:f895:336d%3
IPv4 Address							192.168.1.131
Subnet Mask							255.255.255.0
Default Gateway							192.168.1.1

2.Enter into set up interface, choose net. Choose "DHCP" (automatically get the IP address) or

"static"(type in IP address manually).

Tips: if you want to type in IP address manually, make sure the IP address is in the same internet section with the service, and won't confused with other IP in the LAN.

CHISON	ł		Setup						×
General Measu	urement Comment	BodyMark	ExamMode	Report	DICOM	Net	Sys	tem	
Settings	Storage Net Printer								
Net									
O DHCP	Static								
lp address	192.168.30167								
Netmask	255.255.255.0								
Gateway	192.168.301								
lp address	192.168.301		Test						
							8	×.	Ŷ
			Default	0		Cancel	En		ġ,

3.Choose net "storage" interface, type in service name, IP address, user name, password and the name of shared files, click "add" to add a network storage, you can choose the export route. As shown in the picture.

CHISON	Setu	p	8
General Measurement Comm	nent BodyMark ExamM	lode Report DICOM	Net System
Settings Storage Net	Printer		
Configure Service			
Service Name		P Address	Ping
User Name		Password	
Shared Dir			
		Add	Clear
Service List			
Service Name	IP Address	Shared Dir	Default
lala	192.168.30.93	share	
SS	192.168.30.171		
li	192.168.30.249		yes
Delete			
Delete			
			📓 🖓 🦞
			Cancel En 🖵 🖑

	Ex	portInfo	×
UDisk	Service Name	serer	•
DicomPrint DicomStorage	Image Type	💿 BMP 🔘 JPG 📃 Single Frame wi	th Raw data
Network Storage		🗹 Cine transform to AVI	
Folder Name		2016-09-12-0591_Kim	
		Export	

Ping: test if the IP is connected or not.

Clear: clear all the IP address, user name, password and names of shared file

Update: update the content to the chosen item.

Delete: delete the chosen service item.

Default: set the chosen item as the default net route.

NOTE: you can add multi-numbers of network storage service to realize the transmission among multi systems.

NOTE: If windows turn off password protected sharing, then in set up, you can type in user name and password arbitrarily.

Net Printer Setup

1. First press [Setup] key to enter into the Setup interface, then select Net and Net Printer.

CHISON		Setup				×
General Measurement Comm	nent BodyMark	ExamMode R	eport DICC	OM Net	Syste	m
Settings Storage Net	Printer					
Configure Net Printer						
Name		IP Address				
PPD File						
			Add	Clear		
Name	IP Address			Default		
QQ	192.168.30.51	No				
	Import PPD Files					
					8	₽ \$
		Default		Cancel		9

2. Configure Net Printer

Operation:

a) Input name and IP Address which has installed printer equipment.

b) Choose [Ping] to confirm whether the net connect normally.

c) Select PPD File and press **[Add]**, the net connecting information can be saved on the Printer List.

3. Finally choose the information on the Printer List to finish the printing according to the requirement.

ΕN

CHISON



CHISON Medical Technologies Co., Ltd.

No.3 Changjiang South Road, Xinwu District, Wuxi, 214028 Jiangsu, P.R. China