

User's manual

Head immobilizers and ankle immobilizers

General information	page	1
Warnings	page	1
Description of product	page	1
Operating instructions	page	2
Maintenance and cleaning	page	3
Spare parts	page	3



1



2



3



4



5



6




7



8



9



10



11



12



13



14

1. GENERAL INFORMATION

1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

1.2 Conservation of the instruction manual

The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

1.3 Symbols used



General or specific warning



See instructions for use



Lot number



Product code



Product compliant with specifications of the Directive 93/42/CEE



Single use

1.4 Servicing request

For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on tel. 0039 0521 541111, fax 0039 0521 541222, e-mail info@spencer.it or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY.

1.5 Demolition

Follow the current regulations.

1.6 Labelling

The serial number as indicated below can be found on each appliance and must not be removed or covered. In order to facilitate assistance please indicate or communicate the lot number (LOT) on the label.

2. WARNINGS



2.1 General warnings

- Before carrying out any kind of operation on the appliance, the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Regularly check the appliance.
- In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and the safety, the appliance must be immediately removed from service.

- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.

- The appliance must not in any way be tampered with. In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself.

- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.

- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.

- Handle with care.



2.2 Specific warnings

- The product must be used by trained personnel only.

- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.

- When the device is being used, the assistance of qualified staff must be guaranteed.

- The head immobiliser or ankle immobiliser should not be exposed nor get in contact with heat sources or flammable agents.

- The device must be applied to the patient by at least two operators.

- Always check the integrity of all the parts of the device before use.

- Replace the belts immediately if worn out or damaged (if present).

- Select the fixation points of the belts accurately (if present).

- Position and regulate the device so not to obstruct any operations the rescuers may need to carry out and the use of the rescue appliances.

2.3 Contraindications and side effects


The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

3. PRODUCT DESCRIPTION

3.1 Intended use

The Spencer head immobiliser and ankle immobiliser have been designed to increase the stability of the cranium and the spine and likewise the distal areas of the lower limbs during movement and transport of patients on the spine board. Their function is to eliminate the results caused any movement that may occur and in any direction with a guarantee of stability even if subjected to stress by vertical components (perpendicular to the fixing board). Designed to immobilise the head or the ankles when moving a patient who has suspected spinal injuries. The Spencer head immobiliser and ankle immobiliser, if used in combination with the spine board, restraint belts and rigid cervical collar, offer efficient immobilising which permits extended transport or transport over uneven surfaces even of patients in a critical state.

3.2 Main components (fig. 1 - 2 page I)

 Nylon belts to optimise stability of device

- B** Side support cushions with ear holes
- C** Semi rigid base for application of cushions with Velcro® system
- D** Removable cushion for use as paediatric head immobiliser
- E** Restraint belts for use with spine board
- F** Shaped surface for restraint of the ankles

3.3 Models

SH00300A
 Fix Head Advanced universal head immobiliser, orange
 SH00310A
 Fix Head Advanced universal head immobiliser, black
 SH00201A
 Spencer Contour Anatomic universal head immobiliser, yellow/black
 SH00200A
 Spencer Contour Anatomic universal head immobiliser, blue
 SH00111C
 Super Blue Compact universal head immobiliser, yellow
 SH00112C
 Super Blue Compact universal head immobiliser, orange
 SH00117C
 Super Blue Compact universal head immobiliser, green
 SH00110C
 Super Blue Compact universal head immobiliser, blue
 SH00104C
 Super Blue Compact universal head immobiliser, black
 SH00240A
 Tango Fix Integrated adult/paediatric head immobiliser
 SH00250A
 U-Fix Adult head immobiliser
 SH00260A
 Pedi Fix Paediatric head immobiliser
 SH00262A
 Pedi Go Paediatric head immobiliser, only cushions, for Baby
 Go and Pedi Loc spine boards
 ST02605A
 Pedi Roll Paediatric head immobiliser
 SH00160A
 751 Inflatable disposable head immobiliser
 SH00130C
 755 Universal disposable head immobiliser
 SH00120C
 756 Stabilized disposable head immobiliser
 SH00150A
 FXA Universal ankle immobiliser
 SH00151A
 FXA Pro Anatomic universal ankle immobiliser

3.4 Technical data (chart 3.4 page 4)

4. OPERATING INSTRUCTIONS

4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself. Keep the original packaging for use in case of any further transport.

Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client.

The appliance must be stored in a dry place free from humidity. During storage take care not to put heavy materials onto the device. In no way and under no circumstances should the device be considered as a work top.

4.2 Preparation

On receipt of the product:

- remove the packaging and display the material so that all components are visible
- check that all the components/pieces on the accompanying list are present

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage.

Before putting the device in service make the following careful checks:

- general integrity (absence of cuts, holes, abrasions)
- general functionality
- Velcro® and belts tightness (if present)

If everything is regular, the device is ready to be used.

4.3 Functioning

4.3.1 Application of head immobiliser Fix Head, Spencer Contour, Super Blue, Tango Fix, U-Fix, Pedi Fix (fig. 3 - 4 - 5 - 6 page II)

The product must be used by trained personnel only and must be applied to the patient by at least two operators.

Attach the device to the spine board before use using the belts supplied (E fig. 1 page I), making sure that the side of the spine board with Velcro® is topside. Check the position and the correct fixing of the belts before use with the patient. Load the patient on the spine board maintaining the correct alignment. Block the patient's head by placing the special anatomical cushions symmetrically on each side and without applying any pressure. Fix the two belts which restrain the head (A fig. 1 page I) by passing them over the anatomical cushions and uniting them with the buckles on the opposite side.

Proceed and attach the belts, adjust the tension and block them using the Velcro® system. Proceed in the same way with both belts making sure to block the chin strap first and then the forehead strap. Before transporting the patient it is necessary to proceed with the application of the special restraints for the spine board which will block the patients shoulders so as to avoid compression of the vertebrae in the neck area.

4.3.2 Double use adult/paediatric head immobiliser: Spencer Contour and Tango Fix

Spencer Contour and Tango Fix head immobilisers can be used both for adult and paediatric patients.

They have a removable cushion in the centre of the base applied with Velcro® (fig. 7 page II). If using the head immobiliser with an adult patient, the cushion must be in position in the centre between the two cushions. When in use for a paediatric patient the cushion must be removed. This will also permit the two anatomical cushions to be positioned closer and the Velcro® permits fixing in the correct position.

The Tango Fix head immobiliser is characterised by the special "U" shape of the cushions. The cushion for immobilising the paediatric patient is situated inside the main cushion. Therefore in a confined space you have both adult and paediatric head immobiliser which will permit the correct immobilising for both paediatric and adult patient (fig. 8 page II).

4.3.3 Application of the Pedi Go, Pedi Roll, 751, 755, 756 head immobiliser

The device must be applied only by trained operators. At least two operators are essential for the application.

The Pedi Go, Pedi Roll and 751 head immobilisers are extremely simple to use and to apply. When the patient is in the correct position on the spine board, maintain the correct spinal alignment and proceed to immobilise the head by positioning the cushions at the sides of the head. The Pedi Go and Pedi Roll should then be fixed to the spine board and to the patient's head using the special belts (fig. 9 page II).

The 751 head immobiliser can be fixed using 50 mm surgical tape (fig. 10 page II).

To immobilise the patient with the 755 and 756 an immobiliser base must first be placed on the spine board. Once the patient has been loaded and correctly aligned, with the patient's head positioned on the head immobiliser base proceed to block in position by lifting the extension on the sides until they are positioned on the sides of the head but do not exert any pressure. Fixing of the device should be carried out using the purpose made straps which block the device and the head in position and then maintain them in position on the spine board (fig. 11 - 12 page II).

4.3.4 Application of FXA, FXA Pro ankle immobiliser (fig. 13 - 14 page II)

The device must be applied only by trained operators. At least two operators are essential for the application. Assemble the device on the spine board before use using the purpose made belts supplied (E fig. 1 page I). Load the patient in the correct alignment. Continue to position the ankles in the hollows on the cushion. Block the ankles in a stable, correctly aligned position using the top fixing wrap (fig. 13 - 14 page II).

4.4 Troubleshooting (chart 4.4 page 4)

5. MAINTENANCE AND CLEANING

5.1 Cleaning

Failure to carry out cleaning operations can cause the risk of cross infection.

Do not use aggressive substances or solvents of any type. The device can be washed using delicate detergents or water and soap. At the end always dry the product with a soft and clean cloth. Using high pressure water will raise the risk of lesions. Ensure you have dried the device perfectly so to avoid the presence of residual detergents, which could deteriorate or compromise the integrity and duration of the device.

5.2 Maintenance

5.2.1 Precautionary Maintenance

The person responsible for every day maintenance can only substitute the spare parts indicated on paragraph 6 "Spare Parts". All other substitutions or repairs can be carried out only by the manufacturer or by a centre authorised by the manufacturer. For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device does not require programmed servicing.

5.2.2 Special servicing



Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years. The life span can be expanded only following a general revision of the product that must be carried out by the manufacturer or by a centre authorised by the manufacturer.

6. SPARE PARTS

SH00303A

Base for orange Fix Head

SH00302A

Cushion for orange Fix Head

SH00301A

Kit head/chin straps for Fix Head

SH00242A

Base for Spencer Contour, Tango Fix, U-Fix and Pedi Fix

SH00214A

Cushion for yellow/black Spencer Contour

SH00212A

Cushion for blue Spencer Contour

SH00244A

Kit head/chin straps for Spencer Contour, Super Blue, Tango Fix, U-Fix, Pedi Fix, Pedi Go and Pedi Roll

SH00243A

Base for Super Blue

SH00100C

Cushion for yellow Super Blue

SH00101C

Cushion for orange Super Blue

SH00102C

Cushion for green Super Blue

SH00108C

Cushion for blue Super Blue

SH00103C

Cushion for black Super Blue

SH00241A

Adult/paediatric cushion for Tango Fix

SH00251A

Adult cushion for U-Fix

SH00261A

Paediatric cushion for Pedi Go

Head immobilizers and ankle immobilizers		Head immobilizers				Ankle immobilizers			
Code	Fix Head SH00300A	Spencer Contour SH00200A	Super Blue SH00110C	Tango Fix SH00240A	Pedi Fix SH00260A	Pedi Roll ST02605A	FXA SH00150A	FXA Pro SH00151A	
Tidal weight (kg)	2,3	1,5	1	2,5	1,1	0,092	0,4	0,376	
X-ray compatible	●	●	●	●	●	●	●	●	
Immobilization directions	►►►►	►►►►	►►►►	►►►►	►►►►	►►►►	►►►►	►►►►	
Floatation	●								
Compatibility with spine board	All	All	All	All	Pedi Loc / Baby Go	Pedi Loc	All	All	
Compatibility with scoop stretcher									
Fixation straps	2	2	2	2	2		1	1	
Type of fixation	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps				
Reusable	●	●	●			●	●	●	
Base									
Dimensions (mm)	375 x 10 x h275	390 x 15 x h260	384 x 21 x h290	445 x 28 x h 280	445 x 28 x h 280		345 x 120 x h90		
Weight (kg)	1,2	0,5	0,3	0,5	0,5		0,4		
Type	Semi rigid collection	Removable layer	Stuffed	Removable layer	Removable layer				
Material	Plastic compound	Expanded PE / Spentex®	Expanded PE / Spentex®	Polyethylene / rubber	Polyethylene / rubber		Closed cell expanded PE	Plastic compound	
Manufacturing system	Injection	Injection	Injection	Blow moulding / injection	Blow moulding / injection	Hand made	Hand made	Injection	
Colour	Orange	Blue	Blue	Black / yellow	Black / yellow	Blue	Orange	Black	
Head blocks									
Dimensions (mm)	225 x 90 x h150	235 x 130 x h160	245 x 80 x h160	230 x 125 x h160	190 x 70 x h123	160x80xh80			
Weight (kg)	1,2	1	0,5	1,2	0,6	0,092			
Structure		Stuffed pillows	Stuffed pillows	Rigid with soft pillows	Rigid with soft pillows				
Material	Plastic compound	Plastic compound	Plastic compound	Polyethylene / rubber	Polyethylene / rubber	Vinyl covered			
Manufacturing system	Injection	Injection	Injection	Blow moulding / injection	Blow moulding / injection	Hand made			
Auricular inspection	●	●	●	●	●	●			
Colour	Orange	Blue	Blue	Black / yellow	Black / yellow	Blue	Orange	Black	
Functioning temperature (°C)	from -20 to +40	from -20 to +40	from -20 to +40	from -20 to +60	from -20 to +60	from -20 to +40	from -20 to +40	from -20 to +40	
Storage temperature (°C)	from -20 to +40	from -20 to +40	from -20 to +40	from -20 to +60	from -20 to +60	from -20 to +40	from -20 to +40	from -20 to +40	

4.4

PROBLEM	CAUSE	REMEDY
Excessive mobility of the head/ankles	The basis presents an excessive mobility in relation to the spine board it is strapped to	Tighten the belts at the basis around the spine board
	Loose belts	Check the correct closure and tension of the fastening belts
Lesions to the device	Anatomic supports positioned wrong	Ensure that the hollow part of the anatomic support is in contact with the head/ankles
	Wear and tear, improper use	Put immediately the device out of service and contact the service centre



We declare that the appliances conform with the Directive 93/42/CEE "Medical Devices".

Si dichiara che i dispositivi sono conformi alla Direttiva 93/42/CEE "Dispositivi Medici".

Se declara que los dispositivos son conformes a la Directiva 93/42/CEE "Dispositivos Médicos".

Nous déclarons que les dispositifs sont conformes à la Directive 93/42/CEE "Dispositifs Médicaux".

ISO 9001	Quality management systems - Requirements
EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
EN 980	Graphical symbols for use in the labelling of medical devices
CEI EN ISO 14971	Medical devices - Application of risk management to medical devices

Guarantee of Quality System for the production and the final control of the products certified by the notifying body TÜV SÜD Product Service GmbH.

Sistema di Garanzia di Qualità per la produzione ed il controllo finale dei prodotti certificato dall'organismo notificato TÜV SÜD Product Service GmbH.

Sistema de Garantía de Calidad para la producción y el control final de los productos certificado por el organismo notificado TÜV SÜD Product Service GmbH.

Système de Garantie de Qualité pour la production et le contrôle final des produits certifié par l'organisme TÜV SÜD Product Service GmbH.

Warning

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.

Avvertenza

Le informazioni contenute in questo documento sono soggette a modifica senza preavviso e sono da intendersi come impegno da parte della Spencer Italia S.r.l. con riserva di modifiche. I prodotti Spencer vengono esportati in molti paesi nei quali non valgono sempre identiche regole. Per questa ragione possono esserci differenze fra quanto qui descritto ed i prodotti consegnati. Spencer lavora costantemente al perfezionamento di tutti i tipi e modelli dei prodotti venduti. Contiamo perciò sulla Vostra comprensione se dovessimo riservarci la facoltà di apportare in qualsiasi momento modifiche alla fornitura nella forma, equipaggiamento, allestimento e tecnica rispetto a quanto qui convenuto.

© Copyright Spencer Italia S.r.l.

Tutti i diritti sono riservati. Nessuna parte del documento può essere fotocopiata, riprodotta o tradotta in un'altra lingua senza previo consenso scritto della Spencer Italia S.r.l.

Advertencia

Las informaciones presente en este documento están sometidas a modificación sin aviso y se deben considerar como compromiso por parte de Spencer Italia S.r.l. con discreción de modificaciones. Los productos Spencer son exportados en muchos países en los cuales no valen siempre idénticas reglas. Por ésta razón pueden haber diferencias entre cuanto aquí descrito y los productos efectivamente entregados. Spencer trabaja constantemente en el perfeccionamiento de todos los tipos y modelos de sus productos. Contamos con su comprensión si debemos reservarnos la facultad de hacer modificaciones en cualquier momento al abastecimiento en la forma, equipamiento, instalación y técnica respecto a cuanto aquí convenido.

© Copyright Spencer Italia S.r.l.

Todos los derechos son reservados. Ninguna parte del documento puede ser fotocopiada, reproducida o traducida en otro idioma sin previa aprobación escrita de parte de Spencer Italia S.r.l.

Réserveres

Les informations contenues dans ce document peuvent être modifiées sans aucun préavis et doivent s'entendre comme engagement de la part de Spencer Italia S.r.l. sous réserve de modification. Les produits Spencer sont exportés dans beaucoup de pays où ne sont pas toujours valables les mêmes normes. Pour cette raison il peuvent y avoir des différences parmi ce que l'on a ici décrit et les produits livrés. Spencer opère constamment au perfectionnement de tous les types et modèles de produits vendus. Nous comptons donc sur votre compréhension si l'on doit se réserver la faculté d'apporter à n'importe quel moment des changements à la fourniture dans la forme, équipement, aménagement et technique par rapport à ce qui est ici convenu.

© Copyright Spencer Italia S.r.l.

Tous les droits sont réservés. Aucune partie du document peut être photocopiée, reproduite ou traduite sans aucun pré consentement écrit par Spencer Italia S.r.l.

