



GIMA

PROFESSIONAL MEDICAL PRODUCTS

***RESERVOIR MONOUSO
OXYGEN RESERVOIR***

ModelREF **34257** Oxygen Reservoir 1000 ml with valve - single useREF **34258** Oxygen Reservoir 2500 ml with valve - single useREF **34275** Oxygen Reservoir bag - 2500ml**Accessories**REF **34279** Reservoir Valve**Product Description**

The device designed for use in a breathing circuit as a reservoir for breathing gas and to assist, control, or monitor a patient's ventilation. This product is intended for use by qualified medical personnel.

Usage instruction












1. Eye check for clear airflow passage prior to use.
2. Connect the breathing bag to breathing circuit.
3. Make sure the connection is compatible and fit firmly.

Warning

1. The product should be placed in a dry and clean place and sealed conditions.
2. The operator should be trained and know the skill in CPR resuscitation.

Precautions

- Disposable product-single use only
- The patient should be closely monitored by qualified medical personnel when using this device.
- This product is hygienically packed.
- Read this instruction manual carefully before use.
- Store in clean dry conditions. Away from heat and light.

Symbols					
	Medical Device complies with Directive 93/42/EEC		Product code		Expiration date
	Caution: read instructions (warnings) carefully		Manufacturer		Keep in a cool, dry place
	Consult instructions for use		Date of manufacture		Keep away from sunlight
	Disposable device, do not re-use		Lot number		

REF **34257 - 34258 - 34275 - 34279**

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China

