# GIMA

**PULSOXIMETRO**  
***PULSE OXIMETER***  
**PULSOXYMÈTRE**  
***PULSOXIMETRO***

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**ATTENTION:** *The operators must carefully read and completely understand the present manual before using the product.*

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



## FEATURES

### Product description

COMDEK MD-600P portable pulse oximeter (cod. 34290) is an instrument with high accuracy and reliable performance.

This hand-held unit is easy to operate with the functions of adjustable alarm settings for oxygen saturation and pulse rate.

The product presents the following features:

- high accuracy
- quick assessment of Oxygen saturation and Pulse Rate
- SpO<sub>2</sub> and Pulse Rate display
- adjustable SpO<sub>2</sub> and Pulse alarm settings
- clear numerical LED display
- compact, light weight and easy to operate
- reliable quality and excellent sales service

### Specifications

Oxygen Saturation Range: 0% to 100%

Pulse Rate Range: 30 bpm to 250 bpm

Display:

Oxygen saturation, Pulse Rate, Pot, Prb, Pulse search indicator & Low battery indicator

SpO<sub>2</sub> Accuracy: ±2% between 100% to 80%  
±3% between 79% to 65%  
below 65%, not specified

Pulse Accuracy: ±1% of full scale

Default Alarm Settings:

- SpO<sub>2</sub> Upper Limit:: 100%
- SpO<sub>2</sub> Lower Limit : 85%
- Pulse Rate Upper Limit: 150 bpm
- Pulse Rate Lower Limit: 30 bpm

Power Requirements:

- DC 6 Volts by 4 pieces AA size batteries
- AC AC adapter output DC 6 V, 600 mA
- rechargeables AA batteries (optional)

Operating Temperature:

- 0°C to 55°C

SpO<sub>2</sub> Sensors:

Finger clip sensor & wrap probe (optional)

Dimensions:

172 mm x 90 mm x 36 mm

Weight:

240g

## Field of application

This is a pulse oximeter for the 21st century: reliable, portable and affordable. It is designed to the fields of hospitals, clinical surgery and home care applications.

## PRESCRIPTIONS

### Accident prevention

The product shall be used by qualified personnel only.

*Do not use the equipment in case it is damaged. Apply to your retailer.*

*Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.*

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

## TRANSPORT

### Transport procedures

It is not necessary to follow particular transport procedures. It is however recommended to carefully handle it and avoid shocks. The equipment is delivered in an imitation leather case allowing a comfortable shoulder transport.

### Unpacking

*Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.*

*They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).*

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

## INSTALLATION

### Installation and connections

The pulse oximeter is delivered with the following components: a practical imitation leather case, 4 "AA" alkaline batteries, AC adapter with 6 V output, 600 mA, a finger sensor for heart beat survey. To feed the equipment with 4 "AA" batteries, open the apposite compartment on the back and insert the batteries. Pay attention not to invert the direction of poles (see the picture inside the battery compartment) and then replace the cover.

To feed the equipment with an AC adapter, insert the adapter jack into the corresponding socket on the equipment base and plug in the adapter.



*With the exception of rechargeable batteries, remove common or alkaline batteries when using the AC adapter for charging.*

Insert the finger sensor connector into the pulse oximeter (when the sensor is not inserted, the display views the following writing: **Prb**). It is sufficient to make the grooves present on the connector and on the socket match.

## OPERATION

### Instrument description and symbols

- Select** Pushbutton for the selection of SpO<sub>2</sub> Max and Min values and PULSE/min.
-  Pushbutton to increase SpO<sub>2</sub> Max and Min values and PULSE/min.
-  Pushbutton to decrease SpO<sub>2</sub> Max and Min values and PULSE/min.
-  Turn on/off audio signals
-  Indicator for disabled audio signals.
-  Run-down battery indicator. When the battery charge is low, the battery indicator lights up and gives out two beeps. In this case it is necessary to replace the batteries. At this point the run-down battery indicator turns off.
- Prb** Sensor not connected with the equipment
- Pot** No patient connected
- Lo** Low signal

### Operating instructions

To turn the equipment on, press the green ON/OFF pushbutton. To confirm power-up, an acoustic signal is released and **PWR** pilot light turned on. The right display views the following writing: **Pot** (no patient connected).

Insert the finger of the patient into the examination sensor. It is possible to use any finger, however it is advisable to select the finger suitable to the sensor dimensions. (The PULSE OXIMETER can be fitted with different sensors, follow the instructions).

*In order to perform an accurate detection of the value, the patient shall keep still.*

Wait 6 seconds and the sensor starts up measuring. The left display shows the value of blood oxygen saturation rate (SpO<sub>2</sub> %) and the right display shows the heart rate (Pulse/min).

COMDEK PULSE OXIMETER is provided with a sound signal warning that the values detected by the sensor are higher or lower than the ones stored in the equipment.

Default alarm settings are:

SpO <sub>2</sub>	Max	100%
SpO <sub>2</sub>	Min	50%
Heart rate	Max	250 bpm
Heart rate	Min	30 bpm

To change these values, proceed as follow:

• **SpO 2 alarm max. value**

Press **Select** once, the red led (SpO 2 Min) lights up and a value blinks on the left display. This value represents the **maximum** value for oxygen saturation (SpO 2 ).

To change this value, press key ▲ to increase it and key ▼ to reduce it.

• **SpO 2 alarm min. value**

Press **Select** two times, the red led (SpO 2 Min) lights up and a value blinks on the right display. This value represents the **minimum** value for oxygen saturation (SpO 2 ).

To change this value, press key ▲ to increase it and key ▼ to reduce it.

• **Heart rate alarm max. value**

Press **Select** three times, the red led (Pulse Max) lights up and a value blinks on the left display. This value represents the **maximum** value for heart rate (Pulse/min). To change this value, press key ▲ to increase it and key ▼ to reduce it.

• **Heart rate alarm min. value**

Press **Select** four times, the red led (Pulse Max) lights up and a value blinks on the right display. This value represents the **minimum** value for heart rate (Pulse/min). To change this value, press key ▲ to increase it and key ▼ to reduce it.

When it stops blinking, the equipment returns to the original conditions.

The new alarm settings will be automatically saved.

When the batteries run down or when the equipment is disconnected from power supply, the original alarm settings are restored.

During measurement, the pulse oximeter display shows an intermittent scale of 10 leds indicating the search for beating. As seconds go by, the scale grows longer and becomes stable. When the signal is good, the search bar should show at least 4 leds.

If the loud speaker is on, the equipment gives out intermittent audio signals. Press  pushbutton to stop the sound and the red led  on the display lights up. Press the pushbutton again and the sound starts up again.

The alarm keeps on sounding for a minute and then it stops.

## Switching off

After the measurement has been taken, turn the equipment off by means of ON/OFF pushbutton. A beep is given out and then the equipment stops.

## MAINTENANCE

### General maintenance instructions

The product has been conceived to last forever; it is resistant to most part of chemicals except for methyl ethyl ketones (M.E.K.), diluents and nail polish removers.

*To clean the instrument, make use of a cloth dampened with disinfectant.*

## GUARANTEE

### General guarantee conditions

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, doppler probes, SpO<sub>2</sub> sensors, lamps, batteries, electrodes, handpieces, washers and packings, resistances, etc.

No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

*This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.*

The guarantee will be invalidated if the registration number has been removed, cancelled or altered.

Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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---

**Disposal:** The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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