## HEINE NT300° Ladegerät HEINE NT300° Charger

HEINE Optotechnik GmbH & Co. KG  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49 (0) 81 52 / 38 - 0  
Fax +49 (0) 81 52 / 38 - 202  
E-Mail: info@heine.com · www.heine.com  
med 0809 1 / 3.12

## Instructions

### HEINE NT300° charger

Before using the HEINE NT300° read these instructions carefully and keep them in a safe place for future reference.

#### Safety Information

##### Symbols

In accordance with local regulations this product should be disposed of as an electronic device separately.

Double insulated equipment of Protection Class II

Follow the safety information.

Medical Electrical Equipment  
WITH RESPECT TO ELECTRICAL SHOCK FIRE, AND  
MECHANICAL HAZARDS ONLY IN ACCORDANCE  
WITH UL60601-1/CAN/CSA C22.2 No. 601.1  
28KY

"UL certification is valid for 120VAC/60Hz and BETA L  
and F.O. Laryngoscope handles"

**Range of applications:** The HEINE NT300° is designed exclusively for charging HEINE rechargeable handles.

**Please note: The performance of this instrument can only be guaranteed if genuine HEINE instruments and spare parts are used.**

**General safety:** Under no circumstances should the mains transformer housing or charger be opened. Repairs may only be carried out by qualified personnel (exception: changing the rechargeable battery).

**Safe installation:** The charger must be placed on a solid surface. Do not place the NT300 on soft furniture, cloth or carpet. Do not insert any objects in any opening in the housing.

**Safety in use:** The NT300 may only be used at the mains voltage specified on the identification label. Unplug the unit from the mains if not in use over long periods. If liquid enters the housing, ensure that the unit is tested by authorized personnel immediately. The NT300 must not be used in rooms where an explosion hazard is present.

The following handles are not compatible with the NT300 charger:  
BETA TR, BETA R, BETA cord handle.

#### Installation and use

– The NT300 fulfils the requirements regarding the electromagnetic compatibility. However, its function might be affected by electrical and magnetic fields, such as those which occur in connection with radiological devices or mobile telephones.

– Plug the mains transformer into the mains and the charger is ready for use.

– NT300 ensures that batteries are kept in good condition irrespective of handle size.

– When replacing the handles in the charging compartments, ensure that the internal contacts in the compartment are clean and that the instruments are switched off.

– Instruments placed next to each other must not touch each other.

– If different-sized handles are used, it is possible to fit adaptors so that each handle is held in place correctly. The NT300 must be disconnected from the mains supply when the adaptors are being fitted.

– Ideally, the rechargeable handle should always be replaced and stored in the charger. Handles cannot be overcharged.

– The rechargeable battery can only be charged when fitted to the handle.

– If the second charging compartment is not in use it can be sealed off by the dust cap which is supplied with the unit.

– The pulsating light around the charging compartment shows that current is flowing and that the battery is being recharged.

– When the light is permanently lit, the charging sequence is completed. The battery in the handle is fully-charged and is ready to supply maximum power.

– If the light is switched off, this indicates that the charging compartment is empty or the unit is defective.

#### Cleaning and Maintenance

– We recommend that external surfaces of the NT300 be wiped clean with a disinfectant cleaning wipe. This should only be carried out by qualified personnel.

– Any disinfectant used must be approved by the disinfectant manufacturer for use on plastics. The manufacturer's instructions must be followed.

– The external surfaces may also be wiped clean with a cloth moistened with water or with a neutral pH or slightly alkaline cleansing agent.

– Isopropyl alcohol may only be used at concentrations up to 70%.

– Spray disinfection and soaking are not permitted and may cause damage.

– The article can be neither machine-cleaned, nor thermal disinfected or sterilised. These processes lead to irreparable damage to the product.

#### Technical specification

Mains supply:	100 – 240 Vac / 50 – 60 Hz for Europe 120 Vac / 50 – 60Hz for USA
Current consumption:	max. 30 VA
Fuse:	Integral overload protection
Charging current:	Optimised charging cycle CC-CV for Li-ion BETA L, maximum charging current 1,5 A deltaU for all other HEINE batteries except Li-ion BETA L, maximum charging current 0,5 A
Charging time:	Li-ion BETA L → ca. 2h Other HEINE rechargeable batteries → up to 4h
Operation:	continuous
Classification:	Type BF
Safety class:	double-insulated
Protection class:	IP 20
Dimensions:	142 x 66 x 61 mm
Weight:	830 g
Ambient temperature:	+10 ° to +40 °C

#### Applicable standards:

**IEC 60601-1** (Electrical medical products – Part 1: General conditions for safety including the main performance features)

**UL 60601-1** Standard For Safety For Medical Electrical Equipment, Part1: General requirements for Safety, First Edition, Revision date 2006-04-26)

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### Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries).We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Any modification of a HEINE product with parts or additional parts which do not conform to the original HEINE specification will invalidate the warranty for the correct function of the product and further invalidate any warranty claims which result from such a change or modification. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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