

**Concentrated instrument cleaner
with proteolytic enzymes****CAUTION****Properties**

GERMOCID ENZIMATICO is a proteolytic enzymatic detergent for cleaning and cleansing of medical and dental instruments. It acts effectively on organic waste, even dry ones, and stubborn dirt: it can be used on glass, plastic and metal materials. FOR PROFESSIONAL USE, suitable for manual washing and also for use in ultrasonic devices.

Instruction to use: Dilute GERMOCID ENZIMATICO at the rate of 0.6-10% in a liter of tap water, preferably warm, and mix thoroughly. Before soaking the instruments in the prepared solution, it is recommended to decontaminate them with a specific product, then rinse them using cold water and immerse them in GERMOCID ENZIMATICO solution for 10 minutes. In the presence of dry residue or stubborn dirt, extend the immersion to 1 hour and possibly use a toothbrush.

Composition: non-ionic surfactants 15-30%; anionic surfactants 5-15%; enzymes; perfume (limonene, citral), preservative (benzisothiazolinone).

Warnings: H318 Causes serious eye damage. H315 Causes skin irritation. EUH208 Contains: CITRAL May cause an allergic reaction. P264 Wash hands thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305-P351-P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a POISON CENTER/doctor. Contains: Polymer alkoxylated fatty alcohol. For professional use. Store in a cool, dry place and away from heat sources. After opening, if properly stored, the product retains the validity stated on the label.

Detergent in accordance with EC Regulation n° 648/2004

Validity: 3 years

Lot and expiry date: see packaging

EC MEDICAL DEVICE

Manufacturer:

GERMO S.p.A. - Via Giotto 19/21 - 20032 Cormano (MI) Tel. 02 663 01 938

www.germo.it - info@germodis.com

MADE IN ITALY