GIMA

PROFESSIONAL MEDICAL PRODUCTS

LACCI EMOSTATICI MONOUSO
DISPOSABLE TOURNIQUETS
GARROTS JETABLES
TORNIQUETES DESECHABLES
GARROTES DESCARTÁVEIS
EINWEG-STAU BINDE
JEDNORAZOWE OPASKI UCISKOWE
STASSLANGAR FÖR ENGÅNGSBRUK
ΤΟΥΡΝΙΚΕ ΜΙΑΣ ΧΡΗΣΗΣ
EGYSZER HASZNÁLATOS ÉRSZORÍTÓK
ÛHEKORDSEÐ ŽGUTID
TOURNIQUETS VOOR EENMALIG GEBRUIK
ÛHEKORDSEÐ ŽGUTID

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REF JS-TQ001 (GIMA 37510-37511-37512)

Ningbo Jenius Polymer Sci-Tech Co. Ltd
 N° 2112, Ningheng South Road, Hengxi Town, Yinzhou District,
 Ningbo, Zhejiang, China (Mainland)
 Made In China

MedNet GmbH Borkstrasse 10,
 48163 Münster, Germany

Gima S.p.A.
 Via Marconi, 1 - 20060 Gessate (MI) Italy
 gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Model JS-TQ001, size 450mm*25mm*0.635mm

This product is made of medical polymer thermoplastic elastomer, milky white or color, long flat type,

The new design point is used to pack the box, which is flexible and can be extracted continuously.

This product is suitable for medical institutions in routine treatment and treatment of infusion, blood transfusion, blood transfusion, hemostatic one-time use;

Emergency hemostasis in case of body bleeding or bleeding from field snake and insect bite.

Contraindications, matters needing attention, warnings and reminders:

1. This product is disposable and shall be destroyed after use;
2. If the package is damaged, it is strictly prohibited to use.
3. After use, please cover the suction cover to avoid dust and bacterial contamination;
4. The production date is shown on the side of the packaging box or the certificate of qualification;

Warnings and reminder:

A tourniquet can block the flow of blood, and tying it too long can seriously damage tissue, even lead to limb necrosis.

A tourniquet should only be used to bind the limbs, not the head, neck, or trunk.

Don't cover it with anything else, and don't cover the tourniquet tied to the limb.

If you have to leave the patient alone in the camp, over his head, with a pen or lipstick, write down the timing of the bundle and the location of the kit.

Check blood circulation:

After the tourniquet is bandaged, the toe or finger should be checked frequently to see if there is a dark purple spot at the end of the tourniquet, and whether the body temperature at the end of the tourniquet has dropped.

If there is such a phenomenon, should loosen the cloth band, otherwise for a long time, will cause tissue necrosis.

For the same reason, a tourniquet should not be used for long periods of time to tie a limb, unless it is necessary to do so during an arterial ligation. In general, direct pressure should be tried first, followed by pressure points until direct pressure is maintained on the wound.

If the wound is at the end of the limb, elevate it.

Storage conditions:

the packed tourniquet should be stored in the relative humidity of no more than 80%, without corrosive gas and well ventilated

the packaged tourniquet shall be valid for 3 years from the date of production under the conditions of storage, transportation and storage

Symbols:

	Keep in a cool, dry place
	Manufacturer
	Lot number
	Medical Device compliant with Regulation (EU) 2017/745
	Medical Device
	Keep away from sunlight
	Date of manufacture
	Product code
	Consult instructions for use
	Authorized representative in the European community
	Disposable device, do not re-use
	Expiration date
	Imported by

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies