



GIMA

PROFESSIONAL MEDICAL PRODUCTS

TIRALATTE MANUALE MANUAL BREAST PUMP TIRE-LAIT MANUEL MANUAL DEL SACALECHES BOMBA TIRA-LEITE MANUAL MANUELLE MUTTERMILCHPUMPE

REF 41700/LD - 101




JOYTECH HEALTHCARE CO.LTD.
No.365, Wuzhou Road,
Yuhang Economic Development Zone,
Hangzhou city, 311100 Zhejiang, China
Made in China



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537, Hamburg, Germany

Importato da / Imported by / Importé par /
Importado por / Importado por / Importiert von:
Gima S.p.A. Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Index

Product Description	14
Safety Notice.....	14
Unit Illustration	15
Installation Instructions	16
Product Assembly	16
Pumping	19
Size of Breastshield.....	20
Cleaning	21
Trouble Shooting	22
Waste Treatment	22
Maintenance	22
Specifications	23
Explanation of symbols	23

Product Description

Thank you for purchase and use LD-101 manual breast pump. The manual breast pump is used to express and collect breast milk from lactating mothers.

Please read this manual carefully before using stimulation and expression mode.

This is a single use product that should not be shared between moms.

Intended use

1. This breast pump is for household use only. As the pump is compact and discrete to use, you can take it with you anywhere, allowing you to express milk at your own convenience and maintain your milk supply.
2. If your breasts are engorged (painful or swollen), you can express a small amount of milk before or between feeds to ease the pain and to help your baby latch on more easily.
3. If you are separated from your baby and wish to continue to breastfeed when you are reunited, you should express your milk regularly to stimulate your milk supply.

Safety Notice



Caution

1. Use the product only for its intended use as described in this manual.
2. Please check that all the parts of breast pump is complete before using.
3. Do not use attachments not recommended by the manufacturer.
4. Never use the breast pump while you are pregnant, as pumping can induce labor.
5. Inspect all appropriate pump components before each use.
6. Clean and sanitize all parts that come in contact with your breast and breastmilk prior to first use.
7. Wash all parts that come in contact with your breast and breastmilk after every use.
8. For hygienic reasons, this product is intended for use by a single user.
9. Do not allow children or pets to play with accessories.

Unit Illustration

- 1. 1 × Elastic diaphragm
- 2. 1 × Stem
- 3. 1 × Handle
- 4. 1 × Seal seat
- 5. 4 × Seal hoops
- 6. 1 × Shock pad

- 7. 1 × Breastshield
- 8. 1 × Pump body
- 9. 2 × Valves
- 10. 4 × White membranes
- 11. 2 × Breastmilk bottles
- 12. 1 × Bottle stand

- 13. 1 × Dome cap
- 14. 1 × Nipple
- 15. 1 × Screw ring
- 16. 2 × Lids

Installation Instructions

Product Assembly


Note: Make sure you have cleaned and optionally disinfected the appropriate parts of the breast pump.

After cleaning, follow these steps to assemble your collection units:

Step 1:


Insert the breastshield into the funnel-shaped section of the pump body. Two breastshield sizes are available: 24mm and 27mm.

You can find the size on the breastshield itself.




Step 2:


Please install the white membrane on valve smoothly.
Make sure the white membrane won't curl up.

**Step 3:**

Insert the valve into the pump body from underneath. Push the valve in as far as possible.

**Step 4:**

Screw the pump body clockwise onto the bottle until it is securely fixed.



Step 5:

Insert the stem into the hole of the elastic diaphragm. Push the stem in as far as possible. Place the diaphragm with stem into the seal seat, make sure the diaphragm edge and seal seat edge close fitting. If there's no "O" shaped sealing ring on the plug at the bottom of the sealing seat, please put in on.




Figure 5-1




Figure 5-2




Figure 5-3


Step 6:

Attach the handle on diaphragm with stem by hooking the U- shape hole in the groove of the stem. Push down the handle onto the seal seat until its well placed.




Step 7:

Insert the seal seat in the hole of the connector and screw the seal clockwise until it is securely fixed.

**Pumping**

Relax in a comfortable chair (you may wish to use cushions to support your back).


Press the assembled pump body against your breast. Make sure that your nipple is centered. Hold the breastshield onto your breast with your thumb and index finger. Support your breast with the palm of your hand. Gently start to press down the handle.



Note: Do not worry if your milk does not flow immediately. Relax and continue pumping. Repositioning the breast pump on your breast from time to time can help stimulate milk flow. Try to express at another time during the day. If the process becomes very uncomfortable or painful, stop using the pump and consult your breastfeeding advisor.

Close the bottle with a lid after pumping. Disassemble and clean the parts that contact breast and milk.




 Do not continue pumping for more than 5 minutes at a time if you do not succeed in expressing any milk.


Use bottle holder to prevent the bottle from tipping over. Do not fill the bottle too full to prevent overfilling and spillage.

Size of Breastshield

Make sure that your nipple fits correctly in the breast shield. Not too tight, it is important that the nipple can move freely, while expressing breast milk.

The breastshield should not be too large, because this can cause pain or less effective expression. Two breast shield sizes are available: 24mm and 27mm. There's one size breast shield come with the product.



 If there is any pain or discomfort during pumping, please consider to choose a larger or smaller size.

Cleaning

Caution

1. Clean and sanitize the parts below before using your pump for the first time.
2. Only use drinking-quality tap or bottled water for cleaning.
3. Take apart and wash all parts that come in contact with the breast and breastmilk immediately before and after use to avoid dry up of milk residues and to prevent growth of bacteria.
4. Please disassemble all parts which require clean to make sure cleaning thoroughly.
5. Please place parts on a clean environment to avoid being polluted.

Cleaning before use

1. Take apart and wash/sanitize all parts that come in contact with the breast and breastmilk.
2. Put all separated parts into pot. Fill in the pot with enough drinking-quality tap or bottled water to cover all parts.
3. Bring the water to boil. Place parts in boiling water for 5 minutes.
4. Allow water to cool and gently remove the parts from the water. Be careful and don't scald your skin.
5. Place the parts neatly on a clean paper towel or in a clean drying rack and allow them to air dry. Avoid using cloth towels to dry the parts because they can carry germs and bacteria that are harmful to your baby.

- Breastshiled
- Values
- White membranes
- Nipple
- Breastmilk bottle and lid



Cleaning after use

1. Take apart and wash all parts that come in contact with the breast and breastmilk.
2. Rinse in cool water all separated parts that came in contact with breast and breastmilk in order to remove breastmilk residue.
3. Place the parts neatly on a clean paper towel or in a clean drying rack and allow them to air dry.

IMPORTANT

The flanges and milk storage containers can also be washed on upper rack of a standard, household dishwasher. DO NOT use solvents or abrasives. Allow to air dry on a clean surface. DO NOT put parts in Steam Bottle Sterilizers or Microwave Sterilizer Bags.

Trouble Shooting

Abnormal phenomenon	Cause analysis	Solutions
Low or No Suction	Connection points is loose	Inspect all connection points to ensure attachment is secure
	The valve is chipped or cracked or there's holes or tears on membrane is	Replace valve and/or membrane prior to pumping
	Breastshield size is not fit	Replace a appropriate size breastshield

Waste Treatment

When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center. Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product is free of hazardous materials.

Maintenance

1. Avoid dropping, slamming, or throwing the unit.



2. Avoid extreme temperatures. Do not expose unit directly under sunshine.














3. Do not use the unit for other purpose.

Specifications

Product Description	Manual Breast Pump	
Model	LD-101	
Vacuum (approx.)	Approx. -250mmHg	
Breastshield size	24mm	
Unit Weight	approx. 159g(5.61oz)	
Unit Dimensions	205 x 186 x 90mm (L x W x H)	
Operating Environment	Temperature	5°C~40°C (41°F~104°F)
	Humidity	15~93%RH
	Pressure	700hPa~1060hPa
Storage Environment	Temperature	-20°C~55°C (-4°F~131°F)
	Humidity	15~93%RH

Specifications are subject to change without notice.

Explanation of symbols

	Caution: read instructions (warnings) carefully		Authorized representative in the European community		Medical Device complies with Directive 93/42/EEC
	Manufacturer		Keep away from sunlight		Serial number
	Product code		Lot number		Consult instructions for use
	Temperature limit		Humidity limit		

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.