TUTORE PER CAVIGLIA
ANKLE BRACE
ATELLE POUR LA CHEVILLE
ORTESIS DE TOBILLO
IMOBILIZADOR DE TORNOZELO
KNÖCHELORTHESE
ΝΑΡΘΗΚΑΣ ΑΣΤΡΑΓΑΛΟΥ
ORTEZA STAWU SKOKOWEGO
ΚΟΤΝÍΚΟVÁ ORTÉZA
FOTLEDSBANDAGE
NILKKATUKI
OPORNICA ZA GLEŽENJ
ČLENKOVÁ ORTÉZA
ORTEŽĀ PENTRU GLEZNA
ENKELSPALK
STEZNIK ZA NOŽNI ZGLOB
ΒΟΚΑΜΕΡΕVÍΤΌ
ANKELBØJLE
ΟΡΤΕΖΑ ΖΑ ΓΛΕΖΕΗ
ČIURNOS JTVARAS
POTĪTES BALSTS
ΠΑΗΚΛΥΟΥΡΤΟΟΣ

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REF AN013 (GIMA 42050 – 42051 – 42052 - 42053)

XIAMEN ORTOCARE SPORTS HEALTH CO., LTD.
Unit 218-2, New Times Building, No.1 Huli Avenue, Xiamen, China
(Fujian) Pilot Free Trade Zone
MADE IN CHINA

SUNGO Cert GmbH
Harffstr. 47,40591 Düsseldorf, Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

and attach to the corresponding hook & loop closure on the outer side of the ankle.

5. Next, extend the inner strap from the dorsum of the foot and pass it under the outer side of the foot. Finally, apply an upward tension and attach to the corresponding hook & loop closure on the inner side of the ankle.

6. Fasten the top strap around the shin with comfortable pressure.

Cleaning Methods / Maintenance:

- Close all hook and loop closures before washing.
- Hand wash in cold water with mild soap, gentle cycle..
- Do not dry clean or tumble dry and iron; do not use bleach.
- Rinse well and air dry.
- Do not expose this brace to any heat source.

Precaution and Warning:

Please consult with your physician if any allergic symptom occurs. Do not use more than 24 hours of continuous.

Storage:

Store the product in a cool dry place away from sunlight and heat and preferably in the original packaging.

Shelf Life:

5 years

Disposal:

Please dispose the product after use to comply with local regulation.

Notice to the user and/or patient:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

ENGLISH

Name:

Ankle Brace

Intended use:

Correct postural and medical immobilization Helps immobilize and provide support on ankle joint.

Usage:

1. Detach all straps and loosen laces.
2. Place foot into brace and ensure heel is in aperture.
3. Fasten the laces by hook and loop closures.
4. Extend the outer strap from the dorsum of the foot and pass it under the inner side of the foot. Finally, apply an upward tension

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