

SPLINT PER DITA
FINGER SPLINT
ATELLE POUR LE DOIGT
FÉRULA DE DEDO
TALA DE DEDO
FINGERSCHIENE
ΝΑΡΘΗΚΑΣ ΔΑΚΤΥΛΩΝ
SZYNA NA PALEC
PRSTOVÁ DLAHA
FINGERSKENA
SORMITUKI
PRSTNA OPORNICA
ATELĂ PENTRU DEGET
PRSTOVÁ DLAHA
VINGERSPALK
UDLAGA ZA PRST
UJJSÍN
FINGERSKINNE
ШИНА ЗА ПРЪСТИ
PIRŠTO JTVARAS
PIRKSTU ŠINA
SÖRMELAHAS

جبيرة الاصبع

REF FS004 (GIMA 42102 – 42103 – 42014 – 42105 – 42106)



XIAMEN ORTOCARE SPORTS HEALTH CO., LTD.
Unit 218-2, New Times Building, No.1 Huli Avenue, Xiamen, China
(Fujian) Pilot Free Trade Zone
MADE IN CHINA



SUNGO Cert GmbH
Harffstr. 47,40591 Düsseldorf, Germany



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Product name:

Finger Splint

**Indications :**

Protects fingertip and nail bed injuries. Protects against Mallet finger as well.

CONTRAINDICATIONS:

Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity. The determination of when to apply the brace rests strictly at the discretion of the treating physician.

This product was designed to compliment the variety of medical treatments common to the above afflictions. The frequency and duration of use should be determined by your prescribing health care professional.



Features & Benefits:

1. Used to immobilize and hold the distal finger in hyperextension
2. Supports finger in desired alignment
3. Protects finger from further injury


Indications:

1. Mallet finger
2. Distal finger fracture, dislocation and sprain

assorted Finger Splints
No. 1, 2, 3, 4, 5, 5,5, 6



Instructions for use:




Notice to the user and/or patient:


Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established




GIMA WARRANTY TERMS








The Gima 12-month standard B2B warranty applies.

Indice dei simboli - Symbol index - Index des symboles - Índice de símbolos- índice de símbolo - Symbolindex - Ευρετήριο συμβόλων - Indeks symboli - Index symbolů - Symbol index - Symboli-indeksi - Indeks simbolov - Index symbolov - Index de simbol - Symbol index - Indeks simbola - Szimbólum index - Symbolindeks - Индекс на символа - Simbolių rodyklė - Simbolu rādītājs - Sūmbolite indeks

	<p> IT - Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB - Caution: read instructions (warnings) carefully FR - Attention: lisez attentivement les instructions (avertissements) ES - Precaución: lea las instrucciones (advertencias) cuidadosamente PT - Cuidado: leia as instruções (avisos) cuidadosamente DE - Achtung: Anweisungen (Warnungen) sorgfältig lesen GR - Προσοχή: διαβάστε προσεκτικά τις οδηγίες (ενοτάσεις) PL - Ostrzeżenie - Zobacz instrukcję obsługi CZ - Pozor: Pečlivě si přečtěte a dodržujte pokyny (varování) k použití SE - Varsamhet: läs anvisningarna (varningar) noga FI - Huomio: Lue käyttöohjeet (varoitukset) ja noudata niitä huolellisesti SI - Pozor: Preberite in skrbno sledite navodilom (opozorilom) za uporabo SK - Pozor: Pozorne si prečítajte a dodržiavajte pokyny na použitie (výstrahy) RO - Atenție: Citiți și respectați cu atenție instrucțiunile (avertismentele) de utilizare NL - Opgelet: Lees en volg aandachtig de gebruiksaanwijzing (waarschuwingen) HR - Pozor: Pročitajte i pažljivo slijedite upute (upozorenja) za upotrebu HU - Figyelem: Figyelmesen olvassa el és kövesse a használati utasításokat (figyelmeztetéseket) DK - Forsigtig: Læs instruktioner (advarsler) omhyggeligt BG - Forsigtig: Læs instruktioner (advarsler) omhyggeligt LT - Dėmesio: perskaitykite ir atidžiai laikykitės naudojimo instrukcijų (įspėjimų). LV - Uzmanību: Ilasiet un uzmanīgi ievērojiet lietošanas instrukcijas (brīdinājumus) EE - Tähelepanu! Lugege kasutusjuhised (hoiatused) läbi ja järgige neid hoolikalt AA - ال حذوق اللات ع ليم اللات ح ذري راجع ن اية </p>
---	--

	<p>IT - Leggere le istruzioni per l'uso GB - Consult instructions for use FR - Consulter les instructions d'utilisation ES - Consultar las instrucciones de uso PT - Consulte as instruções de uso DE - Gebrauchsanweisung beachten GR - Διαβάστε προσεκτικά τις οδηγίες χρήσης PL - Przeczytaj instrukcję użytkowania CZ - Přečtěte si návod k použití SE - Läs bruksanvisningen FI - Lue käyttöohjeet SI - Preberite navodila za uporabo SK - Prečítajte si návod na použitie RO - Citiți instrucțiunile de utilizare NL - Lees de gebruiksaanwijzing HR - Pročitajte upute za uporabu HU - Olvassa el a használati utasításokat DK - Se brugsvejledningen BG - Se brugsvejledningen LT - Perskaitykite naudojimo instrukcijas LV - Izlasiet lietošanas instrukcijas EE - Lugege kasutusjuhendit</p> <p>AA - اقرأ بدقة وحرص تعليمات الاستخدام</p>		
<table border="1" data-bbox="69 922 140 954"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	<p>IT - Rappresentante autorizzato nella Comunità europea GB - Authorized representative in the European community FR - Représentant autorisé dans la Communauté européenne ES - Representante autorizado en la Comunidad Europea PT - Representante autorizado na União Europeia DE - Autorisierter Vertreter in der EG GR - Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση PL - Upoważniony przedstawiciel we Wspólnocie Europejskiej CZ - Zplnomocnený zástupce v Evropském společenství SE - Auktoriserad representant i Europeiska gemenskapen FI - Valtuutettu edustaja Euroopan yhteisössä SI - Pooblaščenī zastopnik za Evropsko skupnost SK - Splnomocnený zástupca v Európskom spoločenstve RO - Représentant autorizat pe teritoriul Comunității Europene NL - Bevoegde vertegenwoordiger in de Europese Gemeenschap HR - Ovlašteni predstavnik u Europskoj zajednici HU - Meghatalmazott képviselő az Európai Közösségekben DK - Autoriseret repræsentant i det Europæiske Fællesskab BG - Autoriseret repræsentant i det Europæiske Fællesskab LT - Įgaliotasis atstovas Europos bendrijoje LV - Pilnvarotais pārstāvis Eiropas Kopienā EE - Volitatud esindaja Euroopa Ühenduses</p> <p>AA - ممثل معتمد في الاتحاد الأوروبي</p>
EC	REP		

	<p>IT - Dispositivo medico conforme al regolamento (UE) 2017/745 GB - Medical Device compliant with Regulation (EU) 2017/745 FR - Dispositif médical conforme au règlement (UE) 2017/745 ES - Producto sanitario conforme con el reglamento (UE) 2017/745 PT - Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 DE - Medizinprodukt im Sinne der Verordnung (EU) 2017/745 GR - Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 PL - Wyrób medyczny zgodny z Rozporządzenie (UE) 2017/745 CZ - Zdravotnický prostředek v souladu s nařízením (EU) č. 2017/745 SE - Den medicintekniska produkten överensstämmer med förordning (EU) 2017/745 FI - Lääkinnällinen laite, joka vastaa asetusta (EU) 2017/745 SI - Medicinski pripomoček, skladen z uredbó (EU) 2017/745 SK - Zdravotnícka pomôcka v súlade s nariadením (EÚ) 2017/745 RO - Dispozitiv medical conform regulamentului (UE) 2017/745 NL - Medisch hulpmiddel in overeenstemming met verordening (EU) 2017/745 HR - Medicinski proizvod skladan propisu (EU) 2017/745 HU - A 2017/745/EU rendeltetnek megfelelő orvostechnikai eszköz DK - Medicinsk udstyr i overensstemmelse med forordning (EU) 2017/745 BG - Medicinsk udstyr i overensstemmelse med forordning (EU) 2017/745 LT - Medicinos prietaisais, atitinkantis reglamentą (ES) 2017/745 LV - Medicīniska ierīce, kas atbilst Regulai (ES) 2017/745 EE - Määrusele (EL) 2017/745 vastav meditsiiniiseade (UE) 2017/745</p> <p>AA - جهاز طبي يتوافق مع التوجيه</p>
	<p>IT - Data di scadenza GB - Expiration date FR - Date d'échéance ES - Fecha de caducidad PT - Data de validade DE - Ablaufdatum GR - Ημερομηνία λήξεως PL - Data ważności CZ - Datum ukončení platnosti SE - Utgångsdatum FI - Viimeinen voimassaolopäivä SI - Rok uporabnosti SK - Dátum expirácie RO - Valabil până la data de NL - Vervaldatum HR - Datum isteka HU - Lejárati dátum DK - Udløbsdato BG - Udløbsdato LT - Galiojimo laikas LV - Derīguma termiņš EE - Aegumiskuupäev</p> <p>AA - تاريخ انتهاء الصلاحية</p>
	<p>IT - Codice prodotto GB - Product code FR - Code produit ES - Código producto PT - Código produto DE - Erzeugniscode GR - Κωδικός προϊόντος PL - Numer katalogowy CZ - Kód výrobku SE - Produktkod FI - Tuotekoodi SI - Koda izdelka SK - Kód výrobku RO - Cod produs NL - Productcode HR - Šifra proizvoda HU - Termékkód DK - Produktkode BG - Produktkode LT - Prekės kodas LV - Produkta kods EE - Toote kood</p> <p>AA - كود المنتج</p>

	<p>IT - Dispositivo medico GB - Medical Device FR - Dispositif médical ES - Producto sanitario PT - Dispositivo médico DE - Medizinprodukt GR - Ιατροτεχνολογικό προϊόν PL - Wyrób medyczny CZ - Zdravotnický prostředek SE - Medicinteknisk produkt FI - Lääkinnällinen laite SI - Medicinski pripomoček SK - Zdravotnícka pomôcka RO - Dispozitiv medical NL - Medisch hulpmiddel HR - Medicinski uređaj HU - Orvostechnikai eszköz DK - Medicinsk udstyr BG - Medicinsk udstyr LT - Medicininis prietaisai LV - Medicīniskā ierīce EE - Meditsiiniline seade</p> <p>جهاز طبي - AA</p>	 <p>IT - Importato da GB - Imported by FR - Importé par ES - Importado por PT - Importado por DE - Eingeführt von GR - Εισαγωγή από PL - Importowane przez CZ - Dovezeno uživatelem SE - Importerad av FI - Tuoja SI - Uvozil SK - Dovážal RO - Importat de NL - Geïmporteerd door HR - Uvezeno od strane HU - Importáltá DK - Importeret af BG - Importeret af LT - Importavo LV - Importēja EE - Importija</p> <p>مستورد عن طريق - AA</p>
	<p>IT - Numero di lotto GB - Lot number FR - Numéro de lot ES - Número de lote PT - Número de lote DE - Chargennummer GR - Αριθμός παρτίδας PL - Kod partii CZ - Číslo šarže SE - Satsnummer FI - Eränumero SI - Številka partije SK - Číslo šarže RO - Număr de lot NL - Partijnummer HR - Broj serije HU - Tételszám DK - Batchnummer BG - Batchnummer LT - Partijos numeris LV - Partijas numurs EE - Partii number</p> <p>رقم الدفعة - AA</p>	 <p>IT - Identificatore univoco del dispositivo GB - Unique device identifier FR - Identifiant unique de l'appareil ES - Identificador de dispositivo único PT - Identificador exclusivo do dispositivo DE - Unique Device Identifier (Eindeutige Kennung des Geräts) GR - Μοναδικό αναγνωριστικό συσκευής PL - Unikalny identyfikator urządzenia CZ - Jediněčný identifikátor zařízení SE - Unik identifierare för enheten FI - Laitteen yksilöllinen tunniste SI - Enolični identifikator naprave SK - Jediněčný identifikátor zariadenia RO - Identificatorul unic al dispozitivului NL - Unieke identificatie van het apparaat HR - Jedinstveni identifikator uređaja HU - Az eszköz egyedi azonosítója DK - Unik identifikator for enheden BG - Уникален идентификатор на устройството LT - Unikalus įrenginio identifikatorius LV - Unikāls ierīces identifikators EE - Seadme kordumatu identifikaator</p> <p>معرف فريد للجهاز - AA</p>
	<p>IT - Non sterile GB - Non-sterile FR - Pas stérile ES - No estéril PT - Não estéril DE - nicht steril GR - όχι αποστειρωμένο PL - Nie sterylne CZ - Nesterilní SE - Ej steril FI - Ei-steriili SI - Ni sterilno SK - Nesterilné RO - Nesteril NL - Niet steriel HR - Nije sterilno HU - Nem steril DK - Ikke-steril BG - Нестерилен LT - Ne sterilus LV - Nav sterilis EE - Mittesteriilne</p> <p>ليس معقم - AA</p>	 <p>IT - Conservare in luogo fresco ed asciutto GB - Keep in a cool, dry place FR - À conserver dans un endroit frais et sec ES - Conservar en un lugar fresco y seco PT - Armazenar em local fresco e seco DE - An einem kühlen und trockenen Ort lagern GR - Διατηρείται σε δροσερό και στεγνό περιβάλλον PL - Przechowywać w suchym miejscu CZ - Skladujte na větraném a suchém místě SE - Förvara på svalt och torrt ställe FI - Säilytä kuivassa ja viileässä SI - Hraniti na suhem in hladnem mestu SK - Skladujte na chladnom a suchom mieste RO - A se păstra într-un loc răcoros și uscat NL - Koel en droog opslaan HR - Čuvati na hladnom i suhom mjestu HU - Száraz, hűvös helyen tárolandó DK - Opbevares køligt og tørt BG - Opbevares køligt og tørt LT - Laikyti vėsioje ir sausoje vietoje LV - Uzglabāt vēsā, sausā vietā EE - Hoida jahedas ja kuivas kohas</p> <p>يحفظ في مكان بارد وجاف - AA</p>
	<p>IT - Data di fabbricazione GB - Date of manufacture FR - Date de fabrication ES - Fecha de fabricación PT - Data de fabrico DE - Herstellungsdatum GR - Ημερομηνία παραγωγής PL - Data produkcji CZ - Datum výroby SE - Tillverkningsdatum FI - Valmistuspäivämäärä SI - Datum proizvodnje SK - Dátum výroby RO - Data fabricației NL - Productiedatum HR - Datum proizvodnje HU - Gyártás dátuma DK - Fabrikationsdato BG - Fabrikationsdato LT - Pagaminimo data LV - Izgatavošanas datums EE - Valmistamise kuupäev</p> <p>تاريخ التصنيع - AA</p>	<p>IT - Fabbricante GB - Manufacturer FR - Fabricant ES - Fabricante PT - Fabricante DE - Hersteller GR - Παραγωγός PL - Producent CZ - Výrobce SE - Tillverkare FI - Valmistaja SI - Proizvajalec SK - Výrobca RO - Producător NL - Fabrikant HR - Proizvođač HU - Gyártó DK - Fabrikant BG - Fabrikant LT - Gamintojas LV - Ražotājs EE - Tootja</p> <p>الشركة المصنعة - AA</p>



IT - Conservare al riparo dalla luce solare **GB**
 - Keep away from sunlight **FR** - À conserver à
 l'abri de la lumière du soleil **ES** - Conservar al
 amparo de la luz solar **PT** - Guardar ao abrigo
 da luz solar **DE** - Vor Sonneneinstrahlung
 geschützt lagern **GR** - Κρατήστε το μακριά
 από ηλιακή ακτινοβολία **PL** - Przechowywać
 z dala od światła słonecznego **CZ** - Skladujte
 mimo sluneční světlo **SE** - Skyddas från
 solljus **FI** - Säilytä auringonvalolta suojaassa
SI - Hraniti zaščiteno pred sončno svetlobo
SK - Skladujte mimo slnečného svetla **RO**
 - A se păstra ferit de razele soarelui **NL**
 - Afgeschermd van zonlicht opslaan **HR** -
 Čuvati zaštićeno od sunčeve svjetlosti **HU** -
 Napfénytől védve tárolandó **DK** - Må ikke
 udsættes for sollys **BG** - Må ikke udsættes for
 sollys **LT** - Saugoti nuo saulės spindulių **LV** -
 Uzglabāt prom no saules gaismas **EE** - Hoida
 eemal päikesevalgusest
AA - يحفظ بعيدًا عن ضوء الشمس