Name: Bath Chair

Model: KDB-793, KDB-794, KDB-795, KDB-796, KDB-797, KDB-798, KDB-799, KDB-799B, KDB-787, KDB-788

Intended Use: The Bath Chair is used as bathing tools for people with disabilities, patients and the elderly who are difficult to maintain balance while showering.

Operator: The device can be operated by the user self, the nursing staff or the caregiver.

Patient populations: People with disabilities, patients and the elderly.

Contraindications: Supervise children, when they use the bath chair.

Instruction for Use

1. Before use, adjust the bath seat to a suitable height, so that four feet are on the ground smoothly.
2. Check whether the backrest is installed firmly and whether there is a swing during use. Please sit in the middle of the shower chair during use.
3. The height of the bath seat can be adjusted. When adjusting the height, first push down the marbles, and then pull it down. Just adjust to the height suitable for your own use.
4. The chair should not be used to push the user while seated.
5. Before use, please make sure that the adjusting screw has been tightened and the button screw has been ejected from the slot
6. Wash with clean water and air dry after use
7. Do not use this product on uneven pavement
8. keep away from fire.

Maintenance

1. Routinely clean and disinfect once a day by medical staff or users, and thoroughly disinfect the chair when changing patients.
2. Do not contact with strong acids and alkalis to prevent corrosion.

Cautions

Please do not support on one side when getting up or sitting down, to prevent the Bath Chair from being overturned and injured.

Use it with care on uneven ground. Please clean and disinfect the bath chair in time.

This product is safe to use with a load of 125kg. Be careful when using it in excess of this weight, or it can be customized by the company.

Storage:The product should be stored in normal temperature and humidity and keep dry.

Disposal: Disposal restrictions depend upon the contamination incurred during use and are subject to national or local legislation.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is locate.

Indice dei simboli - Index of symbols - Index des symboles - Índice de símbolos - índice de símbolo - Symbolindex - Ευρετήριο συμβόλων - Indeks symboli - Index symbolů - Symbol index - Symboli-indeksi - Indeks simbolov - Index symbolov - Index de simbol - Symbol index - Indeks simbola - Szimbólum index - Symbolindex - Индекс на символа - Simbolių rodyklė - Simbolu rādītājs

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GIMA 42902

Zhongshan Kangdebao (KDB) Rehabilitation Equipment Co., Ltd.
Card 1, No. 5, 24# Longcheng Road, Dongsheng Town, Zhongshan City,
Guangdon Province, China

REF KDB-797

SUNGO Europe B.V.
Fascianatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com