X36-1 LETTINO DA VISITA ELETTRICO
X36-1 ELECTRIC EXAMINATION TABLE
X36-1 DIVAN D'EXAMEN ÉLECTRIQUE
X36-1 CAMILLA DE EXPLORACIÓN ELÉCTRICA
X36-1 ELEKTRISCHE UNTERSUCHUNGSLIEGE
X36-1 MESA DE EXAME ELÉTRICA
X36-1 ΗΛΕΚΤΡΙΚΟ ΕΞΕΤΑΣΤΙΚΟ ΚΡΕΒΑΤΙ
X36-1 ЕЛЕКТРИЧЕСКА МАСА ЗА ПРЕГЛЕД
X36-1 ELEKTRICKÝ VYŠETŘOVACÍ STŮL
X36-1 ELEKTRISK UNDERSØGELSESELEJE
X36-1 ELEKTRILINE LÄBIVAATUSLAUD

X36-1 SÄHKÖINEN TUTKIMUSPÖYTÄ
X36-1 ELEKTROMOS VIZSGÁLÓASZTAL
X36-1 ELEKTRINĚ APŽIŮROS KUŠETĚ
X36-1 ELEKTRISKAIS IZMEKLĒŠANAS GALDS
X36-1 ELEKTRISCHE ONDERZOEKSTAFEL
X36-1 ELEKTRYCZNY STÓŁ DO BADAŃ
X36-1 CANAPEA ELECTRICĂ DE EXAMINARE
X36-1 ELEKTRICKÝ VYŠETROVACÍ STŮL
X36-1 ELEKTRIČNA PREGLEDNA MIZA
X36-1 ELEKTRISKT UNDERSÖKNINGSBORD

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Jiangsu Saikang Medical Equipment Co., Ltd
 No. 35 Lehong Road, Modern Agriculture Demonstration Park,
 Zhangjiagang City, Jiangsu Province, China
 Made in China


X36-1

SUNGO Europe B.V.
 Fascinat Boulevard 522, Unit 1.7,
 2909VA Capelleaan den IJssel, The Netherlands

Gima S.p.A.
 Via Marconi, 1 - 20060 Gessate (MI) Italy
 gima@gimaitaly.com - export@gimaitaly.com
 www.gimaitaly.com

ENGLISH

INTENDED USE

An examination table is integral to patient comfort while suffering from injury, trauma or shock.

WARNING

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.


1. Use method

Table Size	Length	mm	1960±10
	Width		680±10
	The height of the table above the ground		470-910 ±10
Weight	Maximum patient weight	kg	180
Weight	Safe working load	kg	200

1.1 Height Adjustment

The whole examination table can be controlled by hand operator

2. Electric examination table-maintenance

1. In order to use the examination table safely, it is necessary to perform regular safety inspections on the examination table. It is recommended to check it every six months to ensure that the connections are not loose and all operations are normal.
2. When the examination table frame ages and it reaches a certain servicing time, which the metal part of the examination table and plastic parts of the stainless steel basin frame can be recycle.
3. Avoid scratching the panel with sharp-angle appliances or knives during use, and clean it frequently to keep it clean and dry.
4. If the panel is accidentally stained with stains, it is recommended to clean it. Do not use alkaline or corrosive chemicals to clean the examination table, which will cause the stainless steel surface to rust.

3. Examination table-cleaning

1. When cleaning, please wring out the cloth soaked in neutral detergent diluted with water and wipe it, then wring out the cloth soaked in water to wipe off the remaining detergent ingredients, and finally use a dry wipe dry with a cloth.
2. Do not use volatile substances (thinner, volatile agent, gasoline, etc.), which may cause chemical reactions and damage the examination table.
3. When using a disinfectant to clean, be sure to dilute it according to its specified concentration before use. Depending on the composition of the disinfectant, it may corrode metal parts, resin parts, etc., causing undesirable phenomena such as discoloration and deformation. Therefore, the content of the disinfectant is recommended as follows:
0.05 - 0.2% ammonium chloride disinfectant
0.05 - 0.2% Chlorinated phenyl disinfectant 0.05% diclofenac ethane solution 0.05 - 0.2% sodium hypochlorite disinfectant
Do not use smoking sterilizers, autoclaves, and do not use methyl (phenol) to clean the headboard and foot of the examination table. It may cause corrosion, discoloration and deterioration.
Note: When cleaning the examination table or changing the examination table clothes, please pay attention to the corners, edges and screws of the frame to prevent scratches.

4. Packaging, transport, storage

1. The packaging of examination table is carried out according to the contract or product standard.
2. Avoid rushing out, violent vibration, and protection from sunlight and rain during the transportation of the examination table.
3. The examination table should be stored:
a) Ambient temperature: -40°C ~ 70°C.
b) Relative humidity: ≤ 95%.
c) Atmospheric pressure: 500hPa ~ 1060 hPa.

5. After-sales service

1. Please keep the files that comes with the machine and invoices of this product properly, and you need to present these files when the company performs warranty and maintenance for the product.
2. If there is any problem in the process of using, please contact our company in time, so that our company can provide you with accurate and fast technical support and maintenance services in a timely manner.
3. From the date of sale, if the product is broken or damaged due to the correct installation and use according to the specification, the product will enjoy one-year free warranty and lifelong maintenance service with the "certificate" or invoice.
4. From the date of purchase, if it is indeed damaged or does not work normally due to quality problems within one year, the company will provide to the users a free repairing of the product.
5. Lifetime service from manufacturer: JIANGSU SAIKANG MEDICAL EQUIPMENT CO LTD.

EMC

- 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- 2) * Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) * Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity		
The V8v8c is intended for use in the following electromagnetic environment. The customer or user of the V8v8c should make sure it is used in such an environment		
Emission test	Compliance	Electromagnetic Environment Guidance

RF Emission to CISPR 11	Group 1	The bed uses RF energy for its internal function
RF Emission to CISPR 11	Class B	The bed is intended for use in all types of establishment including residential and similar user that are directly connected to a public supply network that also serves buildings used for residential purposes
Harmonics according to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker acc. to IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity		
The V8v8c is intended for use within healthcare facility environment. Users or customers must use bed in this environment		
Immunity test	IEC 60601-1-2 EN 60601-1-2 Test level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)..... IEC 61000-4-2 /EN 61000-4-2	Contact discharge: +8kV Air discharge: +2kV, +4kV, +8kV, +15kV	Contact discharge: +8kV Air discharge: +2kV, +4kV, +8kV, +15kV Floor should be made of wood and concrete or be tiled with ceramic tiles. If the floor is covered with synthetic flooring material, the relative air humidity must be at least 30% Can be used when higher ESD levels are present
Electrical fast transient/burst IEC 61000-4-4 /EN 61000-4-4	+2kV Power supply line +1kV Input/Output lines	+2kV Power supply line The quality of the supply voltage should be equivalent to that of a typical business or hospital environment.
Surge IEC 61000-4-5 /EN 61000-4-5	Line - line: +0.5kV, +1.0kV Line - ground: +0.5kV, +1.0kV, +2.0kV	Line - line: +0.5kV, +1.0kV Line - ground: +0.5kV, +1.0kV, +2.0kV The quality of the supply voltage should be equivalent to that of a typical business or hospital environment.
Power Frequency Magnetic Fields IEC 61000-4-8 /EN 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz Magnetic fields with a network frequency should be equivalent to those to be found in a typical business or hospital environment.
Voltage Dips IEC 61000-4- 11 /EN 61000-4-11	0% UT: 0.5 cycles Phase: 00, 450, 900, 1350, 1800, 2250, 2700, 3150	0% UT: 0.5 cycles Phase: 00, 450, 900, 1350, 1800, 2250,
	0% UT: 1 cycle 70% UT: 25/30 cycles Single phase: 00	0% UT: 1 cycle 70% UT: 25/30 cycles Single phase: 00 0% UT ; 300 cycle The quality of the supply voltage should be equivalent to that of a typical business or hospital environment.
	0% UT ; 300 cycle	If the person using the bed requires that the bed functions must continue despite any interruptions in the energy supply, it is recommended that the bed be connected to an uninterruptible electricity supply or a battery
Short interruptions IEC 61000-4-11 /EN 61000-4-11	0% UT: 250/300 cycles	0% UT: 250/300 cycles The quality of the supply voltage should be equivalent to that of a typical business or hospital environment. If the person using the bed requires that the bed functions must continue despite any interruptions in the energy supply, it is recommended that the bed be connected to an uninterruptible electricity supply or a battery
Note: UT is the AC Power voltage before applying in the test level		

Guidance and Manufacturer's Declaration - Electromagnetic Immunity							
The V8v8c is intended for use in an electromagnetic environment that used wireless RF transceiver. Electromagnetic interference can be prevented by maintaining the minimum distance between RF transceiver and V8v8c depend on the maximum output and frequency of the communication device as recommended below.							
Immunity test	IEC 60601-1-2 /EN 60601-1-2 Test level	Compliance level (V/m)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)

Enclosure port to RF wireless communication equipment IEC 61000-4-3 /EN 61000-4-3	27V/m 385MHz	27	380-390	TETRA 400			
	28V/m 450MHz	28	430-470	GMRS 460, FRS 460			
	9V/m 710MHz	9	704-787	LTE Band 13, 17			
	9V/m 745MHz						
	9V/m 780MHz						
	28V/m 810MHz	28	800-960	GSM 800/900, TETRA 800, iDen 820, CDMA 850 LTE Band 5			
	28V/m 870MHz						
	28V/m 930MHz						
	28V/m 1720MHz	28	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4,25; UMT S			
	28V/m 1845MHz						
	28V/m 1970MHz						
	28V/m 2450MHz	28	2400-2570	Bluetooth, WLAN, 802.11 b/g/n RFID 2450, LTE band 7			
	9V/m 5240MHz	9	5100-5800	WLAN 802.11 a/n			
	9V/m 5500MHz						
	9V/m 5785MHz						

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The V8v8c is intended for use in the following electromagnetic environment. The customer or user of the V8v8c should make sure it is used in such an environment.			
Immunity test	IEC 60601-1-2 /EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment - Guidance
RF Conductivity IEC 61000-4-6 EN 61000-4-6	3V 0.15MHz-80MHz	3V 0.15MHz-80MHz	Field intensity from the fixed RF transmitter determined by a field investigation of the magnetic field must be lower than the compliance level of the respective frequency band. Near devices with the following symbol, interference could occur. 
	6V in ISM bands between 0.15MHz- 80MHz 80%AM at 1kHz	6V in ISM bands between 0.15MHz- 80MHz 80%AM at 1kHz	
Radioactive RF electromagnetic field IEC 61000-4-3 /EN 61000-4-3	3V/M 80MHz-2.7GHz 80%AM at 1kHz	3V/M 80MHz-2.7GHz 80%AM at 1kHz	

a. The intensity of a magnetic field from fixed transmitters such as a radio (mobile/radio) telephone base station, land mobile radio, amateur radio, AM and FM radio broadcasting, and TV broadcasting cannot be theoretically predicted accurately. To determine the electromagnetic environment for a fixed RF transmitter, consider investigating the electromagnetic field in the location where the bed is used. If the intensity of the magnetic field in the location where the V8v8c is used exceeds the above RF compliance level, make sure to monitor if the V8v8c operates appropriately. If abnormal movement is found, take additional measures as needed, such as changing the direction or position of the V8v8c.

b. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

GIMA WARRANTY TERMS

THE GIMA 12-MONTH STANDARD B2B WARRANTY APPLIES.

Simboli - Symbols - Symboles - Symbole - Símbolos - Símbolos - Σύμβολα - Символи - Symboly - Symboler - Symboldid - Symbolit - Simboli - Szi mbólumok - Simboliai - Simboli - Symboler - Symbolen - Symbolika - Simboluri - Symboly - Simboli - Symboler

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