## Automatic Upper Arm Blood Pressure Monitor

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**Read Instruction manual ① and ② before use.**

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## 1. Introduction

Thank you for purchasing the OMRON Automatic Upper Arm Blood Pressure Monitor. This blood pressure monitor uses the oscillometric method of blood pressure measurement. This means this monitor detects your blood movement through your brachial artery and converts the movements into a digital reading.

### 1.1 Safety Instructions

This instruction manual provides you with important information about the OMRON Automatic Upper Arm Blood Pressure Monitor. To ensure the safe and proper use of this monitor, READ and UNDERSTAND all of the safety and operating instructions. **If you do not understand these instructions or have any questions, contact your OMRON retail outlet or distributor before attempting to use this monitor. For specific information about your own blood pressure, consult with your physician.**

### 1.2 Intended Use

This device is a digital monitor intended for use in measuring blood pressure and pulse rate in adult patient population. The device detects the appearance of irregular heartbeats during measurement and indicates this via a symbol with the measurement results. It is mainly designed for general household use.

### 1.3 Receiving and Inspection

Remove this monitor from the packaging and inspect for damage. If this monitor is damaged, DO NOT USE and consult with your OMRON retail outlet or distributor.

## 2. Important Safety Information

Read the Important Safety Information in this instruction manual before using this monitor. Follow this instruction manual thoroughly for your safety.

Keep for future reference. For specific information about your own blood pressure, CONSULT WITH YOUR PHYSICIAN.

### 2.1 Warning

**Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.**

- DO NOT use this monitor on infants, toddlers, children or persons who cannot express themselves.
- DO NOT adjust medication based on readings from this blood pressure monitor. Take medication as prescribed by your physician. ONLY a physician is qualified to diagnose and treat high blood pressure.
- DO NOT use this monitor on an injured arm or an arm under medical treatment.
- DO NOT apply the arm cuff on your arm while on an intravenous drip or blood transfusion.
- DO NOT use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- DO NOT use this monitor in oxygen rich environments or near flammable gas.
- Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation; arterial sclerosis; poor perfusion; diabetes; pregnancy; pre-eclampsia or renal disease. NOTE that any of these conditions in addition to patient motion, trembling, or shivering may affect the measurement reading.
- NEVER diagnose or treat yourself based on your readings. ALWAYS consult with your physician.
- To help avoid strangulation, keep the air tube and AC adapter cable away from infants, toddlers and children.
- This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers and children.

### Data Transmission

- This product emits radio frequencies (RF) in the 2.4 GHz band. DO NOT use this product in locations where RF is restricted, such as on an aircraft or in hospitals. Turn off the **Bluetooth®** feature in this monitor, remove batteries and/or unplug the AC adapter when in RF restricted areas.

## AC Adapter (optional accessory) Handling and Usage

- DO NOT use the AC adapter if this monitor or the AC adapter cable is damaged. If this monitor or the cable is damaged, turn off the power and unplug the AC adapter immediately.
- Plug the AC adapter into the appropriate voltage outlet. DO NOT use in a multi-outlet plug.
- NEVER plug in or unplug the AC adapter from the electric outlet with wet hands.
- DO NOT disassemble or attempt to repair the AC adapter.

## Battery Handling and Usage

- Keep batteries out of the reach of infants, toddlers and children.

### 2.2 Caution

**Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property.**

- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.
- Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an arteriovenous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury.
- Consult with your physician before using this monitor if you have had a mastectomy.
- Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders as cuff inflation can cause bruising.
- DO NOT take measurements more often than necessary because bruising, due to blood flow interference, may occur.
- ONLY inflate the arm cuff when it is applied on your upper arm.
- Remove the arm cuff if it does not start deflating during a measurement.
- DO NOT use this monitor for any purpose other than measuring blood pressure.
- During measurement, make sure that no mobile device or any other electrical device that emit electromagnetic fields is within 30 cm of this monitor. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.

- DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.
- DO NOT use in a location where there is moisture or a risk of water splashing this monitor. This may damage this monitor.
- DO NOT use this monitor in a moving vehicle such as in a car or on an aircraft.
- DO NOT drop or subject this monitor to strong shocks or vibrations.
- DO NOT use this monitor in places with high or low humidity or high or low temperatures. Refer to section 6.
- During measurement, observe the arm to ensure that the monitor is not causing prolonged impairment to blood circulation.
- DO NOT use this monitor in high-use environments such as medical clinics or physician offices.
- DO NOT use this monitor with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation and/or cause an inaccurate reading.
- Avoid bathing, drinking alcohol or caffeine, smoking, exercising and eating for at least 30 minutes before taking a measurement.
- Rest for at least 5 minutes before taking a measurement.
- Remove tight-fitting or thick clothing from your arm while taking a measurement.
- Remain still and DO NOT talk while taking a measurement.
- ONLY use the arm cuff on persons whose arm circumference is within the specified range of the cuff.
- Ensure that this monitor has acclimated to room temperature before taking a measurement. Taking a measurement after an extreme temperature change could lead to an inaccurate reading. OMRON recommends waiting for approximately 2 hours for the monitor to warm up or cool down when the monitor is used in an environment within the temperature specified as operating conditions after it is stored either at the maximum or at the minimum storage temperature. For additional information on operating and storage/transport temperature, refer to section 6.
- DO NOT use this monitor after the durable period has ended. Refer to section 6.
- DO NOT crease the arm cuff or the air tube excessively.
- DO NOT fold or kink the air tube while taking a measurement. This may cause an injury by interrupting blood flow.
- To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.

- ONLY use the AC adapter, arm cuff, batteries and accessories specified for this monitor. Use of unsupported AC adapters, arm cuffs and batteries may damage and/or may be hazardous to this monitor.
- ONLY use the approved arm cuff for this monitor. Use of other arm cuffs may result in incorrect readings.
- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied. NOTE: refer to "If your systolic pressure is more than 210 mmHg" in section 13 of instruction manual (2) for additional information.
- Read and follow the "Correct Disposal of This Product" in section 7 when disposing of the device and any used accessories or optional parts.

### Data Transmission

- DO NOT replace batteries or unplug the AC adapter while your readings are being transferred to your smart device. This may result in incorrect operation of this monitor and failure to transfer your blood pressure data.

### AC Adapter (optional accessory) Handling and Usage

- Fully insert the AC adapter into the outlet.
- When unplugging the AC adapter from the outlet, be sure to safely pull from the AC adapter. DO NOT pull from the AC adapter cable.
- When handling the AC adapter cable:
  - Do not damage it. / Do not break it. / Do not tamper with it.
  - DO NOT pinch it. / Do not forcibly bend or pull it. / Do not twist it.
  - DO NOT use it if it is gathered in a bundle.
  - DO NOT place it under heavy objects.
- Wipe any dust off of the AC adapter.
- Unplug the AC adapter when not in use.
- Unplug the AC adapter before cleaning this monitor.

### Battery Handling and Usage

- DO NOT insert batteries with their polarities incorrectly aligned.
- ONLY use 4 "AA" alkaline or manganese batteries with this monitor. DO NOT use other types of batteries. DO NOT use new and used batteries together. DO NOT use different brands of batteries together.
- Remove batteries if this monitor will not be used for a long period of time.

- If battery fluid should get in your eyes, immediately rinse with plenty of clean water. Consult with your physician immediately.
- If battery fluid should get on your skin, wash your skin immediately with plenty of clean, lukewarm water. If irritation, injury or pain persists, consult with your physician.
- DO NOT use batteries after their expiration date.
- Periodically check batteries to ensure they are in good working condition.

## 2.3 General Precautions

- To stop a measurement, press the [START/STOP] button while taking a measurement.
- When you take a measurement on the right arm, the air tube should be at the side of your elbow. Be careful not to rest your arm on the air tube.

- Blood pressure may differ between the right and left arm, and may result in a different measurement value. Always use the same arm for measurements. If the values between both arms differ substantially, check with your physician on which arm to use for your measurements.
- When using an optional AC adapter, make sure not to place your monitor in a location where it is difficult to plug and unplug the AC adapter.

### Battery Handling and Usage


- Disposal of used batteries should be carried out in accordance with local regulations.
- The supplied batteries may have a shorter life span than new batteries.

Remember to have a record of your blood pressure and pulse readings for your physician. A single measurement does not provide an accurate indication of your true blood pressure. Please use the Blood Pressure Diary to keep records of several readings over a certain period of time. To download PDF files of the diary, visit [www.omron-healthcare.com](http://www.omron-healthcare.com).

### 3. Error Messages and Troubleshooting

If any of the below problems occur during measurement, check to make sure that no other electrical device is within 30 cm. If the problem persists, please refer to the table below.

Display/Problem	Possible Cause	Solution
 appears or the arm cuff does not inflate.	The [START/STOP] button was pressed while the arm cuff is not applied.	Press the [START/STOP] button again to turn the monitor off. After inserting the air plug securely and applying the arm cuff correctly, press the [START/STOP] button.
	Air plug is not completely plugged into the monitor.	Insert the air plug securely.
	The arm cuff is not applied correctly.	Apply the arm cuff correctly, then take another measurement. Refer to section 7 of instruction manual (2).
	Air is leaking from the arm cuff.	Replace the arm cuff to the new one. Refer to section 14 of instruction manual (2).
 appears or a measurement cannot be completed after the arm cuff inflates.	You move or talk during a measurement and the arm cuff does not inflate sufficiently.	Remain still and do not talk during a measurement. If "E2" appears repeatedly, inflate the arm cuff manually until the systolic pressure is 30 to 40 mmHg above your previous readings. Refer to section 13 of instruction manual (2).
	Due to the systolic pressure is above 210 mmHg, a measurement cannot be taken.	
 appears	The arm cuff is inflated exceeding the maximum allowable pressure.	Do not touch the arm cuff and/or bend the air tube while taking a measurement. If inflating the arm cuff manually, refer to section 13 of instruction manual (2).
 appears	You move or talk during a measurement. Vibrations disrupt a measurement.	Remain still and do not talk during a measurement.
 appears	The pulse rate is not detected correctly.	Apply the arm cuff correctly, then take another measurement. Refer to section 7 of instruction manual (2). Remain still and sit correctly during a measurement.
 appears		If the "  " symbol continues to appear, we recommend you to consult with your physician.
 does not flash during a measurement.		

Display/Problem	Possible Cause	Solution
 appears	The monitor has malfunctioned.	Press the [START/STOP] button again. If “Er” still appears, contact your OMRON retail outlet or distributor.
 appears	The monitor cannot connect to a smart device or transmit data correctly.	Follow the instructions shown in the “OMRON connect” app. If the “Err” symbol still appears after checking the app, contact your OMRON retail outlet or distributor.
 flashes	The monitor is waiting for pairing with the smart device.	Refer to section 5 of instruction manual (2) for pairing your monitor with your smart device, or press [START/STOP] button to cancel pairing and turn your monitor off.
 flashes	The monitor is ready to transfer your readings to the smart device.	Open the “OMRON connect” app to transfer your readings.
 flashes	More than 48 readings are not transferred. The date and time is not set.	Pair or transfer your readings to the “OMRON connect” app so you can keep them in memory in the app, and this error symbol disappears.
 appears	There are 60 readings in memory to be transferred	
 flashes	Batteries are low.	Replacing all 4 batteries with new ones is recommended. Refer to section 4 of instruction manual (2).
 appears or the monitor is turned off unexpectedly during a measurement	Batteries are depleted.	Immediately replace all 4 batteries with new ones. Refer to section 4 of instruction manual (2).
Nothing appears on the display of the monitor.	Battery polarities are not properly aligned.	Check the battery installation for proper placement. Refer to section 4 of instruction manual (2).
Readings appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, and/or how you apply the arm cuff, may affect your blood pressure. Review section 2 of instruction manual (2).	

Display/Problem	Possible Cause	Solution
Any other communication issue occurs.	Follow the instructions shown in the smart device, or visit the "Help" section in the "OMRON connect" app for further help. If the problem still persists, contact your OMRON retail outlet or distributor.	
Any other problem occurs.	Press the [START/STOP] button to turn the monitor off, then press it again to take a measurement. If the problem continues, remove all batteries and wait for 30 seconds. Then re-install batteries. If the problem still persists, contact your OMRON retail outlet or distributor.	

## 4. Limited Warranty

Thank you for buying an OMRON product. This product is constructed of high quality materials and great care has been taken in its manufacturing. It is designed to give you every satisfaction, provided that it is properly operated and maintained as described in the instruction manual.

This product is warranted by OMRON for a period of 5 years after the date of purchase. The proper construction, workmanship and materials of this product is warranted by OMRON. During this period of warranty OMRON will, without charge for labour or parts, repair or replace the defect product or any defective parts.

The warranty does not cover any of the following:

- A. Transport costs and risks of transport.
- B. Costs for repairs and / or defects resulting from repairs done by unauthorised persons.
- C. Periodic check-ups and maintenance.
- D. Failure or wear of optional parts or other attachments other than the main device itself, unless explicitly warranted above.
- E. Costs arising due to non-acceptance of a claim (those will be charged for).
- F. Damages of any kind including personal caused accidentally or from misuse.
- G. Calibration service is not included within the warranty.
- H. Optional parts have a one (1) year warranty from date of purchase. Optional parts include, but are not limited to the following items: cuff and cuff tube.

Should warranty service be required please apply to the dealer whom the product was purchased from or an authorised OMRON distributor. For the address refer to the product packaging / literature or to your specialised retailer. If you have difficulties in finding OMRON customer services, contact us for information:

[www.omron-healthcare.com](http://www.omron-healthcare.com)

Repair or replacement under the warranty does not give rise to any extension or renewal of the warranty period.

The warranty will be granted only if the complete product is returned together with the original invoice / cash ticket issued to the consumer by the retailer.

## 5. Maintenance

### 5.1 Maintenance

To protect your monitor from damage, follow the directions below:

Changes or modifications not approved by the manufacturer will void the user warranty.

#### **Caution**

---

DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

---

### 5.2 Storage

- Keep your monitor in the storage case when not in use.
  1. Remove the arm cuff from the monitor.

#### **Caution**

---

To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.

---

2. Gently fold the air tube into the arm cuff. Note: Do not bend or crease the air tube excessively.
  3. Place your monitor and other components in the storage case.
- Store your monitor and other components in a clean, safe location.
  - Do not store your monitor and other components:
    - If your monitor and other components are wet.
    - In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
    - In locations exposed to vibrations or shocks.

### 5.3 Cleaning

- Do not use any abrasive or volatile cleaners.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to clean your monitor and arm cuff, and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to clean your monitor and arm cuff or other components.

### 5.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

## 6. Specifications

Product description	Automatic Upper Arm Blood Pressure Monitor		
Product category	Electronic Sphygmomanometers		
Model (code)	M4 Intelli IT (HEM-7155T-EBK) X4 Smart (HEM-7155T-ESL)	Display	LCD digital display
Cuff pressure range	0 to 299 mmHg	Pulse measurement range	40 to 180 beats / min.
Blood pressure measurement range	SYS: 60 to 260 mmHg / DIA: 40 to 215 mmHg		
Accuracy	Pressure: $\pm 3$ mmHg / Pulse: $\pm 5\%$ of display reading		
Inflation	Automatic by electric pump	Deflation	Automatic pressure release valve
Measurement method	Oscillometric method	Transmission method	<b>Bluetooth</b> ® Low Energy
Wireless communication	Frequency range: 2.4 GHz (2400 - 2483.5 MHz) / Modulation: GFSK Effective radiated power: < 20 dBm		
Operating mode	Continuous operation		
IP classification	Monitor: IP20 / Optional AC adapter: IP21 (HHP-CM01) or IP22 (HHP-BFH01)		
Rating	DC6 V 4.0 W	Applied part	Type BF (arm cuff)
Power source	4 "AA" batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50 - 60 Hz 0.12 - 0.065 A)		
Battery life	Approximately 1000 measurements (using new alkaline batteries)		
Durable period (Service life)	Monitor: 5 years / Cuff: 5 years / Optional AC adapter: 5 years		
Operating conditions	+10 to +40 °C / 15 to 90% RH (non-condensing) / 800 to 1060 hPa		
Storage / Transport conditions	-20 to +60 °C / 10 to 90% RH (non-condensing)		
Contents	Monitor, arm cuff (HEM-FL31), 4 "AA" batteries, storage case, Instruction Manual ① and ②, setup instructions		
Protection against electric shock	Internally powered ME equipment (when using only batteries) Class II ME equipment (optional AC adapter)		
Weight	Monitor: approximately 340 g (not including batteries) / Arm cuff: approximately 163 g		
Dimensions (approximately value)	Monitor: 105 mm (W) × 85 mm (H) × 152 mm (L) Arm cuff: 145 mm × 532 mm (air tube: 750 mm)		
Memory	Stores up to 60 readings per user		

## Note

- These specifications are subject to change without notice.
- This monitor is clinically investigated according to the requirements of EN ISO 81060-2:2014 and complies with EN ISO 81060-2:2014 and EN ISO 81060-2:2019 + A1:2020. In the clinical validation study, K5 was used on 85 subjects for determination of diastolic blood pressure.
- This device has been validated for use on pregnant and pre-eclampsia patients according to the Modified European Society of Hypertension Protocol\*.
- This device has been validated for use on diabetic (Type II) population\*\*.
- IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. This monitor and optional AC adapter are protected against solid foreign objects of 12.5 mm diameter and greater such as a finger. The optional AC adapter HHP-CM01 is protected against vertically falling water drops which may cause issues during a normal operation. The optional AC adapter HHP-BFH01 is protected against oblique falling water drops which may cause issues during a normal operation.

\*Topouchian J et al. Vascular Health and Risk Management 2018;14 189–197

\*\* Chahine M.N. et al. Medical Devices: Evidence and Research 2018;11 11–20

## About a wireless communication interference

This product operates in an unlicensed ISM band at 2.4 GHz. In the event this product is used near other wireless devices such as microwave and wireless LAN, which operate on the same frequency band as this product, there is a possibility that interference may occur. If interference occurs, stop the operation of the other devices or relocate this product away from other wireless devices before attempting to use it.

## 7. Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

This marking shown on the product or its literature, indicates that it should not be disposed of, with other household wastes at the end of its working life.

To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this product from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can return this item for environmentally safe recycling.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial waste for disposal.

## 8. Important Information Regarding Electromagnetic Compatibility (EMC)

HEM-7155T-EBK and HEM-7155T-ESL conform to EN60601-1-2:2015 Electromagnetic Compatibility (EMC) standard.

Further documentation in accordance with this EMC standard is available at [www.omron-healthcare.com](http://www.omron-healthcare.com)

Refer to the EMC information for HEM-7155T-EBK and HEM-7155T-ESL on the website.

## 9. Guidance and Manufacturer's Declaration

- This blood pressure monitor is designed according to the European Standard EN1060, Non-invasive sphygmomanometers Part 1: General Requirements and Part 3: Supplementary requirements for electromechanical blood pressure measuring systems.
- Hereby, OMRON HEALTHCARE Co., Ltd., declares that the radio equipments type HEM-7155T-EBK and HEM-7155T-ESL are complied with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: [www.omron-healthcare.com](http://www.omron-healthcare.com)
- This OMRON product is produced under the strict quality system of OMRON HEALTHCARE Co., Ltd., Japan. The Core component for OMRON blood pressure monitors, which is the Pressure Sensor, is produced in Japan.
- Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

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### Applied part - Type BF Degree of protection against electric shock (leakage current)

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### Class II equipment. Protection against electric shock

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**IP XX**

### Ingress protection degree provided by IEC 60529

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The **Bluetooth®** word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by OMRON HEALTHCARE Co., Ltd. is under license. Other trademarks and trade names are those of their respective owners. App Store and App Store logo are service marks of Apple Inc., registered in the U.S and other countries. Google Play and Google Play logo are trademarks of Google LLC.

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
Intellisense

IM2-HEM-7155T-E-05-01/2022  
2870441-5D

## Read instruction manual ① and ② before use.


- FR Lire le mode d'emploi ① et ② avant l'utilisation.
- DE Lesen Sie vor der Verwendung die Gebrauchsanweisung ① und ②.
- IT Leggere il manuale di istruzioni ① e ② prima dell'uso.
- ES Lea el manual de instrucciones ① y ② antes del uso.
- NL Lees de gebruiksaanwijzing ① en ② voor gebruik.
- RU Прочтите руководство по эксплуатации ① и ② перед использованием.
- TR Kullanmadan önce, kullanım kılavuzu ① ve ②'yi okuyun.

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

## The date and time will automatically be set when your monitor is paired with the app.

- FR La date et l'heure seront automatiquement réglées lorsque votre tensiromètre sera jumelé avec l'application.
- DE Wenn Ihr Messgerät mit der App verknüpft ist, werden Datum und Uhrzeit automatisch eingestellt.
- ES La fecha y la hora se ajustarán automáticamente cuando el monitor se sincronice con la aplicación.
- NL De datum en tijd worden automatisch ingesteld wanneer uw bloeddrukmeter is gekoppeld aan de app.
- RU Дата и время установятся автоматически после синхронизации с приложением.
- TR Ölçüm cihazınız uygulamayla eşleştirildiğinde tarih ve saat otomatik olarak ayarlanır.

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



## If your monitor is paired with your smart device, date and time is set automatically. When you need to set them manually, set year > month > day > hour > minute.

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
## A Tube side of the cuff should be 1 - 2 cm above the inside elbow.


- FR Le côté tuyau du brassard doit être positionné 1 à 2 cm au-dessus de l'intérieur du coude.
- DE Das Manschettenstück mit dem Schlauch muss 1 bis 2 cm oberhalb des Ellbogens liegen.
- IT Il lato del bracciale con il tubo deve trovarsi al di sopra dell'interno del gomito, a una distanza di circa 1 o 2 cm.
- ES El lado del tubo del manguito deberá quedar 1 o 2 cm por encima de la parte interna del codo.
- NL De kant met de slang van de manchet moet 1 - 2 cm boven de binnenkant van de elleboog liggen.
- RU Край манжеты с трубкой должен находиться на 1-2 см выше локтевого сгиба.
- TR Kolluğun boru tarafı, dirsek içinin 1 ila 2 cm üstünde olmalıdır.

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## B Make sure that air tube is on the inside of your arm and wrap the cuff securely so it can no longer slip round.

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
## If taking measurements on the right arm, refer to:

- FR Pour la prise de mesures au bras droit, voir :
- DE Bei Messungen am rechten Arm siehe:
- IT Se la misurazione viene eseguita al braccio destro, fare riferimento a:
- ES Si va a realizar mediciones en el brazo derecho, consulte:
- NL Als u metingen aan de rechterarm uitvoert, raadpleeg dan:
- RU При измерении на правой руке см:
- TR Sağ koldan ölçüm yapıyorsanız aşağıdakilere bakın:

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
## 1 Sit comfortably with your back and arm supported.

- FR S'asseoir confortablement de manière à ce que le dos et le bras soient bien soutenus.
- DE Bequem hinsetzen, mit dem Rücken anlehnen und den Arm auf eine Unterlage legen.
- IT Sedere comodamente con la schiena e il braccio ben sostenuti.
- ES Siéntese cómodamente de modo que la espalda y el brazo estén bien apoyados.
- NL Neem een comfortabele zitpositie in, met uw rug en arm ondersteund.
- RU Сядьте удобно, чтобы спина и рука опирались на что-либо.
- TR Sirtınız ve kolunuz desteklenerek rahat bir şekilde oturun.

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## 2 Place the arm cuff at the same level as your heart.

- FR Le brassard doit se trouver au même niveau que votre cœur.
- DE Die Manschette auf Herzhöhe platzieren.
- IT Posizionare il bracciale allo stesso livello del cuore.
- ES Coloque el manguito al mismo nivel que el corazón.
- NL Plaats de armmanchet op hetzelfde niveau als uw hart.
- RU Манжета должна находиться на уровне сердца.
- TR Kolluğu kalbinizle aynı düzeye getirin.

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## 3 Keep feet flat, legs uncrossed, remain still and do not talk.


- FR Garder les pieds à plat, les jambes non croisées, ne pas bouger et ne pas parler.
- DE Füße flach auf den Boden stellen, die Beine nicht kreuzen, still halten und nicht sprechen.
- IT Tenere i piedi ben poggiati, non incrociare le gambe, rimanere fermi e non parlare.
- ES Mantenga los pies planos, las piernas sin cruzar, quédese quieto y no hable.
- NL Houd de voeten plat, benen niet gekruist, zit stil en praat niet.
- RU Поставьте ноги ровно, не перекрещивая их, не двигайтесь и не разговаривайте.
- TR Ayaklarınızı düz, bacaklarınızı açık tutun, hareketsiz kalın ve konuşmayın.

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## Switching user ID enables you to save readings for 2 people.

- FR Changer d'ID Utilisateur permet d'enregistrer les résultats pour 2 personnes.
- DE Mit verschiedenen Benutzer-IDs lassen sich Messwerte für 2 Personen speichern.
- IT Cambiare ID utente permette di registrare i risultati relativi a 2 persone.
- ES Cambiar de ID de usuario le permite guardar las lecturas de 2 personas.
- NL Door het wisselen van gebruikers-ID kunt u metingen voor 2 personen opslaan.
- RU Переключение идентификатора пользователя позволяет сохранять результаты для 2-x людей.
- TR Kullanıcı kimlikleri arasında geçiş yapma, 2 kişi için ölçüm değerlerini kaydetmenizi sağlar.

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**When the [START/STOP] button is pressed, the measurement is taken and saved automatically. Open the app to transfer the reading.**

**FR** Après une pression sur le bouton [START/STOP], la mesure démarre et le résultat est enregistré automatiquement. Ouvrir l'application pour transférer le résultat.  
**DE** Beim Drücken der Taste [START/STOP] erfolgt die Messung. Die Messwerte werden automatisch gespeichert. Öffnen Sie die App, um die Messwerte zu übertragen.  
**IT** Quando si preme il pulsante [START/STOP], la misurazione viene eseguita e salvata automaticamente. Aprire la app per trasferire i risultati.  
**ES** Al pulsar el botón [START/STOP], se obtiene la medición, que se guarda automáticamente. Abra la aplicación para transferir la lectura.  
**NL** Wanneer op de knop [START/STOP] wordt gedrukt, wordt de meting automatisch uitgevoerd en opgeslagen. Open de app om de meting te versturen.  
**RU** После того, как кнопка [START/STOP] нажата, измерение выполняется и сохраняется автоматически. Откройте приложение, чтобы передать результаты.  
**TR** [START/STOP] düğmesine basıldığında ölçüm yapılır ve otomatik olarak kaydedilir. Ölçüm değerini aktarmak için uygulamayı açın.

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
**If your systolic pressure is more than 210 mmHg:**

After the arm cuff starts to inflate, press and hold the [START/STOP] button until the monitor inflates 30 to 40 mmHg higher than your expected systolic pressure.

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**The guest mode can be used to take a single measurement for another user. No readings are stored in the memory.**

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**2 Appears when an irregular rhythm\*\* is detected during a measurement. If it appears repeatedly, OMRON recommends to consult your physician.**

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**3 Appears when your body moves during a measurement. Remove the arm cuff, wait 2-3 minutes and try again.**

**FR** S'affiche si vous bougez pendant une mesure. Retirez le brassard, attendez 2 à 3 minutes et essayez à nouveau.  
**DE** Wird angezeigt, wenn Sie sich während der Messung bewegen. Nehmen Sie die Manschette ab, warten Sie 2–3 Minuten und versuchen Sie es erneut.  
**IT** Viene visualizzato se l'utilizzatore si muove durante la misurazione. Rimuovere il bracciale, attendere 2-3 minuti e riprovare.  
**ES** Aparece cuando se mueve el cuerpo durante una medición. Retire el manguito, espere unos 2 o 3 minutos e inténtelo de nuevo.  
**NL** Verschijnt wanneer uw lichaam tijdens een meting beweegt. Verwijder de armband, wacht 2-3 minuten en probeer het opnieuw.  
**RU** Отображается при движении тела во время измерения. Снимите манжету, подождите 2-3 минуты и попробуйте еще раз.  
**TR** Bir ölçüm esnasında vücudunuz hareket ettiğinde görünür. Kolluğu çıkarın, 2-3 dakika bekleyip tekrar deneyin.  
**AR** تظهر عند تحرك جسمك في أثناء عملية القياس. قم بإزالة الشريط الضاغط للذراع وانتظر من ٢ إلى ٣ دقائق المحاولة مرة أخرى.

**NL** Metingen bekijken  
**RU** Проверка результатов  
**TR** Ölçüm Değerlerini Kontrol Etme

**4 Cuff is tight enough.**

**FR** Le brassard est suffisamment serré.  
**DE** Manschette sitzt ausreichend straff.  
**IT** Il bracciale è stretto a sufficienza.  
**ES** El manguito está lo suficientemente prieto.  
**NL** Manchet zit strak genoeg.  
**RU** Манжета затянута достаточно туго.  
**TR** Kolluk yeterince sikirdi.  
**AR** الشريط الضاغط محكم بما يكفي.

**4 Apply cuff again MORE TIGHTLY.**

**FR** Poser le brassard en le serrant davantage.  
**DE** Manschette STRAFFER ziehen.  
**IT** Applicare di nuovo il bracciale STRINGENDOLO DI PIÙ.  
**ES** Vuelva a poner el manguito MÁS PRIETO.  
**NL** Breng de manchet STRAKKER aan.  
**RU** Наложите манжету еще раз БОЛЕЕ ТУГО.  
**TR** Kolluğu tekrar, DAHA SIKI bir şekilde takın.  
**AR** قبل استخدام وظائف الذاكرة، حدد معرف المستخدم الخاص بك.

**\*The high blood pressure definition is based on the 2018 ESH/ESC Guidelines.**


**FR** \*La définition de l'hypertension est basée sur les recommandations ESH/ESC 2018.  
**DE** \*Die Definition für Bluthochdruck basiert auf den 2018 ESH/ESC Richtlinien.  
**IT** \*La definizione di alta pressione arteriosa si basa sulle linee guida 2018 di ESH/ESC.  
**ES** \*La definición de presión arterial alta se basa en las guías 2018 ESH/ESC.  
**NL** \*De definitie van hoge bloeddruk is gebaseerd op de ESH/ESC-richtlijnen uit 2018.  
**RU** \*Определение высокого артериального давления основано на Рекомендациях 2018 ESH/ESC.  
**TR** \*Yüksek kan basıncı tanımı, 2018 ESH/ESC Kilavuzlarını temel alır.

**\*\*An irregular heartbeat rhythm is defined as a rhythm that is 25% less or 25% more than the average rhythm detected during a measurement.**

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
**Before using memory functions, select your user ID.**

**FR** Avant d'utiliser les fonctions de mémoire, sélectionner votre ID Utilisateur.  
**DE** Wählen Sie vor dem Verwenden der Speicherfunktionen Ihre Benutzer-ID aus.  
**IT** Prima di utilizzare le funzioni di memoria selezionare il proprio ID utente.  
**ES** Antes de usar las funciones de memoria, seleccione su ID de usuario.  
**NL** Selecteer uw gebruikers-ID voordat u geheugenfuncties gebruikt.  
**RU** Перед использованием функций памяти выберите идентификатор пользователя.  
**TR** Hafıza fonksiyonunu kullanmadan önce kullanıcı kimliğinizi seçin.  
**AR** قبل استخدام وظائف الذاكرة، حدد معرف المستخدم الخاص بك.

**12.1 Readings Stored in Memory**

**FR** Mesures stockées en mémoire  
**DE** Gespeicherte Messungen  
**IT** Risultati conservati in memoria  
**ES** Lecturas guardadas en la memoria  
**NL** Meetwaarden opgeslagen in het geheugen  
**RU** Сохранение результатов в памяти  
**TR** Hafızada Saklanan Ölçüm Değerleri  
**AR** القراءات المخزنة في الذاكرة

**12.2 Average of the Latest 2 or 3 Readings Taken within a 10 Minute Span**

**FR** Moyenne des 2 ou 3 dernières mesures prises en l'espace de 10 minutes  
**DE** Mittelwert der letzten 2 oder 3 Messwerte, die in einem 10-Minuten-Zeitraum erfasst wurden  
**IT** Media degli ultimi 2 o 3 risultati ottenuti nell'arco di 10 minuti  
**ES** Promedio de las 2 o 3 últimas lecturas realizadas en un intervalo de 10 minutos  
**NL** Gemiddelde van de laatste 2 of 3 metingen uitgevoerd binnen een tijdsperiode van 10 minuten  
**RU** Среднее значение из последних 2-х или 3-х результатов полученных в течение 10 минут  
**TR** 10 Dakika Aralıkta Alınan Son 2 veya 3 Ölçüm Değerinin Ortalaması  
**AR** متوسط آخر قراءتين أو ٣ قراءات تم الحصول عليها في مدة زمنية تبلغ ١٠ دقائق

**12.3 Deleting All Readings**

**FR** Suppression de toutes les mesures  
**DE** Löschen aller Messwerte  
**IT** Cancellazione di tutti i risultati  
**ES** Eliminación de todas las lecturas  
**NL** Alle metingen verwijderen  
**RU** Удаление всех результатов  
**TR** Tüm Ölçüm Değerlerini Silme  
**AR** حذف جميع القراءات

**13 Other Settings**  
**FR** Autres réglages  
**DE** Weitere Einstellungen  
**IT** Altre impostazioni  
**ES** Otros ajustes  
**NL** Andere instellingen  
**RU** Другие настройки  
**TR** Diğer Ayarlar  
**AR** الإعدادات الأخرى

**13.1 Disabling/Enabling Bluetooth**


**FR** Désactivation/activation de Bluetooth  
**DE** Deaktivieren/Aktivieren der Bluetooth-Funktion  
**IT** Disattivazione/attivazione della funzione Bluetooth  
**ES** Activar o desactivar el Bluetooth  
**NL** Bluetooth uitschakelen/inschakelen  
**RU** Выключение/включение Bluetooth  
**TR** Bluetooth'u Devre Dışı Brakma/Etkinleştirme  
**AR** تعطيل/تفعيل خاصية Bluetooth

**Bluetooth is enabled by default.**


**FR** Bluetooth est activé par défaut.  
**DE** Die Bluetooth-Funktion ist standardmäßig aktiviert.  
**IT** La funzione Bluetooth è attiva per impostazione predefinita.  
**ES** El Bluetooth está activado por defecto.  
**NL** Bluetooth is standaard ingeschakeld.  
**RU** Bluetooth включен по умолчанию.  
**TR** Bluetooth varsayılan olarak etkindir.  
**AR** يتم تمكين خاصية Bluetooth بشكل افتراضي.

**13.2 Restoring to the Default Settings**

**FR** Réinitialisation aux réglages par défaut  
**DE** Wiederherstellen der Standardinstellungen  
**IT** Ripristino delle impostazioni predefinite  
**ES** Restablecimiento a los ajustes de fábrica  
**NL** De standaardinstellingen herstellen  
**RU** Восстановление настроек по умолчанию  
**TR** Varsayılan Ayarları Geri Yükleme  
**AR** استعادة الإعدادات الافتراضية

**14 Optional Medical Accessories**  
**FR** Accessoires médicaux optionnels  
**DE** Optionales medizinisches Zubehör  
**IT** Accessori medicali opzionali  
**ES** Accesorios médicos opcionales  
**NL** Optionele medische accessoires  
**RU** Дополнительно принадлежности  
**TR** Opsiyonel Tıbbi Aksesuarlar  
**AR** الملحقات الطبية الاختيارية

**Do not throw the air plug away. The air plug can be applicable to the optional cuff.**

**FR** Ne pas jeter la prise de gonflage. La prise de gonflage peut être utilisée pour le brassard en option.  
**DE** Entsorgen Sie den Luftschauchstecker nicht. Der Luftschauchstecker wird für die optionale Manschette verwendet.  
**IT** Non gettare via l'attacco del tubo dell'aria. L'attacco del tubo dell'aria può essere applicato al bracciale opzionale.  
**NL** No tire el conector para tubo de aire. El conector para tubo de aire puede ser utilizado con el manguito opcional.  
**RU** Goodi de plug van de luchtslang kan worden gebruikt op de optionele manchet.  
**TR** He vıyabırsayınca vıduşlu ştekere. On mekte pıdodıyık k dıplonlelınen mınketıe. Hava tırasını atmayın. Hava tırasını isteğe bağılı kolluğa uygulayınabılır.  
**AR** لا تفسد الليف مع الضيف مع الليف الضاغط الاختياري.

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