


ambu.com
Tél. +45 72 25 20 00
Denmark
DK-2750 Ballerup
Ballerupbæken 13
Ambu A/S
EU Importer:

Medico Electronics Ltd.
Bld. Saint-Michel 65/6,
1040 Brussels, Belgium
Made in India
[AN ISO 13485 COMPANY]
www.medicoelectronics.com
Plot 142A/1, 1, 2, 7, 28 & 29, Noida Special Economic Zone,
Noida, Uttar Pradesh - 201305, India



ECG Electrodes

Instructions for use

These instructions for use are intended for the following Ambu White Sensor electrodes:
WS, WSRT, WP25, WP3A, 0215M, 0315M, 0415M, 0715M, 2742*, 2742RT*, 2837*, 3345C/RT*, 3351*,
RT*, 3552*, 3552/RT*, 3552DUAL/RT*, 4200, 4300, 4242, 4440, 4500M, 4500M-H, 4530, 4535M, 4540,
4560M, 4570M, 4831/2, 4841P, 5050, 6060, 7841P, CMM1.
Ambu Blue Sensor electrodes: BS-TAB-2300.

Instructions for use

Инструкции за употреба

Návod k použití

Brugervejledning

Bedienungsanleitung

Οδηγίες χρήσης

Instrucciones de uso

Kasutusjuhend

Käyttöohjeet

Mode d'emploi

Upute za upotrebu

Használati útmutató

Istruzioni per l'uso

使用説明書

Naudojimo instrukcija

Lietošanas pamācība

Gebruiksaanwijzing

Bruksanvisning

Instrukcja użycia

Instruções de utilização

Instruções de utilizare

Инструкция по применению

Návod na použití

Navodila za uporabo

Bruksanvisning

Kullanım talimatları

使用说明



Instructions for use

ECG electrodes

1.0 Intended Use

The electrode is applied to the surface of the body to transmit the electrical signal at the body surface to a recording and monitoring device that produces an electrocardiogram or vectorcardiogram.

1.2 Warnings and Cautions

Warnings

1. Surface electrodes should be applied only to intact, clean skin (e.g., not over open wounds, lacerations, infected, or inflamed areas).

2. For single use only. Reuse of the product may cause contamination, leading to infections.

Cautions

- During surgical procedures electrodes should be placed as far as possible from the electro-surgical area to minimize unwanted RF current flow which could lead to burns.
- Electrodes should only be used by or in consultation with a health care provider familiar with their proper placement and use.
- Do not soak, rinse, or sterilize the electrodes as these procedures may leave harmful residues or cause malfunction of the electrodes.
- To minimize the risk of skin irritation, when replacing

1.3 Potential Adverse Events

Potential adverse events in relation to the use of surface electrodes are skin reactions such as skin burns, skin redness, skin itching and small blisters.

2.0 Explanation of Symbols Used

Symbols for the Ambu White/Blue Sensor electrodes	Indication	Symbols for the Ambu White/Blue Sensor electrodes	Indication	Symbols for the Ambu White/Blue Sensor electrodes	Indication
	Medical Device		Wet gel		Cloth
	Radiolucent		Polymer		Foam
	Solid gel		Electrical Safety Type II Applied Part		
	Magnetic Resonance Conditional: Static magnetic field of 1.5 Tesla and 3 Tesla only. Maximum spatial gradient field of 25,000 Gauss/cm or 250 Tesla/m. Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning. * The Maximum Spatial Gradient Field of 12,900 Gauss/cm or 129 Tesla/m.				

3.0 Directions for Use

- If needed, clean application areas with a mild soap and water and dry thoroughly.
- For adults only: Clip off hair on excessively hairy areas if necessary. Only gently shave if skin irritation is too high. The correct placement of electrodes is typically determined by the hospital physician and/or local protocol procedures. Protective liners may have sharp corners, and care should be taken when removing electrodes from the pouch. Pull the electrode away from the protective liner and place it on the skin. Dispose of the protective liner immediately after removal of electrodes.
- When using electrodes with connectors, the connector should be positioned in the direction so that leadwires run in the same direction towards the equipment.
- Place the electrode on the skin by gently pressing around the adhesive edge with a circular smooth and firm motion.
- Connect the leadwire from the ECG monitoring device to the electrode.
- Electrode Removal: Gently lift the tab or the edge of the electrode, maintaining a low profile, and slowly peel off the electrode from the skin, holding and supporting down the newly exposed skin as the electrode is removed. Keep the electrode close to the skin surface when pulling back over itself. Avoid exposing the electrode at 45 degree angle to reduce the risk of mechanical trauma.
- Self-adhesive electrodes should be replaced if they no longer stick firmly to the skin.
- Gel residue can be removed using paper towels or water.
- Electrodes in an opened pouch may dry out, therefore do not open pouch until necessary.
- The electrodes are disposable and must be disposed of according to local hospital procedures.

4.0 Duration of application (time in hours)

Item number	Gel System	Adhesive System	Sensor Material	Patient Indication	Options: Application of use	Duration of application (hours)
BS-TAB-2300	Solid Gel	Solid Gel	Ag/AgCl	Adult	Resting	1/2
WS	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
WP5RT	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
WP25	Solid Gel	Acrylic	Ag/AgCl	Pediatric	Resting/Long term	48
WP3D	Solid Gel	Acrylic	Ag/AgCl	Adult/Pediatric	Resting/Long term	48
0215M	Solid Gel	Solid Gel	Ag/AgCl	Pediatric	Resting	1/2
0315M	Solid Gel	Solid Gel	Ag/AgCl	Adult	Resting	1/2
0415M	Solid Gel	Solid Gel	Ag/AgCl	Adult	Resting	1/2
0715M	Solid Gel	Solid Gel	Ag/AgCl	Adult	Resting	1/2
2742	Wet gel	Acrylic	Ag/AgCl	Adult/Pediatric	Resting	24
2742RT	Wet gel	Acrylic	Ag/AgCl	Adult	Long term**	48
2837	Wet gel	Acrylic	Ag/AgCl	Adult	Resting	1/2
3345C/RT	Wet gel	Acrylic	Ag/AgCl	Adult	Resting/Long term	48
3351/RT	Solid Gel	Acrylic	Ag/AgCl	Adult	Stress-Resting/Monitoring	48
3552	Wet gel	Acrylic	Ag/AgCl	Adult	Resting	24
3552/RT	Wet gel	Acrylic	Ag/AgCl	Adult	Long term**	48
3552DUAL/RT	Solid Gel	Acrylic	Ag/AgCl	Adult	Resting/Long term	48
4200	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
4300	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
4242	Wet gel	Acrylic	Ag/AgCl	Adult	Stress/Long term	48
4440	Wet gel	Acrylic	Ag/AgCl	Adult	Long term	48
4500M	Wet gel	Acrylic	Ag/AgCl	Adult	Short term	24
4500M-H	Wet gel	Acrylic	Ag/AgCl	Adult	Short term	24
4530	Wet gel	Acrylic	Ag/AgCl	Adult	Resting	24
4535M	Wet gel	Acrylic	Ag/AgCl	Adult	Short term	24
4540	Wet gel	Acrylic	Ag/AgCl	Adult	Resting	24
4560M	Wet gel	Acrylic	Ag/AgCl	Adult	Stress test	24
4570M	Wet gel	Acrylic	Ag/AgCl	Pediatric	Short term	24
4831/2	Solid Gel	Acrylic	Ag/AgCl	Adult/Pediatric	Short term	24
4841P	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
5050	Wet gel	Acrylic	Ag/AgCl	Adult	Stress/Long term	48
6060	Wet gel	Acrylic	Ag/AgCl	Adult	Stress/Long term	48
7841P	Solid Gel	Acrylic	Ag/AgCl	Adult	Short term	24
CMM1	Solid Gel	Acrylic	Ag/AgCl	Adult	Long term	48