



PROFESSIONAL MEDICAL PRODUCTS

REF 80210/80211
DIGITAL VET SCALE

SPECIFICATION:

- Precision measuring system
- Tare function
- Overload indicator "Err"
- Low battery indicator
- with LCD wall-mount detachable display

TECHNICAL PARAMETERS:

Voltage: 9 V alkaline battery and mains adapter plug in (included)

Max. measurable weight: 300 kg

Weighing increments: 100 g

Min. measurable weight: 2 kg

Platform: 600x800 mm

Weighing units: kg/lb/st

Working/Storage temperature: 10-40°C/-10-50°C

IMPORTANT SAFETY INSTRUCTIONS:

- Make sure you follow the marked polarity of battery when inserting and replacing it.
- If the appliance will be out of use for a longer time, remove the battery.
- Do not leave the unit unattended when in use.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- Use the product only for the intended purpose specified in the Operating Manual.
- Use the unit only on a stable horizontal surface.
- Before cleaning, assembly, or disassembly of the cover, switch the unit off and remove the battery!
- Do not use abrasive or chemically aggressive substances for cleaning the unit.
- Do not use the unit if it does not operate properly, or if it has been dropped, damaged or exposed to any liquids. Take the unit to an authorized service center for testing and repair.
- The unit is suitable for domestic use only. It is not designed for commercial operation.
- Keep the unit away from heat sources, such as radiators, ovens etc. Do not expose the unit to direct sunlight or moisture.
- Do not touch the unit with wet or damp hands.
- The product is a fragile measuring device. Prevent it from crashing or falling.
- Do not repair the appliance yourself. Contact an authorized service center.

FAILURE TO FOLLOW THE MANUFACTURER'S INSTRUCTIONS MAY LEAD TO REFUSAL OF WARRANTY REPAIR!

OPERATING BUTTONS:

1. The scale can be switched on with the ON/OFF/Tare button.
2. Then "0.0" is displayed, the scale is ready for measuring, put weight on it.
3. After measuring when you press the "Hold" button once, the weight value of the display will freeze and it will display the weight continuously, and store in the memory.
4. After weighing by pressing the ON/OFF /Tare button twice the scale switches off, then the information from the memory will be cleared.
5. If you press the ON/OFF/Tare button just once, then the display will go back to "0.0".
6. Then put weight on the scale. The measured value will be displayed. With pressing the "D" button we can get the difference between the actual weight and the previous weight in the memory.


FURTHER SPECIFICATIONS:

- The scale is able to measure continually without power off.
- Memory for the comparison of the actual and the previous weighing.
- The adapter plug-in can be found on the box of display.
- Buttons are functional with sound-signal
- Recommended working environment:
10- 40 Celsius degree.
The most accurate measuring
result: on + 22 Celsius degree.
- The prevailing weight value is able to see continually on the display.

ATTENTION!

- Avoid injuries on the mains adapter plug in.
- Don't push the scale aside when weight is on it.

DISPLAY:



OTHER

INFORMATION ON THE DISPLAY:

If the "Lo" warning is displayed after switching on (Fig. 4), or there is nothing on the display, the battery is dead and has to be replaced.



CHANGE BATTERY:

Please open the battery holder at the backside of the scale and carefully take out the battery by pressing the small silver colour holding plate to get the battery out of the battery holder. Replace the low battery to a new battery. The low battery should be thrown to the appropriate place.

BATTERIES ARE NOT INCLUDED IN WARRANTY!

OVERLOAD WARNING:

When the scale is overloaded, the following is displayed:



CLEANING AND MAINTENANCE:

The surface of the unit may only be cleaned using a wet cloth. Do not use any detergents or hard objects as they may damage the unit's surface!

Never clean the unit under running water, do not rinse it, and do not submerge it in water!

REPAIR AND MAINTENANCE:

Any extensive maintenance or repair that requires access to the inner parts of the product must be performed by an expert service centre.

ENVIRONMENTAL PROTECTION:

- Packaging materials and obsolete appliances should be recycled.
- The transport box may be disposed of as sorted waste.
- Polyethylene bags shall be handed over for recycling.



DISPOSAL

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel travelling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or

with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from.

Products sent to GIMA will be rejected.

Importato da / Imported by/Importé de / Importado de:

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

www.gimaitaly.com

Made in Hungary

M 80210 -80211 - EN - Rev. 2.02.2020