Test Marker Cardiaco Combo a dischetto per sangue Intero, plasma o siero
Step Cardiac Combo Marker test disk for whole blood, plasma or serum
Test Marqueur Cardiaque Combo sur cassette sur sang total, plasma ou sérum
Marker Test Cardiaco Combo de disco para sangre entera, plasma o suero

MANUALE D'USO
OPERATOR'S MANUAL
MANUEL D'UTILIZATION
MANUAL DE USO

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- 1, Andersen PT, Miller-Petersen J, Henneberg EW, et al, "IpermioglibiHypermyoglobinemia After Successful Arterial Embolectomy," *Surgery*, 1987, 102:25-31.
- 2, Cummins B, Russell GJ, Chandler ST, et al, "Uptake of Radioiodinated Cardiac Specific Troponin I Antibodies in Myocardial Infarction," *Cardiovasc Res*, 1990, 24(4):317-27.
- 3, Delanghe J, Chapelle JP, el Allaf M, et al, "Quantitative Turbidimetric Assay for Determining Myoglobin Evaluated," *Ann Clin Biochem*, 1991, 28(Pt 5):474-9.
- 4, Donnelly R and Hillis WS, "Cardiac Troponin T," *Lancet*, 1993, 341(8842): 410-1.
- 5, Gibler WB, Gibler CD, Weinshenker E, et al, "Myoglobin as an Early Indicator of Acute Myocardial Infarction," *Ann Emerg Med*, 1987, 16:851-6.
- 6, Mair J, Smidt J, Artner-Dworzak E, et al, "Rapid Diagnosis of Myocardial Infarction by Immunoturbidimetric Myoglobin Measurement," *Lancet*, 1991, 337(8753):1343-4.
- 7, Parmacek MS and Leiden JM, "Structure, Function, and Regulation of Troponin C," *Circulation*, 1991, 84(3):991-1003.
- 8, "Troponin T and Myocardial Damage," *Lancet*, 1991, 338(8758):23-4, (editorial).
- 9, Serrano S, Chueca P, Carrasco E, et al, "Predictive Value of Myoglobin in Early Diagnosis of Acute Mtocardial Infarction," *Ann Emerg Med*, 1990, 19(8):953.
- 10, Silva DP Jr, Landt Y, Porter SE, et al, "Development and Application of Monoclonal Antibodies to Human Cardiac Myoglobin in A Rapid Fluorescence Immunoassay," *Clin Chem*, 1991, 37(8):1356-64.
- 11, Vrenna L, Castaldo AM, Castaldo P, et al, "Comparison Between Nephelometric and RIA Methods for Serum Myoglobin, and Efficiency of Myoglobin Assay for Early Diagnosis of Myocardial Infarction," *Clin Chem*, 1992, 38(5):789-90. Marker Test Cardíaco Combo de disco

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GIMA S.p.A.