ENGLISH Instructions for Use

2800 series Filters/2800 series Filters with attachments

Intended use: Spirometer pre filters are used in pulmonary function tests for trapping bacteria and viruses and other particulates, reducing the risk of cross-contamination between the patient and the machine.

Filters can be supplied with accessories, such as nose clip, mouthpiece and bite grip.

Bidirectional Filter.

Intended use of nose clip (A508): A nose clip is a device designed to hold the nostrils closed to prevent the air from entering, or escaping during the spirometry testing.

Intended use of mouthpiece (A571): A device designed for partial insertion into a patient's mouth (in between the lips) to facilitate access to the respiratory system. It is typically attached, either directly or through other tubes/adaptors, to diagnostic and/or therapeutic respiratory devices (e.g., ventilator, incentive spirometer, breath analyser, breath collector).

Intended use of bite grip (A539): A mouthpiece to be used in conjunction with the spiroguard filter used in pulmonary lung function testing.

Cautions

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in harm to the patient.

Any use of the medical device requires full understanding and strict observation of all portions of these Instructions for Use. The medical device may only be used for the purpose specified under "Intended Use". Observe all CAUTION statements throughout this manual and all statements on medical device labels. The manufacturer disclaims any liability for patient injury caused by noncompliance with these statements.

The product may only be used on the patient side.

The medical device IFU is not available for each individual filter. Only one copy of the Instructions for Use is included in the distribution package and should therefore be kept in an accessible location for users.


Medical device is only to be used by healthcare professionals.

User has the responsibility to report any serious incident that has occurred in relation to the medical device, to the manufacturer, EC REP and the competent authority of the Member State in which the user and/or the patient is established.

To avoid contamination and soiling, the product should remain packaged until ready to be used.

Do not use the product if the packaging and/or product is damaged.

If product was kept outside the stated storage temperature, product functionality cannot be guaranteed. Storage Temperature: +5°C to 40°C.

Store the product at room temperature, in dry place, free of dust and not exposed to direct sunlight.

Handle filter with care. Abusive handling can damage the filter.

For Your Safety and that of Your Patients, follow the Instruction for Use of the machine!

Before installation, check that all system components are free of obstructions and of foreign bodies. Otherwise, use is limited, or faulty operation is possible.

Inappropriate connection to other devices can result in a patient hazard.

Medical device products have been designed, tested and manufactured exclusively for disposable use and for a period of use not exceeding 24 hours following the first opening of the primary packaging.

Do not use after expiry date stated on primary packaging.

Manufacturing date

Expiry date

Monitor filter performance during the use. If any changes to filter performance or appearance are observed, change the filter to new one.

Always keep spare filter ready to exchange with the one in use, if any of the below are observed:

- Increase in resistance;
- Filter is visibly soiled or damaged or if it rattles;
- Excessive condensation;
- Associated alarm.

This product is for single patient use only.

The products must not be reused, cleaned or sterilised / not to be sterilised after use in order to re-use. Exceptions from this rule will be stated within the Instruction For Use for the particular products code if applicable. Reuse, cleaning or sterilization may lead to failure of the products, risk infection and injury to the patient.

This product is latex and PHT free.

Batch code symbol.

Prescription only.

Following use, the products must be disposed off in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.

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REF	2800/01XXXX Spiroguard	2800/02XXXX Spiroguard	2800/03XXXX Spiroguard	2800/10XXXX Spiroguard	2800/11XXXX Spiroguard Integral Mouthpiece	2800/15XXXX Spiroguard	2800/17XXXX Spiroguard	2800/21XXXX Spiroguard Integral Mouthpiece	2800/22XXXX Spiroguard	2800/23XXXX Spiroguard with Integral Mouthpiece	2800/24XXXX Spiroguard	2800/25XXXX Spiroguard with Integral Mouthpiece	2800/26XXXX Spiroguard	2800/27XXXX Spiroguard with Integral Mouthpiece	2800/30XXXX Spiroguard with Integral Mouthpiece	2800/31XXXX Spiroguard Integral Mouthpiece	2800/32XXXX Spiroguard	2800/72XXXX Compact Electrostatic Spirometry Filter	2800/729XXXX Compact Electrostatic Spirometry Filter
Product Code; Code du produit; Artikelnummer; Codice prodotto; Código de producto; Productcode; Código de produto; Produktkod; Produktkode; Produktkode; Tuotekoodi; Kod produktu; Код на продукта; Sifra proizvoda; Kód výrobku; Τοοτε κωδ; Κυβωδός προϊόντος; Termékkód; Cód Tárge; Izstrādājuma kods; Gaminio kodas; Kодиcт tal-prodot; Cod produs; Kód výrobku; Koda izdelka.																			
Product description; Description du produit; Produktbeschreibung; Descripción producto; Descripción de producto; Productomschrijving; Descrição do produto; Produktbeskrivning; Produktbeskrivelse; Produktbeskrivelse; Tuotteen kuvaus; Opis produktu; Описание на продукта; Opis proizvoda; Opis výrobku; Τοοτε κειδελος; Περιγραφή προϊόντος; Termékleírás; Tuarisc ar Thárige; Izstrādājuma apraksts; Gaminio aprásmas; Deskripcija tal-prodot; Descrierea produsului; Popis výrobku; Opis izdelka.	Inlet/Outlet Connectors: - OD 34mm / ID 30.1mm - Machine Side; - OD 29.22mm / ID 26.14mm - Patient Side. Approx. Dimensions: 96.7mm diameter x 79mm height. Weight: 37g (approx.).	Inlet/Outlet Connectors: - OD 34mm / ID 28.2mm - Machine Side; - OD 29.22mm / ID 26.7mm - Patient Side. Approx. Dimensions: 97mm diameter x 79mm height. Weight: 37g (approx.).	Inlet/Outlet Connectors: - OD 34mm / ID 30.1mm - Machine Side; - OD 29.22mm / ID 26.14mm - Patient Side. Approx. Dimensions: 96.7mm diameter x 79.2mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD34mm / ID30.1mm - Machine Side; - OD 30.1mm / ID 26.2mm - Patient Side. Approx. Dimensions: 97mm diameter x 79mm height. Weight: 37g (approx.).	Inlet/Outlet Connectors: - OD30.65mm / ID25.9mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 86mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD30.65mm / ID25.9mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 86mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD29.2mm / ID26.7mm - Patient Side; - OD 29.4mm / ID 25.9mm - Machine Side. Approx. Dimensions: 96.8mm diameter x 78.8mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD34mm / ID29.25mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 97mm diameter x 93mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD31.2mm / ID26.09mm - Patient Side. Approx. Dimensions: 97mm diameter x 85mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD 48.36mm / ID 45.35mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 79.2mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD 48.36mm / ID 45.35mm - Machine Side; - OD 31.2mm / ID 29.09mm - Patient Side. Approx. Dimensions: 96.8mm diameter x 71.4mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD34.75mm / ID27.8mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 86mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD34.75mm / ID27.8mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 92.9mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD34mm / ID29.25mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.8mm diameter x 92.9mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD29.3mm / ID27.2mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 97mm diameter x 92mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD 29.4mm ID 25.9mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 96.7mm diameter x 86mm height. Weight: 35gm (approx.).	Inlet/Outlet Connectors: - OD 29.2mm ID 26.7mm - Patient Side; - OD 34.0mm ID 29.25mm - Machine Side. Approx. Dimensions: 96.8mm diameter x 84.9mm height. Weight: 35gm (approx.).	Inlet/Outlet Connectors: - OD 34mm / ID 29.25mm - Machine Side; - Integral Mouthpiece - Patient Side. Approx. Dimensions: 68.5mm diameter x 92.6mm height. Weight: 37gm (approx.).	Inlet/Outlet Connectors: - OD 25.9mm / ID 24.9mm approx. - Machine Side; Integral Bite Style Mouthpiece OD 25.2mm approx. / ID 23.8mm approx. - Patient Side. Approx. Dimensions: 68.6mm diameter x 94mm height. Weight: 21gm (approx.). Bidirectional Filter.
NON STERILE	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.	Yes; Oui; Ja; Si; Si; Ja; Sim; Ja; Ja; Joo; Tak; DA; Da; Ano; Jah; Nai; Igen; Tá; Já; Taip; Iva; Da; Ano; Da.

Ambient Conditions: Conditions ambiantes : Umgebungsbedingungen : Condizioni ambientali : Condiciones ambientales : Omgevingscondities : Condições ambientais : Omgivningsförhållanden : Omgivelsesbetingelser : Krav til omgivelsene : Ympäristöolosuhteet : Warunki otoczenia : Условия на околната среда : Увjeti okoline : Okolni podminky : Keskkonnatingimused : Συνθήκες περιβάλλοντος : Környezeti feltételek : Coínollacha Comhthreallacha : Apkárthige apstáikí : Aplinkos sąlygos : Kundizjonijiet Ambientali : Condiții ambientale : Okolité podmienky : Okolje.

Temperature during storage; Temperatura de stocaj; Temperatur bei Lagerung; Temperatura durante l'immagazzinamento; Temperatura durante el almacenamiento; Temperatur tijdens opslag; Temperatura em armazenamento; Temperatur under förvaring; Temperatur under opbevaring; Temperatur under oppbevaring; Lämpötila varastoinnin aikana; Temperatura podczas składowania; Температура по време на съхранение; Temperatura tijekom skladištenja; Teplota během skladování; Temperatur ladustamise ajal; Οερυχοκρασία κατά την αποθήκευση; Hőmérséklet tárolás közben; Teocht le linn na stórála; Temperatura uzglabšanas laikā; Temperatura sandēlavimo metu; Temperatura waqt il-hazna; Temperatura în timpul depozitării; Teplota počas skladovania; Temperatura med skladištenjem.	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)	5 to 40 °C (41-104 °F)

Performance Data: Caractéristiques ventilatoires : Leistungsdaten : Dati di rendimento : Características de rendimiento : Karakteristieke waarden : Dados de desempenho : Prestandauppgifter : Ydelse : Kapasitetsdata : Suorituskykytiedot : Parametry sprężu; Дани за ефективността; Podaci o izvedbi; Údaje o výkonu; Toimivusandmed; Δεδομένα απόδοσης; Teļestitmnē adati; Sonrai Feidhmíochta; Veiktspējas dati; Pasirodymo data; Dejta tal-Prestazzjoni; Date de performanță; Údaje o výkonnosti; Podatki o uspešnosti.

Internal Volume; Volume interne; Innenvolumen; Volumen interno; Volume interno; Intern volume; Volume interno; Innenvolum; Indre volum; Sisältävolum; Objętość wewnętrzna; Вътрешен обем; Unutarnji volumen; Vnitřní svazek; Sisemine maht; Εσωτερικός όγκος; Belső kötet; Imleabhar Innheánach; Iekšējais tilpums; Vidinis tūris; Volum Intern; Volum intern; Vnútrotný objem; Notranja glasnost.	72ml	76ml	73ml	74ml	74ml	65ml	67ml	77ml	74ml	68ml	64ml	75ml	72ml	77ml	73ml	70ml	49ml	42ml

Resistance; Résistance; Widerstand; Resistencia; Resistencia; Weerstand; Resistência; Resistans; Modstand; Motstand; Vastus; Опор; Спротивля; Опротив; Odpor; Vastupanu; Αντίσταση; Ellenállás; Frittaloch; Pretestība; Paspiresināmas; Resistēnca; Resistēnā; Odpor; Odpornost.	Max.36.3Pa @ 30L/min; Max.70.4Pa @ 60L/min; Max.118.8Pa @ 90L/min.	Max.33Pa @ 30L/min; Max.60.5Pa @ 60L/min; Max.94.6Pa @ 90L/min.	Max.36.3Pa @ 30L/min; Max.70.4Pa @ 60L/min; Max.118.8Pa @ 90L/min.	Max.37.4Pa @ 30L/min; Max.71.5Pa @ 60L/min; Max.119.9Pa @ 90L/min.	Max.29.7Pa @ 30L/min; Max.63.8Pa @ 60L/min; Max.104.5Pa @ 90L/min.	Max.38.5Pa @ 30L/min; Max.72.6Pa @ 60L/min; Max.113.3Pa @ 90L/min.	Max.33Pa @ 30L/min; Max.57.2Pa @ 60L/min; Max.89.1Pa @ 90L/min.	Max.35.2Pa @ 30L/min; Max.75.9Pa @ 60L/min; Max.127.6Pa @ 90L/min.	Max.35.2Pa @ 30L/min; Max.75.9Pa @ 60L/min; Max.132Pa @ 90L/min.	Max.33Pa @ 30L/min; Max.68.2Pa @ 60L/min; Max.112.2Pa @ 90L/min.	Max.35.2Pa @ 30L/min; Max.63.8Pa @ 60L/min; Max.108.9Pa @ 90L/min.	Max.33Pa @ 30L/min; Max.59.4Pa @ 60L/min; Max.91.3Pa @ 90L/min.	Max.52.8Pa @ 30L/min; Max.96.8Pa @ 60L/min; Max.156.2Pa @ 90L/min.	Max.29.7Pa @ 30L/min; Max.63.8Pa @ 60L/min; Max.104.5Pa @ 90L/min.	Max.37.4Pa @ 30L/min; Max.73.7Pa @ 60L/min; Max.121Pa @ 90L/min.	Max.37.4Pa @ 30L/min; Max.70.4Pa @ 60L/min; Max.113.3Pa @ 90L/min.	Max.35.2Pa @ 30L/min; Max.78.1Pa @ 60L/min; Max.93.5Pa @ 90L/min.	Max.79Pa @ 30L/min; Max.160Pa @ 60L/min; Max.259Pa @ 90L/min.	Max.72.6Pa @ 30L/min; Max.155.1Pa @ 60L/min; Max.244.2Pa @ 90L/min.

Filter Efficiency @ 30L/min; Efficacité du filtre à 30L / min; Filterleistung bei 30 l / min; Efficacia del filtro @ 30L / min; Eficacia del filtro a 30 l / min; Filtrerefficiēnā @ 30L / min; Eficácia do filtro @ 30L / min; Filtreringseffektivitet @ 30L / min; Filtreringseffektivitet @ 30L / min; Filter Effektivitet @ 30L / min; Suodattimen tehokkuus @ 30L / min; Wydajność filtra @ 30 l / min; Ефективност на филтъра @ 30L / min; Учинковитост филтра @ 30L / min; Учинност филтру pH 30 l / min; Filtri efektiivsus @ 30L / min; Αποδοτικότητα φίλτρου @ 30L / min; Szűrőhatékonyság @ 30L / perc; Είεαhtacht Scagaire @ 30L / nóm; Filtra efektivitate @ 30L / min; Filro efektyvumas @ 30L / min; Efficjenza tal-filtra @ 30L / min; Eficencia filtrului @ 30L / min; Učinost filtra @ 30L / min; Učinkovitost filtra @ 30L / min.	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 99%	Min. 94%

Bacterial Efficiency; Rétenion bactérienne; Bakterienretention; Ritenzione di batteri; Retención de bacterias; Bacterieretentie; Retenção bacteriana; Bakterieretention; Bakterieretention; Bakterieretensjon; Bakterien suodatus; Zatrzymywanie bakterii; Ефективност на бактериите; Učinkovitost bakterija; Bakteriální účinnost; Bakterite efektiivsus; Βακτηριακή αποτοκρομακότητα; Bakteriális hatékonyság; Eifeachtúlacht Baictéarach; Bakteriū efektiivitate; Bakterijū efektyvumas; Efficjenza Batterjalj; Eficiență bacteriană; Bakteriálna účinnost; Učinkovitost bakterij.	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%

Viral Efficiency; Rétenion virale; Virenretention; Ritenzione di virus; Retención de virus; Virusretentie; Retenção de vírus; Virusretention; Virusretention; Virusreten; Virussten suodatus; Zatrzymywanie wirusów; Βιρυχνα ефективност; Virusna učinkovitost; Virová účinnost; Virussik efektiivsus; Ixrti αποτοκρομακότητα; Virushatékonyaság; Eifeachtúlacht Vireasach; Virusū efektiivitate; Virusinis efektyvumas; Efficjenza Viralj; Eficiență virală; Virusová účinnost; Virusna učinkovitost.	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%	Min.99.999%

Material Data; Matériau; Material; Dati materiale; Material; Material; Materiale; Materiale; Materiali; Material; Дани за материала; Podaci o materijalu; Materijalové údaje; Materiaised andmed; Δεδομένα υλικού; Anyagi adatok; Sonrai Ábhar; Materijalni dati; Esminiai duomenys; Dejta Materjalj; Date materiale; Materijalové údaje; Materialni podatki.

Housing; Boliter; Gehäuse; Alloggiamento; Carcasa; Behuizing; Estrutura; Kabinet; Hus; Hus; Pidkie; Obudowa; Корпус; Kucište; Kryt; Ümbris; Пери́фери́я; Burkloch; Títhloch; Májokj; Bústas; Djar; Carcasa; Kryt; Ohšje.	HIPS	HIPS	HIPS	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP	PP

Filter Media; Matière du filtre : Filter-Medium : Materiale del filtro : Filtro : Filter medium : Material do filtro : Filtermaterial : Filtermateriale : Filtermedium : Suodattimen väliaine : Medium filtrujące; Филтромаща среда; Filtrarsi medij; Filtrační prosředí; Filter; Матра дйрђишњ; Szűrőközeg; Meán Scagaire; Filtrēt multivdi; Filtruoti laikmena; Filtra I-Midja; Medi de filtrare; Filtračné médium; Filtrini medij.	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber	Electrostatic Blended Synthetic Fiber