ITALIANO

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ENGLISH

Uric Acid Test Strips

Use only with GIMACARE Multi-Functional Monitoring System.

Warnings

- For *in vitro* diagnostic use (outside of the body) only.
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this manual and your GIMACARE Multi-Functional Monitoring System Owner's Manual before you begin. Only use GIMACARE Uric Acid Test Strips with GIMACARE Multi-Functional Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure, or those who are in shock.
- For patients with impaired peripheral circulation, collection of capillary blood from the approved sample sites is not advised as the results may not be a true reflection of the physiological uric acid level.
- This product is not intended for the diagnosis of diseases. Never change your medication without consulting your physician first.
- Keep test strips and lancets away from young children. If swallowed, consult a doctor immediately.

Intended Use

GIMACARE Uric Acid Test Strips, when used together with GIMACARE Multi-Functional Monitoring System, allow you to measure uric acid levels by yourself at home or by healthcare professionals. Users should only obtain fresh capillary whole blood samples from the fingertips. This system should not be used for diagnosis or screening of diseases.

Limitations

- Hematocrit: The hematocrit level is limited to between 20% and 60%. Please ask your healthcare professional if you do not know your hematocrit level.
- In vitro* acetaminophen up to 2 mg/dL, ascorbic acid up to 4 mg/dL, levo-dopa up to 1 mg/dL, bilirubin (unconjugated) up to 2 mg/dL and methyl-dopa up to 0.3 mg/dL showed no interference.
- Metabolites: high concentration of acetaminophen, bilirubin, hemoglobin, methyl-dopa and ascorbic acid may affect the test results. The user is suggested to perform a uric acid test in the morning after fasting for 12 hours.
- Altitude Effects: Altitudes up to 11,500ft (3,500m) do not affect test results

Storage and Handling

- Δ Do not use the test strips if they have expired.
- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it. (For strip vial only)
- Store the test strips in a cool, dry place between 2°C and 30°C (35,6°F and 86°F) and between 10% and 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers. (For strip vial only)
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial or individual foil packet. Close the vial immediately after taking out a strip. (For strip vial only)
- Keep the vial closed at all times. (For strip vial only)
- Do not bend, cut, or alter the test strip.

Strip Appearance

- Absorbent Hole**
Apply blood sample here. The blood will be automatically absorbed.
- Confirmation Window**
This is where you confirm if enough blood has been applied to the absorbent hole in the strip.
- Test Strip Handle**
Hold this part to insert the test strip into the slot.
- Contact Bars**
Insert this end of the test strip into your meter. Push it in firmly until it will not go any further.

Calibration

Calibrate the meter every time you begin to use a new package of test strips by setting the meter with the correct code. Test results may be inaccurate if the code number displayed on the meter does not match the number printed on your test strip vial label/packet.

Code Chip

- Insert the code chip with the meter switched off. Wait until a number and "UA" appears on display.
 - Remove the code chip. The display will show "OFF", and then the meter will automatically switch off.
- Checking the Code Number**
- Make sure that the number and "UA" displayed on the meter matches the number on your test strip vial label/packet before you proceed. If the numbers match, you can proceed with the test. If they do not match, please stop testing and insert the correct code chip. If the problem persists, contact customer service for help.

Testing Your Uric Acid

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTING. ALWAYS REFER TO THE OWNER'S MANUAL AND LANCET INSERT ON HOW TO COLLECT A BLOOD SAMPLE.

- Insert the test strip fully into the slot of your meter until it will not go any further. When the strip is fully inserted, the meter will perform several self-checks.
- Collect a blood sample with the test strip. Wipe off the first drop of blood with a clean cotton swab. Make sure there is a sufficient quantity of blood in order to provide accurate test results. Apply the blood drop to the absorbent hole of the test strip, and wait until the confirmation window is fully filled. The meter will start counting down. Do NOT apply a smeared blood sample.
- After a few seconds, the meter will display your blood uric acid level. The last reading will be automatically saved in the meter. The meter will turn off automatically after the test strip is removed. Make sure to dispose the used lancet and test strip carefully.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations. Please refer to your Owner's Manual for further information.

Reading Your Result

Your uric acid readings deliver plasma equivalent results and are displayed in milligrams of uric acid per deciliter of blood (mg/dL) or micromoles of uric acid of blood (μmol/L). The measurement unit is based on the actual meter you used.

Reference values

Male	3.5 to 7.2 mg/dL (208 to 428 μmol/L)
Female	2.6 to 6 mg/dL (155 to 357 μmol/L)

Source: National Kidney Foundation (2014).

Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
- Check the expiration date of the test strips.
- Check the performance of your meter and test strip with the control solutions
- Δ If your test results are significantly different from what you expect, or in unusually high or low levels, please repeat the test with a new test strip or contact your healthcare professional.

Expiry Date Reminder

For your convenience, the expiry date reminder will activate and notify you the

number of days remaining until the strip's expiry date shown on the vial label or on the foil packet. The count down begins from 30 days to 1 day, which will be shown at the center of the display screen. When you see the date reminder, please use the remaining test strips before they expire.

Δ The error message E-2 will appear in the following situations:

- The test strip is expired;
- The code chip is expired; or
- On the initial set-up, the date has been set incorrectly on the meter.

If the error message E-2 appears, please repeat the test with a new lot of test strip to get accurate results.

Chemical Components:

Uric acid catalyst $\geq 0,15$ mg

Non-reactive ingredients $\geq 0,05$ mg

Quality Control Testing

Our control solutions contain a known amount of uric acid that will react with the test strips. If you are concerned about the meter or test strips are not working properly, you can check the performance of meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial or on the test strip package. Please refer to the Owner's Manual for the step-by-step quality control test instructions.

Δ LThe reference range of the control solutions may vary with each new vial or package of test strips.

Make sure you check the range on the label of your current vial or on the current package

Additional Information

Always wear gloves and follow your local biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh capillary whole blood samples only.

Sample Size: 0.5 μL

Reaction Time: 15 seconds

System Measurement Range: 3 to 20 mg/dL (178 to 1189 μmol/L).

Hematocrit Range: 20% to 60%

Accuracy

Accuracy describes how well the readings from a testing system (meter and test strips) agree with the readings from an internationally accepted reference system (laboratory analyzer) and is performed according to an internationally recognized standard. The prepared samples were then tested on this meter and compared to the readings from the Cobas C311 analyzer.

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Symbol Information

Symbol	Referent	Symbol	Referent
	<i>In vitro</i> diagnostic medical device		Humidity limitation
	Consult instructions for use or consult electronic instructions for use		Manufacturer
	Temperature limit		Model number
	Use-by date		Do not re-use
	Batch code		Unique device identifier
	Caution		CE mark
	Authorized representative in the European Community		

	MedNet EC-REP GmbH Borkstraße 10 48163 Münster, Germany		TaiDoc Technology Corporation B1-7F, No.127, Wugong 2nd Rd., Wugu Dist., 24888 New Taipei City, Taiwan www.taidoc.com
--	--	--	---

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