***Test Marker Cardiaco Combo a dischetto per sangue Intero, plasma o siero***  
***Step Cardiac Combo Marker test disk for whole blood, plasma or serum***  
***Test Marqueur Cardiaque Combo sur cassette sur sang total, plasma ou sérum***  
***Marker Test Cardiaco Combo de disco para sangre entera, plasma o suero***

**MANUALE D'USO**  
**OPERATOR'S MANUAL**  
**MANUEL D'UTILIZATION**  
**MANUAL DE USO**

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### BIBLIOGRAPHIE

- 1, Andersen PT, Miller-Petersen J, Henneberg EW, et al, "IpermiogliHypermyoglobinemia After Successful Arterial Embolectomy," *Surgery*, 1987, 102:25-31.
- 2, Cummins B, Russell GJ, Chandler ST, et al, "Uptake of Radioiodinated Cardiac Specific Troponin I Antibodies in Myocardial Infarction," *Cardiovasc Res*, 1990, 24(4):317-27.
- 3, Delanghe J, Chapelle JP, el Allaf M, et al, "Quantitative Turbidimetric Assay for Determining Myoglobin Evaluated," *Ann Clin Biochem*, 1991, 28(Pt 5):474-9.
- 4, Donnelly R and Hillis WS, "Cardiac Troponin T," *Lancet*, 1993, 341(8842): 410-1.
- 5, Gibler WB, Gibler CD, Weinschenker E, et al, "Myoglobin as an Early Indicator of Acute Myocardial Infarction," *Ann Emerg Med*, 1987, 16:851-6.
- 6, Mair J, Smidt J, Artner-Dworzak E, et al, "Rapid Diagnosis of Myocardial Infarction by Immunoturbidimetric Myoglobin Measurement," *Lancet*, 1991, 337(8753):1343-4.
- 7, Parmacek MS and Leiden JM, "Structure, Function, and Regulation of Troponin C," *Circulation*, 1991, 84(3):991-1003.
- 8, "Troponin T and Myocardial Damage," *Lancet*, 1991, 338(8758):23-4, (editorial).
- 9, Serrano S, Chueca P, Carrasco E, et al, "Predictive Value of Myoglobin in Early Diagnosis of Acute Myocardial Infarction," *Ann Emerg Med*, 1990, 19(8):953.
- 10, Silva DP Jr, Landt Y, Porter SE, et al, "Development and Application of Monoclonal Antibodies to Human Cardiac Myoglobin in A Rapid Fluorescence Immunoassay," *Clin Chem*, 1991, 37(8):1356-64.
- 11, Vrenna L, Castaldo AM, Castaldo P, et al, "Comparison Between Nephelometric and RIA Methods for Serum Myoglobin, and Efficiency of Myoglobin Assay for Early Diagnosis of Myocardial Infarction," *Clin Chem*, 1992, 38(5):789-90.

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GIMA S.p.A.